
USER MANUAL E-connect S+

E-connect S+

E-connect S+

P/N: IFU-6035268

Version: 02

Issued: 2026.03.26

Size: 160mm*92mm

English Instruction	Page 2
Deutsch Anleitung	Page 33
FrançaisInstructions	Page 69
Italiano Istruzioni	Page 104
Español Instrucciones	Page 139
Português instruções	Page 173
Polski Instrukcje	Page 206
Română instruciuni	Page 241
Nederlands instructies	Page 274
İngilizce Talimatlar	Page 308

Content

1. Scope of E-connect S+	3
1.1 Parts Identification	3
1.2 Components	3
1.3 Accessories	3
2. Symbols used	4
3. Before Use	4
3.1 Scope of application	4
3.2 Contraindications	5
4. Installing the E-connect S+	6
4.1 Installation of the contra angle	6
4.2 Install the file	6
4.3 Installation of the Insulating Sleeve	6
4.4 Connecting measuring wire	7
4.5 Connecting charging base	7
5. Use Interface	7
5.1 Panel key	7
5.2 Screen display	8
5.3 Terms and definition	8
6. Setting	9
6.1 Set memory mode	9
6.2 Setting parameters	9
6.3 Preset programs	11
6.4 Advanced setting	13
6.5 Parameter logic	14
7. Operation	16
7.1 Charge	16
7.2 Motor Operation	17
7.3 Apex operation and not suitable condition	18
8. Cleaning, Disinfection and Sterilization	20
8.1 Foreword	20
8.2 General recommendations	20
8.3 Autoclavable Components	21
8.4 Disinfection Components	24
9. Error Indication	24
10. Troubleshooting	24
11. Technical Data	26
12. EMC Tables	27
13. Statement	31

1.Scope of E-connect S+

1.1 Parts Identification

(Fig 1.1) Note: This product does not contain root canal file

- ① Contra Angle ② Motor Handpiece ③ Charging Base ④ Insulating Sleeve ⑤ File Clip (2 pcs) ⑥ Lip Hook (2 pcs) ⑦ Measuring Wire ⑧ Adapter ⑨ Spray Nozzle ⑩ Apex Tester ⑪ Handpiece Base(optional)

1.2 Components

Motor Handpiece (1pc) Fig 1.2a Part No: 6051153	Charging Base (1pc) Fig 1.2b Part No: 6051075	Contra Angle (1pc) Part No: 6036010 Fig 1.2c
Measuring Wire (1pc) Fig 1.2d Part No: 6015015	Lip Hook (2pcs) Fig 1.2e Part No: 6072002	File Clip (2pcs) Fig 1.2f Part No: 6151036


For different regions, there are several different adapter options to be selected as follows.

Standard	Adapter	Power plug
European standard	Adapter (1pc) Fig 1.2g Part No: 6016021	/
American standard	Adapter (1pc) Fig 1.2h Part No: 6516003	American standard power plug (1pc) Fig 1.2i Part No: 6016011
Multi-standard	Adapter (1pc) Fig 1.2j Part No: 6516003	British standard power plug (1pc) Fig 1.2k Part No: 6016009
		Australian standard power plug (1pc) Fig 1.2l Part No: 6016010
		Argentina standard power plug (1pc) Fig 1.2m Part No: 6016014

1.3 Accessories

Apex Tester (1pc) Fig 1.3a Part No: 6151005	Spray Nozzle (1pc) Fig 1.3b Part No: 6051108	Insulating Sleeve (1pc) Fig 1.3c Part No: 6004027	Handpiece Base(optional) Fig 1.3d Part No: 6005002
--	---	---	--

2. Symbols used

	General warning sign		Keep dry
	Caution		CE marking
	Serial number		Direct current
	Catalogue number		Consult instructions for use
	Batch code		Type B applied part
	Medical device		Humidity limitation
	Authorized representative in the European Community		Sterilizable in a steam sterilizer (autoclave) at the temperature specified
	Manufacturer		Temperature limitation
	Country of manufacture + Date of manufacture		Atmospheric pressure limitation
	Class II equipment		Manufacturer's Logo
	Washer-disinfector for thermal disinfection		Dispose of in accordance with the WEEE directive

3. Before Use

3.1 Scope of application

E-connect S+ is a cordless endodontic treatment motorized handpiece with root canal measuring capability. It can be used to enlarge canals while monitoring the position of the file tip inside the canal. It can be used as a low speed motorized handpiece and device for measuring canal length.

This device must only be used in hospital environments, clinics or dental offices by qualified dental personnel and not used in the oxygen-rich environment.

3.2 Contraindications

The integrated apex locator of the E-connect S+ is contraindicated in cases where patient/user carry medical implants such as pace makers or cochlear implants etc.

Do not use the device for implants or other non-endodontic dental procedures.

Safety and effectiveness have not been established in pregnant women and children.

Read the following warnings before use:

- The device must not be placed in humid surroundings or anywhere where it can come into contact with any type of liquids.
- Do not expose the device to direct or indirect heat sources. Do not use the device in the presence of free oxygen, anesthetic gas or combustible materials. The device must be operated, used and stored in a safe environment.
- The device requires special precautions with regard to electromagnetic compatibility (EMC) and must be installed and operated in strict compliance with the EMC information. In particular, do not use the device in the vicinity of fluorescent lamps, radio transmitters, remote controls and do not use this system near the active Surgical Equipment and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high. Portable RF communications device (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the E-connect S+, including cables specified by the manufacturer. Otherwise, degradation of the performance of this device could result.
- Please do not charge, use or store this device at high temperature. Please pay attention to the use and storage conditions.
- Gloves and a rubber dam are compulsory during treatment.
- No modification of this equipment is allowed. Never open or repair the device yourself, otherwise, void the warranty.
- If irregularities occur in the device during treatment, switch it off. Contact the local dealer for treatment.
- Please use the original power adapter when charging.
- If liquid flows out of the handpiece, it can be judged as battery leakage. Please stop using immediately and contact the local dealer for treatment.
- Do not dismount the contra angle during the operation of the main engine, otherwise the contra angle and motor gear will be damaged.
- Please use the original contra angle, which gear ratio is 1:1. The contra angle is not field repairable.
- Use continuous file in continuous mode; use reciprocating file in reciprocating mode, and use according to rotation speed, torque and contra angle recommended by the root canal file manufacturer.
- The user or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the

competent authority of the Member State in which the user or patient is established.

- It is forbidden to use non-original parts for the equipment.
- Do not load and use the device for a long time. Otherwise, the temperature of the device will rise, which may cause minor burn to the operator or the patient. (The surface of some applied parts, such as contra angle, will reach maximum to 48°C if loading the device continuously for more than 1 minute. The surface temperature of the motor handpiece will reach maximum to 52°C if loading the device continuously for more than 10 minutes.)
- Do not use the device on patients who are allergic to Nickel.

4. Installing the E-connect S+

4.1 Installation of the contra angle

Make sure 4 pins on contra angle alignment the slots of handpiece, plug them together until a sound of "click" to make sure to install into place. (Fig 4.1a)

The contra angle can be 360 degrees rotated without take off, this make it easy to see the screen during treatment by rotating the contra angle. (Fig 4.1b)

- After connecting the contra angle and handpiece, pull it gently to make sure the connection is good
- Make sure the assembly is connected properly, otherwise might cause unexpected motor reverse, even hurt the patients.
- It is strictly forbidden to start the equipment before the contra angle is plugged in, otherwise it will hurt the operator.

4.2 Install the file

Install: insert the root canal file and rotate it slightly from left to right to ensure that the file needle is aligned with the internal bayonet, and then push it in slightly to complete the installation.

Remove: press and hold the back cover button on the contra angle to release the internal bayonet and gently pull out the root canal file. (Fig 4.2a)

- Check the file before inserting the file. Do not use the damaged file.
- Pull the file gently to make sure that the file is secure in handpiece properly, otherwise it may pop out and hurt the patient.
- Be careful when inserting and removing files to avoid injury to fingers.
- When removing the file, press the button tightly to release the internal bayonet. If the bayonet is not fully released, the bearing will be damaged.
- Make sure the motor is stopped when inserting and removing files.

4.3 Installation of the Insulating Sleeve

Install: assemble according to the figure. (Fig 4.3a)

Remove: pull out the installation of the Insulating Sleeve in the opposite direction of installing. (Fig 4.3b)

- You can also use disposable sleeve (sold separately) instead of insulating sleeve
- Without the insulating sleeve, when performing the apex measurement with handpiece, wear appropriate insulated gloves, and make sure the contra angle does not touch the lips. It is advisable to use a rubber dam when performing such treatments. It's not necessary to connect file clip during motor combine apex function, only during single apex function.

4.4 Connecting measuring wire

If want activity apex measurement function, uncap the USB cover on handpiece, insert measuring wire. (Fig 4.4a)

- It's not necessary to connect file clip during motor combine apex function, only during single apex function. (Fig 4.4b)
- Match colors to connect the lip hook and file clip, if connect lip hook with black slot, apex auto start will have no function. (Fig 4.4c)

4.5 Connecting charging base

Plug the USB of adapter into the charging base, and plug the other end into a power outlet, the Power LED on charging base will light up (green). (Fig 4.5a)

Put the handpiece all the way into the charging base, the charge state will show on the screen. (Fig 4.5b)

- Put the handpiece into the charging base in the right direction, otherwise the handpiece will not be recharged. (Fig 4.5c)
- Only the original adapter could be used.
- Don't position the device where it is difficult to operate the disconnection device

5. Use Interface

5.1 Panel key

(Fig 5.1) ① ● Main switch ② Display screen ③ S Setting key
④ < Decrease key ⑤ > Increase key

Power on

Press ● more than 0.5 seconds to turn on the instrument.

Memory Change

Press < or > to change at standby mode.

Operation Mode Change

Press **S** once at standby mode, press **<** or **>** to switch, then press **S** or wait 5 seconds to confirm.

Parameter Adjustment

Press **S** till target parameters are displayed at standby mode, press **<** or **>** to adjust, then press **•** or wait 5 seconds to confirm.

Preset Program Selection

Long press **S** at standby mode to enter preset program, press **<** or **>** to select deserve program and press **•** to confirm.

Power off

Holding **S** and press **•**.

Advanced Setting

Holding **S** then press **•** to enter advanced setting at power off mode, press **S** till target parameters are displayed, press **<** or **>** to adjust the parameters, then press **•** to confirm.

5.2 Screen display

Standby interface (Fig 5.2a)

①Memory name ②Speed ③Voice ④Residual battery power ⑤Unit of speed (Revolutions Per Minute) ⑥Unit of torque (Newton Centimeter) ⑦Torque ⑧Operation mode ⑨Memory mode number

Working mode interface (Fig 5.2b)

①Rotate direction: Forward ②Set speed ③Set torque limit ④Scale mark of set torque limit ⑤Torque display scale ⑥Realtime torque

Reference point interface (Fig 5.2c)

①Flash bar of apical reverse position ②2mm meter reading (Indicative dimension) ③Apex (Major/Anatomic apical foramen) ④1mm-3mm (Estimate dimension) distance scale from apex (Estimate dimension)

Canal measurement standby interface (Fig 5.2d)

①Memory mode number M0 is apex stand-alone memory

Canal measurement start interface (Fig 5.2e)

①Indication number The number has no represent of actual length, only for indication ②Canal length indicator bar

5.3 Terms and definition

Fwd	Forward (Clockwise rotation)
Rev	Reverse (Counter clockwise rotation)
REC	Reciprocation Be applied to reciprocating file, path file and rotary file protection by setting some special angle
ATC	Adaptive torque control Up to setting torque, the motor will move with reciprocating mode; when torque reduce to normal value, the motor will clockwise rotate
EAL	Electronic apex locator In the mode, the device will work like a stand-alone apex locator
AP	Apex

	Major apical foramen or Anatomic apical foramen
Reference point	During combined length determination, normally apical reverse must active before reaching major apical foramen, setting apical reverse position by change the flash bar
FWD angle	Angle of rotate forward (can be set in REC & ATC mode.)
REV angle	Angle of rotate reverse (can be set in REC & ATC mode.)
Memory mode	Such as M0-M10
Operation mode	Such as Fwd, Rev, REC and ATC.
Separation of instruments	The file used in root canal therapy is broken accidentally.

6.Setting

6.1 Set memory mode

Fig 6.1a	The device has 10 memory modes (M1-M10), press < or > at standby mode to switch, the memory number (①) will change together. Each memory mode includes its own speed (②), operate mode(④) and torque(③). These parameters can be set separately. (Adjust parameters according to chapter 7.2).
Fig 6.1b	M0 is special memory for stand-alone apex locator function (See chapter 7.3 Apex operation and not suitable condition).

6.2 Setting parameters

- All the parameters must be set according to files, make sure all the parameters are expected before starting the motor, otherwise has risk of file broken.

Fig 6.2a	Before starting the motor, check whether the operation mode (①) is correct. If it is not the expected operation mode, press S once on to enter the operation mode selection, and press < or > to switch, then press S once or wait for 5 seconds to confirm.
Fig 6.2b	This device has four built-in operation modes: Fwd, Rev, REC and ATC (Refer to chapter 5.3 Terms and definition for more information).  ● When using Rev mode, continuous slow alarm sound will appear after starting the motor, which is used to

	remind the operator that the motor is rotating counter clockwise.
	<p>Press S several times, the set of the other parameters of current mode will be shown. Including speed, torque and angle. Confirm all the parameters are correct, if any of them is not desired, press < or > to switch</p> <p></p> <ul style="list-style-type: none"> ● The logic of the parameters in different operation modes are not exactly the same. (Refer to chapter 6.5 Parameter logic).
Fig 6.2c	<p>Speed can be set from 50 rpm to 1500 rpm. Press S several times at standby mode till speed is displayed. Press < or > to switch and press S or wait for 5 seconds to confirm.</p> <p></p> <ul style="list-style-type: none"> ● When using REC mode, there are five speed can be set (100rpm ~ 500rpm). (Refer to chapter 6.5 Parameter logic). <p></p> <ul style="list-style-type: none"> ● When using ATC mode, the speed setting will be different. (Refer to chapter 6.5 Parameter logic).
Fig 6.2d	<p>Torque can be set from 0.5N·cm to 5.0N·cm. Press S several times at standby mode till torque is displayed. Press < or > to switch and press S or wait for 5 seconds to confirm.</p> <p></p> <ul style="list-style-type: none"> ● When using REC or ATC mode, or at high speed of Fwd and Rev mode, the torque settings will be different.(Refer to chapter 6.5 Parameter logic).
Fig 6.2e	<p>E-connect S+ integrated apex locator, if the lip hook is connecting with patient's lip, when the endo file entering root canal, the motor will start automatically. Press < or > to shut off this function if not expected, press ● to start and stop the motor.</p> <p></p> <ul style="list-style-type: none"> ● The motor will start automatically if handpiece (without insulating sleeve) or file touch the patient's lip or operator's fingers (without insulating glove), take care to avoid this, the file rotated by motor has risk of injure someone.
Fig 6.2f	<p>When the endo file out of root canal, the motor will not auto stop with default setting, Press < or > to select auto stop "ON" if needed.</p>
Fig 6.2g	<p>Because of integrated apex locator, when the file reaches the reference point, the motor will response according to setting, it can be Reverse, SlowDown, Stop and Off.</p> <p>Press < or > to change.</p>

	<p>Reverse: rotation direction changing till the file upward a little bit by operator, rotation direction will change back again.</p> <p>SlowDown: rotation slowdown when approach the reference point, will reverse if reach.</p> <p>Stop: rotation stop when reach the reference point, upward a little bit and will rotate again.</p> <p>Off: rotating as usual even if reach the reference point.</p>
Fig 6.2h	<p>During combined length determination, normally apical reverse must active before reaching major apical foramen, Press < or > to set apical reverse position by change the flash bar (①), the motor will reverse while reaching the flash bar every time.</p>
Fig 6.2i	<p>Activating in REC and ATC operation mode. forward angle (Clockwise rotation angle) can be adjusted by operator from 30° to 370°, Press < or > to change.</p>
Fig 6.2j	<p>Activating in REC and ATC operation mode. forward angle (Clockwise rotation angle) can be adjusted by operator from 30° to 370°, Press < or > to change.</p>
	<div style="background-color: #cccccc; padding: 5px;">  <ul style="list-style-type: none"> ● The sum of FWD Angle and REV Angle must be greater than 120°, the motor system has closed the angle not needed. For example: if you set FWD Angle 30°, the REV Angle must be setting greater than 90°. </div>

6.3 Preset programs

Fig 6.3a	<p>For convenience, we preset some common file system. Long press S to entry preset program during standby state, the interface will show as left.</p> <p>New File (①) Indicates a new folder. You can use this function to create a mode. Myfile (②) represents the current memory mode. The operator can replace the current memory mode with preset program (③) .press < or > to change, then press ● to confirm.</p>
Fig 6.3b	<p>In new pattern, long press S to change the pattern name, press < or > to select a letter, press S to confirm the pattern name, and select "√"</p> <p>When you need to delete the letter can be quickly deleted by pressing the ● .</p>
Fig 6.3c	<p>In new mode, continue to press S to enter the "Taper" mode. Press < or > to select the number, and press S to confirm.</p>
Fig 6.3d	<p>In the new mode, continue to press S to enter the "Number" mode, press < or > to select the Number, press S to confirm.</p>

Fig 6.3e	In new mode continue to press S , select the operation mode, press < or > to select forward, reverse, reciprocating or adaptive, press S to confirm.
Fig 6.3f	In new mode, continue to press S to enter the Speed Setting mode. Press < or > to select a speed from 50rpm to 1500rpm, and press S to confirm.
Fig 6.3g	In the new mode, continue to press S to enter the "Torque Setting" mode, press < or > , select torque from 0.5N·cm, 0.8n·cm and 1N·cm, press S to confirm.
Fig 6.3h	Continue pressing S in New mode to enter "Color ring" mode. Press < or > to select the color and press S to confirm.
Fig 6.3i	To delete the new mode, hold down S in the corresponding mode. The interface is as shown in the left picture. Select "Yes" and press S to delete
Fig 6.3j	<p>If you select one of the preset programs, such as "OneCurve" (①) , the rotational speed (②) and torque (③) values will be automatically set according to the manufacturer's recommendations for the file.</p>  <p>Protaper®, GATES®, Pro.Glider®, and Wave one® are registered trademarks of Denberg; Mt two®, Flex.Master®, Reciprocating iproc® and R-Pilot® are registered trademarks of VDW Corporation; K3XF®, TF® are registered trademarks of Spoon Corporation; OneG®, OneShape®, OneFlare®, 2Shape® and OneCurve® are registered trademarks of Micro-Mega, Inc.; XPendo.Shaper®, XPendo.Finisher®, iRace®, BT-Race® and BioRace® are registered trademarks of FKG Corporation. E-flex is our file needle</p>
Fig 6.3k	Some preset modes, such as Flex.Master(①), can be set to different numbers and taper of the file needle: in the corresponding preset mode press < or > to select the appropriate number and taper (②), press S to confirm.
Fig 6.3l	<p>When using preset programs, the current use mode code will be changed to the corresponding preset needle number and taper(①), and the operation mode (④), speed (②) and torque (③) will also be automatically set.</p>  <p>● All memory modes (from M1 to M10) can be replaced by preset programs in this way.</p>
Fig 6.3m	When using preset programs, the parameters still can be changed, and the changed parameters will be marked(①)(②)(Refer to chapter 6.2 Set parameter).

	<p>If you need to reset the parameters of the current preset program, just long press S at standby mode to enter preset mode and press S to choose the preset program again (Such as "One Curve"), the parameters will be reset and the mark will disappear.</p> <p>If want to change back to default setting, long press S to entry preset program during standby state, select "OneCurve" and press S to confirm, the default setting will be recalled, and the 4 corners around will disappear.</p> <p>If want to change back to M1(or M2-M10), long press S to entry preset program during standby state, press < or > to select M1(or M2-M10) again, then press S to confirm</p> <div style="background-color: #cccccc; padding: 5px;">  <ul style="list-style-type: none"> ● It is not recommended to change the preset parameters, unless the operator confirms that the parameters need to be changed, otherwise the instrument have risk to be separated. </div>
--	---

6.4 Advanced setting

Fig 6.4a	Hold S then press • for about 0.5 seconds at power off mode will enter advanced setting, the version will display.
Fig 6.4b	Press S again, the "Auto Power Off" time can be change, press < or > to adjust, then press • to confirm. The "Auto Power Off" time can be set from 3-15 minutes.
Fig 6.4c	Press S again, the "Auto Return time" can be change, it means when setting parameters just like speed and torque, the system will back to standby interface if there is no operation in 5 seconds. press < or > to adjust, then press • to confirm. The "Auto Return time" can be set from 3-60 seconds.
Fig 6.4d	Press S again, the "Beeper Volume" can be change, press < or > to adjust, then press • to confirm. The "Beeper Volume" can be set from 0-3.
Fig 6.4e	In the advanced setting mode, press S again to enter the "Backlight Setting" interface, where you can set the backlight intensity. Press < or > to adjust and press S to confirm.
Fig 6.4f	Press S again, the "Habit hand" can be change, press < or > to adjust, then press S to confirm. The right hand and the left hand can be set. After switching to the left hand usage, the display interface will be rotated 180° for the convenience of the left hand operator to observe the display.
Fig 6.4g	In the advanced setting mode, press S again to enter the "root test sensitivity" setting interface. Press < or > to adjust the root test sensitivity and press

	S to confirm.
Fig 6.4h	In the advanced setting mode, press S again to enter the "Language" setting interface. Press < or > to switch between Chinese and English. Press S to confirm.
Fig 6.4i	Press S again, the "Start memory" can be change, it means every time turn power on, which memory mode will appear first. press < or > to adjust, then press ● to confirm. M1 and Last (the memory mode number when you turn power of) can be set.
Fig 6.4j	Press S again, entry "Calibration" function, press < or > to select "ON", press ● to start calibration.  <ul style="list-style-type: none"> • Before calibrating, making sure the original contra angle is installed, and do not install the file. • The torque will not correct if calibration without original contra angle or any load on contra angle chuck, and has risk of file broken. • When it is found that there is torque bar displayed in the screen of the no-load starting motor, it can be solved by recalibrating.
Fig 6.4k	Press S again, entry "Restore setting" function, press < or > to select "ON", press S to start recovering, all the parameters be set by operator will be recovered by default factory setting  <ul style="list-style-type: none"> • After restoring the settings, all most everything will be recovered to factory settings. If necessary, please record the important parameters in other places.

6.5 Parameter logic

The factory default parameters of the ten memory modes are shown in the table below. The parameters can be adjusted as needed.

Parameter	M0	M1	M2	M3	M4	M5	M6	M7	M8	M9	M10
Operation mode	APEX	Fwd	Fwd	REC	REC	ATC	ATC	Rev	Rev	Fwd	Fwd
Speed (rpm)	/	300	400	350	350	350	300	350	500	800	1000
Torque (N·cm)	/	3.0	2.0	N/A	N/A	1.5	1.5	2.5	2.0	1.5	1.0
Fwd angle	/	N/A	N/A	30	40	370	210	N/A	N/A	N/A	N/A
Rev Angle	/	N/A	N/A	150	160	50	50	N/A	N/A	N/A	N/A

The default advanced settings parameters are shown in the following table. The parameters can be adjusted as needed

Auto Power Off	10min
Auto Return	5s
Volume	3

Apex Sensitivity	Mid
language	English
Calibration	OFF

Back Light	4
Habit Hand	Right

Restore Settings	OFF
Startup Memory	M1

Speed (rpm) settings are different at different operation modes. See the table below for details.

Fwd		Rev					REC	ATC
50	60	70	80	90	100	110	100	100
120	150	200	250	280	300	200	200	
350	400	450	500	550	600	300	300	
650	700	750	800	850	900	400	400	
950	1000	1100	1200	1300	500	500		
1400	1500							

Torque (N·cm) settings are different at different operation modes. Even at the same operation mode, the torque settings will be limited by the speed settings. See the table below for detail.

Fwd/Rev(50-400rpm)						Fwd/Rev(450-700rpm)				
0.5	0.8	1.0	1.5	1.8		0.5	0.8	1.0	1.5	
2.0	2.2	2.5	3.0	3.2	3.5	1.8	2.0	2.2	2.5	
4.0	4.5	5.0	R.L			3.0	3.2	3.5	4.0	
Fwd/Rev(750-1200rpm)						Fwd/Rev(1300-1500rpm)				
0.5	0.8	1.0				0.5	0.8	1.0		
1.5	1.8	2.0								
REC(100-400rpm)						REC(400-500rpm)		ATC(100-500rpm)		
3.0	3.2	3.5	4.0			3.0	3.2	0.5	0.8	1.0
4.5	5.0					3.5	4.0	1.5	1.8	2.0
								2.2	2.5	3.0

Fwd angle (°) and **Rev angle (°)** settings are different at different operation modes. See the table below for details.

	Fwd	Rev	REC								ATC
Fwd angle	/		30	35	40	45	50	55	60	65	The same with the left table
		70	75	80	85	90	95	100	105		
		110	115	120	125	130	135	140			
		145	150	155	160	165	170	175			
		180	185	190	195	200	205	210			
		215	220	225	230	235	240	245			
		250	255	260	265	270	275	280			
		285	290	295	300	305	310	315			
		320	325	330	335	340	345	350			
		355	360	365	370						
		Rev angle	/	The same with the front table							

- The sum of forward angle and reverse angle should be bigger than 120°. Any angle settings not follow this rule are disabled. Such as: Fwd angle is 30°, Rev angle can only be set bigger than 90°.

7. Operation

7.1 Charge

Fig 7.1a	The left figure on the top left of the screen shows the remaining battery level. When it turns right as shown left means low battery, please charge in time.
	 <ul style="list-style-type: none">● If the battery level is less than 15%, it must be recharged within 30 days, otherwise the battery will be irretrievably damaged due to low power.● If you do not use this product for a long time, please charge the product at least once a month.
Fig 7.1b	If the battery power is lower than 15%, the speed and torque may be lower than the set value. As shown in the left figure, the low power alarm will appear on the display screen with continuously use, and the device will automatically shut down.
	 <ul style="list-style-type: none">● As the display of the remaining battery level is based on the battery voltage level, when there is a sudden large torque load during the operation, the display may decrease at the same time.
Fig 7.1c	When not using charging base, connecting the adapter directly to the device can also charge the device, and charging mode will display on the screen. Charging base is recommended for charging (Refer to chapter 4.4 Parameter logic).  <ul style="list-style-type: none">● Please use the original power adapter
Fig 7.1d	Charging mode will display on the screen when charging (①). When the battery is full or nearly full, the screen will stop flashing, and will display as figure shows left(②). It takes about 4 hours to fully charge the battery. If the remaining battery level or situation of the battery is different, the time to fully charge will be different. According to the use state of the battery, the battery can be recharged 300-500 times, and then the battery power will be significantly reduced.

- Do not use the equipment for treatment during charging.

- It is forbidden to replace batteries by non-professionals or untrained personnel. If the wrong battery is used or installed incorrectly, the electronic components will be damaged. Product information is printed on the battery surface covered with film.
- The device should be placed where it is easy to charge and operate disconnected device.

7.2 Motor Operation

Fig
7.2a

Press ● at standby mode to start the motor, then the torque bar will display on the screen (**Refer to chapter 5.2 Screen display for more about Torque Bar**).

- Before using in the treatment, please try it out of the mouth to ensure that the function of the device is normal.
- The root canal file may be damaged suddenly when it enters the root canal which is too curved or not in good shape. When the user feels that the root canal is abnormal, please stop using the device immediately and confirm the correct operation parameters and methods.
- Even if the normal parameters are set, due to the metal fatigue of the root canal file, the instrument will be separated. Therefore, when using the root canal file, do not exceed the times recommended by the manufacturer, and replace it in time.
- When the root canal file is subjected to excessive external force, it may break. When using this device, do not apply excessive external force to the root canal file.
- Do not press the back cover of the contra angle during the treatment, otherwise the device will be damaged, and even the flying file will hurt the patient.
- The electromagnetic noise in the surrounding environment may interfere with the normal operation of the device. Please do not completely rely on the automatic control of the device, and always pay attention to the feedback information on the screen.

- When there is any abnormality, please stop using the device. This device is not suitable for all types of root canals. It is recommended to use according to the instructions of root canal file.
- The root canal file is easy to fracture at high speed. Please follow the rotation speed recommended by the manufacturer. Please check the set speed before use.
- Be careful when using the root canal file with materials other than nickel titanium.
- Please use disposable gloves and rubber barrier for treatment.
- After the treatment, please take out the root canal file to avoid

damage to the root canal file.	
Fig 7.2b	When using motor combine apex function, the measure wire must be connecting with motor by USB socket, and white slot connects with lip by lip hook, keep the black slot idle. The reference point bar will show on the screen (more information about reference point bar, please see chapter 5.2 Screen display).
Fig 7.2c	 <ul style="list-style-type: none"> ● We strongly recommend check the function every time before use. ● Touch the lip hook with the file in the contra angle and check that all the bars on the meter on the screen light up, and the motor should be reversed continuously.

7.3 Apex operation and not suitable condition

Fig 7.3a	When using as a standalone apex locator. We suggest put the handpiece on the charging base to get better visual angle. the measure wire must be connecting with motor by USB socket, white slot connects with lip by lip hook, and black slot connect with file clip. the canal length indicator bar will show on the screen (more information about canal length indicator bar, please see chapter 5.2 Screen display).
Fig 7.3b	The reference point can be adjusting when use stand-alone apex function. Press S to active reference point interface during M0 standby state, Press < or > to change reference point by change the flash bar (①), a continuous beeping with appear when reach the reference point.
Fig 7.3c	 <ul style="list-style-type: none"> ● The file clip must hold the file correctly. ● Push the button of file clip in the direction shown by the arrow, clip the holder onto the metal upper part of the file and then release the button. ● The file clip must be nearly vertical with file handle, otherwise will damage the tip of the file holder.
Fig 7.3d	 <ul style="list-style-type: none"> ● We recommend using Apex Tester to check the accuracy of apex locator every time before use. ● Insert the Apex Tester into handpiece USB socket during M0 mode, check the highlight number must between 01-03 (0.3mm-0.8mm on the top)
Fig 7.3e	 <ul style="list-style-type: none"> ● We recommend check the connection of apex locator every time before use. ● Touch the lip hook with metal part of file clip, check all

	the bars on the meter on the screen light up, and “over” flash on the top.
Unsuitable situation of root canals for Electric Measurement Cannot obtain precise measurements if the root canal conditions as below	
Fig 7.3f	Root canal with a large apical foramen The root canal cannot be accurately measured because of the lesion or incomplete development of the apical foramen. The results may show that the length measured is shorter than the actual one.
Fig 7.3g	Root canal blood overflow from the opening If blood spills from the root opening and contacts the gums, it will cause leakage of electricity, which cannot be accurately measured. Wait for the bleeding to stop completely. Clean the root canal and the opening, completely empty the root canal blood, and then measure it. The root canal uses a chemical solution to flow out from the opening If a chemical solution flows out of the root canal, it is impossible to get an accurate measurement. It is important to remove the overflow from the opening.
Fig 7.3h	Broken crown If the crown is broken, a segment of the gingival tissue enters the lumen, and the contact between the gingival tissue and the root file causes electrical leakage, which cannot be accurately measured. In this case, the appropriate material should be used to isolate the gingival tissue.
Fig 7.3i	The crack tooth Leakage through branch of the root canal Broken teeth can cause electrical leakage and cannot be accurately measured. Branch tubes can also cause leakage.
Fig 7.3j	Retreatment canal which was filled with gutta-percha The gutta-percha must be completely removed to eliminate its insulation, then pass a small file all the way through the apical foramen and then put a little saline in the canal, but do not let it overflow the canal opening.
Fig 7.3k	Crown or metal prosthesis that touches gingival tissue Accurate measurement cannot be obtained if the file touches a mental prosthesis that is touching gingival tissue. In this case, widen the opening at the top of the crown so that the file will not touch the mental prosthesis before taking a measurement.
Fig 7.3l	Cutting debris on tooth Pulp inside canal Remove all cutting debris on the tooth. Remove all the pulp inside the canal. Otherwise an accurate measurement cannot be obtained.
Fig 7.3m	Caries touching the gums In this case, electrical leakage through the caries infected area to the gums are impossible to obtain an accurate

	measurement.
Fig 7.3n	Blocked canal The meter will not run if the canal is blocked. Opening the canal all the way to the apical construction to measure it.
Fig 7.3o	Extremely dry canal If the canal is extremely dry, the meter may not work until it is quite close to the apex. In this case, try to moisten the canal with oxydol or saline.
Difference measuring result between Apex locator reading and Radiography Sometimes the reading of the apex locator reading does not correspond to the X-ray image. this does not mean inaccurate of apex locator or X-ray, depending on the angle of the X-ray beam, the root tip may not be displayed correctly. The position of the root tip seems to differ from its true position.	
Fig 7.3p	The X-ray photo shows that the actual apex of the root canal is not the same as the anatomic end. In fact, the apical foramen is located at the coronal end. in this case, X-ray may indicate that the file needle has not reached the apical foramen, even if it has actually reached the apical foramen.

8. Cleaning, Disinfection and Sterilization

8.1 Foreword

For hygiene and sanitary safety purpose, the components (contra angle, lip hook, file clip and insulating sleeve) must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent uses. Comply with your national guidelines, standards and requirements for cleaning, disinfection and sterilization.

Reprocessing procedures have only limited implications to this dental device. The limitation of the numbers of reprocessing procedures is therefore determined by the function / wear of the device. The contra angle, lip hook, file clip and insulating sleeve are verified to be able to withstand 250 times reprocessing cycles. The device should no longer be reused in case of signs of material degradation.

In case of damage the device should be reprocessed before sending back to the manufacturer for repair.

8.2 General recommendations

- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments, where applicable after sterility.
- For your own safety, please wear personal protective device (gloves, safety glasses, etc.).
- Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, and FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer.
- The water quality must meet the requirements of EN 285/EN 13060/EN ISO 17665.

- Thoroughly clean and wash the components before autoclaving.
- Do not lubricate the motor handpiece.
- Do not clean the contra angle with an ultrasonic cleaning device.
- Do not use bleach or chloride disinfectant materials.

8.3 Autoclavable Components

Autoclavable Components	
Contra Angle Fig 1.2c	Lip Hook Fig 1.2e
File Clip Fig 1.2f	Insulating Sleeve Fig 1.3c
 <ul style="list-style-type: none"> ● Only the components above can be autoclaved. ● Before first use and after each use, sterilize the above components. 	
Reprocessing Instructions	
Preparation at the Point of Use	<p>Disconnect the components (contra angle, lip hook, file clip and insulating sleeve) from the handpiece. Refer to "Chapter 4-Installing the E-connect S+" of this manual for disassembly instructions. Remove gross contaminations from the components with cold water (0-40°C) immediately after use. Don't use a fixating detergent or hot water (40-100°C) as this can cause the fixation of residuals which may influence the result of the reprocessing process.</p> <p>Store the instruments in a humid surrounding.</p>
	 <ul style="list-style-type: none"> ● Do not submerge the components or wipe them with any of the following functional water (acidic electrolyzed water, strong alkaline solution, or ozone water), medical agents (glutaral, etc.), or any other special types of water or commercial cleaning liquids. Such liquids may result in metal corrosion and adhesion of the residual medical agents to the components.
Transportation	Safe storage and transportation to the reprocessing area to avoid any damage and contamination to the environment.
Preparation for Decontamination	<p>The devices must be reprocessed in a disassembled state.</p>  <ul style="list-style-type: none"> ● Do not fail to take out the file before cleaning the contra angle. ● Observe suitable personal protective measures.
Pre-Cleaning	Do a manual pre-cleaning, until the components are visually clean. Submerge the components in a cleaning solution and flush the lumens with a water jet pistol with cold tap water for at least 10 seconds. Clean the surfaces with a soft bristol brush.
Cleaning	Regarding cleaning/disinfection, rinsing and drying, it is to distinguish between manual and automated

	<p>reprocessing methods. Preference is to be given to automated reprocessing methods, especially due to the better standardizing potential and industrial safety.</p> <p>Automated Cleaning:</p> <p>Carefully put the components into the washer-disinfector on a tray and set the parameters as follows, then start the program:</p> <ul style="list-style-type: none"> • 4 min pre-washing with cold water (0-40°C); • emptying • 5 min washing with a mild alkaline cleaner(pH value between 7.5 and 8.5) at 55±2°C; • emptying • 3 min neutralising with warm water (40-60°C); • emptying • 5 min intermediate rinsing with warm water (40-60°C); • emptying <p>Note: The automated cleaning processes have been validated by using 0.5% Neodisher MediClean forte (Dr. Weigert) and Rapid-M-320 Washer-disinfector from SHINVA.</p> <div style="background-color: #cccccc; padding: 5px;">  <ul style="list-style-type: none"> ● Use only approved washer-disinfectors according to EN ISO 15883, maintain and calibrate it regularly. ● Follow instructions and observe concentrations given by the manufacturer (see general recommendations). ● Avoid any contact between the contra angle and any instrument, kit, support or container. </div>
Disinfection	<p>Automated Thermal Disinfection in washer/disinfector under consideration of national requirements in regards to A0 value (see EN ISO 15883).</p> <p>A disinfection cycle of 5 min disinfection at 93±2°C has been validated for the device to achieve an A0 value of 3000.</p> <p>After automated cleaning, the instruments should be automated disinfected immediately. A manual disinfection is not recommended.</p>
Drying	<p>Automated Drying: Drying of outside of instrument through drying cycle of washer/disinfector. If needed, additional manual drying can be performed through lint free towel. Insufflate cavities of instruments by using sterile compressed air.</p>
Functional Testing, Maintenance	<p>Visual inspection for cleanliness of the instruments and reassembling. Functional testing according to instructions of use. If necessary, perform reprocessing process again until instrument is visibly clean.</p> <p>Before packaging and autoclaving, make sure that the components have been maintained acc. to manufacturer's instruction.</p>

	<p>Only the contra angle needs to be lubricated. (Fig 8.3)</p>  <ul style="list-style-type: none"> ● Before autoclaving, the contra angle must be lubricated. ● Attaching the spray nozzle to oil can and contra angle, press the oil can button more than 3 seconds, till all the black oil flow out from the head of the contra angle.
Packaging	<p>The devices that require sterilization can be packaged together in one pouch, but it is necessary to ensure that the packaging pouch is large enough and will not be damaged due to excessive volume.</p>  <ul style="list-style-type: none"> ● Check the validity period of pouch given by the manufacturer to determine the shelf life. ● Use pouches which resist to a temperature up to 141°C and in accordance with EN ISO 11607.
Sterilization	<p>Sterilization of instruments by applying a fractionated pre-vacuum steam sterilization process (according to EN 285/EN 13060/EN ISO 17665) under consideration of the respective country requirements. Minimum requirements: 5 min at 134±2°C Drying time: at least 8min</p>  <ul style="list-style-type: none"> ● Use only approved autoclave devices according to EN 13060 or EN 285. ● Use a validated sterilization procedure according to EN ISO 17665. ● Respect the maintenance procedure of the autoclave device given by the manufacturer. ● Use only this recommended sterilization procedure. ● Control the efficiency (packaging integrity, no humidity, color change of sterilization indicators, physicochemical integrators, digital records of cycles parameters). ● Wait for cooling before touching.
Storage	<p>Storage of sterilized instruments in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.</p>  <ul style="list-style-type: none"> ● Sterility cannot be guaranteed if packaging is open, damaged or wet. ● Check the packaging and the contra angle before using it (packaging integrity, no humidity and validity period).
 <ul style="list-style-type: none"> ● The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing 	

a medical device for use. It remains the responsibility of the processor to ensure that the processing, as actually performed using device, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

8.4 Disinfection Components

Disinfection components				
Motor Handpiece Fig 1.2a	Charging Base Fig 1.2b	Adapter Fig 1.2g	Measuring Wire Fig 1.2d	Handpiece Base Fig 1.3d
Wipe all the surfaces with a cloth lightly moistened with Ethanol for Disinfection (Ethanol 70~80 vol%) at least 2min,repeat for 5 times.				
 <ul style="list-style-type: none"> ● Do not use disinfectants other than alcohol for disinfection. ● Do not use excessive alcohol to prevent alcohol from seeping into the parts and damaging the internal parts. ● Disinfect before and after each use. 				

9. Error Indication

Fig 9a	When motor is running in over torque reversing or R.L. torque set mode, if the load exceeds the device limit, this warning will appear on the screen. Press ● to leave this page to return to standby mode.
Fig 9b	When the battery is nearly empty, this warning will appear on the screen. Press ● to leave this page to return to standby mode and please charge in time.

10. Troubleshooting

When trouble occurs, check this chart before contacting your distributor. If none of these are applicable or the trouble is not remedied even after action has been taken, the product may have failed. Contact your distributor.

Problem	Cause	Solution	Reference
Not able to power on	Low battery	Charge the handpiece.	7.1
	Duration of pressing switch is too short.	Press the switch for more than 0.5 seconds.	5.1
The LED light on	Wrong adapter is used.	Please use original adapter	7.1

the charging base does not work	Poor connection.	Please check the connection	4.4
Charge page does not show on the screen with correct charge operation	Hand piece does not exactly put into the base	Please check the connection	4.4
	The thimble on the charging base cannot rebound normally	Remove the foreign matter between the movable part and the fixed part of the charging thimble	/
	The head contact point of the charging thimble on the charging base is dirty	Clean contact points	/
	The charging base is damaged	Connect the power adapter directly to the hand piece for charging and contact the dealer.	4.4
Screen does not display	Screen does not work	Try to power on the device normally to check if there is any beep sound. If so, press the main switch again to watch whether the motor run. Then contact the dealer.	/
Motor does not work	Contra angle stuck	Take off the contra angle and press the switch to see whether the motor is running. If so, please clean the contra angle or contact the dealer to maintain the contra angle.	/
	System protect or hand piece damage.	Check through Error warning	9
Motor can not be stopped	Problems in the internal circuit.	Contact the dealer	/
Motor auto reverse	Torque exceed the set torque	Check the set torque	6.2
	Reverse mode is set	Check the operation mode set	6.2
Motor does not	Torque set is R.L. mode	Check the set torque	6.2

reverse	Torque does not reach the set torque	Check the set torque	6.2
Motor rotation direction changes frequently	Operation mode is REC or ATC	Check the operation mode set	6.2
No beep sound	Beep volume is "0"	Check the beep volume set	6.4
Beep all the time	Operation mode is REC or torque set is R.L.	Check the operation set or torque set	6.2

11. Technical Data

Manufacturer	Changzhou Sifary Medical Technology Co., Ltd.
Model	E-connect S+
Dimensions	20.1cm x 10.6cm x 13.0±1cm (package)
Gross weight	1.13kg±10%
Contra angle	Gear ratio: 1:1 Compatible with rotary and reciprocating instruments, equipped with $\varnothing 2.35$ mm nickel titanium root canal file conforming to ISO 1797:2017, Type 1, Files length 11-31mm.
Motor handpiece	Input: DC 5V/1A
Battery	DC 3.7V/1900mAh
European standard adapter	Model No: UE05LV2-050100SPA Input: AC 100-240 V,50/60Hz,0.2A Output: DC 5V/1A, 5W
Multi-standard adapter	Model No: UES06WOCP-050100SPA Input: AC 100-240 V,50/60Hz,0.2A Output: DC 5V/1A
Charging base	Input: DC 5V/1A Output: DC 5V/1A
Torque range	0.5 ~ 5.0 N·cm
Speed range	50 ~ 1500 rpm
Classification of anti-electric shock types	Class II at charging mode; Internally powered device at running mode.
Applied part	B (Contra angle, File clip, Lip hook, Insulating sleeve)
Operation mode	Non-continuous, duty cycle: ON 5 mins, OFF 5 mins
Ingress protection	IPX0

Operation conditions	Use: in enclosed spaces Ambient temperature: 10°C ~ 40 °C Relative humidity: 30% ~ 75% Atmospheric pressure :70kPa ~ 106kPa
Transport and storage conditions	Ambient temperature: -20 °C ~ +55 °C Relative humidity: 20% ~ 80 % Atmospheric pressure: 70kPa ~ 106kPa

12. EMC Tables

Guidance and manufacturer's declaration - electromagnetic emissions		
The E-connect S+ is intended for use in the electromagnetic environment specified below. The customer or the user of the E-connect S+ should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Professional healthcare facility environment and Home healthcare environment
RF emissions CISPR 11	Class B	Professional healthcare facility environment
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	
 <p>The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.</p>		

Guidance and manufacturer's declaration - electromagnetic immunity			
The E-connect S+ is intended for use in the electromagnetic environment specified below. The customer or the user of the E-connect S+ should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge	+/- 8 kV contact +/- 2 kV, +/- 4	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/-	Floors should be wood, concrete, or ceramic tile. If floors are covered with

(ESD) IEC 61000-4-2	kV, +/- 8 kV, +/- 15 kV air	15 kV air	synthetic material, the relative humidity should be at least 30 %.
Electric al fast transients/bursts IEC 61000-4-4	±2kV 100kHz repetition frequency	±2kV 100kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Line to line: ±0.5kV, ±1kV	Line to line: ±0.5kV, ±1kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips IEC 61000-4-11	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of devices requires continued operation during power mains interruptions, it is recommended that devices be powered from an uninterruptible power supply or a battery
Voltage interruptions IEC 61000-4-11	0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0° 0% UT; 250/300 cycle	0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0° 0% UT; 250/300 cycle	
Rated Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz or 60Hz	30 A/m 50Hz or 60Hz	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: UT: rated voltage(s); E.g., 25/30 cycles mean 25 cycles at 50Hz or 30 cycles at 60Hz			

Guidance and manufacturer's declaration – electromagnetic immunity

The **E-connect S+** is intended for use in the electromagnetic environment specified below. The customer or the user of the **E-connect S+** should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V 0.15 MHz – 80 MHz, 6 V in ISM bands between 0.15 MHz and 80 MHz, 80 % AM at 1 kHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the E-connect S+ , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF EM fields IEC 61000-4-3	3 V/m, 80 MHz – 2,7 GHz, 80 % AM at 1 kHz	3V/m	
Proximity fields from RF wireless communication equipment IEC 61000-4-3	See the RF wireless communication equipment table in "Recommended minimum separation distances"	Complies	Recommended minimum separation distances See the RF wireless communication equipment table in "Recommended minimum separation distances"

Recommended minimum separation distances

Nowadays, many RF wireless equipment has being used in various healthcare locations where medical equipment and/or systems are used. When they are used near medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. The **E-connect S+** has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2020. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and the **E-connect S+** as recommended below.

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380	TETRA 400	Pulse	1.8	0.3	27

	- 390		modulati on 18Hz			
450	430 - 470	GMRS 460 FRS 460	FM \pm 5 kHz deviatio n 1 kHz sine	2	0.3	28
710 745 780	704 - 787	LTE Band 13, 17	Pulse modulati on 217Hz	0.2	0.3	9
810 870 930	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulati on 18Hz	2	0.3	28
1720 1845 1970	170 0- 199 0	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulati on 217Hz	2	0.3	28
2450	240 0- 257 0	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulati on 217Hz	2	0.3	28
5240 5500 5785	510 0- 580 0	WLAN 802.11 a/n	Pulse modulati on 217Hz	0.2	0.3	9

Guidance and manufacturer's declaration – electromagnetic immunity

The **E-connect S+** is intended for use in the electromagnetic environment specified below. The customer or the user of the **E-connect S+** should assure that it is used in such an environment.

Proximity magnetic fields	IEC 61000-4-39 test level	Compliance level	Electromagnetic environment – guidance
Proximity magnetic fields	134.2kHz Pulse modulation 2.1 kHz	65A/m	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Proximity magnetic fields	13.56MHz Pulse modulation 50 kHz	7.5A/m	

- Use of accessories and cables other than those specified or provided by the manufacturer of E-connect S+ could result in increased electromagnetic emissions or decreased electromagnetic immunity of E-connect S+ and result in improper operation.

Cable information:

Cable Name	Cable Length (m)	Shielded or not	Remark
Adapter Cable	1.2	No	/

- Use of E-connect S+ adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, E-connect S+ and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the E-connect S+, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- If the use location is near (e.g. less than 1.5 km from) AM, FM or TV broadcast antennas, before using this equipment, it should be observed to verify that it is operating normally to assure that the equipment remains safe with regard to electromagnetic disturbances throughout the expected service life.

13.Statement

Service Life

The service life of E-connect S+ series products is 3 years. It is recommended that the equipment be checked and repaired at the dealer once a year.

Maintenance

MANUFACTURER will provide circuit diagrams, component part lists, descriptions, calibration instructions to assist to SERVICE PERSONNEL in parts repair. MANUFACTURER needs to carry out technical training and technical support for SERVICE PERSONNEL, so as to realize the correct maintenance of the product. Repairing the product by unauthorized personnel may cause damage to the product and become unusable.

Disposal

The package should be recycled. Metal parts of the device are disposed as scrap metal. Synthetic materials, electrical components, and printed circuit boards are disposed as electrical scrap. The lithium batteries are disposed as special refuse. Please deal with them according to the local environmental protection laws and regulation.

Rights

All rights of modifying the product are reserved to the manufacturer without further notice. The pictures are only for reference. The final interpretation rights belong to Changzhou Sifary Medical

Technology Co., Ltd. The industrial design, inner structure, etc, have claimed for several patents by SIFARY, any copy or fake product must take legal responsibilities.

Changzhou Sifary Medical Technology Co., Ltd.

Add: No. 26 Yandanghe Road, Xinbei District, 213000 Changzhou, Jiangsu, China.

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All rights reserved.