# bestdent

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This package insert is not a reference to surgical techniques. It is designed to assist in using this product.

**DESCRIPTION:** The absorbable, sterile, malleable, and water-insoluble haemostatic sponge is designed for haemostatic use by direct application to bleeding surfaces. It is non-pyrogenic and biocompatible. It is manufactured from highly purified pharmaceutical-grade granules and purified water for use in dental surgical procedures. When implanted in vivo and used in adequate amounts, it is completely absorbed in less than 4 weeks.

**INDICATIONS:** In oral and dental surgery, haemostatic sponges serve as an aid in providing haemostasis. They may be used either dry or moistened, depending on the conditions present during the operation and the surgeon's preference. Sterile isotonic sodium chloride solution (sterile saline) is suitable for use.

**CONTRAINDICATIONS:** The sponge should not be used:

- For closure of skin incisions, as it may interfere with healing of skin edges. This is due to mechanical interposition, not to intrinsic interference with wound healing.

**IN PATIENTS WITH KNOWN ALLERGIES TO PORCINE COLLAGEN (SEE WARNINGS)**
**HANDLING INSTRUCTIONS:**

- The sponge can be used dry or saturated with a sterile isotonic sodium chloride solution (sterile saline solution) and can be cut to the desired size.

**CONTRAINDICATIONS:**

- When applied dry, cut a piece of sponge to the desired size, roll it between your fingers, and compress it slightly so that it fits into the cavity or socket. After insertion into the cavity or socket, slight pressure should be applied with a finger for one to two minutes.

**When used with sterile saline solution, it should first be immersed in the solution and then removed, squeezed between gloved fingers to expel air bubbles, and then placed back into the saline solution until needed.**

The sponge should quickly return to its original size, with a slight increase in thickness and shape in the solution.

- If it does not, it should be removed again and kneaded vigorously until all the air is expelled and it expands back to its original size, with slight increases in thickness and shape when returned to the sterile saline solution. If it used wet, it may be blotted with gauze until just damp before applying it to the bleeding site. It should be held in place by applying moderate pressure with a small gauze pad or gauze until hemostasis is achieved.

**Removal of the pad or gauze should be carried out after moistening it with a few drops of sterile saline solution, in order to minimize the risk of adherence and prevent displacement of the hemostatic sponge, which at that point should contain a stable clot.**

- For additional applications, fresh pieces prepared as described above should be used. Use only the minimum amount necessary to achieve hemostasis. Once hemostasis has been achieved, any excess sponge should be carefully removed.

**WARNING:**

- Patients with a history of allergies to porcine products may be at risk of severe acute hypersensitivity reactions, including anaphylaxis (see **CONTRAINDICATIONS**). If an anaphylactic reaction occurs, administration of the absorbable sponge must be discontinued immediately and any application removed.

**Overpacking with the sponge should be avoided, as it returns to its original volume and may interfere with normal function and/or cause compression necrosis of surrounding tissue and nerve damage. If used for packing, the sponge should be removed after haemostasis.**

**Use with caution in contaminated areas of the body. If signs of infection or abscess develop at the placement site, reoperation may be necessary to remove the infected material and allow drainage.**

- Do not re-sterilize. Do not use if the package is open or damaged. This device is intended, tested, and manufactured for single use only.

**The absorbable haemostatic sponge is not intended to substitute for meticulous surgical technique or proper application of sutures or other conventional haemostatic procedures.**

**PRECAUTIONS:**

- Use only the minimum amount necessary to achieve haemostasis, leaving it in place until bleeding stops and removing excess material afterwards.

The sponge should not be used with methyl methacrylate adhesives.

**It is not recommended for the primary treatment of coagulation disorders.**

- It should not be used in the presence of infection unless accompanied by antibiotics in infected wounds.
- This product should be used only by healthcare professionals trained in procedures and surgical techniques involving absorbable haemostatic sponges.

**The sponge is packaged in a sterile blister. Once the package is opened, the contents are subject to contamination. It is recommended to use the sponge immediately after opening and discard any unused material.**

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