optiprint clara

Germany's favorite 3D Resin

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1. Produktbeschreibung

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EN - instruction manual

1. Product description

optiprint lumina is a light-curing 3D printing resin for the additive manufacturing of custom-made bite splints for various medical indications prescribed by dental professionals. It is suitable for use in devices with a wavelength of 385 nm and 405 nm, unless otherwise stated on the label. The target patient group of this Class IIa medical device is defined as adults and adolescents.

2. Intended purpose

optiprint clara is a 3D printing resin for the additive manufacturing of bite splints.

3. Contraindications

The material should not be used for any purpose other than additive manufacturing of the intended purpose. Do not use the polymer material in the case of allergies to any of the ingredients (contains methacrylate monomers and oligomers).

4. Risk minimization and safety instructions

• Improper use and deviations from the described processing will result in impaired quality and biocompatibility as well as undesirable mechanical properties of the finished molded part.

• Biocompatibility is only guaranteed with proper application (post-polymerization without oxygen). Only pro-cess all molded parts when they are fully polymerized.

• The post-polymerization of the molded parts takes place in a suitable polymerization device (e.g., otoflash G171), see 5.4.

• After completion of the construction process, the molded part should be cleaned with suitable cleaning solution (e.g., isopropanol 99%) in an ultrasonic tank.

• The batch no. must be indicated in any procedure requiring the material to be identified.

• For additively manufactured medical devices made of printed resins, storing in water for 24 hours is recommended.

• Observe all recommended settings of the printer and the light-curing device.

• Read and observe the Safety Data Sheet (SDS) before use.

• For rupture safety, usual material wall thicknesses must be observed when creating the custom-made products (see 5.1)

• Proper personal protective equipment (nitrile gloves, safety glasses, protective clothing) must be worn when handling optiprint products as well as uncured printed molded parts.

• Avoid contact with skin and eyes before post-curing. The optiprint product can irritate eyes and skin.

• In rare cases, allergic reactions to components of optiprint products may occur. In case of accidental contact, follow first aid measures (rinse thoroughly with water and consult a doctor if necessary). See the SDS.

• To be used by qualified professionals only. Keep out of reach of children!

5. Processing steps

5.1 Designing

The minimum thickness must not be less than 0.9 mm.

Occlusal thickness	≥ 0.9 mm
Peripheral thickness	≥ 1.3 mm

Make sure that the construction is positioned at an angle of 45-50° and that it is sufficiently supported by supports. Depending on the printer, a base plate may be required for supported objects.

5.2 Printing

Observe the instructions for use of the printer and the software. Ensure that work is carried out in clean conditions. Contamination on the 3D printer can cause defects on the molded part and damage the tray. Corresponding material parameters for optiprint printing resins can be downloaded from the printer manufacturer's database.

Shake the optiprint bottle before use, and fill the material tray of the printer with sufficient material. Remove any bubbles with a clean instrument/spatula.

If you are using a printer without integrated heating, it is recommended to bring both the printer and the optiprint product to an operating temperature of 30°C/86°F. A cold start should be avoided.

5.3 Cleaning after printing

Due to the high viscosity, the printing resin does not drip completely off the building platform after the end of printing. Strip off the residue of the printing resin with a clean wooden spatula. Then remove the building platform from the printer and detach the molded parts from the building platform. Disconnect the support structures before cleaning.

Recommendation

In a non-heated ultrasonic bath: place cleaning fluid and the molded part in a second, sealable container (5 minutes). Isopropanol (99%) and optiprint clean can be used as cleaning fluids. When using optiprint clean, a 2-minute follow-up cleaning with isopropanol is always required. optiprint clean is not suitable for cleaning surfaces and equipment.

Following this, dry the molded part with compressed air and check the results of the cleaning.

Complete the post-treatment of the molded parts quickly and adhere to the time specifications! In particular, avoid an unnecessarily long bath in the cleaning fluid and long waiting times between steps!

Alternatively, use the corresponding Rapid Shape Wash program with isopropanol (99%).

You can determine cleaning has been sufficient if the surface is matt. Shiny areas require spot post-cleaning with cleaning fluid and a brush. Clean until noshiny areas are visible.

5.4 Finishing off

The final properties and also the final color depend on the post-curing process!

To achieve the desired material properties and biocompatibility as well as the coloration, the fully cleaned and dried molded parts must be post-cured in oxygen-free conditions.

Recommendation

Otoflash G171 flash-light device (NK Optik) with nitrogen flooding.

Ensure you use the Plexiglas tub with a UVB filter! It can be identified by the inscription: NK Optik 360N2.

When post-curing in other light curing devices, make sure that the curing device has a product-specific, predefined setting or a comparably high energy input (200 W).

OTOFASH G171 LIGHT OUTPUT	2 x 3000 flashes (turn after 3000 flashes)
ALTERNATIVE LIGHT OUTPUT (200W, 315-400 NM)	10 minutes

6. Final cleaning

Dental molded parts produced with optiprint can be polished in the conventional way. A final cleaning is performed with detergent and water in a cold ultrasonic bath (5 minutes) until no polishing agent residues can be seen and felt on the molded part. Then rinse the molded part in clear water for 30 seconds until no more foaming from the detergent is visible.

7. Instructions for the laboratory/dentist for delivery of the molded part to the patient

Storing in water for 24 hours is recommended before inserting the molded parts. Disinfection is not necessary as standard.

If you want to ensure hygiene when handing over the medical device to your customer, use the optiprint prevente disinfection bath (NW-Chemie GmbH), which is especially suitable for 3D printed molded parts. Spectrum of activity: Bactericidal incl. TBC, levurocidal, limited virucidal (HIV, HBV, HCV, Sars-CoV-2).

optiprint clara must not be autoclaved. Medical devices made of optiprint clara are thermoflexible. Before placing/removing the mold-ed part on the model or inserting it into the patient, please place it in warm water (approx. 40°C) for 2 minutes. This will make the material adapt even better to the shape of the teeth. The user must be careful not to cause reinfestation with the water bath. Pass on the following handling information to the patient and check the splint during follow-up examinations.

Instructions for the patient

Comfort: becomes flexible due to body temperature. Before inserting the splint, place it in water at about 40°C for 2 minutes. Do not wear during eating.

Cleaning: after wearing, clean with a soft toothbrush under running water. Do not use toothpaste, denture cleaners, cleaning solutions or alcohol-based disinfectants. Storage: keep dry until next use.

Protect the splint from excessive heat and cold.

Avoid direct sunlight

8. Technical data

Please see the separate TDS (Technical Data Sheet).

9. Composition

Acrylic mixture, photoinitiator, inhibitor.

10. Storage

Protect the product from strong light and heat sources, recommended storage temperature: 5°C to 30°C.

Close the bottle after each use. If storage out of light and dust is ensured, the material can remain in the material tray after printing.

Do not use after the expiry date.

11. Disposal

Dispose of liquid printing resin in accordance with official regulations. Must not be disposed of with household waste. Do not discharge into drains.

12. Information

Report all serious incidents (death, serious deterioration of health, serious risk to public health) related to the medical device to the manufacturer (email: mailbox@dentona.de) and the competent authority of the Member State.

If defects in the material occur within the warranty period, the user is only entitled to replacement of the material if the defect is the responsibility of the manufacturer. dentona AG shall not be liable for any loss or damage caused by this material, whether direct or indirect, or special incidental or consequential damages, regardless of the legal basis. dentona AG is liable for direct material damage to the material caused by intent or gross negligence on the part of its legal representatives or executive employees, as well as for personal injury in accordance with the statutory provisions. Any liability for the material and consequential damages resulting from its use are excluded if the user has not followed the specified steps of the procedure.

13. Symbols

 Please follow the instructions for use (also available at www.dentona.de)

 Use by/expiry date

 Protect from sunlight

 Caution: hazardous to health

 Medical device

 Unique Device Identification


 Article number

 Caution: hazardous to the environment

 CE marking with number of the notified body

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