INSTRUCTIONS FOR USE

FIRST STEP

- 1. Place the guide on the implant and keep it in an axial position with the locking key
- 2. Place the Initial extraction drill LEFT Ø 1mm in the handpiece.

IMPORTANT:

- A) Program the surgical device at 1200 rpm in counterclockwise mode.
- B) Use abundant irrigation to cool the drill during milling.

- 3. Insert the Initial extraction drill LEFT Ø 1mm inside the guide, and begin to mill with pressure until the screw engages with the threaded part of the drill. The milling residues have to be sucked. If this operation is not efficacious:

- 4. Place the complementary drill LEFT in the handpiece.

IMPORTANT:

- C) Program the surgical device at 100 rpm in COUNTERCLOCKWISE mode.
- D) Use abundant irrigation to cool the drill during milling.

- 5. Place the complementary drill LEFT inside the guide, and begin the extraction with a pressure until the screw engages with the threaded part of the drill. The residues of this friction have to be sucked.

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COMPONENTS TO USE (not included in the screw removal kit Proclinic®).

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In the hypothetical case that you have not been able to perform the extraction of the fragment of the fractured screw, following the instructions of the first phase, and before to judge the implant as lost, you can follow the steps below. None of the elements described here is included in the screw removal kit Proclinic®, and have to be purchased separately.

- A) For implants with internal metric 1.8mm: This procedure has to be done without the drill's guide. Use the ball drill diameter 1.4 mm, inserted in the handpiece and set it in the milling done previously. Press intermittently, using a maximum of 850 rpm, in the clockwise direction with abundant irrigation, axially, until to remove most of fractured screw, taking care not to deviate from the direction of the original drilling. The residues generated by the cutting have to be removed.
- B) For implants with internal metric 2.0mm: This procedure has to be done without the drill's guide. Use the ball drill diameter 1.6 mm, inserted in the handpiece and set it in the milling done previously. Press intermittently, using a maximum of 850 rpm, in the clockwise direction with abundant irrigation, axially, until to remove most of fractured screw, taking care not to deviate from the direction of the original drilling. The residues generated by the cutting have to be removed.

WARNING: The presence of fragments of screw inside the implant, poor irrigation and failing of cleaning of the thread, or even an excessive force in the insertion of the tap, could cause its breakage.

In both cases, remove any residue of screw inside the implant with a forceps or a probe.

At this point must be used taps for the internal thread of the implant diameter. Follow these instructions about the internal thread of the implant for its proper recovery.


- A) Turn 180 ° clockwise mode.
- B) Turn back 90 ° counterclockwise mode.
- C) Repeat this sequence until, after some repetitions, you will have restored the expected functionality.

WARNINGS:

Reuse of disposable involves a possible deterioration of its characteristics, which implies the risk of tissue infection, surgical failure and / or worsening of the patient's health.

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PRECAUTIONS

Unfavorable load conditions and mechanical stresses could cause fracture of the screw. It is important to determine the cause of the fracture, because their removal and replacement does not solve the problem completely.

INITIAL EXTRACTION DRILL AND COMPLEMENTARY DRILL ARE TO BE CONSIDERED DISPOSABLE.

Store the product in a dry controlled environment, at an average temperature between + 10 ° C and + 40 ° C, like any other surgical instrument.

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RESPONSIBILITY AND WARRANTY

The user must ensure that the product used is suitable for the intended purpose. Neither these instructions or procedures discharge the user from the obligation to use it properly. Please follow the instructions to make the process as efficient as possible. All accessories used with the extraction kit Proclinic® must be original.

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