VOCO

Bifix® Temp

EN Instructions for use

MD EU Medical device

Product description:

Bifix Temp is a dual-curing, composite-based temporary luting material. The balanced translucency of **Bifix Temp** enables highly aesthetic results to be achieved for non-permanent luting of temporary restorations made from crown and bridge materials as well as temporary CAD/CAM-fabricated pieces. The low film thickness ensures precise placement of temporary restorations at all times. The stability of **Bifix Temp** facilitates particularly simple and reliable removal of excess material. Optionally available as **Bifix Temp fluorescent**: Fluorescent when illuminated with UV-A light source for easy detection of **Bifix Temp fluorescent**.

Working time: > 1.5 min (23 °C)
Setting time: < 4 min (37 °C)

Indications:

– Non-permanent luting of temporary restorations

Contraindications:

Bifix Temp contains dimethacrylates and diacrylates, benzoyl peroxide, amines and BHT. **Bifix Temp** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Bifix Temp is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Bifix Temp should only be used by a professionally trained dental practitioner.

Caution:

The surface of composite build-ups and restorations must be coated with a thin layer of insulation (glycerine gel, Vaseline) in order to prevent **Bifix Temp** from luting permanently.

The use of (meth)acrylate-based bonding agents and protective varnishes prior to luting with **Bifix Temp** must be avoided as this can result in a permanent bond.

Use:

1. Preparation:

Ensure that the surface is clean and dry; if necessary it should be cleaned and dried. Follow the instructions for application on composite build-ups and restorations (see **caution** above). Bring material to room temperature before use.

2. Placing the restoration:

Remove the protective cap from the QuickMix syringe. Due to technical reasons, there may be slight differences in the filling levels of the syringe chambers prior to first use. Therefore remove the syringe cap and extract material until both chambers deliver uniformly. Then fit a mixing tip type 9. The notches on the syringe and the mixing tip must line up exactly. Turn the mixing tip 90° clockwise to lock the tip into place. Apply **Bifix Temp** to the dry temporary restoration and/or the prepared tooth surface. Then place the restoration by applying gentle pressure, generating excess material in the process. The working time at room temperature is at least 1.5 min. The attached mixing tip is intended for single use only. After use, store the QuickMix syringe tightly sealed. Check that material flows freely through the openings before reusing. The warranted product characteristics can only be achieved when using the supplied original mixing tips, which are also available as accessories.

3. Removal of excess material:

Any excess material can be removed using the following methods:

- Wait until the material turns elastic (monitor change in consistency as a consequence of intraoral chemical curing) and remove excess material using a suitable instrument.
- Alternatively, the surplus material can be rendered elastic after short curing with a polymerisation light (around 2 s) and then removed using a suitable instrument.
- If not cured, **Bifix Temp** can be removed using a foam pellet or a similar object.

Removal of excess material following brief curing is recommended. The restoration must be held in position during the polymerisation phase. Contact points and interproximal areas must be cleaned with dental floss. There must be no excess material remaining in the sulcus as this can cause gingival irritation.

4. Curing:

Self-curing is completed 4 min after insertion.

The polymerisation time can be controlled by additional light-curing.

For this, each side is cured for a minimum of 10 s using a polymerisation light with an intensity of at least 500 mW/cm².

Warnings, precautionary measures:

- UV light can cause damage. Patients and dental personnel should always take note of the corresponding safety measures. This includes wearing appropriate protective UV goggles, among other things. Do not look into the light source and do not aim the light directly into the eyes. Only use the light at the treatment site.
- **Bifix Temp** is eugenol-free and well-suited for use prior to subsequent adhesive permanent luting as it has no influence on the curing of the materials used.
- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Bifix Temp**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Bifix Temp**.
- If, during the course of treatment, the temporary restoration is to be luted again, the application of a thin layer of, for example, glycerine gel or Vaseline to the restoration prior to luting can make it easier to subsequently remove any **Bifix Temp** residue from the temporary restoration.
- In case of contact with the eyes, immediately rinse with plenty of water and seek medical advice.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Storage:

Storage at 4°C-23°C. Refrigeration is recommended. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Bifix Temp** must be reported to VOCO GmbH and the responsible authority.

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VC 60 AA1055 E1 1124 99 © by VOCO