## TruNatomy®

|       |                                  |     |
|-------|----------------------------------|-----|
| EN    | Instructions For Use .....       | 2   |
| FR    | Instructions d'utilisation ..... | 7   |
| IT    | Istruzioni per l'uso .....       | 13  |
| DE    | Gebrauchsanweisung .....         | 19  |
| ES    | Instrucciones de uso .....       | 25  |
| PT    | Instruções de Utilização .....   | 31  |
| NL    | Gebruiksaanwijzing .....         | 37  |
| SV    | Användarinstruktioner .....      | 43  |
| DA    | Brugsanvisning .....             | 48  |
| NO    | Brukerinstruksjoner .....        | 54  |
| LT    | Naudojimo instrukcija .....      | 59  |
| LV    | Lietošanas norādījumi .....      | 65  |
| ET    | Kasutusjuhised .....             | 71  |
| CS    | Návod k použití .....            | 77  |
| SK    | Návod na použitie .....          | 83  |
| HU    | Használati utasítás .....        | 89  |
| RO    | Instrucțiuni de utilizare .....  | 95  |
| SR    | Uputstva za upotrebu .....       | 101 |
| BS    | Uputstva za upotrebu .....       | 107 |
| CNR   | Uputstva za upotrebu .....       | 113 |
| SQ    | Udhëzimet për përdorim .....     | 119 |
| EL    | Οδηγίες χρήσης .....             | 125 |
| TR    | Kullanım Kılavuzu .....          | 131 |
| PT-BR | Instruções de uso .....          | 136 |
| KO    | 사용 지침 .....                      | 142 |

TruNatomy®

EN

FOR DENTAL USE ONLY  
FOR PROFESSIONAL USE ONLY  
STERILE – STERILIZED BY RADIATION  
LIABLE FOR REPROCESSING

## INSTRUCTIONS FOR USE

### REUSABLE STERILE ENDODONTIC INSTRUMENTS

#### 0) INDICATIONS FOR USE

- Indication for use:  
The products are indicated for treatment of endodontic disease.
- Intended purpose:  
Engine-driven instrument intended for root canal preparation (glide path, shaping and debridement of the root canal).
- Intended User:  
Endodontic instruments are to be used only in a clinical or hospital environment, following good dental practice, by qualified dental professionals such as general practitioners as well as Endo specialists (Endodontist) and Dental Assistants.
- Instruments shall be used in combination with an endodontic motor (manually in severe curvatures).

#### 1) CONTRAINDICATIONS

None.

#### 2) CONTENT, COMPOSITION AND COMPATIBLE DEVICES

TruNatomy® instruments are made of three main components: the working part made of a nickel-titanium alloy, a colored silicone stopper and a plated brass shank with a colored ring (except for Orifice Modifier).

The files are manufactured with a process that results in a file that has a colored appearance.

Due to this proprietary processing, the files may appear slightly curved. This is not a manufacturing defect. While the file can be easily straightened using only your fingers, it is not necessary to straighten the file prior to use. Once inside the canal, the file will follow the root canal anatomy.

The TruNatomy® range comprises following endodontic instruments:

| Code<br>xx-length<br>yyy-file size | Product                        | ISO XXX YYYv<br>XXX-ISO size<br>YYY-taper size<br>v-variable taper                           | Available<br>lengths | Packaging and content                              |
|------------------------------------|--------------------------------|----------------------------------------------------------------------------------------------|----------------------|----------------------------------------------------|
| BST TNMM<br>4xx GLR                | TruNatomy® Glider              | ISO 017 002v                                                                                 | 21/25/31 mm          | Blister with 4 instruments (sterile)               |
| BST TNMM<br>416 OMD                | TruNatomy® Orifice<br>Modifier | ISO 020 008                                                                                  | 16 mm                | Blister with 4 instruments (sterile)               |
| BST TNMM<br>4xx yyy                | TruNatomy® files               | Small: ISO 020 004v<br>Prime: ISO 026 004v<br>Medium: ISO 036<br>003v<br>Large: ISO 046 002v | 21/25/31 mm          | Blister with 4 instruments (sterile)               |
| BST TNMM<br>5xx AST                | TruNatomy®<br>Assortment       | ISO 017 002v<br>ISO 020 008<br>ISO 020 004v<br>ISO 026 004v<br>ISO 036 003v                  | 21/25/31 mm          | Blister with 1 instrument of each ISO<br>(sterile) |
| BST TNMM<br>6xx AST                | TruNatomy®<br>Assortment       | ISO 017 002v<br>ISO 020 008<br>ISO 020 004v<br>ISO 026 004v<br>ISO 036 003v<br>ISO 046 002v  | 21/25/31mm           | Blister with 1 instrument of each ISO<br>(sterile) |
| BST TNMM<br>3xx SQC                | TruNatomy®<br>Sequence         | ISO 017 002v<br>ISO 020 008<br>ISO 026 004v                                                  | 21/25/31 mm          | Blister with 1 instrument of each ISO<br>(sterile) |

- Torque control devices ensure optimal usage.
- Use with endodontic motors in a constant rotation at a speed of 500 rpm.
- Set torque at: 1.5 Ncm.
- Lubricants such as NaOCl, EDTA, ProLube®, Glyde™ shall be used.
- The use of radiographs in combination with an apex locator and a tool for adjusting the silicon stopper to the correct working length is the appropriate method of working length determination.

### 3) WARNINGS

- Strictly follow this Instruction for Use and the General Processing Instructions For Endodontic Products (see section 8) to minimize following risks to the device, the patient and/or the user:
  - Breakage of instruments.
  - Cross-contamination
  - Heat generation due to insufficient lubrication and irrigation.
  - Swallowing of working part of the instrument.
  - Toxic or allergic reactions caused by processing residues.

- Similar to all machine-driven root canal instruments, TruNatomy® instruments should not be used in a root canal with an abrupt apical curvature due to heightened risk of separation. In this case, pre-curved hand files should be used in the apical region.

#### 4) PRECAUTIONS

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Use a rubber dam system during the endodontic procedure.
- For your own safety, wear personal protective equipment (gloves, glasses, mask).
- Inspect the packaging before use and do not use the instruments if the packaging is damaged.
- Do not use the instruments after expiration date.
- Check the instrument before each use for signs of defects such as deformations (bent, unwound), breakage, corrosion, damaged cutting edges, loss of color coding or marking. With these indications the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Before using any instrument, make sure it is well connected to the contra-angle head.
- Check instrument and clean working part frequently during instrumentation, inspecting for signs of distortion, elongation or wear, such as uneven flutes, dull spots. With these indications that the devices are not able to fulfil the intended use with the required safety level, instruments should be discarded.
- Instruments should not be fully immersed in a sodium hypochlorite solution (NaOCl). Only the working part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- Exercise caution in the apical area and in canals that divide, and/or exhibit abrupt curvatures or recurvatures.
- Irrigate abundantly and frequently the canal throughout the procedure and after each instrument used (according to good dental practices).
- Always use minimal apical pressure. Never force the files down the canal.
- When instrument does not easily progress, clean and inspect the cutting flutes, then irrigate, recapitulate with a manual file and reirrigate.
- For shaping extremely curved canals it is safer to use the file only to shape one canal in order to reduce the risk of breakage. Pay attention to the following good practices:
  - Use a new file and discard it after the canal was treated (single canal use).
  - Use manual instead of rotary files.
  - Use small size, flexible or/and NiTi files (this will help avoid canal transportation).
  - Visually inspect the working part for all the defects listed in the former paragraph during use (i.e after each wave).
  - Avoid the standard reaming continual rotational motion and instead use small angle motions (filing motion, watch winding oscillation motion, or balanced force technique) in order to limit the rotational bending fatigue on the instruments and improve their expected life.

## 5) ADVERSE REACTIONS / UNDESIRABLE SIDE EFFECTS

Following undesirable side effects may occur:

- Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

## 6) STORAGE CONDITIONS

Keep products in a dry and clean environment, away from sources of moisture and direct sunlight.

## 7) STEP BY STEP INSTRUCTIONS

- 1) Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root canal.
- 2) Estimate the working length using well-angulated preoperative radiographs.
- 3) Prepare a conservative access cavity sufficient enough to reveal all root canal orifices.
- 4) Scout the coronal 2/3 of canals with an ISO 010 K-File in the presence of lubricant such as ProLube® and Glyde™ FILE PREP and irrigate.
- 5) Followed by a TruNatomy® Orifice Modifier at 500 rpm and 1.50 Ncm. With irrigant in the canal advance TruNatomy® Orifice Modifier in 2-3 gentle amplitudes approximately 2-5 mm in-and-out of the canal. Repeat until the coronal third is shaped. The instrument has 7 mm of cutting flutes, insertion into the canal should not exceed the length of the cutting portion. Irrigate the canal and clean cutting flutes routinely.
- 6) Scout the whole root canal with an ISO 010 K-File, determine working length by using an electronic apex locator (EAL) in combination with radiographs, irrigate and confirm patency.
- 7) With irrigant in the canal, create and confirm a reproducible glide path using a TruNatomy® Glider in 2-3 gentle amplitudes approximately 2-5 mm. Irrigate and repeat until previously confirmed working length with an EAL has been reached.
- 8) ALWAYS begin shaping with TruNatomy® PRIME (500 rpm / 1.5 Ncm) file passively in the presence of sodium hypochlorite with no more than 2-3 gentle amplitudes approximately 2-5 mm in-and-out of the canal. Irrigate and repeat as necessary until full working length is reached. Upon reaching length, remove the file to avoid over-enlarging the apical canal portion.
- 9) Routinely irrigate the canal and clean the instruments cutting flutes of debris upon removal.
- 10) If TruNatomy® PRIME does not progress easily, remove, irrigate, and recapitulate with an ISO 010 K-File to confirm canal patency and move to the TruNatomy® SMALL.
- 11) Inspect cutting flutes routinely upon removal for presence of unwinding and straightening. If deformation is noted, discard and use a new TruNatomy® instrument.
- 12) Advance TruNatomy® SMALL passively in the presence of sodium hypochlorite with no more than 2-3 gentle amplitudes approximately 2-5 mm in-and-out and remove instrument. Irrigate and repeat as necessary to working length in a gentle/passive in-and-out motion and then use TruNatomy® PRIME to working length to optimize the shape.
- 13) If TruNatomy® PRIME is loose at length with no dentinal debris in the apical flutes, continue shaping with TruNatomy® MEDIUM and subsequently with LARGE if needed. Care should be taken to guide each instrument gently to full working length to ensure completeness of shape. Use TruNatomy® LARGE in larger and straighter canals only, such as maxillary central incisors, some palatal or distal canals of molars.

- 14) Following proper cleaning and shaping, confirm the final working length measurement using the last TruNatomy® instrument by hand. Verifying the shape at working length and passive fit ensures that TruNatomy® Conform Fit® gutta-percha points matches the correct apical canal diameter. This is achieved by passively placing the last TruNatomy® instrument in the canal.
  - If the instrument reaches the working length passively, choose and try-in the matching TruNatomy® Conform Fit® gutta-percha points.
  - If the fit is not passive, instrument (as described above) with the last TruNatomy® shaping instrument after reconfirming working length, irrigating and confirming patency.
- 15) Once appropriateness of the shape is confirmed, proceed with 3-D disinfection protocols.
- 16) Use dedicated TruNatomy® paper points to dry the root canals and dedicated TruNatomy® Conform Fit® gutta-percha points to obturate, selecting either AH Plus® or AH Plus® Bioceramic Sealer.

## 8) HYGIENE, DISINFECTION, CLEANING AND STERILIZATION

- Please follow the General Processing Instructions For Endodontic Products (<https://www.dentsplysirona.com/en/service-contact/ifu.html>) for the disinfection, cleaning and sterilization steps.
- Do not re-use silicone stops. Remove and discard silicone stops after each use.
- Products shall be disposed according to local regulations for the safe disposal of sharp and contaminated devices.

## 9) ADDITIONAL INFORMATION

- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
- Sterility cannot be guaranteed if packaging is open, damaged or wet.
- To get a free printed copy of IFU please see section “Get a print copy of IFU” on website <https://www.dentsplysirona.com/en/service-contact/ifu.html>.
- Explanation of non-harmonized symbols for IFUs and labels, see IFU Symbols (<https://www.dentsplysirona.com/en/service-contact/ifu/ifu-symbols.html>).

| Authorised Representative                                                                  |                                                                                                                              | Manufacturer                                                                                                                                                                                                         |  |
|--------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| <br>2797 | <div>EC REP</div> <div>DENTSPLY DeTrey GmbH</div> <div>De-Trey-Straße 1</div> <div>D-78467 Konstanz</div> <div>Germany</div> |  <div>Maillefer Instruments Holding Sàrl</div> <div>Chemin du Verger 3</div> <div>CH-1338 Ballaigues</div> <div>Switzerland</div> |  |
|                                                                                            |                                                                                                                              |  <div><a href="https://www.dentsplysirona.com/ifu">www.dentsplysirona.com/ifu</a></div>                                           |  |
|                                                                                            |                                                                                                                              |                                                                                                                                                                                                                      |  |

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TruNatomy®

EN

FOR DENTAL USE ONLY  
FOR PROFESSIONAL USE ONLY  
SINGLE USE – REPROCESSING PROHIBITED  
STERILE – STERILIZED BY RADIATION

## INSTRUCTIONS FOR USE

### SINGLE USE STERILE ENDODONTIC INSTRUMENTS

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### 3) WARNINGS

- Strictly follow this Instruction for Use to minimize following risks to the TruNatomy® devices, the patient and/or the user:
  - Breakage of instruments.
  - Cross-contamination
  - Heat generation due to insufficient lubrication and irrigation.
  - Swallowing of working part of the instrument.
  - Toxic or allergic reactions caused by processing residues.

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- TruNatomy® instruments are single use devices and cannot be reused. Re-use can increase the risk of cross contamination or breakage.
- Disinfection, cleaning and sterilization is not applicable.
- Products shall be disposed according to local regulations for the safe disposal of sharp and contaminated devices.

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