Hemostatic sponge

Instructions for use

- en Instructions for use
- es Instrucciones de uso
- fr Mode d'emploi
- pt Instruções de uso

Instructions for use

en

Description:

This is a non-toxic, non-allergenic, non-immunogenic, non-pyrogenic and overall bio-compatible product. It's a surgical haemostatic sponge, manufactured from highly purified first extract grade gelatine material for use in various surgical procedures, where traditional haemostatic is difficult or impractical and use of other non-absorbable materials is undesirable. This product absorbs approximately 40-50 times its weight of water/blood and adheres easily to the bleeding site. The uniform porosity guarantees a favourable haemostatic. When implanted in vivo, it is completely absorbed within 3-4 weeks. Proclinic gelatine sponges have a porous structure which activates the thrombocytes at the moment blood comes in contact with the matrix of the sponge. This causes the thrombocytes to release a series of substances which promote their aggregation at the same time as their surfaces change character, thus enabling them to act as a catalyst for the formation of the fibrin.

Indications:

It can be effectively used in Dental Surgery. This product may be used in oral surgery, dry or saturated with thrombin solution to control bleeding in extractions and apicectomies especially in patients with a tendency to bleed extensively. To between gloved fingers to expel air bubbles and then replaced in saline until needed. The sponge should promptly return to its original size, with slight expansion in thickness and shape in the solution. If it does not, it should be removed again and kneaded vigorously until all air is expelled and it does expand to its original size, with slight increases in thickness and shape when returned to the sterile saline. If used wet, it may be blotted to dampness on gauze before application to the bleeding site. It should be held in place with moderate pressure, using a pledge of cotton or small gauze sponge until haemostasis results. Removal of the pledge or gauze is made easier by wetting it with a few drops of sterile saline, to prevent pulling up the sponge which by then should enclose a firm clot. Use of suction applied over the pledge of cotton or gauze to draw blood is unnecessary, as the sponge will draw up sufficient blood by capillary action. The first application of the product will usually control bleeding, but if not, additional applications may be made. For additional applications, fresh pieces should be used, prepared as described above. Use only the minimum amount of sponges, cut to appropriate size, necessary to produce haemostasis. It may be left in place at the bleeding site, when necessary. Since this product causes little more cellular reaction than does the blood clot, the wound may be closed over it. It may be left in place when applied to mucosal surfaces until it liquefies.

Achieve immediate haemostasis in:

- Dentoalveolar surgery and extraction of teeth in patients with risk of haemorrhage
- Regeneration of mucosal tissues removal of mandibular cysts
- Surgical removal of oral tumours and leucoplakias
- Surgical treatment of periapical and dental abscesses
- Maxillo-facial surgery

Directions for use:

Sterile technique should always be used to remove this product from its packaging. Cut to the desired size, a piece, either dry or saturated with sterile, isotonic sodium chloride solution (sterile saline), can be applied with pressure directly to the bleeding site. When applied dry, a single piece should be manually applied to the bleeding site, and held in place with moderate pressure until haemostasis results. When used with sterile saline, it should be first immersed in the solution and then withdrawn, squeezed.

Handling and storage:

This product is packed in sterile blisters with guaranteed absolute sterility. Once the package is opened, contents are subject to contamination. It is recommended that this product be used as soon as the package is opened and unused contents discarded. It should not be reused or re-sterilized as it may lose its physical characteristics as well as sterility.

Storage:

The product should be stored in its original packaging and the outer package must therefore be closed immediately after use. The box should be stored in a clean, dry room at a temperature not more than 30°C.

Precaution:

Extensive packing/application of this product, particularly within bony cavities, should be avoided, since recovering to its initial volume may interfere with normal tissue function and/or could cause possible or eventual compression necrosis of surrounding tissue. The sterility of the product cannot be guaranteed if the package bears evidence of damage or has been.

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OBELIS S.A.
Laboratorios Clarben S.A.
AEGIS LIFESCIENCES PVT. LTD.