
Composite Warmer

(en) Composite Warmer
Instructions for Use

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en ENGLISH**Safety**

CAUTION! Please read these instructions carefully before you set up and start the device. As with all technical devices, this unit will function properly and be safe to operate only if both the general safety regulations and the specific safety instructions laid out in these operating instructions are followed.

1. The device must be operated in accordance with the following instructions. Solventum accepts no liability for any damage arising from use of the device for any other purpose.
2. The device is for indoor use only.
3. Placement or operation of the device on non-level surfaces is not permissible. Set up the device on a dry, level, and non-slip surface.
4. Do not store materials that require cooling next to the device.
5. Do not cover the device to avoid damages to the cover and device electronics.
6. Do not place device in reachable distance to patient.
7. Check the power supply, base station, composite holder and silicone sleeve for damage at regular intervals. If a defect is found, stop operation of the device.
8. If there is any reason to believe that the safety of the device may be compromised, the device must be taken out of operation and labelled appropriately to ensure the device is not inadvertently used by another person. The safety may be compromised, for example, when the device fails to function properly or shows visible signs of damage.
9. Only service centers authorized by Solventum may open the device housing and repair the device.
10. Only adequately labelled original Solventum replacement parts may be used in the repair and maintenance of this device. Any damage resulting from the use of third-party replacement parts voids the liability of Solventum.
11. Keep solvents and sources of intense heat away from the device as they may damage the plastic housing of the base station.
12. When cleaning and disinfecting the device, ensure that no cleaning or disinfecting agent enters the base station since this may lead to electric short circuiting.
13. Clean device only when disconnected and cooled down to protect user from heat.
14. **Caution:** During operation, the cone shaped surface of the base station becomes hot (indicated with the warning symbol “Caution, hot surface”). To prevent burns, avoid direct contact with this area. The composite holder is covered with a silicone sleeve to ensure safe handling. Pay attention to the exposed areas (especially on the sides and on the bottom of the device), as they also become hot and mean a risk of burns. Therefore, the silicone sleeve is also indicated with the warning symbol “Caution, hot surface”.
15. Please report a serious incident occurring in relation to the device to Solventum and the local competent authority (EU) or local regulatory authority.

Product Description

Solventum™ Filtek™ Composite Warmer is a device for warming Filtek dental composites in capsules and flowable syringes indicated for warming up to 70 °C / 158 °F (e.g., Filtek™ Universal, Filtek™ One, Filtek™ Supreme XTE/Ultra, Filtek™ Easy Match, Filtek™ Supreme Flowable, Filtek™ Bulk Fill Flowable) to increase flowability. Filtek Composite Warmer has been specifically designed for use of Filtek composites. The device requires no maintenance.

 The product may only be used when the product labeling is clearly readable. These Instructions for Use should be kept for the duration of product use. For details on all additional products mentioned, please refer to the corresponding Instructions for Use. The instructions for use are available at eFU.Solventum.com.

Intended Purpose

Intended purpose: Warming device for Filtek composites in capsules and flowable syringes used in dentistry.

Intended users: Educated dental professionals who have theoretical and practical knowledge on usage of dental products.

Patient target group: All patients requiring dental treatment, unless the patient's condition limits the use.

Clinical benefit: Increase flowability of Filtek dental composites.

Indications

Warming of Filtek dental composites in capsules and flowable syringes indicated for warming; applicable for all patients requiring composite treatment.

Contraindications



None known.

Potential Undesirable Side Effects and Complications

None known.

Residual Risks

In situations with improper in-office infection control, or if the Instructions for Use are not followed, the following risks may occur: carpopedal spasm, cardiac arrhythmia, infection, laceration or other wounds, thermal burn.

Technical Data**Base station****Composite holder with silicone sleeve**

Operating voltage (power supply):	100 V – 240 V, 50/60 Hz
Power supply:	USB-C (5 V, 3 A)
Relative Humidity:	0% to 80%
Operating temperature:	70 °C max./ 158 °F max.
Ambient temperature:	5 °C to 40 °C/41 °F to 104 °F
Atmospheric pressure:	784 hPa – 1060 hPa
Preheating time of device:	10 min max. at room temperature 23 °C/73 °F
Heating time composite capsules:	2 min (to at least 50 °C/122 °F after preheating)
Heating time flowable syringes:	5 min (to at least 50 °C/122 °F after preheating)
Timer setting:	2 min
Dimensions:	Diameter: 76 mm Height: 76 mm
Weight:	190g (without power supply)
Year of manufacture:	see type plate

Transport and Storage Conditions

Ambient temperature:	– 20 °C to +60 °C/ – 4 °F to +140 °F
Relative humidity:	0% to 80%
Atmospheric pressure:	784 hPa – 1060 hPa

LED Signals

Status LED yellow:	Device is operating/heating
Status LED green:	Device has reached operating temperature
Timer LED blue:	Timer is running
Timer LED red:	Device malfunction

Installation

Option 1

- ▶ Connect USB-C plug to the wall adapter.
- ▶ Plug wall adapter to an adequate electrical socket.

Option 2

- ▶ Connect USB-C plug directly to an adequate USB-C port (e.g., treatment unit, computer, powerbank, etc.).
- ⚠ Do not use a USB-C extension cable or USB adapter, as this affects the performance of the device.

Operation

Switch on device/Preheating

- ▶ Place the composite holder on the heating cone of the base station.
- ▶ Make sure the composite holder is covered with its silicone sleeve to avoid contact with hot surfaces. Never use the device without a silicone sleeve. If sleeve has been disconnected e.g., for reprocessing, reattach the sleeve to the holder using the guiding lines on both prior to heating the device.

- ▶ Attach the silicone sleeve to the composite holder, letting it slightly protrude.

- ▶ Make sure the sleeve is attached correctly by checking the fit.
- ▶ Press operation button once to turn on device and start preheating. Yellow LED light indicates preheating of device. Required time for preheating: appr. 10 min. at room temperature 23 °C / 73 °F.
- ▶ The change to green LED indicates the device has come to its operating temperature. Composite can be placed now.
- ⚠ When placing a fresh or cooled down holder on the warm base station it will take appr. 25 sec until the status LED switches to yellow indicating heating/reheating.

Inserting Filtek composite capsules into composite holder

- ▶ After required preheating time (indicated by green status LED), insert capsules in the outlets provided.
- ▶ Ensure the capsules are securely positioned with the unopened cap facing downwards and outside. Insert closed capsules only. In case of contamination with material, saliva or blood follow instructions for reprocessing the device.
- ⚠ Each capsule outlet is marked with its own symbol on the top of the silicone sleeve to better distinguish different capsules (e.g., shades).

Timer setting/Capsule heating

- ▶ Press operating button shortly to start timer. Blue LED light confirms that the timer is set. The blue LED disappears after the required capsule heating time of 2 min.
- ▶ If capsules are warmed without use of the timer function, make sure to keep the required heating time (2 min).
- ⚠ Please note that the timer can only be set when the device has reached its operating temperature indicated by the green status LED.

Inserting Filtek flowable composite syringe into composite holder

- ▶ After required preheating time (indicated by green status LED) insert flowable syringe carefully in provided outlet.

- ▶ Ensure the syringe with attached dispensing tip is securely positioned with tip facing downwards. Make sure not to bend the tip.
- ▶ Only insert syringes with attached dispensing tip.
- ⚠ In case of contamination with material, saliva or blood, please follow instructions for reprocessing the device.

Flowable syringe heating

- ▶ Allow the flowable syringe to heat for at least 5 min.
 - 6 capsules and 1 flowable syringe can be warmed simultaneously. Please note that the syringe requires a longer heating time to reach the required temperature (in total 5 min for the syringe, 2 min for the capsules).
- ⚠ It is recommended to warm only the amount of composite required for immediate use.

Removing capsules/flowable syringe

- ▶ Remove the capsules/syringe from the holder after the respective required heating time and apply the material immediately.
- ▶ Remove the flowable syringe by the syringe shaft or finger plate. Do not remove the flowable syringe by pulling the plunger as this may lead to air inclusions in the composite material.

Optionally the composite holder can be fully removed from the base station and brought to the point of use/patient tray where the material is taken out and applied. Once the composite holder is removed from the base station, it preserves heat (at least 50 °C / 122 °F) for up to 6 min depending on the ambient conditions. Longer procedures require reheating of material. Remove composite from the device if not immediately used.

Power off

To switch off the device, keep operation button pressed until all LED lights disappear.

Errors

Error	Possible cause ▶ Remedy
Status LED does not appear after power button is pushed. Device is not heating.	Power failure. ▶ Check both the power cable and the power source.
	Defect of device electronics. ▶ Have device repaired/ exchanged.
Status LED does not appear after power button is pushed. Device is heating.	Defect of device electronics. ▶ Have device repaired/ exchanged.
Status LED does not turn green after required preheating time.	Defect of device electronics. ▶ Have device repaired/ exchanged.
Timer LED does not flash after button is pushed while status LED shows green light.	Defect of device electronics. ▶ Have device repaired/ exchanged.
	Operating has been pressed too long. ▶ Press operating button again shortly.
Timer LED does not disappear after 2 min capsule heating time.	Defect of device electronics. ▶ Have device repaired/ exchanged.
Red LED flashes.	Device malfunction. ▶ Have device repaired/ exchanged.
Device is not heating (enough)/ device gets too hot.	Device malfunction. ▶ Have device repaired/ exchanged.

Notes

- Filtek Composite Warmer is a multiple patient, multiple use device.
- Always check Filtek instructions for use for warming indications, heating times and instructions.
- Do not use this device for dental materials from a category other than the above listed Filtek composites.
- Filtek composite warmer is intended for brief warming of dental composites prior to application. Longer warming beyond the treatment period may damage the material and should therefore be avoided.
- Remove composites from the device if not immediately used.
- Discard composites if maximum heating time is exceeded. Consult Filtek instructions for use for more information.

Cleaning and disinfection of base station

Warnings

Clean plastic housing and heating cone of the base station after it has cooled down and while it is disconnected from power source. Do not immerse. Cleaning and disinfecting solutions must not enter the device. Do not use solvents or scouring agents, as they may damage the device components.

Make sure that the chosen cleaning and disinfectant agents do not contain any of the following materials:

- Oxidation agents (e.g., hydrogen peroxide)
- Oils
- Glutaraldehydes

Use only validated methods for cleaning and disinfection.

Use solely cleaning and disinfectant agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). CaviWipes™ have been tested for use as a cleaning and disinfectant agent. Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.

General Basics

The following procedure should be applied after use or before any subsequent use. The device must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

Manual cleaning by wiping

- ▶ Clean the device after or before every application with ready-to-use disinfectant wipes (e.g., CaviWipes™ based on the active substances alcohols and quaternary compounds) or equivalent cleaning agents. Follow the cleaning agent manufacturer's directions for use.
- ▶ Thoroughly clean the surface of the device with two ready-to-use CaviWipes™ for at least 30 seconds or until no more dirt is visible.
- ▶ Special care must be exercised in cleaning of the gaps and transitions.

Manual disinfection by wiping (Intermediate level disinfection)

- ▶ Disinfect the device with ready-to-use disinfectant wipes (e.g., CaviWipes™ based on the active substances alcohols and quaternary compounds) or an equivalent disinfection agent. Follow the disinfectant agent manufacturer's directions for use, observing in particular the contact times.
- ▶ Disinfect the surface of the device for three minutes using at least two new disinfectant wipes. Keep all surfaces moist during the three-minute period.
- ▶ Special care must be exercised in disinfection of the gaps and transitions.

Inspection, Maintenance, Testing

- ▶ Examine the device for damage, discoloration, and contamination before every use.
- ▶ Do not use damaged devices under any circumstances. If the device is not visibly clean, repeat the cleaning and disinfection procedure.

Cleaning and disinfection of Composite holder and Silicone Sleeve

Warnings

For reprocessing, the silicone sleeve must be disconnected from the composite holder. Only clean when cooled down to avoid touching hot surfaces.

Observe the following requirements when selecting the cleaning and disinfectants agents used:

- The following agents should not be included: oxidizing agents (e.g. hydrogen peroxides), solvents and oils.
- The pH value of the cleaning solution should be between pH 7.5 and 11
- The concentrations, temperatures and exposure times specified by the manufacturer of the cleaning agents and disinfectants as well as the instructions for rinsing must be strictly followed.

Use solely solutions that have been freshly prepared with water that has following microbiological quality. The water should be free of facultative pathogenic microorganisms and should not exceed 100 CFU/ml of microorganisms.

- Use only a soft, clean and lint-free cloth for drying.
- Never use metal brushes or steel wool for cleaning. Always use a clean, soft brush for cleaning.
- Only use validated methods for cleaning and disinfection.
- Only use cleaning agents and disinfectant agents that have been tested for effectiveness and compatibility with the medical device used (e.g. VAH approval, FDA clearance).
- Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.

Limitations on reprocessing

The composite holder and silicone sleeve can be reprocessed without restriction if handled with care, provided the following instructions are followed and no damage is visible. The composite holder and silicone sleeve must be carefully checked for damage before each reuse.

General basics

The medical device is supplied non-sterile and must be reprocessed in accordance with these instructions before first use and before each subsequent use.

Preparation for decontamination

Reprocess the medical device as quickly as possible after use.

Reprocessing

The reprocessing of the composite holder and silicone sleeve include the steps of pre-treatment, cleaning, and disinfection. After these steps, reprocessing is complete. **Optionally** a sterilization can be performed after pre-treatment, cleaning, and disinfection.

Cleaning and disinfection should, if possible, always be done automatically (washer-disinfector). If automatic reprocessing is not possible, manual cleaning and disinfection can be carried out in an immersion bath. The sole use of manual methods is not recommended as they are less effective and reproducible.

The reprocessing procedure must be carried out before the first use and after each subsequent use.

Pre-treatment

Pre-treatment must be carried out before any manual or automated cleaning and disinfection. Gross contaminations must be removed from the product immediately after use (within a maximum of 2 h).

- ▶ Remove all visible impurities manually under running water (tap water quality, temperature 22 °C – 44 °C/ 71 °F – 111 °F).
- ▶ Thoroughly clean the recesses in the composite holder under running water for at least 1 minute using a clean, soft and dry disposable brush.

- ▶ Rinse the surface of the composite holder and the silicone sleeve under running water for at least 1 minute using a clean, soft cloth.

Check the composite holder and silicone sleeve for visible contaminations and repeat the pre-cleaning process if necessary.

Automated cleaning and thermal disinfection

Automated washer-disinfectors that meet the following requirements are suitable for use:

- The device has been approved for effectiveness (e.g. DGHM or FDA approval/clearance/registration or CE labelling in accordance with ISO 15883).
- If available, use a tested program for thermal disinfection (A0 value > 3000 or at least 5 min at 90 °C / 194 °F or A0 value > 600 or for

older devices at least 1 min at 90 °C / 194 °F). The program used is suitable for composite holders and silicone sleeves and contains sufficient rinsing cycles (in the case of chemical disinfection, consider the risk of disinfectant residues on the instruments).

- The water used for rinsing is demineralized water. The microbiological quality of the water used should be free of facultative pathogenic microorganisms and should not exceed 100 CFU/ml of microorganisms.
- The air used for drying is filtered (oil-free, low in microbiological and particulate contamination).
- The washer-disinfector is regularly maintained and checked.

Processing in a washer-disinfector

- After pre-treatment of the composite holder and silicone sleeve, place the medical devices with the inner surface facing downwards in the wire or sieve basket of the washer-disinfector. Ensure that the instruments do not touch each other in the washer-disinfector.
- For automated cleaning, use an alkaline cleaning agent (e.g. Neodisher® MediClean dental, Weigert) at a concentration of 0.2% (see product data sheet) or an equivalent agent.
- Rinse with demineralized water.
- Start the cleaning program (Vario-TD from Miele G7836 CD or equivalent program).

Phase	Time	Temperature
Precleaning	5 min	45 °C
Main Cleaning	5 min	50 °C
Rinse	30 min	60 °C
Disinfection	5 min	90 °C
Drying	7 min	90 °C

- Remove the composite holder and silicone sleeve from the washer-disinfector at the end of the program.
- Check the composite holder and silicone sleeve for visible impurities. If any impurities are visible, repeat the cycle in the washer-disinfector.
- Check the composite holder and silicone sleeve again and store at room temperature, protected from contamination and dry after completion of the thermal disinfection.

Manual cleaning and disinfection

Only in case that an automated washer-disinfector is not available manual cleaning can be performed with ultrasonic bath followed by the manual disinfection.

Manual cleaning with ultrasonic bath

- Place the composite holder and silicone sleeve in the ultrasonic cleaning cassette after pre-treatment. Ensure the composite holder and silicone sleeve or instruments do not touch each other. The medical device must be completely immersed in the solution
- Fill the ultrasonic bath with the orthophthalaldehyde-based cleaning agent Cidezyme® in a 0.8% solution (8 ml/l) or an equivalent agent. The water used should be at least drinking water quality and between 25 – 35 °C.
- Start the cleaning program.
- Set the ultrasonic treatment to 1 min.
- Remove the composite holder and silicone cover from the ultrasonic bath and rinse thoroughly under running water for at least 5 minutes. Rinse all recesses with water, using a clean, soft brush or cloth if necessary, to remove all residues of Cidezyme cleaner.
- Check the composite holder and silicone sleeve for visible impurities. If contamination is still present, repeat the cycle in the ultrasonic bath.
- Check the composite holder and silicone sleeve again and continue with the disinfection in a disinfection bath.

Manual disinfection by immersion

- Fill the immersion bath with undiluted or ready-to-use Cidex® OPA disinfectant or an equivalent agent.
- Place the cleaned and checked composite holder and the silicone sleeve in the disinfection bath at 20 °C for the specified exposure time of 5 minutes. Make sure that the medical devices are sufficiently covered with liquid and that the instruments do not touch each other in the immersion bath.
- After the exposure time, remove the composite holder and silicone sleeve from the disinfectant bath and rinse thoroughly with water at least 5 times (for at least 1 min). Rinse all recesses with water, if necessary using a clean, soft brush or cloth to remove all disinfectant residues.
- Blow off/blow out the composite holder and silicone sleeve, dry with filtered compressed air or leave to dry in a clean place at room temperature.
- Store the composite holder and silicone sleeve in a clean place as soon as possible after removal, if necessary after additional drying.

Sterilization (Optional Step)

Optionally a sterilization can be performed after pre-treatment, cleaning, and disinfection.

Packaging for Sterilization

Package composite holder and silicone sleeve in single wrapped sterilization pouches (one-way packaging), meeting following requirements:

- Each part is packaged separately
- Packaging meets ANSI AAMI ISO 11607 (for USA: FDA-Clearance)
- Suitable for steam sterilization (temperature resistance up to at least 142 °C / 288 °F, sufficient steam permeability)
- Sufficient protection of the composite holder and the silicone sleeve or sterilization pouches against mechanical damage

Sterilization

A fractionated vacuum procedure (at least 3 vacuum stages) or a gravity procedure (with adequate product drying) must be used for the sterilization¹. The use of the less effective gravity procedure is permissible solely if the fractionated vacuum procedure is unavailable.

The user is solely responsible for selection of the parameters required for the gravity procedure; the user must therefore calculate or validate these parameters in accordance with the loading configuration. The drying time actually required is also directly dependent on parameters that are in the sole responsibility of the user (loading configuration and density, sterilizer condition, etc.) and must consequently be determined by the user.

- Steam sterilizer in conformity with EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA-Clearance)
- Validated according to ISO 17665 (valid IQ and OQ as well as and product-specific performance assessment (PQ))
- Maximum sterilization temperature 135 °C (287 °F; including tolerance according to ISO 17665)

Package composite holder and silicone sleeve separately in single wrapped sterilization pouches (e.g., Steriking, Wipak) and treat by steam sterilization in autoclave at following parameters:

Overview of the Sterilization Procedures Used:

Method	Exposure time at 132 °C (270 °F) Validation at 132 °C / 270 °F	Exposure time at 134 °C (273 °F) Validation at 132 °C / 270 °F	Drying time
Fractional prevacuum Product: composite holder and silicone sleeve, wrapped	5 min ²	5 min	20 min

Method	Sterilization time at 121 °C (250 °F) Validation at 121°C	Sterilization time at 134 °C (273 °F) Validation at 121°C	Drying time
Gravity displacement¹ Product: composite holder and silicone sleeve, wrapped	20 min	20 min	20 min

¹ The gravity procedure is not applicable for processing in the European Union.

² The sterilization temperature of 132 °C (270 °F) is not applicable for processing in the European Union

- Immediate-use steam sterilization is generally not permitted.
- Sterilization of unwrapped instruments is generally not permitted.
- Do not use hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization or plasma sterilization.

Storage

- After reprocessing, the composite holder and silicone cover must be stored in a dry place, protected from light and dust-free.

Control, maintenance, inspection

- Check the medical device for damage, discoloration and contamination before each use.
- Do not use damaged devices under any circumstances.
- There are symbols on the top of the base station to ensure safe and proper use. The rating plate of the device is located on the underside of the base station. There is also a UDI label on the silicone sleeve. Before each use, check that the symbols, the rating plate, and the UDI labelling are intact and clearly legible.
- If the medical device is not visibly clean, repeat the reprocessing procedure.

Additional Information

The information above is intended to enable the re-processing of the **Filtek Composite Warmer** and complies to the standard ISO 17664-2. The national guidelines also apply.

In addition, all applicable national legal regulations and regulations on hygiene relating to your practice or the hospital apply.

Return of Old Electric and Electronic Equipment for Disposal Collection

Users of electric and electronic equipment are required to collect their old equipment separately from other waste in accordance with the regulations of the specific country. Old electric and electronic equipment must not be disposed of with unsorted household waste. This separate collection is a prerequisite for recycling and reprocessing as an important method for preserving environmental resources. Devices should be disposed of in a safe manner and in compliance with local, state/province, and federal law.

Return and Collection Systems

When your Filtek Composite Warmer is no longer usable, do not dispose the device with household waste.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Symbol Glossary

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer.

Reference Number and Symbol Title	Symbol	Description of Symbol
Swiss Authorized Representative		Indicates the authorised representative in Switzerland.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.1.7 Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.3.8 Humidity limitation		Indicates the range of humidity to which the medical device can be safely exposed.
ISO 15223-1 5.3.9 Atmospheric pressure limitation		Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
ISO 15223-1 5.4.3 Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use.
CE Mark		Indicates conformity to all applicable European Union Medical Device Regulations and Directives.
ISO 15223-1 5.7.7 Medical Device		Indicates the item is a medical device.
ISO 15223-1 5.7.10 Unique device identifier		Indicates a carrier that contains unique device identifier information.
Caution, hot surface		To indicate product can be hot and should not be touched without taking care.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)

Reference Number and Symbol Title	Symbol	Description of Symbol
Non-corrugated fiberboard	21 PAP	Indicates product packaging is made of non-corrugated fiberboard. Official Journal of the EC; Commission Decision (97/129/EC)
ISO 15223-1 5.1.8 Importer		Indicates the entity importing the medical device into the locale.
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle. Directive 2012/19 EC on waste electrical and electronic equipment (WEEE).

For more information, see eIFU.Solventum.com.

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© Solventum 2025. Solventum, the S logo and Filtek are trademarks of Solventum or its affiliates.
Solventum Germany GmbH
U.S.A. 1-800-634-2249.