Adhese® 2

[en] Instructions for Use

Dual-curing, single-component dental adhesive

Rx ONLY

CE 0123

Date information prepared:
2025-05-28 / WW / Rev. 1
web

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

ivoclar

English

1 Intended use

Intended purpose

Bonding of dental restorations to dental hard tissue, conditioning of Tetric CAD restorations, sealing of dentin lesions.

Patient target group

Patients with permanent teeth

Intended users / Special training

- Dentists
- No special training required

Application

For dental use only.

Description

Adhese® 2 is a dual-curing single-component dental adhesive for enamel and dentin that is compatible with all etching techniques (self-etch, selective-enamel-etch and etch & rinse techniques). The areas of application include adhesive cementation of indirect restorations and direct restorative procedures.

Adhese 2 can either be light-cured for 10 s at a light intensity of 500 to 1,300 mW/cm² or for 5 s at a light intensity of 1,800 to 2,200 mW/cm². In the restoration of Class I and II cavities, Adhese 2 can be light-cured in only 3 s from the occlusal aspect at a light intensity of 2,700 to 3,300 mW/cm².

Adhese 2 is available in the VivaPen®, in bottles or in Free Stand® Single Dose units.

Adhese 2 in the VivaPen

- The VivaPen DC (dual-curing) cannula (blue) is coated with the co-initiator required for the self-curing reaction. By activating the VivaPen push button, the flocked cannula tip is saturated with adhesive liquid and mixes with the co-initiator.
- The VivaPen LC cannula (black) does not contain any co-initiator and is only used in combination with light-curing procedures.

Adhese 2 in the bottle

- As applicators for the bottle, commercial applicators or special Adhese 2 DC applicators are available. The Adhese 2 DC applicator is coated with the co-initiator required for the self-curing reaction. When using a DC (dual-curing) procedure, make sure that you are using the correct applicator (Adhese 2 DC applicator). When the applicator is dipped into the Adhese 2 liquid, the co-initiator mixes with the adhesive. When using a LC (light-curing) procedure, an applicator without co-initiator can be used.
- The applicators are available in two different versions:
Regular: for cavities, crown preparations, etc.
Small/Endo: for root canals and micro-preparations

Adhese 2 in Free Stand Single Dose units

The Free Stand Single Dose delivery form already contains the required applicator coated with co-initiator. The Single Dose units are available in two versions:

- Regular: for cavities, crown preparations, etc.
- Small/Endo: for root canals and micro-preparations

Indications

Missing tooth structure in anterior and posterior teeth, partial edentulism in the anterior and posterior region

Areas of application:

- Pure self-curing of Adhese 2 only in combination with VivaPen DC cannulas and DC applicators**
 - Adhesive cementation of indirect restorations and endodontic posts using Variolink® Esthetic DC in combination with all etching techniques.
 - Adhesive cementation of indirect restorations and endodontic posts using dual-curing luting composites in combination with the etch & rinse technique.
- Light-curing of Adhese 2 at 500 to 1300 mW/cm² (10 seconds) or at 1800 to 2200 mW/cm² (5 seconds)**
 - Adhesive cementation of indirect restorations with light- and dual-curing luting composites
 - Directly placed build-up fillings and core build-ups made with light- and dual-curing composites.
 - Directly placed light-curing composite restorations.
 - Repair of fractured composite restorations.
 - Desensitization of hypersensitive cervical areas (Note: Provide additional cooling with an air syringe when applying the product.)
 - Sealing of prepared tooth surfaces before temporary/permanent cementation of indirect restorations
- Light-curing of Adhese 2 using Bluephase® PowerCure in the 3sCure curing mode**
 - Only restorations in the posterior region of permanent dentition (Class I and II, including the replacement of individual cusps) when light-cured from the occlusal aspect.

Contraindication

- The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.
- Direct pulp capping

Limitations of use

- Do not use Adhese 2 as a primer for ceramic restorative materials. In such cases, suitable ceramic primers must be used.
- Adhese 2 must not be used if the stipulated working technique cannot be employed.

Processing restrictions

- ⊗ Do not reuse. The Free Stand Single Dose units, cannulas, applicators and sleeves are not intended to be reprocessed or reused.

Side effects

In rare cases, components of Adhese 2 may lead to sensitization. In such cases, Adhese 2 should no longer be used.

Interactions

- Peroxide-based bleaching agents and disinfectants may inhibit the polymerization of adhesives. Such materials must not be used ≤ 14 days prior to the application of Adhese 2.
- Phenolic substances (e.g. eugenol) may inhibit the polymerization of Adhese 2.
- Before any treatment with Adhese 2, completely rinse off all iron- and aluminium-based styptic agents with water.
- Fluoridation agents must not be used prior to the application of Adhese 2.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

HEMA, tricyclic urethane dimethacrylate, MDP, ethanol, water, D3MA, methacrylate-modified polyacrylic acid, silicon dioxide, KOH, camphor-quinone, ethyl p-dimethyl aminobenzoate, 2-dimethyl aminoethyl methacrylate

2 Application

2.1 Procedure for indirect and direct restorations

2.1.1 Isolation

Adequate relative or absolute isolation is required.

2.1.2 Pulp protection / cavity liner

The 3sCure curing mode must not be used for restorations in areas close to the pulp. In very deep cavities, areas close to the pulp must be selectively coated with a calcium hydroxide liner and subsequently

covered with a pressure-resistant cement (e.g. a glass ionomer cement).

2.1.3 Conditioning with phosphoric acid gel

The bond to enamel can be improved by selectively etching the enamel or by applying the "etch & rinse" technique. Unprepared enamel surfaces must always be conditioned with phosphoric acid gel.

 Please observe the instructions for use for the phosphoric acid gel.

a. Selective enamel etching

Apply phosphoric acid gel onto the enamel and allow it to react for 15–30 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with oil- and water-free compressed air until the etched enamel surfaces appear chalky white.

b. Etch & rinse technique

Apply phosphoric acid gel onto the prepared enamel first, and then onto the dentin. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with oil- and water-free compressed air until the etched enamel surfaces appear chalky white.

2.1.4 Handling of the VivaPen, bottle and Free Stand Single Dose unit

a. VivaPen

- Remove the cap of the VivaPen by turning it beyond the point where resistance is felt and simultaneously pulling it off.
- Remove the seal before using the DC cannula.
- Attach the required VivaPen cannula (LC or DC) by snapping it into place. The attached cannula can be rotated to the desired position up to the point where resistance is felt. Caution: Rotating the cannula beyond the point of resistance will detach the cannula!
- Make sure to use the VivaPen cannula (LC = black, DC = blue) required for the respective procedure.
- For hygienic reasons, always use a VivaPen protective sleeve by sliding the Pen with the cannula into the protective sleeve until the cannula punctures through the sleeve.
- Remove the paper separation strip from the protective sleeve before using the VivaPen.
- The use of the protective sleeve also minimizes the risk of a cannula coming loose.
- The cannula tube can be bent up to 90°.
Caution: Bending the cannula several times or bending it at the cannula hub may break the cannula. Excessive bending may result in blockage of the cannula tube.
-  Do not reuse. A new cannula and protective sleeve must be used for each patient.
- Press the push button firmly down and wait until the flocked cannula tip is saturated with Adhese 2. Repeat if required. The flocked cannula tip is saturated when the yellow colour of Adhese 2 is visible.
- Remove the protective sleeve and cannula after having used the VivaPen. The VivaPen cannula is removed by turning it beyond the point of resistance and simultaneously pulling it off.
- Close the VivaPen with the original cap after use.

Application notes

- When using the VivaPen for the first time, several clicks may be required to wet the cannula tip due to the engineering of the VivaPen. This does not influence the indicated fill-level or the function of the VivaPen.
- When you use the black LC cannula in a DC procedure, the adhesive will not self-cure. This means that no reliable bond between the restoration and the tooth structure will be achieved.
- In the case of narrow and long root canal preparations, please observe the length of the DC cannula.
- The VivaPen DC (dual-curing) cannula is intended for the treatment of one cavity only. For further uses, a new DC cannula must be used.
- The coated flocked DC cannula tip must not come in contact with the VivaPen opening, as there is a risk of polymerization of the adhesive in the VivaPen.
- Following activation of the VivaPen with the DC cannula in place, the adhesive should be applied within 120 seconds.
- Protect the activated DC applicator/DC cannula from light until the adhesive is applied.
- Do not immerse the VivaPen in disinfectant solutions or autoclave it.
- Avoid any contact of the reusable VivaPen not covered with a protective sleeve with the mouth of the patient.
- Dispose of the VivaPen in case of contamination.

- If the cannula has been contaminated (e.g. the cannula has dropped on the floor), it must be replaced by a new one.
- Cannulas of other manufacturers are not compatible with the VivaPen and therefore must not be used.
- The white indicator of the transparent fill-level gauge shows the remaining amount of adhesive in the VivaPen. If the white indicator is no longer displayed, the VivaPen is empty.

b. Bottle

- Dispense the desired amount of Adhese 2 (max. 1 drop when using a DC procedure) into a mixing well and apply it with a disposable applicator (LC or DC applicator).
- When using a DC procedure, mix the adhesive with the co-initiator by rotating the DC applicator for approx. 5 seconds. Once the adhesive and co-initiator are thoroughly mixed, the flocked applicator tip will have a uniformly yellow appearance. If not, repeat the mixing procedure.
- Protect the adhesive from light.
- Carefully close the bottle immediately after each use.
-  Do not reuse. A new disposable applicator must be used for each application.

Application notes

- If an applicator not coated with co-initiator is used in a DC procedure, the adhesive will not self-cure. This means that no reliable bond between the restoration and the tooth structure will be achieved.
- In the case of narrow and long root canal preparations, please observe the length of the DC Endo applicator.
- The coated DC applicator must not come in contact with the dropper opening of the bottle, as there is a risk of polymerization of the adhesive in the bottle.
- The DC applicator, which is coated with co-initiator, is intended for single use only. For further uses, a new DC applicator must be used.
- Do not use the bottle intraorally. For hygienic reasons, suitable application aids should be used.
- Following activation (mixing of the adhesive with the co-initiator), the adhesive can be applied for approx. 120 seconds.
- Protect the activated DC applicator/DC cannula from light until the adhesive is applied.

c. Activation of the Free Stand Single Dose unit

- Remove the seal.
- Take the blue DC applicator from its holder and dip it into the container with the adhesive.
- Then mix the adhesive with the co-initiator by rotating the DC applicator for approx. 5 seconds. Once the adhesive and co-initiator are thoroughly mixed, the flocked applicator tip will have a uniformly yellow appearance. If not, repeat the mixing procedure.
- Do not dip the DC applicator into the adhesive more than twice.
-  Do not reuse. A new Free Stand Single Dose unit must be used for each application.

Application notes

- Always use the DC applicator (blue) included in the delivery form.
- If the DC applicator included in the delivery form is not used in DC procedures, the adhesive will not cure. This means that no reliable bond between the restoration and the tooth structure will be achieved.
- In the case of narrow and long root canal preparations, please observe the length of the DC Endo applicator.
- Following activation, the adhesive can be applied for approx. 120 seconds.

2.1.5 Application of the adhesive

- Starting with the enamel, completely coat the tooth surfaces to be treated with Adhese 2.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Adhese 2 with oil- and water-free compressed air until a glossy, immobile film layer results.

Important information:

- Avoid pooling, since this can compromise the fitting accuracy of the permanent restoration.
- Exposure to intensive light should be avoided during application.

- Prevent any contamination of the applicator/cannula with blood, saliva, sulcus fluid or water during application. In the event of contamination, the cavity must be rinsed again and the procedure must be repeated with a new applicator / cannula.

2.1.6 Light-curing the adhesive

Self-curing of Adhese 2

Under the following conditions, Adhese 2 can be used in the purely self-curing mode in combination with suitable dual-curing core build-up and luting composites.

Application / Indication	Curing mode of the luting or core build-up composite	Etching technique	
		Etch & rinse	Self-etch Selective enamel etching
Adhesive cementation of indirect restorations and endodontic posts	Dual-curing	✓	✓
	Self-curing	✓	Only in combination with Variolink Esthetic DC

Light-curing of Adhese 2

Under the following conditions, Adhese 2 must be light-cured:

Application / Indication	Curing mode of the core build-up or restorative composite material	Etching technique	
		Etch & rinse	Self-etch Selective enamel etching
Direct restorations, indirect translucent single tooth restorations and core build-ups	Light- or dual-curing	✓	✓
Sealing, desensitization, immediate dentin sealing	–	✓	✓

Adhese 2 can either be light-cured for 10 s at a light intensity of 500 to 1,300 mW/cm² or for 5 s at a light intensity of 1,800 to 2,200 mW/cm². In the fabrication of Class I and II restorations, Adhese 2 can be light-cured in only 3 s from the occlusal aspect at a light intensity of 2,700 to 3,300 mW/cm².

If the diameter of the light guide tip is smaller than the restoration diameter, perform an overlapping polymerization in order to ensure that all areas of the restoration are covered.

 The instructions for use of the curing light must be observed.

Exposure times and light intensities are listed in the table below.

Light intensity mW/cm ²	Exposure time
500 – 900	10 s
1,000 – 1,300	10 s
1,800 – 2,200	5 s
2,700 – 3,300	3 s

Safety note

- Avoid direct exposure of the gingiva, mucous membrane or skin to the light emitted by the curing light.
- The 3sCure mode must not be used in the case of caries profunda and very deep cavities.

2.1.7 Application of the restorative or luting composite

-  Proceed according to the instructions for use of the restorative and/or luting composite used.

2.2 Procedure for intraoral repairs of fractured composite fillings/ composite restorations with light-curing materials

- Roughen the surface to be repaired using diamond finishers and thoroughly clean with water spray and dry with oil- and water-free compressed air.
- Proceed according to the direct application for Adhese 2.

2.3 Procedure for the conditioning of Tetric® CAD restorations

-  Proceed according to the instructions for use for Tetric CAD.

3 Safety information

- Adhese 2 is irritant. Avoid contact with the skin, mucous membrane and eyes.
- After skin contact, immediately rinse with copious amounts of water.
- If the material should come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician/ ophthalmologist.
- In rare cases, contact with the skin may lead to sensitization towards the ingredients.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylate.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 76152082ABOND010HK

Warnings

Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).

Warnings for phosphoric acid gel:

Phosphoric acid gel is corrosive. Avoid contact with the eyes, mucous membrane and skin (protective goggles are recommended for both patients and operators). If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician / ophthalmologist.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Failure of the adhesive bond
- Postoperative sensitivity

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- If the indicated storage temperature is exceeded, proper functioning of the product cannot be ensured.
- Do not use the product after the indicated date of expiration.
- Expiration date: See note on packages.
- Keep the foil bag of the VivaPen DC cannulas closed until use. After opening, store cannulas in the foil bag or plastic box with closed lid, and use within six months.
- Store DC applicators exclusively in the plastic box provided. Keep the plastic box closed.

Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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