D. Devices

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en_ Biological incubator Instructions for use

DECLARATION OF CONFORMITY

This declaration is issued to confirm that the device complies with the requirements of the following directives, in accordance with the ISO/IEC 17050-1 standard:

- Directive 2014/30/EU on Electromagnetic Compatibility (EMC)
- Directive 2014/35/EU (Low Voltage Directive) concerning electrical equipment designed for use within certain voltage limits

WARRANTY TERMS

1. Duration:
The product is covered by a 12-month warranty.
2. Exclusions:
The warranty does not cover:
 - Errors resulting from lack of scheduled maintenance or due to negligence or misuse by the user.

- Periodic checks and maintenance.
 - Repair or replacement of parts subject to wear and tear, unless proven defective at the time of delivery.
 - Defects caused by personnel or during transport.
 - Defects or damage caused by improper use or operation inconsistent with the instruction manual, or inconsistent with the intended use.
 - Defects or damage caused by water or air contamination, or by random events such as fire, flood, lightning strike, power surges, etc.
 - Defects or damage caused by the use of detergents, disinfectants, or sterilization methods not clearly described in the operation and maintenance manual.
 - Natural discoloration of plastic components.
3. Limitations:
The warranty entitles the owner to the repair or replacement of defective components free of charge.

The right to replace the entire device is excluded where components integrated into CHIKANG MEDICAL products come from other companies and are covered by their own warranty certificates.

In such cases, the conditions and limitations stated in those certificates shall apply.

Service repairs do not extend or renew the original warranty period.

4. Claims:

Submitting a complaint about the quality of the delivered equipment does not entitle the buyer to delay or withhold payments.

5. Disputes:

In the event of any dispute concerning the warranty or its interpretation, the competent court for the registered office of MMG shall have jurisdiction, regardless of the place of sale.

6. Deviations:

Any deviations from the warranty conditions outlined in this section shall not entitle the buyer to additional rights and will be considered only in exceptional cases.

7. Other Provisions:

In all matters not regulated above, the applicable provisions of the Civil Code shall apply.

GENERAL SAFETY RULES

1. Ensure that the device is connected to the power supply using the correct voltage, as indicated on the nameplate.
2. Do not remove the nameplate.
3. Use dry cloths to clean the device.
4. Before carrying out any repair, disconnect the power cord from the socket.
5. Only use original spare parts.

DEVICE DESCRIPTION

Is a device designed for incubating biological indicators at 56°C. The incubator maintains a constant temperature throughout the incubation period, following a startup time of approximately 3 to 5 minutes. The incubation duration depends on the type of ampoule tests used.

Once the test is initiated by pressing the START button, the device keeps the temperature stable until it is turned off.

The incubator (see Fig. 1) includes:

- A metal frame
- A plastic housing (Fig. 1, item B)
- An aluminum heater (Fig. 1, item C)
- Three compartments for biological tests
- An electronic board for automatic temperature regulation

At the back of the device, there is a socket for connecting the supplied power adapter (Fig. 1, item D).

On the front of the device:

- A START button (Fig. 1, item E)
- A warning LED indicator labeled STAND BY/ALARM (Fig. 1, item A)

On the top of the device:


A transparent insulation lid (Fig. 1, item F)

DEVICE COMPONENTS (Fig. 1)

Item Descriptions:

- A – STAND BY/ALARM warning LED
- B – Plastic housing
- C – Heating element with three test slots
- D – Power socket
- E – START button
- F – Insulation lid

Fig. 1


BIOLOGICAL INDICATOR (Fig. 2)

The biological indicator (Fig. 2) is a thermoplastic container that holds a paper strip inoculated with bacterial spores and a small glass capsule containing the culture medium.

Components of the Biological Indicator:

- G – Nut
- H – Filter
- I – Label
- J – Strip with TYPE-1 chemical indicator
- K – Glass capsule
- L – Culture medium
- M – Paper strip with bacterial spores

Fig. 2

TECHNICAL CHARACTERISTICS

- Operating Temperature Range: 17°C – 27°C
- Device Dimensions (W × H × D): 130 × 85 × 90 mm
- Net Weight: 0.8 kg
- Frame Material: Fe360B steel
- Housing Material: ABS plastic
- Heater: Anodized aluminum block
- Package Dimensions: 220 × 170 × 110 mm
- Supply Voltage: 10–14 V
- Number of Phases: 1
- Power Frequency: 50 Hz
- Power Consumption: 46 W
- Protection Level: PE
- Insulation Class: 1
- Cable Length:
 - L = 1200 mm (2 × 0.75 mm² + ground wire)
 - L2 = 1400 mm
- Maximum Current in Heating Phase: 4.5 A
- Average Current over 24h Heating Cycle: 2 A
- Standby Current: 0.1 A
- Maximum Power Consumption: 60 W
- Standby Power Consumption: 1.2 W

- Operating Temperature: 56°C
- Temperature Stabilization Accuracy: $\pm 1^\circ\text{C}$

INSTALLATION

The biological incubator must be placed on a flat, stable surface. Due to the characteristics of the device, the following precautions must be observed:

- Do not use the device in locations where there is a risk of explosion.
- Remove the device from its packaging and inspect it. If any damage is detected, contact the distributor or manufacturer immediately.
- Place the device on a flat surface, away from heat sources such as ovens, washing machines, or cookers.
- Position the device so that the warning light is clearly visible to the user.
- First, connect the power cord to a 230V electrical outlet. Then, plug the 12V cord into the device's power socket (see Figure 1).

Important: Use only the power adapter provided in the set (12V, 5A). Using any other power supply may damage the module and void the warranty.

DEVICE OPERATION OVERVIEW

A biological test incubator is a device intended for use in medical and/or dental applications.

- The manufacturer is not responsible for any damage or injury caused by improper use or lack of user training.
 - Operating temperature range: $+17^\circ\text{C}$ to $+27^\circ\text{C}$
- There is an internal diode that signals a device malfunction or a temporary power loss.

A thermistor inside the device continuously monitors and regulates the temperature during operation.

INSTRUCTIONS FOR USE

1. Connect the 230V power cord to the main socket, then connect the 12V power cord to the device's power socket (see Fig. 1 – item D).
2. The warning LED (see Fig. 1 – item A) will light up red.
3. Press the Start/Standby button (see Fig. 1 – item E). The LED will begin flashing green, indicating that the heating cycle has started.

Once the device reaches the operating temperature of 56°C, the LED will emit a steady green light.

The electronic control board maintains a constant temperature in the aluminum heating block throughout the incubation period by regulating the heater.

In case of a power failure:

- The device automatically resumes temperature control once power is restored.
- If the temperature remains above 52°C, the device will resume normal operation.
- If the temperature drops below 52°C, the incubator will reheat the block to the correct temperature.
- If the temperature cannot be restored, the LED will turn red, indicating a fault.

TEST VALIDITY

At the end of the incubation period (refer to the instructions specific to the ampoule test used), if the LED indicator lights up **green**, the test is **valid**.

If the LED is **red** at the end of the incubation time, the test is **invalid**.

Safety Warnings:

- Hot surface: Do not touch the heating block while the device is operating.
- Keep out of the reach of children.
- Do not touch metal surfaces while the device is in operation.

INSTRUCTIONS FOR USE

1. Identify the appropriate biological indicator to be used.
2. Insert the biological indicator into a representative package for sterilization, as recommended by the autoclave manufacturer.
3. Run a complete sterilization cycle.
4. After sterilization, remove the indicator from the autoclave and allow it to cool.
5. Let the indicator rest for at least 10 minutes to cool properly.
6. Check for discoloration on the chemical indicator strip inside the biological indicator (see Fig. 2).
7. Crush the glass capsule (see Fig. 2) inside the biological indicator using pliers or by pressing it against a hard surface. This should be done immediately before placing the indicator in the incubator.
8. The indicator must be inserted into the incubator within one hour after sterilization.
9. Insert the indicator into the designated slots of the aluminum heating block.
10. Start the incubation process by pressing the START button. After approximately 10 minutes, verify that the LED has changed from flashing green to steady green.
11. If this change does not occur, refer to the maintenance and troubleshooting section.
12. Incubate the biological indicator for the duration specified in the indicator's manual.
13. For detailed procedures, please refer to the instructions for use included inside the biological indicator packaging.

UTILIZATION

- Directive on the Disposal of Electrical and Electronic Equipment (WEEE):
In accordance with Directive 2002/96/EC, this symbol indicates that the product must not be disposed of as municipal waste at the end of its service life.
- The product must be taken to a company specialized in the disposal of electrical and electronic devices, or returned to the supplier.
- The device owner is responsible for ensuring proper disposal.

MAINTENANCE AND TROUBLESHOOTING

The device does not require any specific maintenance, except for periodic cleaning of the heating element. This is necessary to prevent the build-up of deposits inside the biological indicator sockets.

When cleaning the incubator, it is sufficient to follow basic hygiene practices as recommended for similar devices:

- Switch off the device and disconnect it from the power supply.
- If the device was recently used, allow it to cool down completely.
- Wipe the surfaces with a soft cloth moistened with soapy water, then dry thoroughly with a clean paper towel.
- Do not use harsh detergents that may permanently alter the appearance of plastic parts.
- Do not spray chemicals directly on the device, as this may cause internal electrical damage.

Important:

- Do not attempt to maintain or repair the interior of the device.
- Disassembling the device will cause permanent damage and void the warranty.

ERROR SIGNALING

The device uses LED flashing patterns to indicate different types of errors:

1. Damaged Heater:
 - 3 seconds of red flashes
 - 3 seconds off

- Repeats continuously
- 2. Incorrect Power Adapter:
 - 4 short red flashes
 - Then 2 short red flashes
 - Repeats continuously
- 3. Heater Temperature Drop Below 52°C (due to power failure):
 - Continuous red light (steady)
- 4. Damaged NTC Sensor:
 - 1 second red flash
 - 1 second off
 - Repeats continuously

Resetting the Device to Clear Errors

To reset the device and clear the error signals:
Disconnect the power supply, wait 3 seconds, and reconnect.

Power Failure Recovery

The system constantly monitors both supply voltage and temperature. If a temporary power outage occurs without allowing the heater to cool down significantly:

- The device will automatically return to heating mode once power is restored.
- The error indicator will be cleared automatically, and the LED will show a steady green light, confirming that the device is operating correctly.

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PROCLINIC S.A.U.
LISTER BIOMEDICAL CO.,LTD.