VOCO

V-Print® splint 2.0

EN

Instructions for use

MD

EU Medical device

CE

0482

Product description:
V-Print splint 2.0 is a light-curing acrylic for the generative production of therapeutic splints as well as auxiliary parts and functional parts for use in dental diagnostics. Thanks to its special formula, **V-Print splint 2.0** displays particularly high impact resistance.

- Indications:**
- Therapeutic splints
 - Auxiliary parts and functional parts for diagnostics
 - Bleaching splints (home bleaching)

Contraindications:
V-Print splint 2.0 contains methacrylate and phosphine oxide.
V-Print splint 2.0 should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
V-Print splint 2.0 is suitable for use in all patients without any age or gender restrictions.

Performance features:
The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
V-Print splint 2.0 should only be used by a professionally trained dental practitioner.

Hardware and software requirements:

CAD software¹ dental scanner	Software for the planning and design of therapeutic splints as well as splints, auxiliary parts and functional parts for use in dental diagnostics. The software and dental scanner must satisfy local and current medical device specifications and allow for issuance of the patient-specific design as an STL data set.
CAM software	Software for preparation of the print order. The part will not be modified during this process. Structures that facilitate the 3D printing are simply created. For example: - Autodesk Netfabb version 2020 or later for SoliFlex 3D printing.

¹The designation **Software as Medical Device SaMD** includes standalone (autonomous) software that is a medical device (MD) and not part of one.

Manufacturing equipment	For example: VOCO SoliFlex 170 VOCO SoliFlex 350 VOCO SoliFlex 650 VOCO SoliFlex 170 HD
Post-curing devices	For example: Otofash G171

See also: accompanying list of resources or www.voco.dental/3dprintingpartners
All manuals and/or operating instructions for the respective programmes, and for device, materials and/or parts manufacturers, which are required for the manufacturing process, must be observed.

Clarify ahead of time whether the programmes, devices and/or objects that you intend to use have been designed and approved for the corresponding applications.

CAUTION: Non-authorised changes to the process equipment, parameters, or software could result in the **V-Print splint 2.0** end object not satisfying specifications.

User:
Preparation:
For an indication-appropriate CAD construction, the following design conditions must be observed:

- Minimum wall thickness 1,5 mm / Minimum wall thickness in case of an occlusal construction of adjusted splints: 2-3 mm
- Round off internal edges on the outside.

Prepare a print job using CAM-software.

Processing:

- Notes:**
- Use separate material containers and cleaning baths for each printing material so as to rule out the possibility of cross-contamination.
 - Protect uncured **V-Print splint 2.0** and not yet post-cured objects from exposure to light immediately in order to avoid accidental polymerisation.

The materials container should be filled immediately before the start of the printing process. It is important to ensure that the material is free of bubbles to the extent possible, and filled to the fill level mark.
Start the print job observing the parameters that you previously selected.
Once the printing process has ended, a dripping time of approximately 10 minutes is recommended. Next, carefully detach the printed objects from the build platform.
In the following steps, the printed objects will need to be cleaned, dried and post-exposed, in order to guarantee the required product characteristics. A detailed explanation of the steps outlined above can be found under **Post-processing**.
Post-process the printed object without delay so as to avoid distortion of its shape and ensure the required precision.

Recommendation: After completing your work, transfer the remaining material from the material tray to the original container (use a stainless steel sieve if necessary). This serves to check the material tray and also enables optimum storage of the print material.

Post-processing:
Cleaning

For cleaning purposes, use isopropanol (purity ≥ 98%) as a cleaning solution, in a cleaning device. An unheated ultrasonic bath or an unheated stirring bath may be used as a cleaning device.

The printed objects must be cleaned in two, or optionally in three steps. Position the unclean printed objects inside the cleaning bath so that any openings point downwards. Use tweezers or appropriate lowering baskets for transferring the printed objects into the detergent solution.
Please ensure that the printed objects do not come into contact with one another during cleaning.

	Ultrasonic bath	Stirring bath
Pre-cleaning (optional)	Carefully pre-clean the printed objects by submerging them several times in a beaker with isopropanol.	
Preliminary cleaning*	3 minutes – may be used multiple times	3 minutes – may be used multiple times
Final cleaning	2 minutes – fresh cleaning bath	2 minutes – fresh cleaning bath

*Note: The bath's cleaning efficacy decreases with increased use. Resin residue on the surface could indicate that the cleaning efficacy of the bath has diminished, or that the parts came into contact with one another. When the cleaning efficacy decreases, the respective bath must be replaced.
Next, the printed objects must be dried carefully using compressed air. If there is any resin residue on the printed object after the final cleaning, or if residue escapes from the undercuts when drying, the printed object can be briefly immersed once again in the final cleaning bath. Next, repeat the drying process.

Preparation for post-exposure:
Obstructive support structures should be removed before the post-exposure process, carefully and without exerting pressure, using a rotary instrument and as close to the printed object as possible. Use a suction device. Remove any remaining resin dust carefully using compressed air. Then rinse printed objects for a few seconds with fresh isopropanol and remove any dust using a fine, clean brush if necessary. Then dry the printed objects again carefully using compressed air.

Post-exposure:
Conduct the post-exposure a minimum of 15 minutes after the most recent contact with isopropanol.
It is important to ensure that the printed objects do not overlap or contact each other, as post-exposure would be negatively affected by the shadows that are cast.

Post-exposure can be conducted using the following devices:

Post-exposure device	Programme	
For example: Xenon photoflash unit Otofash G171 (optional with N2)	2 x 2000 flashes	After 2000 flashes, observe a cooling phase of at least 2 minutes with open lid. Next, turn over and light-cure with another 2000 flashes.

See also: accompanying list of resources

Finishing:
Generally work with low contact pressure and reduced speed. This guarantees consistent results and reduces processing marks.
In order to sand the support stubs use a rubber polisher or a fine-toothed carbide bur, for example. The bur can also be used for additional finishing of special structures.
In order to achieve form-fit grinding, e.g. between the support base and the print object, it is recommended to process the surface in the corresponding area with coarser or finer rubber polishers. A corresponding result can also be achieved with sandpaper, possibly with different grit sizes.
In order to obtain a high-gloss finish, polish the object first with a pumice stone. Then, thoroughly remove any pumice dust under running water using a brush. Finally, polish the object to a high-gloss using a buffer and high-gloss polishing paste, without applying excessive pressure.

Final cleaning:
Clean the object thoroughly. Coarse residue can be removed with the steam jet cleaner. The arising deformation is reversible. The final cleaning can be conducted by briefly placing the object in an unheated ultrasonic water bath. In order to remove oily or fatty contaminants, a surfactant solution may be used in place of water.

Disinfection:
Objects manufactured from **V-Print splint 2.0** may be disinfected using alcohol- or aldehyde-based disinfectants (e.g. ethanol (≥ 70 %), MD 520 by Dürr, Cavex ImpreSafe by Cavex). Observe the manufacturer's Instructions for use.

Warnings, precautionary measures:

- Only use **V-Print splint 2.0** intraorally in a fully cured state. Pay attention to the finishing process.
- Contact between uncured **V-Print splint 2.0** and the skin/mucous membranes and eyes can cause irritation and should be avoided.
- The wearing of protective clothing is recommended. Furthermore, it is important to ensure that no vapours and/or dusts are inhaled. The wearing of a suitable mask and/or the use of suction devices is recommended.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for the intended purpose.

Storage:
Storage at **15°C–28°C**. Seal the bottle immediately after use in order to avoid exposure to light and accidental polymerisation. Do not use after the expiry date.

Disposal:
Dispose of the product in accordance with local regulations.

Reporting obligation:
Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition, and a serious risk to public health that arises or could have arisen in association with the use of **V-Print splint 2.0** must be reported to VOCO GmbH and to the responsible authority.

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VC 60 AA6124 E1 0525 99 © by VOCO