# optiprint® interim

Germany's favorite 3D Resin

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### 1. Produktbeschreibung

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## EN - Instructions for use

### 1. Product description

optiprint interim is a light-curing 3D printing resin for the additive manufacturing of removable denture bases, interim partial dentures and their other custom-made resin parts for various medical indications, prescribed by dental professionals. It is suitable for use in devices of 385 nm and 405 nm wavelength. The target patient group for this class I medical device is defined as adults and adolescents.

### 2. Indication

optiprint interim is a 3D printing resin for the additive manufacturing of interim partial dentures and as well as their small parts for up to 30 days permanent use.

### 3. Contraindication

The material should not be used for any purpose other than the additive manufacturing of the indication. Do not use the polymerised material in the case of allergies to any of the ingredients (contains methacrylate monomers and oligomers).

### 4. Risk minimization and safety instructions

Improper use and deviations from the described processing will result in impaired quality and biocompatibility as well as undesirable mechanical properties of the finished printed part. Biocompatibility is only guaranteed with proper application (post-polymerization without oxygen). Only process all printed parts once they are fully polymerized. The post-polymerization of the printed parts takes place in a suitable polymerization device (e.g., Otoflash G171, see 5.4). After completion of the construction process, the printed part should be cleaned with suitable cleaning solution (e.g., isopropanol 99 %) in an ultrasonic bath. The batch no. must be indicated in any procedure requiring a material identification. For additively manufactured medical devices made of print resins, storing in water for 24 hours is recommended. Observe all recommended settings of the printer and the light-curing device. Read and observe the Safety Data Sheet (SDS) before use. For rupture safety, usual material wall thicknesses must be observed when creating the custom-made products (see 5.1). Proper personal protective equipment (nitrile gloves, safety glasses, protective clothing) must be worn when handling optiprint products as well as uncured printed parts. Avoid contact with skin and eyes before post-processing. The optiprint product can irritate eyes and skin. In rare cases, allergic reactions to components of optiprint products may occur. In case of accidental contact, follow the first aid measures (rinse thoroughly with water and consult a doctor if necessary, see the SDS). To be used by qualified professionals only. Keep out of reach of children!

### 5. Processing

#### 5.1 Designing

The minimum thickness should not be less than 1.5 mm. Make sure that the construction is aligned in a most possible horizontal position and that it is sufficiently supported by an appropriate support structure. To improve stability, a perforated base plate should also be selected.

#### 5.2 Printing

Observe the instructions for use of the printer and the software. Ensure that work is carried out in clean conditions. Contamination on the 3D printer can cause defects on the printed part and damage the tray. The relevant material parameters for optiprint resins can be obtained from the dentona homepage or from the device manufacturer. Shake the optiprint bottle before use, and fill the material tray of the printer with sufficient material. Remove any bubbles with a clean instrument/spatula.

If you are using a printer without integrated heating, it is recommended to heat both the printer and the optiprint material to an operating temperature of 30 °C/86 °F. A cold start should be avoided.

#### 5.3 Cleaning after printing

Due to its high viscosity, the printing resin does not drip completely from the build platform after printing. Wipe off the residues of the printing resin with a clean wooden spatula. Then remove the build platform from the printer and detach the printed parts from the platform. Do not separate the support structures and the base plate before cleaning/post-curing to avoid distortion.

#### 5.4 Recommendation

In a non-heated ultrasonic bath: place the cleaning liquid and the printed part in a second, sealable container (5 minutes pre-cleaning and 5 minutes post-cleaning). Isopropanol (99 %) should be used as the cleaning liquid. Following this, dry the part with compressed air and check the cleaning results. Complete the post-treatment of the part quickly and adhere to the time specifications! Avoid an unnecessarily long bath in the cleaning liquid and long waiting times between the steps! Alternatively, use the corresponding Rapid Shape Wash program with isopropanol (99%). You can recognize a sufficient cleaning if the surface is matt. Shiny areas require spot post-cleaning with cleaning liquid and a brush. Clean until no shiny areas are visible any more.

#### 5.5 Finishing off

The final properties and also the final color depend on the post-curing process! To avoid unwanted distortions, it is recommended to perform the post-curing on the model or to leave the support structure including the base plate on the object for post-curing. To achieve the desired material properties and biocompatibility as well as the coloration, the fully cleaned and dried printed parts must be post-cured in oxygen-free conditions.

#### 5.6 Recommendation

Otoflash G171 flash-light device (NK Optik) with nitrogen flooding. Ensure you use the Plexiglas tub with a UVB filter! It can be identified by the inscription: NK Optik 360N2. Alternatively, use the corresponding Rapid Shape Cure program with vacuum. When post-curing in other light curing devices, make sure that they offer a product-specific, predefined setting or a comparably high energy input (200 W).

LIGHT OUTPUT IN THE OTOFLASH G171	2 x 3000 flashes (turn after 3000 flashes)
ALTERNATIVE LIGHT OUTPUT (200W, 315-400 NM)	10 minutes

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