Ivotion Base Print

[en] Instructions for Use
Denture base material for stereolithography

Rx ONLY

CE 0123

Manufacturer:
Ivoclar Vivadent AG
Bendererstrasse 2
9494 Schaan/Liechtenstein
www.ivoclar.com

Date information prepared:
2024-12-12 / Rev. 1

ivoclar

English

1 Intended use

Intended purpose
Fabrication of bases for removable complete dentures

Patient target group
Edentulous adult patients

Intended users / Special training

- Dental laboratory technicians (fabrication of restorations in the dental laboratory)
- Denturists (fabrication of restorations in the dental laboratory)
- No special training required.

Use
For dental use only.

Description
Ivotion Base Print is a light-curing resin material for stereolithography 3D printing. It is suitable for the fabrication of permanent denture bases in removable complete denture prosthetics.

 The material can be processed with approved 3D printing systems.

Technical data

	Specifications in accordance with ISO 20795-1 Type 2* and Type 4*	Typical material properties
Flexural strength	≥ 65 MPa	72 MPa
Flexural modulus	≥ 2000 MPa	2250 MPa
Water absorption	≤ 32 µg/mm ³	15.7 µg/mm ³
Water solubility	≤ 8.0 µg/mm ³	2.3 µg/mm ³
MMA residual monomer content	≤ 2.2 %	≤ 0.1 %

* Currently, there is no standard that defines the physical requirements for denture base materials processed with 3D printing systems. The product complies with the requirements of ISO 20795-1 for light-activated denture base resins. The water solubility value fulfils the specification for autopolymerizable materials defined in ISO standard 20795-1.

Indications
Complete edentulism

Areas of application:
Denture bases in removable complete denture prosthetics

Contraindications
The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use
Unpolymerized material must not come into contact with the patient's mouth.
Adhering to the following framework conditions is essential to ensure the successful processing of this product:

- Observe the required minimum thickness (2 mm).
- Only process in approved 3D printing systems.
- Staining / layering with SR Nexco.

 For guidance regarding processing, please consult the Instructions for Use of the product.

- The final denture is for use by one patient only.

System requirements
Ivotion Base Print can be processed in approved 3D printing systems. A current overview of approved systems can be found on www.ivoclar.com.

Side effects
In individual cases, allergic reactions to (meth)acrylate materials have been reported.

Interactions
There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition
UDMA, tricyclodecane methanol acrylate, copolymer, fumed silica, photoinitiator TPO-L

2 Application
The use of the product in conjunction with approved 3D printing systems is explained in the respective instructions for use. A current overview of approved systems can be found at www.ivoclar.com. The description below is based on the use of the material in conjunction with the PrograPrint System.

2.1 Preparation of the printing procedure
Ivotion Base Print is suitable for processing in PrograPrint devices. Before a new bottle is inserted into the PrograPrint Cartridge, shake the material in the bottle for two minutes. As a result of the shaking, foam may form. To use the resin, remove the screw cap of the bottle and snap on the dispensing valve included in the delivery. Make sure to correctly align the valve. The valve must be in the "closed" position. The bottle can now be inserted into the PrograPrint Cartridge. In case of foam formation, the material must be kept protected from light for one hour before it can be used. If the printer is not used for several weeks, components of the material may settle in the PrograPrint Pool. In this case, the material must be manually mixed.

 For details about the insertion of the cartridge and the operation of your PrograPrint printer, please refer to the Operating Instructions of the device.

To ensure the traceability of the material batch, the batch number of each new resin container must be separately documented. If the material polymerizes prematurely, it must be replaced completely.

2.2 Cleaning
After the printing process is complete, objects made from Ivotion Base Print must be cleaned and post-cured promptly. For this purpose, the objects are left on the PrograPrint Stage.

- If a considerable amount of uncured material adheres to the surface of the object, we recommend leaving the object to drain in the printer for several minutes. During this process, the material must not be exposed to ambient light.
- The cleaning unit PrograPrint Clean is designed for cleaning printed objects, using isopropanol cleaning fluid (≥ 99%). Cleaning should be performed immediately after the printing process in order to ensure optimum material properties.

 For detailed information on the procedure, please refer to the PrograPrint Clean Operating Instructions.

Use the following parameters for cleaning objects in the PrograPrint Clean:

Rough cleaning	Fine cleaning
2 min	2 min
850 rpm	

- After cleaning, the objects must be dried with oil-free air at low pressure from a distance of 20–30 cm. If shiny or insufficiently cleaned spots are still recognizable on the object, it must be finely recleaned with isopropanol from a spray bottle and blown dry again before post-curing.

2.3 Option 1: Bonding with 3D-printed teeth prior to joint post-curing

Check the fit of the teeth in the denture base. Correct any interfering contacts on the tooth itself. Use cross-cut or diamond-coated tungsten carbide burs for adjustments. Before joint post-curing the denture base and teeth, use uncured Ivotion Base Print material to bond the two components. Remove any excess material with a brush.

Also refer to the instructions for use of the tooth material.

2.4 Post-curing

- Post-cure the object in the post-curing unit PrograPrint Cure using a program that is especially suitable for the Ivotion Base Print material.

 For detailed information on the procedure, please refer to the PrograPrint Cure Operating Instructions.

The parameters are stored as programs on the PrograPrint Cure or can be obtained through an update.

- Once the printed objects have been successfully cleaned and post-cured, they may be removed from the build platform. For this purpose, the tools included in the PrograPrint printer delivery can be used. Detach the objects including the support structures from the PrograPrint Stage with the help of these tools.

2.5 Option 2: Bonding with denture teeth after post-curing the denture base

Check the fit of the teeth in the denture base. If there are occlusal interferences, make adjustments to the tooth. Use cross-cut or diamond-coated tungsten carbide burs. After post-curing the denture base, use an appropriate bonding system (i.e. Ivotion Bond).

 For processing information refer to the Instructions for Use.

2.6 Finishing

- Detach the support structures from the object using dental laboratory tools.
- Any remaining residues of the support structures should be removed using cross-cut tungsten carbide burs. These burs can also be used for making shape adjustments. When finishing the objects, pay particular attention to the margins of tooth cavities, interdental spaces and basal surfaces. Make sure to wear protective goggles during finishing.
- Please note that the improper use of grinding tools, polishing brushes and steam jets can lead to overheating of the object, which may result in damage to the material.

2.7 Polishing

- Polish the object in the usual manner using customary brushes and polishing agents suitable for denture base materials. As a result of the additive manufacturing process, layer lines can form in some areas of the denture base. These layer lines can be removed by polishing.
- After polishing, clean the object under running water using a brush or with a steam cleaner or in an ultrasonic bath.

2.8 Repair / Relining

A denture fabricated with Ivotion Base Print can be repaired or relined using the conventional self-curing denture base material ProBase Cold.

 When using ProBase Cold for repairs or adjustments, follow the respective instructions in the Instructions for Use of the product.

2.9 Application notes

- Failure to observe the stipulated processing instructions may lead to failure.
- Avoid overheating of the material when finishing the objects.
- Protect the material from direct exposure to light once you have removed it from the packaging.
- Take care not to mix or confuse different materials.
- The products are intended to be used only by professionals. The fabrication in the dental laboratory must be carried out according to standard norms accepted by the profession.
- If materials are mixed, their biocompatibility may be compromised.
- During the printing process, foam may form. However, this does not affect the qualitative properties of the product in any way.
- The shading of Ivotion Base Print is based on the Ivoclar shade standards for denture base materials. If required, shade samples may be printed using the material to create a shade guide. The samples should be stored in a place protected from sunlight. The respective digital printing data can be found at www.ivoclar.com.
- Typically, white discolouration may occur when stereolithography resins are used in combination with a solvent. In order to prevent these discolourations, allow the object to dry for 10 minutes after cleaning, avoiding exposure to ambient light. If white discolouration persists, additional cleaning with isopropanol in an ultrasonic bath (2 minutes) is recommended. Alternatively, polish the affected surface with glass beads (100 µm) at a pressure of 0.5 bar after post-curing.

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).

Warnings

- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- Do not inhale grinding dust during finishing – use suction equipment and a face mask.
- The processing of the material can result in sharp edges, so that there is a risk of injury.
- Irritating to eyes, skin and respiratory system. Make sure to wear appropriate personal protective equipment (e.g. protective goggles, gloves, etc.).
- Avoid contact of the skin with uncured material and cleaning fluid. Commercially available medical gloves do not provide protection against the sensitizing effect of (meth)acrylates.
- If the material should come in contact with the eyes, rinse with copious amounts of water and consult a physician, if necessary.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

- Liquid, uncured material as well as cleaning fluid must not get into the sewage system, surface water or groundwater. Empty material bottles can be disposed of with household waste once they have been thoroughly cleaned. Uncured material must not be disposed of with household waste.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Fracture of the denture may lead to ingestion of material.
- There is a risk of injury due to fragments.
- Denture stomatitis may occur.

4 Shelf life and storage

- Storage temperature <30 °C (<86 °F)
- Do not use the product after the indicated expiry date.
- Expiry date: See note on the packages. Once the bottle has been opened, the shelf life of the material is six months.
- Store material at room temperature in a dark, well-ventilated place.
- Store opened material bottles in the closed PrograPrint Cartridge.
- To avoid premature curing, minimize exposure of liquid material to daylight as much as possible.

Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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