UTS 3D

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Date information prepared:
2020-08-26/740988/WW/Rev. 0

ivoclar
vivadent[®]
passion vision innovation

Views of the Universal Transferbow System / List of Parts

1.0 Universal transferbow (complete)

- 1.1 Earpiece mounting
- 1.2 Earpiece rod screw
- 1.3 Side arms
- 1.4 Side arm support
- 1.5 Scale
- 1.6 Multifunctional mounting
- 1.7 Width setting screw
- 1.8 Registration joint connection
- 1.9 Connection for nosepiece
- 1.10 Nosepiece holder screw
- 1.11 Registration joint connection screw
- 1.12 Thread for screw 3.1

2.0 Individually adjustable nosepiece (complete)

- 2.1 Nosepiece holder
- 2.2 Nosepiece rod
- 2.3 Nosepad holder
- 2.4 Nosepiece rod screw

3.0 Reference indicator

- 3.1 Reference indicator screw

4.0 Axis pins

- 4.1 Thread
- 4.2 Adapter pin for articulator

5.0 Intercondylar distance (S,M,L)**6.0 Earpieces (complete)**

- 6.1 Earpiece rod
- 6.2 Threaded openings for axis pins
- 6.3 Earpiece
- 6.4 Threaded pin for earpieces
- 6.5 Earpiece support

7.0 3D registration joint

- 7.1 Connector
- 7.2 Tommy screw
- 7.3 Bite fork connector
- 7.4 Bite fork connector screw
- 7.5 3D bite fork connector

8.0 3D bite fork

- 8.1 Connector

9.0 Earpiece, black, 'large'



- 9.1 Threaded connector

10.0 Nosepads

- 10.1 Nosepads, green (medium)
- 10.2 Nosepads, blue (soft)

11.0 3D bite fork "S"

- 11.1 3D bite fork connector "S"

3D registration joint (7.0)

3D bite fork (8.0)
3D bite fork "S" (11.0)

Universal transferbow (complete) (1.0)

Earpiece (complete) (6.0)

Nosepiece (complete) (2.0)

Nosepads

Reference indicator (3.0)

Earpiece (L) (9.0)

Preface

Dear Customer,

Thank you for having purchased the UTS 3D transferbow.
Inappropriate use may damage the equipment and result in personal harm.
Please observe the relevant notes and read the Operating Instructions.
We hope you enjoy working with the UTS 3D transferbow.

1. Intended use

Intended purpose

Bite registration

Use

For dental use only!

Description

The UTS 3D system is a transfer system. It enables the transfer of the spatial position of a patient's individual Bonwill triangle to the joint elements of the articulator. With its adjustable nosepiece, the UTS 3D system is also suitable for model transfer according to Camper's plane (CP) or the Frankfort horizontal (FH). The adjustable earpieces also make the UTS 3D suitable for model transfer according to specific coordinates.

Indications

Determining the position of the occlusal plane in relation to reference planes on the skull for the transfer to an articulator.

Contraindications

If the device is used according to the instructions, no contraindications are known to date.

Side effects

There are no known side effects to date.

Composition

The following components of the product are intended for intraoral use:

- Position 8.0, 3D bite fork, aluminium AlMg1, colourless anodized

Signs and symbols

The signs and symbols in these Operating Instructions facilitate the finding of important points and have the following meanings:

Risks and dangers

Contraindication

Important information


Observe the Operating Instructions (Failure to observe the Operating Instructions may expose the patient or operators to a risk.)

Target group

Dentists, dental technicians, dental staff

2. Application

The UTS 3D system is a transfer system. It is a fact that some jaws are larger and some are smaller. The Bonwill triangle of various patients may differ up to 50 mm, depending on the size of the skull.


Individual, spatial model orientation by means of a transferbow takes these parameters into account and helps avoid occlusion errors.

The UTS 3D system permits the spatial transfer of a patient's individual Bonwill triangle to the joint elements of the articulator.

With its adjustable nosepiece, the UTS 3D system is suitable for model transfer according to Camper's plane (CP) or the Frankfort horizontal (FH). The adjustable earpieces also make the UTS 3D suitable for model transfer according to specific coordinates.


The Stratos 200 permits working according to these reference planes (CP and FH). For this purpose, two different registration joint holders are available (CP and FH). The registration joint holders are used according to the reference plane selected. Thanks to the registration joint holder of the Stratos 200, the transferbow no longer has to be transported to the dental laboratory. Only the 3D registration joint and the bite fork have to be sent to the lab.

e.g. CP version

The UTS 3D system may also be directly connected to the articulator by means of the vertically adjustable type 2 support pins. For this purpose, however, the entire transferbow including the 3D registration joint and bite fork has to be sent to the dental laboratory.

2.1. Installation and initial start-up

2.1.1. Unpacking and checking the delivery contents

Please make sure that the contents are complete.


- 1 Universal transferbow, complete
- 1 Nosepiece, complete
- 1 Reference indicator
- 1 3D registration joint
- 1 3D bite fork
- 2 Earpieces with support, complete
- 5 Nosepads, green (medium)
- 5 Nosepads, blue (soft)
- 1 Operating Instructions

2.1.2. Assembly and initial start-up

The UTS is already partially assembled.


Mount earpieces:

Mount the desired earpieces (6.3) (M or L) on the threaded pins for earpieces (6.4). Loosen the earpiece rod screw (1.2) and mount the entire assembly on the earpiece mounting (1.1). Make sure the earpieces are correctly positioned.


Attach nosepiece with nosepads:

Place the desired nosepad (10.1 or 10.2) in the nosepad holder (2.3).


Insert the nosepiece rod (2.2) into the nosepiece holder (2.1).

The nosepiece holder (2.1) may now be attached to the connection for the nosepiece (1.9) using the nosepiece holder screw (1.10) of the multifunctional mounting (1.6).


Reference indicator:

Mount the reference indicator (3.0) on the thread for screw 3.1 (1.12) by means of the reference indicator screw (3.1).

3D registration joint with bite with 3D bite fork:

Insert the connector (8.1) of the 3D bite fork (8.0) into the bite fork connector (7.5) and tighten the bite fork connector screw (7.4). Insert the connector (7.1) into the 3D registration joint connection (1.8) and secure the 3D registration joint (7.0) with the 3D registration joint connection screw (1.11).

Axis pins:

The axis pins (4.0) are mounted in the laboratory if the UTS 3D is directly connected to the articulator.


When using the UTS 3D on the patient, the axis pins (4.0) must not be used in order to prevent injury.


2.2. Handling, operation

2.2.1. Earpieces


The earpiece rod (6.1) may be adjusted by loosening the earpiece rod screw (1.2).


Please note that the earpiece support (6.5) is touching the earpiece mounting (1.1) if the CP or FH versions are used.

2.2.2. Nosepiece


The height of the nosepiece holder may be adjusted to the patient's individual reference point by loosening the nosepiece holder screw (1.10).


By loosening the nosepiece rod screw (2.4) the nosepiece rod (2.2) may be adjusted horizontally to fit the patient.

2.2.3. Reference indicator


The reference indicator (3.0) may be adjusted to the patient's individual reference point by loosening the reference indicator screw (3.1).


2.2.4. 3D registration joint


The 3D bite fork (8.0) is attached to the 3D bite fork connection (7.5) by means of the 3D bite fork connector screw (7.4).

The bite fork (older version) can be attached to the connector (7.3), as the connector is equipped with the corresponding adapter on the opposite side. For this purpose, the connector screw (7.4) has to be changed to the opposite side.


The 3D registration joint (7.0) is secured on the multifunctional mounting (1.6) with the registration joint connection screw (1.11).

The tommy screw (7.2) secures the 3D registration joint.



2.2.5. Adjusting the width

The width of the transferbow can be adjusted by loosening the width setting screw (1.7). If the minimum width is not sufficient, the black earpieces 'L' (9.0) can be used. In this way, an additional 20 mm of space is available.

2.2.6. Additional height in the area of the nosepiece

If the nosepiece rod (2.2) is rotated by 180°, an additional 10 mm is gained vertically.

2.3. Practical use on the patient

2.3.1. Preparing the bite fork

Coat the bite fork with sufficient hard wax (warmed) or silicone. Press wax or silicone on the maxillary tooth row. Make sure that clear occlusal support on both sides is provided.

2.3.2. Reference plane

2.3.2.1. CP version

Camper's plane runs from the spina nasalis to the lower edge of the external auditory meatus. Mark the spina nasalis.

2.3.2.2. FH version

The Frankfurt horizontal runs from the suborbital point to the upper outer edge of the auditory meatus. Mark the orbital point.

2.3.3. Positioning the transferbow

The patient should be in an upright position. Loosen the width setting screw (1.7) and completely open the transferbow. Introduce the earpieces (6.3) into the external auditory meatus by slowly pushing the side arms together. To relieve some of the weight, the patient may support the UTS 3D by holding the side arms with both hands.

Do not apply too much pressure when positioning the earpieces to prevent injury to the auditory meatus.

Next, tighten the width setting screw (1.7) to set the width.

2.3.4. Positioning the nosepiece


With the 'nosepiece' screw open (2.4), exert slight pressure with the thumb to position the nosepiece rod (2.2) on the glabella so that the bow sits firmly.

2.3.5. Aligning the transferbow

Align the transferbow by vertical adjustment of the nosepiece, pointing with the reference indicator (3.0) to the desired reference plane (CP or FH).

CP version = spina nasalis

FH version = orbital point


When manipulating the reference indicator (3.0), make sure to observe a certain distance to the eye of the patient in order to prevent injury.

2.3.6. Other versions

By means of the vertically adjustable earpieces, the UTS 3D transferbow is also suitable for model transfer according to specific coordinates. Please refer to the literature for further details on this subject.

2.3.7. Securing the 3D registration joint

Check again if the transferbow is tightly positioned with regard to the relevant reference points. Attach the 3D bite fork (8.0) to the open 3D registration joint (7.0) and position the 3D bite fork on the mandibular tooth row. Additional support for the 3D bite fork may be provided by cotton rolls. The 3D registration joint (7.0) is now secured at the 3D registration joint connection (1.8).

Tighten the tommy screw (7.2) for the rotary joint. The registration procedure is now complete and the corresponding parameters are recorded in the 3D registration joint.


2.3.8. Removing the transferbow

Loosen the width setting screw (1.7) for the nosepiece (1.9). Pull the side arms outward and instruct the patient to open the mouth. Now the whole transferbow including the registration can be removed.

2.3.9. Removing the 3D registration joint

Loosen the registration joint connection screw (1.11) and remove the 3D registration joint (7.0).



The transferbow may be used in combination with the Gnathometer M or the Centric Tray in the manner described above.

Important information regarding the UTS 3D registration joint

Please note that the tommy screw of the UTS 3D registration joint must always be attached and secured towards the right as shown in the picture. If this point is not observed, problems may occur at the registration joint holder (CP or FH version) during transfer to the articulator.

3. Maintenance, cleaning, diagnosis

This chapter describes the user maintenance and cleaning procedures for UTS 3D transferbow. Only those tasks are listed that may be performed by dental professionals and dental staff. All other tasks must be performed by qualified service personnel at a certified Ivoclar Vivadent Service Centre.

Monitoring and maintenance

The time for these maintenance procedures depends on the frequency of use and the working habits of the user. For that reason, the recommended times are only approximations.

Cleaning

Item	Frequency	Cleaning material
Screws that are contaminated with wax or plaster	If required	Rinse with hot water and dry with a dry cloth
Screws that are contaminated with silicone	If required	Clean with a cloth
3D bite fork that has come into contact with saliva or blood	After every use	Sterilization in the autoclave at 134 °C for 5.5 min. using the prevacuum method. Please observe the valid standards and regulations..
3D registration joint	After every use	Clean with a spray disinfectant
UTS 3D transferbow	After every use	Clean with a spray disinfectant
Nosepads	After every use	Use new nosepad (nosepads are disposable)
Earpieces	After every use	Autoclaving device. Please observe the corresponding instructions

Maintenance

Avoid any contact with strong acids and solvents (e.g. MMA, acetone) to prevent damage of the surfaces.

Sterilize and/or disinfect the components before the first use and after every use and/or after contact with saliva or blood.

- Please observe the cleaning and maintenance instructions.
- Use only disinfectants suitable for stainless steel and light metals.
- Avoid any contact with strong acids, brines or solvents.
- Clean in an ultrasonic bath, with water, steam jets or sterilization devices.
- The 3D registration joint must not be cleaned in an ultrasonic bath, with water, steam jets or sterilization devices.

4. What if ...

This chapter will help you to recognize malfunctions and take appropriate measures or, if possible, to perform some minor repairs.

Technical malfunction

Description	Instructions for users	Measures
Patient has pressure sores and complains about pain of the nose	Is the nosepad in place?	Use nosepad.
Patient complains about pain in the ear	Is the width of the transferbow too narrow?	Do not close the transferbow too tightly.
The reference to the corresponding plane (CP or FH) has not been correctly recorded	Were the earpieces correctly mounted so that they came into contact with the earpiece support?	When recording according to the CP and FH method, the earpieces must always come into contact with the earpiece support.
The surface of the transferbow is damaged or partially corroded	Was the surface cleaned with acid or a solvent?	Do not use any acids or solvents to clean the transferbow.
The surface of the transferbow is scratched and damaged	The surface was scratched with a sharp instrument	Do not scratch the surface with a sharp object.

Repair

Repairs may only be carried out by a certified Ivoclar Vivadent Service Centre. Please refer to the addresses of the Service Centres on the back of these Operating Instructions.

If repairs during the warranty period are not carried out by a certified Service Centre, the warranty will expire immediately. Please also refer to the corresponding warranty regulations.

5. Safety information

In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclarvivadent.com and your responsible competent authority.

The current Operating Instructions are available in the download section of the Ivoclar Vivadent AG website (www.ivoclarvivadent.com)

Disposal information

The product must be disposed of according to the corresponding national legal requirements.

Risks and dangers

Risk of injury

When manipulating the UTS 3D while it is attached to the patient, make sure that the reference indicator does not injure the eye.

Do not apply too much pressure when positioning the earpieces in the auditory meatus to prevent injury.

The UTS 3D must only be used for the purpose described in Chapter 2. Further instructions to assure proper use of the UTS 3D transferbow:

- The instructions, regulations, and notes in these Operating Instructions must be observed.
- The UTS 3D transferbow must be properly maintained (see Chapter 3).

6. Shelf life and storage

This product does not require any special storage conditions.

7. Product specifications

Technical data

Weight:

Nosepiece: 25 g
3D registration joint: 94 g
UTS 3D basic bow 160 g

Materials:

Metal parts: aluminium, steel
Side arms: carbon composite material
Earpieces: DELRIN resin

Adjustment possibilities:

Adjustable width:
– with earpiece Large: 105 mm to 175 mm
– with earpiece Short: 85 mm to 155 mm
Adjustable width: 87 mm to 160 mm
Adjustable nosepiece height: 85 mm
Adjustable nosepiece length: 22 mm
Adjustable earpiece height: 22 mm
Adjustable reference indicator height: 66 mm

8. Additional information

Keep out of the reach of children.

The device has been developed solely for use in dentistry. Start-up and operation should be carried out strictly according to the Operating Instructions. Liability cannot be accepted for damages resulting from misuse or failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the apparatus for its suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.