Termo Resin Proclinic

Instructions for use

proclinic

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EN Instructions for use.

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SK Návod na použitie.

SL Navodila za uporabo.

ET Kasutusjuhend.

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SV Bruksanvisning.

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INSTRUCTIONS FOR USE EN

Material/intended purpose: Prosthesis base resin, hot polymer based on methyl methacrylate in powder and liquid form. Classification according to ISO 20795-1 Type 1 Class 1 and according to Medical Devices Directive 93/42/EEC Annex IX Class IIa for removable and fixed dentures. Only for intended use by qualified personnel (dental laboratory/dental practice).

Indications:

- Complete upper and lower dentures using the packing and pressing technique

Contraindications: Patient contact with the unpolymerised acrylic or its individual components is contraindicated. No restrictions on the patient group. However, use should be avoided in cases of proven allergy to one or more components of the product. No known interactions.

Brief information:

Mixing ratio powder/liquid [g]	10:4
Saturation phase approx.	10 min*
Plastic modelling phase approx.	20 min*
Polymerisation at 95°C	20 min*

*The times refer to a material and room temperature of 23°C. Higher temperatures shorten these times, lower temperatures extend them.

Properties:

- Easy processing
- High fit accuracy and guarantee of high-quality results
- Excellent mechanical values
- Optimal mucous membrane compatibility due to low residual monomer content
- Physiological colours and absolute colour stability
- Cadmium free
- Tested and certified biocompatibility
- Phthalate and BPA free

Termo Resin does not adhere to highly cross-linked, synthetic teeth. Synthetic teeth must be pretreated (roughened) and/or provided with mechanical retention.

Recommended mixing ratio: 10 g of powder with 4 g of liquid. Free dosage is possible, but the physical material properties can only be guaranteed if the components are weighed.

Preparation & Isolation: Scald the cuvette cleanly and then evaporate to remove wax residue. The plaster surfaces are thinly coated with alginate isolation twice, making sure that the first layer has not yet dried during the second layer. Important: The alginate isolation must be fully dried before first contact with plastic in order to prevent any white discolouration on the base surfaces.

Processing: Place the liquid in the mixing container (made of PE, PP or silicone) and quickly sprinkle in the appropriate amount of powder. With free dosage, sprinkle in as much powder as is necessary to completely bind the existing liquid. Then mix thoroughly with a wide spatula (made of PE, PP or metal). Make sure there are no bubbles. After about 10 minutes of saturation time, the packing time begins (the paste should no longer be sticky). The plastic paste is then kneaded homogeneously. The paste can be processed within another 20 minutes (end of packing time).

Press: Place the excess plastic paste into the cooled, lukewarm cuvette. Close the cuvette slowly and continuously increase the pressure on the cuvettes. Maintain pressing pressure for 10 minutes after the cuvette is closed. Test pressing is not necessary due to the good plastic consistency of the paste.

Recommended polymerisation: Place the cuvette with clamp in boiling water under polymerisation condition for 20 minutes. Of course, all other polymerisation techniques are also possible with this hot polymer.

After completion: The denture should be stored in water until incorporation. In principle, 24 hours of storage in water is recommended in order to further reduce the residual monomer content (0.8 % according to ISO 20795-1) and thus prevent possible taste disturbances.

Storage and safety instructions:

The material should no longer be used after the expiry date has passed. Do not store above 25°C (77°F). Avoid direct sunlight. Always close the container tightly after use. Keep away from heat/sparks/open flames/hot surfaces. Do not smoke. Dispose of contents/container in accordance with local regulations.
Contains methyl methacrylate; 1,4-butanediol dimethacrylate; dibenzoyl peroxide. Highly flammable liquid and vapour, possible burns. Causes skin irritation. Cause cause allergic skin reactions. Can irritate the airways. Avoid breathing dust/smoke/gas/mist/vapour/aerosol. Wear protective gloves/protective clothing/eye protection/face protection. In case of contact with skin: wash with plenty of water and soap. Please read the additional information in the safety data sheet.
Finished medical device: Taste disorders or injuries possible due to product breakage. Contact with alcohol, coffee, tobacco smoke or certain medications can cause discolouration.

Disposal: Disposal must comply with official regulations. Do not dispose of individual components in household waste or sewers. However, completely polymerised material can be disposed of with household waste.

Reporting: Serious incidents directly related to the device must be reported by the user/patient to the manufacturer and the responsible authority of the country in which they occurred, stating the batch and article number.

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