Auto Resin Proclinic

Instructions for use

proclinic

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EN Instructions for use.

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SK Návod na použitie.

SL Navodila za uporabo.

ET Kasutusjuhend.

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SV Bruksanvisning.

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EN INSTRUCTIONS FOR USE

Material/intended purpose: Prosthesis base resin, autopolymer based on methyl methacrylate in powder and liquid form. Classification as per ISO 20795-1 Type 2 Class 1 and as per Medical Devices Directive 93/42/EEC Annex IX Class IIa for removable and fixed dentures. Only for intended use by qualified personnel (dental laboratory/dental practice).

Indications:

- Completion of model cast prostheses
- Complete upper and lower jaw prostheses using the casting process
- Partial or total refining (indirect application)
- Additions and repairs

Contraindications: Patient contact with the unpolymerised resin or its individual components is contraindicated. No restrictions on the patient group. However, it should not be used in case of proven allergy to one or more components of the product. No known interactions.

Brief information:	
Mixing ratio powder/liquid [g]	10:7
Saturation phase approx.	30 s*
Casting phase approx.	3 min*
Plastic modelling phase approx.	3 min*
Polymerisation (2-3 bar, 45-50 °C) approx.	7 min**

* The times refer to a material and room temperature of 23°C. Higher temperatures shorten these times, lower temperatures lengthen them.
** The 7 min specification refers to a previous working time of up to 8 min. The time in the pressure vessel should last up to 15 minutes after mixing.

Properties:

- Universal application
- Average processing width
- Very good fluidity
- Good mucous membrane tolerance
- Absolute colour stability thanks to a catalyst system without tertiary amines
- Cadmium free
- Tested and certified biocompatibility
- Phthalate and BPA free

Auto Resin does not adhere to highly cross-linked, synthetic teeth. Synthetic teeth must be pretreated (roughened) and/or provided with mechanical retention.

Recommended mixing ratio: 10 g of powder with 7 g of liquid. Free dosage is possible, but the physical material properties can only be guaranteed if the components are weighed.

Preparation: With model cast prostheses, the teeth set in wax are fixed with a plaster or silicone pre-wall. The pre-wall, the model framework and the teeth are carefully cleaned of wax residues by scaling and evaporating. The plaster model (e.g. type 3 plaster) is soaked until no more air bubbles appear.

Isolation: The model and the plaster preliminary walls are coated twice with a thin layer of alginate isolation, making sure that the first layer has not yet dried during the second coating. Important: The alginate isolation must be fully dried before first contact with plastic in order to prevent any white discolouration on the base surfaces.

Processing: Place the liquid in the mixing container (made of PE, PP or silicone) and quickly sprinkle in the appropriate amount of powder. With free dosage, sprinkle in as much powder as is necessary to completely bind the existing liquid. Then mix thoroughly with a wide spatula (made of PE, PP or metal). Make sure there are no bubbles. After about 30 seconds of saturation time, the mixture can be poured for about 2 to 3 minutes. During the casting phase, the plastic is poured into the pre-wall. At the beginning of the plastic modelling phase, the plastic becomes stable, no longer flows out of the pre-wall and can be modelled. The plastic modelling phase lasts about 3 minutes.

Recommended polymerisation: After a maximum of 8 minutes, polymerisation is performed in a polymerisation pressure device at a pressure of 2-3 bar and a water bath temperature of 45 to 50°C for 7 minutes (up to 15 minutes).

After completion: The prosthesis should be stored in water until incorporation. In principle, 24 hours of storage in water is recommended in order to further reduce the residual monomer content (2,4% according to ISO 20795-1) and thus prevent possible taste disturbances.

Storage and safety instructions: The material should no longer be used after the expiry date has passed. Do not store above 25°C (77°F). Avoid direct sunlight. Always close the container tightly after use. Keep away from heat/sparks/open flames/hot surfaces. Do not smoke. Dispose of contents/container in accordance with local regulations. Contains methyl methacrylate; 1,4-butanediol dimethacrylate; dibenzoyl peroxide. Highly flammable liquid and vapour, possible burns. Causes skin irritation. Cause cause allergic skin reactions. Can irritate the airways. Avoid breathing dust/moek/gas/mist/vapour/aerosol. Wear protective gloves/protective clothing/eye protection/face protection. In case of contact with skin: wash with plenty of water and soap. Please read the additional information in the safety data sheet.

Finished medical device: Taste disorders o injuries possible due to product breakage. Contact with alcohol, coffee, tobacco smoke or certain medications can cause discolouration.

Disposal: Disposal must comply with official regulations. Do not dispose of individual components in household waste or sewers. However, completely polymerised plastic can be disposed of with household waste.

Reporting: Serious incidents directly related to the device must be reported by the user/patient to the manufacturer and the responsible authority of the country in which they occurred, stating the batch and article number.

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