Orto Resin Proclinic

Instructions for use

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EN Instructions for use.

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SK Návod na použitie.

SL Navodila za uporabo.

ET Kasutusjuhend.

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SV Bruksanvisning.

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INSTRUCTIONS FOR USE EN

Material/intended purpose: Orthodontic resin, autopolymer based on methyl methacrylate in the form of powder and liquid. Classification according to ISO 20795-2: Type 1 and according to Regulation (EU) 2017/745 Annex VIII, Chapter III: Class IIa for removable and fixed orthodontic appliances. Only for intended use by qualified personnel (dental laboratory/dental practice).

Indications:

- Orthodontic appliances using the mixing method
- Orthodontic appliances using the spray/scatter method
- Production of dental splints (bite ramp, snoring, etc.)

Contraindications: Patient contact with the unpolymerised material or its individual components is contraindicated. No restrictions on the patient group. However, it should not be used in case of proven allergy to one or more components of the product. No known interactions.

Brief information:	
Mixing ratio powder/liquid [g]	10:5
Saturation phase approx.	45 s*
Plastic modelling phase approx.	10 min*
Polymerisation (2-3 bar, 40-45°C) approx.	15 min*

* The times refer to a material and room temperature of 23°C. Higher temperatures shorten these times, lower temperatures lengthen them.

Properties:

- Universal application
- Comfortable processing width
- Good stability
- Safe and fast hardening
- Cadmium free
- Tested and certified biocompatibility
- Phthalate and BPA free

Recommended mixing ratio: 10 g of powder with 5 g of liquid. Free dosage is possible, but the physical material properties can only be guaranteed if the components are weighed.

Recommended equipment: To mix the components, it is recommended to use monomer-resistant materials - cups, spatulas, spray bottles and nozzles made of PE, PP, silicone or inert metals such as stainless steel. A white, type 3 dental plaster is recommended as a model material for the scattering method.

Preparation: The plaster model must be soaked until no more air bubbles emerge from the plaster, so that no air is pressed out of the plaster into the plastic during pressure polymerisation and bubbles are formed.

Isolation: The model and the plaster preliminary walls are coated twice with a thin layer of alginate isolation, making sure that the first layer has not yet dried during the second coating. Important: The alginate isolation must be fully dried before first contact with plastic in order to prevent any white discolouration on the base surfaces.

Processing - Mixing method: Place the liquid in the mixing container and quickly sprinkle in the appropriate amount of powder. With free dosage, sprinkle in as much powder as is necessary to completely bind the existing liquid. Then mix thoroughly with a wide spatula. Make sure there are no bubbles. After about 30-45 seconds of saturation time, the paste can be processed.

Processing - Spray/scatter method: Once the alginate isolation has dried thoroughly, the monomer is applied to the plaster model either from a dosing bottle or with a brush and then immediately sprayed with the polymer. This process must be repeated until the desired layer thickness has been achieved. It is particularly important to note that only as much liquid can be applied as can be absorbed by the powder. Before polymerising the plate, the surface must be wetted again with the liquid and smoothed.

Recommended polymerisation: The polymerisation takes place in a polymerisation pressure device at a pressure of 2 to 3 bar and a water bath temperature of around 40 - 45°C for around 15 minutes.

After completion: The prosthesis should be stored in water until incorporation. In principle, 24 hours of storage in water is recommended in order to further reduce the residual monomer content (3.4 % according to ISO 20795-2) and thus prevent possible taste disturbances.

Storage and safety instructions: The material should no longer be used after the expiry date has passed. Do not store above 25°C (77°F). Avoid direct sunlight. Always close the container tightly after use. Keep away from heat/sparks/open flames/hot surfaces. Do not smoke. Dispose of contents/container in accordance with local regulations. Contains methyl methacrylate; 1,4-butanediol dimethacrylate; dibenzoyl peroxide. Highly flammable liquid and vapour, possible burns. Causes skin irritation. Cause cause allergic skin reactions. Can irritate the airways. Avoid breathing dust/smoke/gas/mist/vapour/aerosol. Wear protective gloves/protective clothing/eye protection/face protection. In case of contact with skin: wash with plenty of water and soap. Please read the additional information on the safety data sheet.

Finished medical device: Taste disorders or injuries possible due to product breakage. Contact with alcohol, coffee, tobacco smoke or certain medications can cause discolouration.
Disposal: Disposal must comply with official regulations. Do not dispose of individual components in household waste or sewers. However, completely polymerised material can be disposed of with household waste.

Reporting: Serious incidents directly related to the device must be reported by the user/patient to the manufacturer and the responsible authority of the country in which they occurred, stating the batch and article number.

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