Provi Resin Proclinic

Instructions for use

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EN Instructions for use.

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SV Bruksanvisning.

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Material/intended purpose: Tooth-colored resin for prosthetics, autopolymer based on methyl methacrylate in the form of powder and liquid. Classification according to Regulation (EU) 2017/445 Annex VIII, Chapter III, Class IIa - invasive product for the production of individual crowns and bridges (veneers) for removable and fixed dentures. Only for intended use by qualified personnel (dental laboratory/dental practice).

Indications:

- Indirect extraoral fabrication of crowns and bridges
- Dental supplementation in the prosthetic application area

Contraindications: Patient contact with the unpolymersed material or its individual components is contraindicated. No restrictions on the patient group. However, it should not be used in case of proven allergy to one or more components of the product. No known interactions.

Brief information:

Mixing ratio powder/liquid [g]	10:5
Saturation phase approx.	60 s*
Casting/flow phase approx.	2 min*
Kneading phase approx.	5 min*
Polymerisation at room temperature approx.	10-14 min*

*The times refer to a material and room temperature of 23°C. Higher temperatures shorten these times, lower temperatures lengthen them.

Properties:

- Easy processing
- Good mechanical values
- Good mucous membrane tolerance
- Absolute colour stability thanks to a catalyst system without tertiary amines
- Natural colouring
- Cadmium free
- Tested and certified biocompatibility
- Phthalate and BPA free

Recommended mixing ratio: 10 g of powder with 5 g of liquid. Free dosage is possible, but the physical material properties can only be guaranteed if the components are weighed.

Mixing: Place the liquid in the mixing vessel (made of PE, PP or silicone) and quickly sprinkle in the appropriate amount of powder. With free dosing, sprinkle in as much powder as is necessary to completely bind the liquid present. Then mix thoroughly with a wide spatula (made of PE, PP or metal). Make sure there are no bubbles.

Production of a crown: Before preparation, take a situation impression with alginate or kneadable silicone and keep it until the work is finished. Prepare the tooth and make a plaster model (e.g. type 4 plaster) and insulate it with alginate insulation. An exhibition of the dump made of wax is possible. Introduce the mixed plastic into the situation impression during its plastic phase. Establish the original occlusion between the die model and the situation impression and allow the plastic to harden. The dentures can then be removed.

Production of a bridge: When making bridges, the gaps are closed using placeholders (e.g. plastic teeth or wax) before the impression is taken. Plastic teeth should be chosen in which grooves can be ground mesially and distally in order to achieve a clamping effect. If several teeth are missing, the plastic teeth must first be set up in wax and inserted into the gap between the teeth. This simulates a closed row of teeth in the impression. After the abutment teeth have been prepared, the dentures are created as described under "Manufacture of a crown". In the posterior region, the transverse fixation of the crown can be carried out by cutting a groove into the situation impression. The workflow after taking the situational impression also corresponds to the production of individual crowns. Material can be subsequently applied to the bar to fix occlusal contacts. Allow the plastic to harden. The dentures can then be removed.

Use in prothetics: The veneering material can also be used to repair or replace chipped or missing teeth on a prosthesis. The procedure is similar to the fabrication of crowns/bridges. A situational impression can be made with kneadable silicone and instead of the plaster model, the prosthesis to be repaired is placed against the impression and the plastic is very carefully dosed into the impression.

Recommended polymerisation: Allow to fully harden outside the pressure vessel in the impression for approx. 10 – 14 minutes. Or after about 5-7 minutes, polymerise in a polymerisation pressure device at a pressure of 2-3 bar and a water bath temperature of 45 to 50°C for about 10 minutes.

After completion: After you have created enough space for the fixing cement, the dentures can be ground, adjusted and polished. Before the incorporation, the dentures should be washed with soap and water. After drying, the dentures can e.g. B. can be used with eugenol-free fastening material. The temporary restoration should be stored in water until incorporation. In principle, 24 hours of storage in water is recommended in order to further reduce the residual monomer content (2.6 % according to ISO 20795-1) and thus prevent possible taste disturbances.

Storage and safety instructions: The material should no longer be used after the expiry date has passed. Do

not store above 25°C (77°F). Avoid direct sunlight. Always close the container tightly after use. Keep away from heat/sparks/open flames/hot surfaces. Do not smoke. Dispose of contents/container in accordance with local regulations. Contains methyl methacrylate; HEMATMDI; 1,4-butanediol dimethacrylate; dibenzoyl peroxide. Highly flammable liquid and vapour, possible burns. Causes skin irritation. Cause cause allergic skin reactions. Can irritate the airways. Avoid breathing dust / smoke / gas / mist / vapour / aerosol. Wear protective gloves / protective clothing / eye protection / face protection. In case of contact with skin: wash with plenty of water and soap. Please read the additional information in the safety data sheet.

Finished medical device: Taste disorders or injuries possible due to product breakage. Contact with alcohol, coffee, tobacco smoke or certain medications can cause discolouration.

Disposal: Disposal must comply with official regulations. Do not dispose of individual components in household waste or sewers. However, completely polymerised material can be disposed of with household waste.

Reporting: Serious incidents directly related to the device must be reported by the user/patient to the manufacturer and the responsible authority of the country in which they occurred, stating the batch and article number.

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