ENGLISH Instructions For Use

Warning: Protective eyewear for patients and staff is recommended to prevent accidental chemically related eye burns and severe eye irritations, when using the Etching Gel System. Avoid Etching Gel contact with oral soft tissue, eyes and skin. If accidental contact occurs, flush immediately with large amounts of water. See a physician if irritation persists. For eye contact, also consult physician after flushing with water.

Precaution: Avoid contact with eyes, skin and oral soft tissue.

First Aid:

Eye contact: Immediately flush with plenty of water. See a physician if irritation persists.

Skin contact: Wash affected area with soap and water. See a physician if irritation persists.

Oral Soft Tissue Contact: Immediately flush oral soft tissue with large amounts of water. If irritation persists, contact a physician.

Description

The Unitek™ Etching Gel Syringe Delivery System is comprised of a syringe having a luer lock end and interlocking, disposable dispensing tips. The syringe contains 35% phosphoric acid by weight. This system is indicated for intraoral delivery of acid gel for enamel etching in the acid etch technique.

Prepare Delivery System

1. Remove cap from Etching Gel syringe (Figure 1) and **SAVE**.

Figure 1

Twist a disposable dispensing tip securely onto the syringe (Figure 2).


Figure 2

A new dispensing tip should be used for each patient. Holding the tip away from the patient and any staff, extrude a small amount of etching gel onto a dispensing pad or a 2 x 2 gauze to assure that the delivery system is not clogged. If clogged, remove the dispensing tip and extrude a small amount of etching gel directly from the syringe. Remove any visible plug if present, from the syringe opening. Replace dispensing tip and again extrude etching gel. If clog remains, discard dispensing tip and replace with a new one.

Etching Procedure

1. Isolate teeth using the Dry Field System or a combination of retractors, dri-angles and cotton rolls.

2. Prophy teeth with an oil-free pumice or paste. Rinse with water.

3. **Enamel Etching:** Using the delivery system, gently extrude etching gel directly onto enamel surfaces to be etched. Leave gel in place for 15 seconds.

a. If desired, the etching gel may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.

b. If a liquid etchant is desired, the etching gel may be dispensed into a dappen dish and stirred to increase its fluidity.

4. **Rinse Etched Surfaces:** Thoroughly rinse for 15 seconds. Remove rinse water with suction. Do not allow patient to rinse. If saliva contacts the etched surfaces, rinse, re-etch for 5 seconds and rinse. Re-isolate the etched surfaces if using cotton rolls.

5. **Dry Etched Enamel:** Thoroughly dry etched surfaces. Air should be oil and water free. The dry etched surfaces should appear frosty white. If not, repeat steps, five, and six.

DO NOT ALLOW THE ETCHED SURFACES TO BE CONTAMINATED BY ANYTHING. If contamination occurs, rinse, re-etch, rinse and dry as above.

6. Immediately proceed with the bonding procedure.

Storage and Use

1. **Delivery System storage:** Remove used dispensing tip (Figure 3) and discard.

Figure 3

Twist storage cap on. (Figure 4).

Figure 4


Storage of the delivery syringe with a used dispensing tip or without the storage cap will allow drying of the etching gel and subsequent clogging of the system.

Replace storage cap with a new dispensing tip at next use.

Note: The ability of the Etching Gel or Liquid to properly change the surface of the enamel varies in each case. Enamel surfaces and their reaction to phosphoric acid treatment vary from patient to patient, from tooth to tooth in a patient, and even from area to area on a patient's tooth. However, the recommended procedure outlined in these instructions, should produce satisfactory results in the majority of cases.

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