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
English

1. Product description

The optiprint devices are light curing 3D printing resins for the production of dental custom-made parts by commercial users. They are suitable in devices with 385 nm and 405 nm wavelength, unless otherwise stated on the label.

The target patient group for medical devices is defined as adults and adolescents.

2. Indication

3D printing resin for the additive production of:

optiprint clara	bite splints	 MD	 C € 0044
optiprint lumina	temporary crowns and bridges	 MD	 C € 0044
optiprint laviva	denture bases	 MD	 C € 0044
optiprint tray 385 / optiprint tray HR	individual impression trays	 MD	 C €
optiprint guide 385 / optiprint guide HR/ optiprint guide 405	drill guides	 MD	 C €
optiprint I-B-T 385	transfer templates for the bracket positioning	 MD	 C €
optiprint mock-up	aesthetic try-ons	 MD	 C €
optiprint gingiva 385	gingiva masks		
optiprint match	scannable dental form parts		
optiprint model	dental models		
optiprint model align	aligner models		

3. Contraindication

The material should not be used for any purposes other than the additive manufacturing of the specified dental parts. Do not use the polymerized material if you are allergic to any of the ingredients (contains methacrylate monomers and oligomers).

4. Risk minimization and safety instructions

- Improper use and deviations from the described processing will lead to an impairment of quality and biocompatibility as well as undesirable mechanical properties of the final part.
- Biocompatibility is only guaranteed if used properly (light curing with nitrogen flooding)! All printed parts should only be processed once fully polymerized.
- The light curing of the printed parts takes place in a suitable light curing unit (e.g. otoflash G171), see table in point 5.3. Finishing.
- After completion of the construction process, the printed part should be cleaned with a suitable cleaning solution (e.g. isopropanol 99%, optiprint clean) in an ultrasonic bath.
- The LOT No. must be specified for each process that requires identification of the material.
- For additively manufactured medical products, storage in water for 24 hours is recommended.

- Observe all recommended settings for the printer and the light curing device.

- Read and observe the safety data sheet (SDS) before use.

- For the safety against breakage, observe the usual material wall thicknesses when creating customized products (see processing guide).

- The correct personal protective equipment (nitrile gloves, protective goggles, protective clothing) must be worn when handling optiprint products and non-cured printed parts.

- Avoid any contact with skin and eyes before light curing. The optiprint product can irritate eyes and skin.

- In rare cases, allergic reactions to components of optiprint products can occur. In the event of accidental contact, follow the “First Aid Measures” (rinse thoroughly with water and consult a doctor if necessary). See separate SDS (Safety Data Sheet).

- Use by qualified personnel only. Keep out of the reach of children!

5. Processing steps

Please see the relevant processing guide

5.1 Printing

Observe the IFU (Instructions for use) of the printer and the software. Pay attention to clean work. Contamination on the 3D printer can cause defects in the printed part and damage the material tray. Corresponding material parameters for optiprint resins can be downloaded in the printer manufacturer’s database. Shake the optiprint bottle before use and fill enough material into the tray. Remove any ocuring bubbles with a clean instrument / spatula. Ensure a processing temperature of 30 ° C.

5.2 Cleaning

Wipe off any residues of the printing resin with a clean spatula. Afterwards remove the build platform from the printer and release the printed parts with a suitable instrument (cutter knife) from the build platform. Then ideally cut off the support structures manually.

Cleaning: In a non-heated ultrasonic bath; add the cleaning liquid and the printed part to a lockable container (5 minutes). Subsequent drying of the printed part with compressed air and controlling of the cleaning result. You can recognize sufficient cleaning by a matt surface. Shiny areas require spot cleaning with cleaning liquid and a brush. Keep cleaning until no more shiny areas occur.

!Complete the post-treatment of the printed parts quickly and adhere to the time specifications. Above all, avoid unnecessarily long baths in the cleaning fluid and long waiting times between the steps. Isopropanol (99%) and optiprint clean can be used as cleaning liquids (when using optiprint clean, a 2-minute post-cleaning with iso-propanol is always necessary. optiprint clean is not suitable for cleaning surfaces and devices).

5.3. Finishing

Light curing in the light curing device (otoflash G171, NK Optik) with nitrogen flooding:

OPTIPRINT PRODUKT	GUIDE 385/405/HR, I-B-T 385, TRAY 385/HR, GINGIVA, MOCK UP, MODEL, MODEL ALIGN, MATCH, LAVIVA, LUMINA	CLARA
LIGHT OUTPUT IN OTOFLASH G171	2x 2000 flashes (turn after 2000 flashes)	2x 3000 flashes (turn after 3000 flashes)
ALTERNATIVE LIGHT OUTPUT (200 W, 315-400 NM)	7 minutes	10 minutes

!The final properties as well as the final color depend on the light curing process.

When light curing in other light curing devices, a comparably high energy input must be ensured (200 W). This depends on the light source used and the exposure time (UVA radiation source (315-400 nm), 10 minutes).

6. Repairing

The dental parts produced with optiprint can be polished and repaired in the conventional way. Please see the relevant processing guide.

7. Delivering

Use the optiprint prevente (dentona AG) disinfectant bath, which is especially suitable for 3d printed parts, to ensure hygiene when handing over the medical device to your customer. Spectrum of activity: bactericidal including TBC, yeasticidal, limited virucidal (HIV, HBV, HCV, Sars-CoV-2).

Advice to the medical device patient (if applicable)

Removable custom-made items should be cleaned under running water after wearing. Don’t use any toothpaste or denture cleaner. Then keep dry until next wear.

8. Technical data

Please see the separate TDS (Technical Data Sheet).

9. Composition

Mixture of acrylate and methacrylate resins, photoinitiators (phosphine oxide), additives, dyes.

10. Storage

Close the bottle carefully after each use. Store optiprint products closed, dry and protected from light at 5-30 ° C. Even a small amount of light can trigger the polymerization. Do not use after the expiry date.

11. Disposal

Dispose of liquid printing resin according to official regulations. Must not be disposed of together with household waste. Do not empty into drains.

Information

dentona AG is certified in accordance with DIN EN ISO 13485 and operates an elaborate quality assurance system. All dentona products are thus free of material and manufacturing defects and always have flawless quality. Our user recommendations are based on a careful determination of practice-related values from our dental laboratory. Achieving the guaranteed quality and usability of the product requires strict compliance with all process steps specified in the product files in a two-tier quality control. The user is responsible for operation and proper use of the material. The medical products are intended for use as custom-made products by dental technicians in accordance with EU regulations.

Report all serious incidents (death, serious deterioration of health, serious risk to public health) that have occurred in connection with the medical device to the manufacturer and the competent authority of the member state.

If defects appear in the material within the warranty period, the user only has a right to replacement of the material. dentona AG is liable for neither losses nor damage caused by this material, irrespective of whether this damage is direct or indirect, and especially collateral or consequential damage, regardless of the legal basis. According to statutory regulations, dentona AG is liable for direct damage to the material based on a deliberate act or gross negligence by its legal representatives or executive staff and for personal damage. Any liability for the material and damage as a consequence of using it are excluded if the user has not observed the specified process steps.

Symbols

 Observe instruction for use (see www.dentona.de)*	 REF Article number
 Use by / Expiry date	 Attention: health risk
 Protect from sunlight	 CE mark
 May cause an allergic skin reaction	 Manufacturer
 Attention: irritant	 Temperature limitation
 Medical device	 Attention: Environmentally hazardous
 LOT LOT number	

Also part of the range optiprint zero: 3D printing resin for the additive manufacturing of cast objects. See the processing guide at www.dentona.de

dentona AG

Otto-Hahn-Straße 27 | 44227 Dortmund

Phone: +49 (0) 231 - 55 56 - 0

Fax: +49 (0) 231 - 55 56 - 30