

3M Science.
Applied to Life.™

3M™ Impregum™ Super Quick Polyether Impression Material

Technical Data Sheet






3M™ Impregum™ Polyether Material. Proven benefits for accuracy you can trust.

With its unique material properties 3M™ Impregum™ Polyether Impression Material has always been the trusted solution for complex situations and implants.


In the past, the longer setting time of polyether has been one of the reasons to choose VPS materials for smaller cases.

New 3M™ Impregum™ Super Quick Polyether Impression Material – the first true superfast polyether - unites the unsurpassed accuracy of polyether with the speed of a VPS.



Superior moisture tolerance.

With its superior moisture tolerance, Impregum Polyether Material displaces blood and saliva right from the first contact with tissue for excellent, void-free detail reproduction.



Unmatched precision.

Excellent flow properties let dentists capture the finest details: Impregum polyether material flows smoothly, even in the toughest spots, producing void-free, high-precision impressions.



Take all the time you need.

Unlike VPS materials, polyether shows consistently great flow behavior. The indicated working time can reliably be used to full capacity – inaccuracies associated with premature setting are avoided. Once the working time is over, setting starts immediately.



Clinical Solutions & Products

3M™ Impregum™ Super Quick
Polyether Impression Material



Portfolio/indications
3M™ Impregum™ Polyether Materials



Related Products:

3M™ Astringent Retraction Paste



3M™ Intra-oral Syringe Purple
3M™ Garant™ Mixing Tips Purple



3M™ Impression Tray



3M™ Pentamix™ 3 Automatic Mixing Unit
3M™ Pentamix™ Lite Automatic Mixing Unit





3M™ Impregum™ Super Quick

Polyether Impression Material

Indications

Recommended for precision impression taking using the monophasic or 1-step technique for smaller cases like single-unit crowns, implants, bridges up to three units, inlays/onlays, veneers.

Features & benefits

Fast like a VPS.

The first superfast polyether with 45 seconds working time and quick 2-minute setting. Specially recommended for single-unit impressions.

100% Polyether. Trusted precision.

A brand new chemistry offers well-known polyether reliability and accuracy.

Patient comfort.

An improved taste together with a short time in the mouth significantly increase patient comfort.

Chemical background


The use of a new initiator required some changes in the material composition. The aziridino-polyether and triglycerides are the beating heart of every polyether impression material. These relatively small, but by no means simple changes led to the desired result: A faster setting reaction and a more neutral taste for increased productivity and more patient comfort.

Composition base/initiator paste +

Ordering information +

Video +

← Back to overview





← Clin

← Chemical background

3M™ Impregum™ Super Quick

Base paste

Aziridino-Polyether

Triglycerides

Fillers

Pigments

Aroma

NEW: Stabilizer

Accelerator

Plasticizers

Initiator paste

NEW: Initiator

NEW: Plasticizers

Fillers

NEW: Pigments



CLOSE

An improved taste and a short time in the mouth significantly increase patient comfort.

← Back to overview



Order code Introductory Kits for 3M™ Pentamix™ 3 Automatic Mixing Unit

- 69383 3M™ Impregum™ Penta™ Super Quick Medium Body**
1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Penta™ Cartridge MB, 10 3M™ Penta™ Mixing Tips – Red, 1 3M™ Penta™ Elastomer Syringe, 1 Polyether Tray Adhesive – 17 ml, 3 3M™ Astringent Retraction Paste capsules
- 69381 3M™ Impregum™ Penta™ Super Quick Medium Body/Light Body**
1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 2 cartridges LB – 50 ml each, 1 3M™ Penta™ Cartridge MB, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple, 5 3M™ Garant™ Intra-oral Tips white, 1 Polyether Tray Adhesive – 17 ml, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Impression Trays
- 69413 3M™ Impregum™ Penta™ Super Quick Heavy Body/Light Body**
1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 2 3M™ Impregum™ Super Quick LB cartridges – 50 ml each, 1 3M™ Impregum™ Penta™ HB Cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit, 10 3M™ Penta™ Mixing Tips – Red, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple, 5 3M™ Garant™ Intra-oral Tips white, 1 3M™ Polyether Tray Adhesive – 17 ml, 3 3M™ Astringent Retraction Paste capsules, 6 3M™ Impression Trays

Order code	Refills	Order code	Accessories
69385	3M™ Impregum™ Penta™ Super Quick Medium Body 2 base pastes – 300 ml each, 2 catalysts – 60 ml each	69387	3M™ Impregum™ Penta™ Super Quick Medium Body Cartridge
69379	3M™ Impregum™ Super Quick Light Body 4 cartridges – 50 ml each, 5 3M™ Intra-oral Syringes Purple, 5 3M™ Garant™ Mixing Tips purple	69401	3M™ Impregum™ Penta™ Super Quick Heavy Body Cartridge
69406	3M™ Impregum™ Penta™ Super Quick Heavy Body 2 base pastes HB – 300 ml each, 2 catalysts HB – 60 ml each	71454	3M™ Garant™ Mixing Tips purple 50 pieces

Order code Value Packs


- 69405 3M™ Impregum™ Penta™ Super Quick Medium Body Value Pack**
6 base pastes MB – 300 ml each, 6 catalyst MB – 60 ml each
- 69412 3M™ Impregum™ Penta™ Super Quick Heavy Body Value Pack**
6 base pastes HB – 300 ml each, 6 catalyst HB – 60 ml each



← Clin

← Video

3M™ Impregum™ Super Quick



The fastest way to polyether precision.

00:02:14

CLOSE ×

An improved taste and a short time in the mouth significantly increase patient comfort.

← Back to overview



← Clinical Solutions & Products

Portfolio overview



	Dispensing system	Viscosity low high	Impression technique	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
Super Quick Setting Materials					
3M™ Impregum™ Penta™ Super Quick (Medium Body)			Monophase or 1-step	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Heavy Body)			1-step	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)			1-step	0:45	2:00
Quick Setting Materials					
3M™ Impregum™ Penta™ H DuoSoft™ Quick (Heavy Body)			1-step	1:00	3:00
3M™ Impregum™ Penta™ Soft Quick (Medium Body)			Monophase	1:00	3:00
3M™ Impregum™ L DuoSoft™ Quick (Light Body)			1-step	1:00	3:00
Regular Setting Materials*					
3M™ Impregum™ Penta™ Soft (Medium Body)			Monophase	1:45	4:15
3M™ Impregum™ Penta™ (Medium Body)			Monophase	1:45	4:15
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body)			1-step	1:45	4:15
3M™ Impregum™ Penta™ L DuoSoft™ (Light Body)			1-step	1:45	4:15
3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)			1-step	1:45	4:15

Recommended indications



Indications	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Veneers	Bridge (> 3 units)	Implant (> 2 units)	Pick-up impression	Functional impression
Super Quick Setting Materials										
3M™ Impregum™ Penta™ Super Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ Super Quick (Medium or Heavy Body) 3M™ Impregum™ Super Quick (Light Body)	1-step	✓	✓	✓	✓					
Quick Setting Materials										
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ Quick (Heavy Body) 3M™ Impregum™ L DuoSoft™ Quick (Light Body)	1-step	✓	✓	✓	✓	✓				
Regular Setting Materials										
3M™ Impregum™ Penta™ Soft (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body) 3M™ Impregum™ Penta™ L DuoSoft™ (Light Body) 3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)	1-step					✓	✓	✓	✓	✓
3M™ Impregum™ Penta™ (Medium Body)	Monophase					✓	✓	✓	✓	✓



← Portfolio overview

3M™ Impregum™ Polyether Impression Material

	Dispensing system	Viscosity low high	Impression technique	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
Super Quick Setting Materials					
3M™ Impregum™ Penta™ Super Quick (Medium Body)			Monophase or 1-step	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Heavy Body)			1-step	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)			1-step	0:45	2:00
Quick Setting Materials					
3M™ Impregum™ Penta™ H DuoSoft™ Quick (Heavy Body)			1-step	1:00	3:00
3M™ Impregum™ Penta™ Soft Quick (Medium Body)			Monophase	1:00	3:00
3M™ Impregum™ L DuoSoft™ Quick (Light Body)			1-step	1:00	3:00
Regular Setting Materials*					
3M™ Impregum™ Penta™ Soft (Medium Body)			Monophase	1:45	4:15
3M™ Impregum™ Penta™ (Medium Body)			Monophase	1:45	4:15
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body)			1-step	1:45	4:15
3M™ Impregum™ Penta™ L DuoSoft™ (Light Body)			1-step	1:45	4:15
3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)			1-step	1:45	4:15

* Setting time from start of mixing 6:00 min.

CLOSE



← Recommended indications

3M™ Impregum™ Polyether Impression Material

Indications	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Veneers	Bridge (> 3 units)	Implant (> 2 units)	Pick-up impression	Functional impression
Super Quick Setting Materials										
3M™ Impregum™ Penta™ Super Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ Super Quick (Medium or Heavy Body) 3M™ Impregum™ Super Quick (Light Body)	1-step	✓	✓	✓	✓					
Quick Setting Materials										
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ Quick (Heavy Body) 3M™ Impregum™ L DuoSoft™ Quick (Light Body)	1-step	✓	✓	✓	✓	✓				
Regular Setting Materials										
3M™ Impregum™ Penta™ Soft (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body) 3M™ Impregum™ Penta™ L DuoSoft™ (Light Body) 3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)	1-step					✓	✓	✓	✓	✓
3M™ Impregum™ Penta™ (Medium Body)	Monophase					✓	✓	✓	✓	✓

CLOSE



3M™ Astringent Retraction Paste

Product description

The faster, gentler way to retract tissue and control bleeding.

Indications

Can be used in a healthy periodontium for indications such as:

- Taking impressions (material-based/digital)
- Preparation of temporary casts
- Preparation of class II and V fillings

Features & benefits

- Extra-fine tip delivers paste right into the sulcus
- Easy access to tight interproximal areas
- Long lasting hemostasis achieved from 15% aluminum chloride paste
- Less risk of tissue trauma and hemorrhage
- Hygienic unit-dose delivery
- Time-saving process and up to 50% faster (vs. cords)
- Gentle on tissue for more patient comfort
- Compatible with most composite dispensers

Ordering information



How to use 3M™
Astringent Retraction Paste





← Clin

← Ordering information

3M™ Astringent Retraction Paste

Order code

Product information

56944

25 capsules 3M™ Astringent Retraction Paste - 0.3 g each

56945

50 capsules 3M™ Astringent Retraction Paste - 0.3 g each

CLOSE

How-to-video 3M™
Astringent Retraction Paste



← Clin

← Video

3M™ Astringent Retraction Paste



How to use 3M™ Astringent Retraction Paste

00:01:13

CLOSE



it into
al areas
ved from

50%

nt comfort
te





3M™ Intra-oral Syringe Purple

Product description

Convenient wash material application with extreme accuracy and less waste.

Indications

Syringing preparations with polyether precision impression materials in accordance with ISO 4823 Type 3 and Type 2, provided the polyether impression materials are intended for syringing.

Features & benefits

- Smooth and easy access to the sulcus and hard-to-reach areas of the mouth
- Easy handling due to small, ergonomic design
- Hygienic single-use syringe:
No disinfection, no messy clean-up
- Less waste of impression material
- Can be (pre)filled directly from the hand dispenser cartridge up to 12 hours in advance

Exclusively designed for use with 3M™ Impregum™ Super Quick Polyether Material only!

3M™ Garant™ Mixing Tips Purple


Shorter and wider mixing tips which optimize the dispensing pressure and allow for more precise syringing.



Ordering information



How to use the 3M™ Intra-oral Syringe Purple



Material waste
3M™ Garant™ Mixing Tips Purple

Material waste
3M™ Intra-oral Syringe Purple



← Ordering information

3M™ Intra-oral Syringe Purple
3M™ Garant™ Mixing Tips Purple

Order code	Product information
71507	3M™ Intra-oral Syringes Purple, 20 pieces
71508	3M™ Intra-oral Syringes Purple, 50 pieces
71454	3M™ Garant™ Mixing Tips purple, 50 pieces
71461N	3M™ Garant™ Intra-oral Tips white, 50 pieces

CLOSE



Intra-oral Syringe Purple

can be prepared directly from the
hand dispenser cartridge up to
12 hours in advance



← Clin

← Video

3M™ Intra-oral Syringe Purple



How to use the 3M™ Intra-oral Syringe Purple

00:00:48

CLOSE ×

12 hours in advance



3M™ Impression Tray

Disposable Impression Tray

Product description

Easy to customize full arch tray that eliminates the need for a tray adhesive and simplifies both preliminary and final impression procedures.

Indications

- Precision impressions
- Implant impressions
- Provisional crown/bridge impressions
- Opposing jaw/study model impressions
- Orthodontic impressions
- Impressions for bleaching trays/mouth guards
- Fabrication of splints

Features & benefits

- Saves time with integrated self-retentive fleece strip and no need for a tray adhesive
- Higher accuracy thanks to optimized tray design
- Hygienic single-use for full infection control
- High patient comfort with design features to reduce gagging

Ordering information





← Clin

← Ordering information

3M™ Impression Tray

Order code	Product information
71615	10 trays, size S, upper jaw
71616	10 trays, size M, upper jaw
71617	10 trays, size L, upper jaw
71618	10 trays, size S, lower jaw
71619	10 trays, size M, lower jaw
71620	10 trays, size L, lower jaw

CLOSE





3M™ Pentamix™ 3 Automatic Mixing Unit 3M™ Pentamix™ Lite Automatic Mixing Unit

Product description

Adds speed and accuracy to the impression procedure through automixing.

Indications

Suitable for VPS and polyether impression materials.

Features & benefits

- Homogeneous, void-free and reproducible mixing quality
- Time-saving procedure at high speed
- Higher productivity
- Easy and intuitive operation
- Hygienic and clean workflow
- Robust, state-of-the-art design



Ordering information





← Ordering information

3M™ Pentamix™ 3 Automatic Mixing Unit

3M™ Pentamix™ Lite Automatic Mixing Unit

Order code	3M™ Pentamix™ 3 Automatic Mixing Unit	Order code	3M™ Pentamix™ Lite Automatic Mixing Unit
77871	3M™ Pentamix™ 3 Automatic Mixing Unit	77903	3M™ Pentamix™ Lite Automatic Mixing Unit – 230 V/50 Hz (INT) incl. 1 Penta™ Cartridge
77949	50 3M™ Penta™ Mixing Tips – Red	77906	3M™ Pentamix™ Lite Automatic Mixing Unit – 230 V/50 Hz (CH) incl. 1 Penta™ Cartridge
77601	3M™ Pentamix™ 3 Wall-Mount Kit	77944	3M™ Penta™ Cartridge for 3M™ Pentamix™ Lite Mixing Unit
69387	3M™ Impregum™ Penta™ Super Quick Medium Body Cartridge	77949	50 3M™ Penta™ Mixing Tips – Red
69401	3M™ Impregum™ Penta™ Super Quick Heavy Body Cartridge		
P3792	3M™ Impregum™ Penta™ Cartridge		
P3775	3M™ Impregum™ Penta™ Soft Quick Cartridge		
P3793	3M™ Impregum™ Penta™ H DuoSoft™ Cartridge		
P3773	3M™ Impregum™ Penta™ H DuoSoft™ Quick Cartridge		
P3795	3M™ Impregum™ Penta™ L DuoSoft™ Cartridge		

CLOSE



Indications & Clinical Application

Operatory guides

Impression tray preparation



Monophase technique



1-step full arch technique



1-step dual-arch technique



Lab handling guidelines for
polyether impressions



Case reports

1-step heavy/light body impression
Replacement of a partial crown



1-step medium/light body impression
Replacement of a partial crown



Dual-arch bite impression
Replacement of a all-ceramic crown



Monophase impression
Replacement of a single-unit crown



Closed implant impression
Production of a crown on implant




Recommended indications 3M™ Impregum™ Polyether Material

3M Science. Applied to Life.™		Recommended Indications 3M™ Impregum™ Polyether Impression Material									
Indications	Impression technique	Single-unit crown	Small implant (+ 2 units)	Majority (+ 2 units)	Small bridge (+ 3 units)	Veneers	Bridge (> 3 units)	Implant (+ 2 units)	Pick-up impression	Functional impression	
Super Quick Setting Materials											
3M™ Impregum™ Penta™ Super Quick (Medium Body)	Monophase										
3M™ Impregum™ Penta™ Super Quick (Medium or Heavy Body)	1-step	✓	✓	✓	✓						
3M™ Impregum™ Super Quick (Light Body)											
Quick Setting Materials											
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	Monophase										
3M™ Impregum™ Penta™ H DuoSoft™ Quick (Heavy Body)	1-step	✓	+	✓	✓						
3M™ Impregum™ L DuoSoft™ Quick (Light Body)											
Regular Setting Materials											
3M™ Impregum™ Penta™ Soft (Medium Body)	Monophase										
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body)	1-step					✓	✓	✓	✓	✓	
3M™ Impregum™ Penta™ L DuoSoft™ (Light Body)											
3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)											
3M™ Impregum™ Penta™ (Medium Body)	Monophase					✓	✓	✓	✓	✓	

3M Deutschland GmbH - ESPE Platz - 82229 Seefeld - Germany - info@3mimpregum.com
3M, DuoSoft®, Garant®, Impregum and Penta are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. © 3M 2016. All rights reserved. 3M.com/Dental

Case reports





← Indications

Recommended indications

Indications	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Veneers	Bridge (> 3 units)	Implant (> 2 units)	Pick-up impression	Functional impression
Super Quick Setting Materials										
3M™ Impregum™ Penta™ Super Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ Super Quick (Medium or Heavy Body) 3M™ Impregum™ Super Quick (Light Body)	1-step	✓	✓	✓	✓					
Quick Setting Materials										
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ Quick (Heavy Body) 3M™ Impregum™ L DuoSoft™ Quick (Light Body)	1-step	✓	✓	✓	✓	✓				
Regular Setting Materials										
3M™ Impregum™ Penta™ Soft (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body) 3M™ Impregum™ Penta™ L DuoSoft™ (Light Body) 3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)	1-step					✓	✓	✓	✓	✓
3M™ Impregum™ Penta™ (Medium Body)	Monophase					✓	✓	✓	✓	✓

CLOSE



← Case report

Indications & clinical application



1-step impression technique using 3M™ Impregum™ Super Quick Polyether Impression Material

00:02:25

CLOSE ×



Indi

Operat

Impres

Monop

1-step f

1-step c

Lab har
polyeth

Case re

1-step h

Replac

1-step r

Replac

Dual-ar

Replac

Monop

Replac

Closed implant impression

Production of a crown on implant



← Indications & Clinical Applications

Operatory guides

Impression tray preparation →

Monophase technique →

1-step full arch technique →

1-step dual-arch technique →

Lab handling guidelines for
polyether impressions →

Case reports

1-step heavy/light body impression
Replacement of a partial crown →

1-step medium/light body impression
Replacement of a partial crown →

Dual-arch bite impression
Replacement of a all-ceramic crown →

Monophase impression
Replacement of a single-unit crown →


Closed implant impression
Production of a crown on implant →

Impression tray preparation


3M™ Impression Tray +




Metal stock tray +



Custom tray +



Dual-arch tray +





← 3M™ Impression Tray

Operatory guides



3M™ Impression Trays are available in 3 sizes for upper and lower jaw. The included retention strip eliminates the need for a tray adhesive. The trays can be customized and adjusted if necessary with grinding or drilling instruments.

CLOSE



← Indi

Operat

Impress

Monop

1-step f

1-step c

Lab har
polyeth

Case re

1-step h

Replac

1-step r

Replac

Dual-ar

Replac

Monop

Replac

Closed implant impression
Production of a crown on implant



← 3M™ Impression Tray

Operatory guides




Due to the design of the 3M™ Impression Tray, less impression material is needed compared to conventional full-arch trays. Directed flow design minimizes flow defects and distal voids for more accurate impressions. Built-in reservoir prevents material overflow to help reduce gagging for better patient comfort.

CLOSE ✕



← Metal stock tray

Operatory guides



Step 1:
Clean and dry the metal tray after try-in.


CLOSE





← Metal stock tray

Operatory guides



Step 2:

Apply a thin homogeneous adhesive layer of uniform red color. Begin at one end of the tray and work towards the other. Let the adhesive fully dry. The required time is dependent on the thickness of the layer. Even for thick layers, two minutes is sufficient.



CLOSE





← Custom tray

Operatory guides



Step 1:

To ensure optimal adhesion remove smear layer and roughen surface.
Optional: Add retention elements.




CLOSE





← Custom tray

Operatory guides



Step 2:
Try-in.




CLOSE





← Custom tray

Operatory guides



Step 3:
Clean the tray after try-in. Thoroughly dry the tray.




CLOSE





← Custom tray

Operatory guides



Step 4:

Apply a thin, homogeneous adhesive layer of uniform red color, beginning from one end of the tray and working towards the other.




CLOSE





← Custom tray

Operatory guides



Step 5:

Let the adhesive fully dry. The required time is dependent on the thickness of the layer and may vary for different custom tray materials. Even for thick layers, five minutes should be sufficient.



CLOSE





← Custom tray

Operatory guides



Step 6:

Let a thin adhesive layer fully dry (at least two minutes) and apply a second thin layer. Let the second layer also dry for at least three minutes.


CLOSE





← Dual-arch tray

Operatory guides



Clean and dry the tray after try-in.
Please follow tray manufacturer's instructions. If required, apply a thin homogeneous adhesive layer of uniform red color at the frame.
In all cases, **Polyether Tray Adhesive must not be applied on the mesh.**

CLOSE





← Indications & Clinical Applications

Operatory guides

Impression tray preparation →

Monophase technique →

1-step full arch technique →

1-step dual-arch technique →

Lab handling guidelines for
polyether impressions →

Case reports

1-step heavy/light body impression
Replacement of a partial crown →

1-step medium/light body impression
Replacement of a partial crown →

Dual-arch bite impression
Replacement of a all-ceramic crown →

Monophase impression
Replacement of a single-unit crown →







Closed implant impression
Production of a crown on implant →

Monophase technique

3M™ Impregum™ Penta™ Super Quick Medium Body Polyether Impression Material

Click on images to enlarge

Gallery →





← Indications

← Monophase technique

Operatory guides

Operatory

Impression

Monophase

1-step for

1-step crown

Lab hand
polyether

Case re

1-step h

Replace

1-step n



Replace

Dual-an

Replace

Monoph

Replace



Step 1:

Insert the 3M™ Impregum™ Penta™ Super Quick cartridge into the 3M™ Pentamix™ 3 Automatic Mixing Unit. With the cartridge loaded, you are ready to mix impression material.




CLOSE





← Monophase technique

Operatory guides



Step 2:
Dispense a small amount of impression material on a mixing pad
to ensure a homogeneous mix.



CLOSE X





← Monophase technique

Operatory guides



Step 3:
Fill the 3M™ Penta™ Elastomer Syringe directly by placing the syringe onto the mixing tip and start dispensing the impression material.



CLOSE





← Monophase technique

Operatory guides



Step 4:

Fill the tray completely while keeping the mixing tip immersed in the material. Syringe the preparation simultaneously.

**Working time of 3M™ Impregum™ Penta™ Super Quick Material:
45 sec.**



CLOSE





← Monophase technique

Operatory guides



Step 5:
Place the filled tray into the mouth.
**Intra-oral setting time of 3M™ Impregum™ Penta™ Super Quick
Material: 2 min.**



CLOSE





← Monophase technique

Operatory guides



Step 6:
After 2 min. the impression can be removed from the mouth.

CLOSE





← Indications & Clinical Applications

Operatory guides

Impression tray preparation →

Monophase technique →

1-step full arch technique →

1-step dual-arch technique →

Lab handling guidelines for
polyether Impressions →

Case reports

1-step heavy/light body impression
Replacement of a partial crown →

1-step medium/light body impression
Replacement of a partial crown →

Dual-arch bite impression
Replacement of a all-ceramic crown →

Monophase impression
Replacement of a single-unit crown →

Closed implant impression
Production of a crown on implant →

1-step full arch technique












3M™ Impregum™ Penta™ Super Quick Heavy or Medium Body Polyether Impression Material

3M™ Impregum™ Super Quick Light Body Polyether Impression Material

3M™ Intra-oral Syringe Purple

Click on images to enlarge

Gallery →





← Indications

Operatory guides

Impression

Monophasic

1-step full arch

1-step dual arch

Lab handling
polyethylene

Case reports

1-step hybrid

Replacement

1-step non-removable

Replacement

Dual-arch



Replacement

Monophasic

Replacement

← 1-step full arch technique

Operatory guides



Step 1:
Insert the 3M™ Impregum™ Penta™ Super Quick Heavy or Medium Body Material cartridge into the 3M™ Pentamix™ 3 Automatic Mixing Unit. With the cartridge loaded, you are ready to mix impression material.


CLOSE ✕



Closed implant impression
Production of a crown on implant

Impression Material




Impression Material →





← 1-step full arch technique

Operatory guides



Step 2:

For the 3M™ Impregum™ Super Quick wash material cartridge either use the 3M™ Garant™ Mixing Tip Purple or the 3M™ Intra-oral Syringe Purple.




CLOSE





← 1-step full arch technique

Operatory guides



Step 3:

Dispense a small amount from the wash material cartridge on a mixing pad. Retain the cap of the cartridge to close it after usage.

CLOSE





← Indications

Operator guides

Impression

Monophasic

1-step full arch

1-step closed

Lab handling

polyethylene

Case re

1-step h

Replace

1-step n

Replace

Dual-arch



Replace

Monophasic

Replace

← 1-step full arch technique

Operator guides



Step 4:

If using the 3M™ Intra-oral Syringe Purple, remove the plunger and place the syringe directly on the wash material cartridge. Pay attention to the index wedge which guides you to the correct position of the syringe on the cartridge.

CLOSE ✕





Closed implant impression
Production of a crown on implant



Session Material




ry →





← 1-step full arch technique

Operatory guides



Step 5:

Load the 3M™ Intra-oral Syringe Purple with the required amount of wash material. Do not exceed the marking line (approx. 1 cm below the mixing tip joint). The loaded syringe can be pre-filled and stored up to 12 hours.



CLOSE





← 1-step full arch technique

Operatory guides



Step 6:

After having removed the 3M™ Intra-oral Syringe Purple from the wash material cartridge, insert the plunger and start pushing the paste no further than to the joint of the mixing tip.


CLOSE





← 1-step full arch technique

Operatory guides



Step 7:
Immediately before wash material application, activate the 3M™ Intra-oral Syringe Purple by bringing the mixing tip into the 180° straight position.

CLOSE ✕





← Indications

Operator guides

Impression

Monophasic

1-step full arch

1-step closed

Lab hard wax
polyethylene

Case re

1-step h

Replace

1-step n

Replace

Dual-ar




Replace

Monop

Replace

← 1-step full arch technique

Operator guides



Step 8:



Dispense a small amount of impression material from both the 3M™ Intra-oral Syringe Purple and the 3M™ Pentamix™ 3 Automatic Mixing Unit on a mixing pad to ensure a homogeneous mix.

CLOSE X



Impression Material




→





← 1-step full arch technique

Operatory guides



Step 9:

Fill the tray completely while keeping the mixing tip immersed in the material. Syringe the preparation simultaneously.

Working time of both 3M™ Impregum™ Super Quick tray and wash material: 45 sec.

CLOSE





← Indications

Operatory

Impression

Monophasic

1-step full

1-step crown

Lab hard
polyether

Case re

1-step h

Replace

1-step n

Replace

Dual-ar



Replace

Monop

Replace

← 1-step full arch technique

Operatory guides



Step 10:

Place the filled tray into the mouth.

**Intra-oral setting time of 3M™ Impregum™ Super Quick Material:
2 min.**


CLOSE ×



Closed implant impression
Production of a crown on implant

ession Material



ry →





← 1-step full arch technique

Operatory guides



Step 11:
After 2 min. the impression can be removed from the mouth.

CLOSE





← Indications & Clinical Applications

Operatory guides

Impression tray preparation



Monophase technique



1-step full arch technique



1-step dual-arch technique



Lab handling guidelines for
polyether impressions



Case reports

1-step heavy/light body impression
Replacement of a partial crown



1-step medium/light body impression
Replacement of a partial crown



Dual-arch bite impression
Replacement of a all-ceramic crown



Monophase impression
Replacement of a single-unit crown



Closed implant impression
Production of a crown on implant



1-step dual-arch technique












3M™ Impregum™ Penta™ Super Quick Heavy or Medium Body Polyether Impression Material

3M™ Impregum™ Super Quick Light Body Polyether Impression Material

3M™ Intra-oral Syringe Purple

Click on images to enlarge

Gallery →





← Indications

Operatory guides

Impression

Monophasic

1-step for

1-step crown

Lab hand
polyether

Case re

1-step h

Replace

1-step n

Replace

Dual-arch



Replace

Monophasic

Replace

← 1-step dual-arch technique

Operatory guides



Step 1:

Insert the 3M™ Impregum™ Penta™ Super Quick Heavy or Medium Body Material cartridge into the 3M™ Pentamix™ 3 Automatic Mixing Unit. With the cartridge loaded, you are ready to mix impression material.


CLOSE



Closed implant impression
Production of a crown on implant

Impression Material




Impression Material





← 1-step dual-arch technique

Operatory guides



Step 2:

For the 3M™ Impregum™ Super Quick wash material cartridge either use the 3M™ Garant™ Mixing Tip Purple or the 3M™ Intra-oral Syringe Purple.




CLOSE





← 1-step dual-arch technique

Operatory guides



Step 3:

Dispense a small amount from the wash material cartridge on a mixing pad. Retain the cap of the cartridge to close it after usage.

CLOSE





← Indications

Operator guides

Impression

Monophasic

1-step fabrication

1-step curing

Lab handling
polyethylene

Case re

1-step h

Replace

1-step n

Replace

Dual-arch



Replace

Monophasic

Replace

← 1-step dual-arch technique

Operator guides



Step 4:

If using the 3M™ Intra-oral Syringe Purple, remove the plunger and place the syringe directly on the wash material cartridge. Pay attention to the index wedge which guides you to the correct position of the syringe on the cartridge.

CLOSE





Closed implant impression
Production of a crown on implant



Impression Material




Operator guides





← 1-step dual-arch technique

Operatory guides



Step 5:

Load the 3M™ Intra-oral Syringe Purple with the required amount of wash material. Do not exceed the marking line (approx. 1 cm below the mixing tip joint). The loaded syringe can be pre-filled and stored up to 12 hours.



CLOSE





← 1-step dual-arch technique

Operatory guides



Step 6:

After having removed the 3M™ Intra-oral Syringe Purple from the wash material cartridge, insert the plunger and start pushing the paste no further than to the joint of the mixing tip.


CLOSE





← 1-step dual-arch technique

Operatory guides



Step 7:

Immediately before wash material application, activate the 3M™ Intra-oral Syringe Purple by bringing the mixing tip into the 180° straight position.

CLOSE





← Indications

Operatory guides

Impression

Monophasic

1-step fabrication

1-step curing

Lab handling
polyethylene

Case reports

1-step handling

Replacement

1-step preparation

Replacement

Dual-arch



Replacement

Monophasic

Replacement

← 1-step dual-arch technique

Operatory guides



Step 8:


Dispense a small amount of impression material from both the 3M™ Intra-oral Syringe Purple and the 3M™ Pentamix™ 3 Automatic Mixing Unit on a mixing pad to ensure a homogeneous mix.

CLOSE X



Impression Material

Impression Material





← Indications

Operatory guides

Impression

Monophasic

1-step fabrication

1-step curing

Lab handling
polyethylene

Case reports

1-step handling

Replacement

1-step handling

Replacement

Dual-arch




Replacement

Monophasic

Replacement

← 1-step dual-arch technique

Operatory guides



Step 9:

Fill both sides of the tray completely while keeping the mixing tip immersed in the material. Syringe the preparation simultaneously.



Working time of both 3M™ Impregum™ Super Quick tray and wash material: 45 sec.

CLOSE ✕



Session Material




ry →





← 1-step dual-arch technique

Operatory guides



Step 10:

Place the filled tray into the mouth.

**Intra-oral setting time of 3M™ Impregum™ Super Quick Material:
2 min.**



CLOSE





← 1-step dual-arch technique

Operatory guides



Step 11:
After 2 min. the impression can be removed from the mouth.

CLOSE





← Indications & Clinical Applications

Operatory guides

Impression tray preparation →

Monophase technique →

1-step full arch technique →

1-step dual-arch technique →

Lab handling guidelines for
polyether impressions →

Case reports

1-step heavy/light body impression
Replacement of a partial crown →

1-step medium/light body impression
Replacement of a partial crown →

Dual-arch bite impression
Replacement of a all-ceramic crown →

Monophase impression
Replacement of a single-unit crown →

Closed implant impression
Production of a crown on implant →

Lab handling guidelines for polyether impressions

Control →

Pouring the cast →

Implant impressions →

Plating →

Disinfection →

Improved removal →

Storage →

Cleaning metal trays →



← Indications & Clinical Applications

Lab handling guidelines for polyether impressions

Operatory guides

Impressions

Monophase

1-step f

1-step c

Lab han

polyeth

Case re

1-step h

Replacement of a partial crown

1-step medium/light body impression

Replacement of a partial crown

Dual-arch bite impression

Replacement of a all-ceramic crown

Monophase impression

Replacement of a single-unit crown

Closed implant impression

Production of a crown on implant



← Lab handling guidelines for polyether impressions

Control

Please inspect the surface quality of the impressions arriving in the lab, (no mushy, melted surfaces). Surfaces should appear detailed and crisp.

CLOSE



← Indications & Clinical Applications

Lab handling guidelines for polyether impressions

Operator guides

Impress

Monop

1-step f

1-step o

Lab han

polyeth

Case re

1-step h

Replac

1-step n

Replacement of a partial crown

Dual-arch bite impression

Replacement of a all-ceramic crown

Monophase impression

Replacement of a single-unit crown

Closed implant impression

Production of a crown on implant



← Lab handling guidelines for polyether impressions

Disinfection

Commercial water-based disinfectants indicated for dental impression materials can be used. Glutaraldehyde solutions are recommended for disinfection. Follow the manufacturer's instructions for use, time, and level of disinfection required. After disinfection, rinse impression well with water and blow dry. The impression should be kept dry when sent to the laboratory. Any temperatures exceeding 60°C/140°F will damage polyether impressions.

CLOSE



← Indications & Clinical Applications

Lab handling guidelines for polyether impressions

Operatory guides

Impress

Monop

1-step f

1-step o

Lab har

polyeth

Case re

1-step h

Replac

1-step medium/light body impression

Replacement of a partial crown

Dual-arch bite impression

Replacement of a all-ceramic crown

Monophase impression

Replacement of a single-unit crown

Closed implant impression

Production of a crown on implant

← Lab handling guidelines for polyether impressions

Pouring the cast

Surfactants (de-bubblizers) are not needed or recommended. Type 3 stone, type 4 or 5 die stone can be used to pour polyether impressions. When pouring polyether impressions with epoxy, or any urethane resin based materials, a separator must be used to prevent adhesion of the materials.

CLOSE ✕





← Indications & Clinical Applications

Lab handling guidelines for polyether impressions

Operator guides

Impress

Monop

1-step f

1-step o

Lab han

polyeth

Case re

1-step h

Replace

1-step n

Replace

Dual-arch bite impression

Replacement of a all-ceramic crown

Monophase impression

Replacement of a single-unit crown

Closed implant impression

Production of a crown on implant

← Lab handling guidelines for polyether impressions

Improved removal

Heating the impression to approx. 40°C/104°F after the cast has set will help in the removal of the cast. Block out or minimize undercuts if information is not needed. Isolating or blocking out the palatal area will help in removal of the model on a maxillary impression with a deep vault. Custom trays are not required during conventional crown and bridge procedures. If you choose to use custom trays, make sure that you have 2 - 3 mm material thickness per 1 mm of undercut. You may section or cut the custom tray if it is “locked-on” the model. Immersing the impression in hot water will soften the adhesive and allow for separation without damage.

CLOSE ✕





← Lab handling guidelines for polyether impressions

Implant impressions

Choosing the gingival mask material in combination with the proper separating material is important:

- a) 3M polyethers (Impregum™, Permadyne™ brands) can be used with common separators based on silicone, wax, or petrolatum (vaseline).
- b) C-silicones (GI-MASK®, Coltène®) can be used without a separator.
- c) Use caution with A-silicones (VPS), they can be used only with a separator that forms a sealed solid waxy layer (Sherasepal-U®, Shera®).

Separating materials can be applied with a brush.

CLOSE

Operate

Impress

Monop

1-step f

1-step c

Lab har

polyeth

Case re

1-step h

Replac

1-step n

Replac

Dual-ar

Replacement of a all-ceramic crown

Monophase impression

Replacement of a single-unit crown

Closed implant impression

Production of a crown on implant





← Indications & Clinical Applications

Lab handling guidelines for polyether impressions

Operatory guides

Impress

Monop

1-step f

1-step o

Lab har

polyeth

Case re

1-step h

Replac

1-step n

Replacement of a partial crown

Dual-arch bite impression

Replacement of a all-ceramic crown

Monophase impression

Replacement of a single-unit crown

Closed implant impression

Production of a crown on implant



← Lab handling guidelines for polyether impressions

Storage

Do not store polyether impressions in sealed plastic bags while still moist. Avoid direct sunlight exposure. Store impressions dry and below 30 °C/86 °F in the dark. Polyether impressions can be poured anytime within 14 days after they have been made if stored properly.

CLOSE



← Indications & Clinical Applications

Lab handling guidelines for polyether impressions

Operatory guides

Impression tray preparation →

Monop...

1-step f...

1-step o...

Lab han...
polyeth...

Case re...

1-step h...

Replace...

1-step medium/light body impression
Replacement of a partial crown →

Dual-arch bite impression
Replacement of a all-ceramic crown →

Monophase impression
Replacement of a single-unit crown →

Closed implant impression
Production of a crown on implant →



← Lab handling guidelines for polyether impressions

Plating

Polyether impressions can be silver-plated, but cannot be copper-plated.


CLOSE ×



← Lab handling guidelines for polyether impressions

Cleaning metal trays

The adhesive layer can be removed from the metal tray with a finger. Peeling under water allows for a fast and easy removal. Alternatively, commercially available agents for cleaning impression trays (e.g. Pluraclean Orange (Pluradent) or Solitine (Kerr)) may be used.



CLOSE



← Indications & Clinical Applications

Operatory guides

Impression tray preparation



Monophase technique



1-step full arch technique



1-step dual-arch technique



Lab handling guidelines for
polyether impressions



Case reports

1-step heavy/light body impression
Replacement of a partial crown



1-step medium/light body impression
Replacement of a partial crown



Dual-arch bite impression
Replacement of a all-ceramic crown



Monophase impression
Replacement of a single-unit crown



Closed implant impression
Production of a crown on implant



Replacement of a partial crown

1-step impression technique using 3M™ Impregum™ Super Quick Polyether Impression Material (Heavy Body/Light Body)

Clinical dentistry and photography by: Dr. Gunnar Reich, Munich, Germany
Lab work by: MDT Nicole Freund, Tutzing, Germany

About the case:

Patient needed a replacement for an insufficient and perforated partial crown on first lower molar #46. Impression was made with 3M™ Impregum™ Penta™ Super Quick Heavy Body and 3M™ Impregum™ Super Quick Light Body Polyether Impression Material in the 1-step technique.

Click on images to enlarge

Gallery





- ← Indications
- Operational
- Impressions
- Monophasic
- 1-step fast
- 1-step dual
- Lab hand
- polyether
- Case reports
- 1-step heavy
- Replacement
- 1-step normal
- Replacement
- Dual-arch
- Replacement
- Monophasic
- Replacement

1-step heavy/light body impression

← Replacement of a partial crown

Case reports



Initial situation: Insufficient and perforated gold alloy partial cast crown on first lower molar #46.

CLOSE



- Closed implant impression
- Production of a crown on implant



1-step heavy/light body impression

Replacement of a partial crown

Case reports



Tooth after minimal preparation with persistent discolored dentin; distal moisture control of subgingival margin achieved by 3M™ Astringent Retraction Paste.

CLOSE





1-step heavy/light body impression

← Replacement of a partial crown

Case reports



Situation prior to impression taking with clean and dry operative field.

CLOSE

← Indications

Operative

Impression

Monophasic

1-step fast

1-step dual

Lab hand
polyether

Case reports

1-step heavy

Replacement

1-step normal

Replacement

Dual-cure

Replacement

Monophasic

Replacement

Closed implant impression


Production of a crown on implant



lick

wn
nta™
polyether

lery →






1-step heavy/light body impression

Replacement of a partial crown

Case reports



Application of 3M™ Impregum™ Super Quick Light Body Polyether Impression Material using 3M™ Intra-oral Syringe Purple enabling easy interproximal access and void-free application of the wash material.

CLOSE

Indic

Operate

Impress

Monop

1-step f

1-step o

Lab har
polyeth

Case re

1-step h

Replac

1-step n

Replac

Dual-ar

Replac

Monop

Replac


Closed implant impression
Production of a crown on implant



ick

wn
nta™
polyether

lery →






1-step heavy/light body impression

← Replacement of a partial crown

Case reports



Detailed view of final impression showing superior detail reproduction of 3M™ Impregum™ Super Quick Polyether Impression Material. The color contrast ensures easy readability.

CLOSE

← Indications

Operational

Impression

Monophasic

1-step fast

1-step dual

Lab hand
polyether

Case reports

1-step heavy

Replacement

1-step normal

Replacement

Dual-arch

Replacement

Monophasic

Replacement

Closed implant impression


Production of a crown on implant



lick

wn
nta™
polyether

lery →






1-step heavy/light body impression

Replacement of a partial crown

Case reports



Final all-ceramic partial crown cemented with 3M™ RelyX™ Unicem 2 Self-Adhesive Universal Resin Cement.

CLOSE





← Indications & Clinical Applications

Operatory guides

[Impression tray preparation](#) →

[Monophase technique](#) →

[1-step full arch technique](#) →

[1-step dual-arch technique](#) →

[Lab handling guidelines for polyether impressions](#) →

Case reports

[1-step heavy/light body impression
Replacement of a partial crown](#) →

[1-step medium/light body impression
Replacement of a partial crown](#) →

[Dual-arch bite impression
Replacement of a all-ceramic crown](#) →

[Monophase impression
Replacement of a single-unit crown](#) →

[Closed implant impression
Production of a crown on implant](#) →

Replacement of a partial crown

1-step impression technique using 3M™ Impregum™ Super Quick Polyether Impression Material


Clinical dentistry and photography by: Dr. Gunnar Reich, Munich, Germany
Lab work by: MDT Nicole Freund, Tutzing, Germany

About the case:

Patient with an insufficient partial crown on first lower molar #36: Impression made with 3M™ Impregum™ Penta™ Super Quick Medium Body and 3M™ Impregum™ Super Quick Light Body Polyether Impression Material in the 1-step impression technique.

Click on images to enlarge

[Gallery](#) →






1-step impression technique

Replacement of a partial crown

Case reports



Initial situation: Insufficient and perforated gold alloy partial cast crown on first lower molar #36.

CLOSE






1-step impression technique

← Replacement of a partial crown

Case reports



Preparation with a partially subgingival preparation line.

CLOSE







1-step impression technique

Replacement of a partial crown

Case reports



Interproximal moisture control of subgingival margin by using 3M™ Astringent Retraction Paste.

CLOSE






1-step impression technique

Replacement of a partial crown

Case reports



Impression was made using the 1-step technique: 3M™ Impregum™ Super Quick Light Body Polyether Impression Material was applied using 3M™ Intra-oral Syringe Purple enabling easy interproximal access and void-free application of the wash material.

CLOSE





1-step impression technique

← Replacement of a partial crown

Case reports



Final 3M™ Impregum™ Super Quick Material impression made with a 3M™ Impression Tray. Due to its self-retentive fleece strip a separate tray adhesive is not required.

CLOSE





1-step impression technique

← Replacement of a partial crown

Case reports



Final glass ceramic restoration cemented with 3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement.

CLOSE





← Indications & Clinical Applications

Operatory guides

[Impression tray preparation](#) →

[Monophase technique](#) →

[1-step full arch technique](#) →

[1-step dual-arch technique](#) →

[Lab handling guidelines for polyether impressions](#) →

Case reports

[1-step heavy/light body impression
Replacement of a partial crown](#) →

[1-step medium/light body impression
Replacement of a partial crown](#) →

[Dual-arch bite impression
Replacement of an all-ceramic crown](#) →

[Monophase impression
Replacement of a single-unit crown](#) →

[Closed implant impression
Production of a crown on implant](#) →

Replacement of an all-ceramic crown

Dual-arch bite impression technique using 3M™ Impregum™ Super Quick Polyether Impression Material


Clinical dentistry and photography by: Dr. Gunnar Reich, Munich, Germany
Lab work by: MDT Nicole Freund, Tutzing, Germany

About the case:

Patient with a chipped all-ceramic crown on second molar: Impression made with 3M™ Impregum™ Super Quick Polyether Impression Material (Medium and Light Body) using the dual-arch bite impression technique.

Click on images to enlarge

[Gallery](#) →






Dual-arch bite impression technique

Replacement of an all-ceramic crown

Case reports



Initial situation: Patient arrived with a largely chipped veneered zirconia crown on second lower molar.

CLOSE






Dual-arch bite impression technique

Replacement of an all-ceramic crown

Case reports



3M™ Astringent Retraction Paste was used for moisture management and slight displacement of gingiva.

CLOSE






Dual-arch bite impression technique

Replacement of an all-ceramic crown

Case reports



For single-unit indications, the dual-arch bite impression technique is a viable alternative to full arch impressions. 3M™ Impregum™ Super Quick Light Body Polyether Impression Material was applied using 3M™ Intra-oral Syringe Purple.

CLOSE





Dual-arch bite impression technique

Replacement of an all-ceramic crown

Case reports



In parallel, the tray was loaded with 3M™ Impregum™ Penta™ Super Quick Medium Body Polyether Impression Material. The impression was done in the dual-arch bite (closed mouth) impression technique.

CLOSE





Dual-arch bite impression technique

Replacement of an all-ceramic crown

Case reports



Final situation: Precisely fitting monolithic 3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia restoration in situ. The crown was cemented with 3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement.

CLOSE





← Indications & Clinical Applications

Operatory guides

[Impression tray preparation](#) →

[Monophase technique](#) →

[1-step full arch technique](#) →

[1-step dual-arch technique](#) →

[Lab handling guidelines for polyether impressions](#) →

Case reports

[1-step heavy/light body impression
Replacement of a partial crown](#) →

[1-step medium/light body impression
Replacement of a partial crown](#) →

[Dual-arch bite impression
Replacement of a all-ceramic crown](#) →

[Monophase impression
Replacement of a single-unit crown](#) →

[Closed implant impression
Production of a crown on implant](#) →

Replacement of a single-unit molar crown

Monophase impression using 3M™ Impregum™ Penta™ Super Quick Polyether Impression Material

Clinical dentistry and photography by: Dr. Gunnar Reich, Munich, Germany
Lab work by: MDT Nicole Freund, Tutzing, Germany

About the case:

Patient with a fractured crown on second upper molar that needed replacement. Impression was made with 3M™ Impregum™ Penta™ Super Quick Medium Body Polyether Impression Material in the monophase technique.

Click on images to enlarge

[Gallery](#) →






Monophase impression

Replacement of a single-unit molar crown

Case reports



Initial situation: Second upper molar needs a crown after endodontic treatment. Insufficient filling on adjacent tooth needs to be renewed.

CLOSE




Closed implant impression
Production of a crown on implant



Monophase impression

Replacement of a single-unit molar crown

Case reports



In order to prepare the situation for impression taking, 3M™ Astringent Retraction Paste was used for moisture management.

CLOSE






Monophase impression

Replacement of a single-unit molar crown

Case reports



Impression was taken in the monophase technique using 3M™ Impregum™ Penta™ Super Quick Medium Body Polyether Impression Material. Syringing of the impression material around preparation with the 3M™ Penta™ Elastomer Syringe while the tray was filled simultaneously.

CLOSE







Monophase impression

Replacement of a single-unit molar crown

Case reports



Final impression made with a single-use 3M™ Impression Tray. Due to its self-retentive fleece strip a separate tray adhesive is not required. The impression shows crisp details and a clear circumferential mapping of the preparation line.

CLOSE





Monophase impression

Replacement of a single-unit molar crown

Case reports



Final situation: Monolithic 3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia crown on the second molar and a renewed filling on the first molar.

CLOSE

← Indications

Operational

Impression

Monophase

1-step filling

1-step crown

Lab handling
polyethylene

Case reports

1-step hand

Replacement

1-step molar

Replacement

Dual-cure

Replacement

Monophase

Replacement

Closed implant impression
Production of a crown on implant



Quick

ment.
n Body

lery →





← Indications & Clinical Applications

Operatory guides

Impression tray preparation →

Monophase technique →

1-step full arch technique →

1-step dual-arch technique →

Lab handling guidelines for
polyether impressions →

Case reports

1-step heavy/light body impression
Replacement of a partial crown →

1-step medium/light body impression
Replacement of a partial crown →

Dual-arch bite impression
Replacement of a all-ceramic crown →

Monophase impression
Replacement of a single-unit crown →

Closed implant impression
Production of a crown on implant →

Production of a crown on implant

Closed implant impression technique using 3M™ Impregum™
Penta™ Super Quick Medium Body Polyether Impression Material


Clinical dentistry and photography by: Dr. Gunnar Reich, Munich, Germany
Lab work by: MDT Nicole Freund, Tutzing, Germany

About the case:

After healing time, patient needed final prosthetics work on implant in first pre-molar regio. Impression made with 3M™ Impregum™ Penta™ Super Quick Medium Body Polyether Impression Material in the closed snap-on implant (monophase) impression technique.

Click on images to enlarge

Gallery →






Closed implant impression

← Production of a crown on implant

Indications & Clinical Applications



Initial situation: Implant (regio first premolar) six months after implantation. Healing cap in situ.



CLOSE X






Closed implant impression

← Production of a crown on implant

Indications & Clinical Applications



Syringing of 3M™ Impregum™ Penta™ Super Quick Medium Body Material around the impression coping with the 3M™ Penta™ Elastomer Syringe.

CLOSE






Closed implant impression

← Production of a crown on implant

Indications & Clinical Applications



Impression coping securely fixed in 3M™ Impregum™ Penta™ Super Quick Medium Body Material. The impression was taken using the monophasic technique.

CLOSE





Closed implant impression

← Production of a crown on implant

Indications & Clinical Applications



Final veneered all-ceramic crown cemented on implant abutment.

CLOSE

← Indi

Operat

Impres

Monop

1-step f

1-step c

Lab har

polyeth

Case re

1-step h

Replac

1-step r

Replac

Dual-ar

Replac

Monop

Replac

Closed implant impression

Production of a crown on implant



erial

pre-
Medium
(base)


ery





Scientific Data

3M™ Impregum™ Super Quick - Testimonial by Dr. C. E. Sabrosa 

Constant flow behavior 

Moisture tolerance 

Super fast setting 

Precision 

Improved taste 

Mechanical properties 



← **3M™ Impregum™ Super Quick**
Testimonial by Dr. C. E. Sabrosa*

Scientific Data



00:02:02

CLOSE

*Associate Professor of the Department of Restorative Dentistry and the Director of the Center for Technology Development at the University of the State of Rio de Janeiro Dental School, Brazil.



← Scientific Data

Constant flow behavior

Consistency Of Light Body Impression
Materials During Working Time



B. Kuppermann et al.

Flowability Throughout Working Time



Dr. C. E. Sabrosa

Constant flow behavior

An impression material should remain flowable throughout the entire working time. If the impression material begins to set during this critical time, an inaccurate impression may result.

Constant flow behavior is even more important in fast setting materials, where every second of the working time counts.





← Scientific Data

Constant flow behavior

Consistency Of Light Body Impression
Materials During Working Time
B. Kuppermann et al.



Flowability Throughout Working Time
Dr. C. E. Sabrosa



Constant flow behavior

Consistency Of Light Body Impression Materials During Working Time *B. Kuppermann, J. Zech, B. Cerny*

Published in: J Dent Res 96 (Spec Iss B):0310 (CED), 2017

In cases where the maximum working time is needed, dentists have to count on the flowability of the impression material even at the end of the working time. Flowability of precision impression materials can be determined by the value of consistency as described in ISO 4823. Objective of the study was a comparison of the consistency at 25 seconds and at the end of the working time according to the manufacturers' instructions.

Only 3M™ Impregum™ Super Quick Polyether Impression Material (named XPE in the study) shows constant consistency over the given working time. All other VPS impression materials tested changed their consistency significantly. The authors stated that the constant flow properties of the polyether material throughout the whole working time may result in an easier and more reliable handling for the clinicians and a better quality of dental impressions.

Materials tested



Method description



Detailed results





← Materials tested

Constant flow behavior

Impression Wash Material (Manufacturer)	Abbreviation
Experimental Polyether Material (3M)	XPE
Aquasil® Ultra XLV Fast Set (Dentsply Sirona)	AQXL
Aquasil® Ultra LV Fast Set (Dentsply Sirona)	AQLV
Take1® Advanced Light Body Fast Set (Kerr)	T1LB
EXAFAST™ NDS Injection Type (GC)	EXA
Honigum Pro Light Fast (DMG)	HON
Flexitime® Light Flow (Heraeus)	FLEX

CLOSE





← Method description

Constant flow behavior



Viscosity test showing the constant flowability of polyether throughout its working time.

00:01:08

CLOSE



ing Time

nt on the Flowability
ncy as
nsistency
cturers'

XPE in
er VPS
thors
ut the
ne






← Detailed results

Constant flow behavior

Consistency of Light Body/Wash materials



line indicates the lower ISO standard consistency limit for light body materials

CLOSE



← Scientific Data

Constant flow behavior

Consistency Of Light Body Impression
Materials During Working Time
B. Kuppermann et al.



Flowability Throughout Working Time
Dr. C. E. Sabrosa



Constant flow behavior

3M™ Impregum™ Super Quick Light Body: Flowability Throughout Working Time

The Shark Fin Test revealed that 3M™ Impregum™ Super Quick Polyether Impression Material shows constant viscosity throughout the whole working time. The indicated working time can actually be used to full capacity and avoids problems associated with premature setting (e.g. inaccuracies at the margins, less accurate reproduction of details).



Method description



Background



Shark Fin Test results





← Method description

Constant flow behavior



Shark Fin Test

00:01:04

CLOSE ×



g time.
problems
accurate





Literature search shows the test in use since 2001.

- Kim M.S., Doherty E.H., Kugel G. “Flow under pressure of four impression materials using shark fin device“. Published in: J Dent Res. 80 AADR Abstracts #624 (2001)
- Richter B., Kuppermann B., Führer C., Klettke T. “Flow Properties of Light Bodied Impression Materials During Working Time“. CED/NOF/IADR 2004 #142.
- Benchimol, Perry R., Kugel G., Hallas M. “Flow of Eight Impression Materials with 2mm slit after 25sec“. IADR 2005 #3083



g time.
problems
curate





← Shark Fin Test results

Constant flow behavior



The Shark Fins are almost the same height, no matter the stage of the working time.

Shark Fin Test results

00:01:33

Courtesy of Dr. Carlos Sabrosa

Associate Professor of the Department of Restorative Dentistry and the Director of the Center for Technology Development at the University of the State of Rio de Janeiro Dental School, Brazil.

CLOSE





← Scientific Data

Moisture tolerance

Material science:
Intrinsic Hydrophilicity



Pairwise Comparison of the Initial
Hydrophilicity of Impression Materials
S. Gaudet et al.



Initial Hydrophilicity of Light Body
Impression Materials
J. Zech et al.



Moisture tolerance

In the unset stage, where the impression material is in touch with the moist oral environment, polyether materials ensure a precise and detailed reproduction of surfaces where it really matters.





← Scientific Data

Moisture tolerance

Material science:
Intrinsic Hydrophilicity



Pairwise Comparison of the Initial
Hydrophilicity of Impression Materials
S. Gaudet et al.



Initial Hydrophilicity of Light Body
Impression Materials
J. Zech et al.




Moisture tolerance

Material science: Intrinsic Hydrophilicity

Polyether is hydrophilic by nature, which is due to its molecular structure. Like water molecules, polyether molecules contain different polar groups which chemically attract and interact with each other (hydrogen bonding). This phenomenon occurs with all bipolar molecules, while VPS materials are intrinsically hydrophobic. Hydrophobic methyl groups hinder interaction of polar water with the siloxane backbone.

What does intrinsic hydrophilicity mean?




Start animation





← What does intrinsic hydrophilicity mean?

Moisture tolerance



What does intrinsic hydrophilicity mean?

00:01:08

CLOSE

← Science

Moisture

Material
Intrinsic

Pairwise
Hydrophilicity
S. Gauc

Initial H
Impress
J. Zech



like water
ically
n occurs
. Hydro-
ckbone.





← Scientific Data

Moisture tolerance

Material science:
Intrinsic Hydrophilicity

Pairwise Comparison of the Initial
Hydrophilicity of Impression Materials
S. Gaudet et al.

Initial Hydrophilicity of Light Body
Impression Materials
J. Zech et al.

Moisture tolerance

Pairwise Comparison of the Initial Hydrophilicity of Impression Materials

S. Gaudet, S. Murali, S. Pagni, G. Kugel

Published in: J Dent Res 96 (Spec Iss A): 1618 (AADR), 2018

When an impression material gets in touch with the moist oral environment, the material needs to ensure a precise and void-free reproduction of the surfaces. In this study, the hydrophilicity of impression materials were directly compared in a side-by-side setup to find out if the new 3M™ Impregum™ Super Quick Polyether Impression Materials are more hydrophilic than VPS materials.

The study confirmed that the new Impregum Super Quick Materials show a superior hydrophilicity compared to all VPS materials tested.

Materials tested

Method description

Results



← Materials tested

Moisture tolerance

Material (Manufacturer)	Type
3M™ Impregum™ Penta™ Super Quick Medium Body (3M) mentioned as experimental material	Medium bodied polyether
Aquasil® Ultra+ Medium Fast Set (Dentsply Sirona)	Medium bodied VPS
Aquasil® Ultra+ Medium Fast Set DECA (Dentsply Sirona)	Medium bodied VPS
Aquasil® Ultra Monophase DECA (Dentsply Sirona)	Medium bodied VPS
EXAFast™ NDS Monophase (GC)	Medium bodied VPS
Flexitime® Monophase Dynamix (Kulzer)	Medium bodied VPS
Honigum Mixstar Mono (DMG)	Medium bodied VPS
Take 1™ Advanced™ Medium Fast Set (Kerr)	Medium bodied VPS
3M™ Impregum™ Super Quick Light Body (3M) mentioned as experimental material	Light bodied polyether
Aquasil® Ultra LV Fast Set (Dentsply Sirona)	Light bodied VPS
Aquasil® Ultra+ LV Fast Set (Dentsply Sirona)	Light bodied VPS
Aquasil® Ultra XLV Fast Set (Dentsply Sirona)	Light bodied VPS
Aquasil® Ultra+ XLV Fast Set (Dentsply Sirona)	Light bodied VPS
EXAFast™ NDS Injection Type (GC)	Light bodied VPS
Flexitime® Light Flow (Heraeus Kulzer)	Light bodied VPS
Panasil® Initial Contact Light (Kettenbach)	Light bodied VPS
Take 1™ Advanced™ Light Body Wash Fast Set (Kerr)	Light bodied VPS


CLOSE





← Method description

Moisture tolerance



Wettability comparison

00:01:04

CLOSE ×

← Scie

Moistur

Material
Intrinsic

Pairwise
Hydroph
S. Gaud

Initial H
Impress
J. Zech



, the
es. In
d in a
yether






← Results

Moisture tolerance

The tested polyether materials showed a larger horizontal spread and thus a higher hydrophilicity than the respective VPS material. In a control measurement with polyether material placed on both glass slides a comparable horizontal water spread was observed. More details with exact values can be found in the publication or the abstract thereof.



The sample video does not show the exact test but shows the behavior which was analyzed.

00:00:32

CLOSE



← Scientific Data

Moisture tolerance

Material science:
Intrinsic Hydrophilicity

Pairwise Comparison of the Initial
Hydrophilicity of Impression Materials
S. Gaudet et al.

Initial Hydrophilicity of Light Body
Impression Materials
J. Zech et al.

Moisture tolerance

Initial Hydrophilicity of Light Body Impression Materials
J. Zech, H. Hoffmann, H. Grupp, B. Cerny

Published in: J Dent Res 96 (Spec Iss A): 1007 (ADDR), 2018

The term hydrophilic is generally used to describe materials with a strong affinity for water or aqueous solutions. This property is crucial for impression materials as they are used in a moist oral environment and have to flow in close proximity to the teeth and gingival tissues in order to capture every important detail. In the study, hydrophilicity was determined by measuring initial water contact angles on surfaces of unset impression materials.

Both polyether materials, the existing 3M™ Impregum™ as well as the new 3M™ Impregum™ Super Quick Polyether Impression Materials showed a significantly higher initial hydrophilicity and thus the lowest contact angles compared to VPS materials.

Materials tested

Method description

Detailed results



← Materials tested

Moisture tolerance

Impression Material (Manufacturer)	Abbreviation
Honigum Pro Light Fast (DMG)	HPL
Panasil® Initial Contact Light (Kettenbach)	PIC
EXAFAST™ NDS Injection Type (GC)	ENI
Aquasil® Ultra LV Fast Set (Dentsply Sirona)	AUL
Aquasil® Ultra+ LV Fast Set (Dentsply Sirona)	AUP
3M™ Impregum™ Garant™ L DuoSoft™ (3M)	IMP
3M™ Impregum™ Super Quick Light Body (3M)	EXP

CLOSE




affinity for
s as they
the teeth
surfaces
3M™
ntly
VPS





← Method description

Moisture tolerance



Contact angle measurement

00:00:37

CLOSE ×



Affinity for
as they
the teeth

surfaces

3M™
ntly
VPS







← Detailed results

Moisture tolerance

Droplet placed on unset material after 20 sec. from start of mixing



Droplet placed on unset material after 45 sec. from start of mixing



CLOSE



affinity for
as they
the teeth
surfaces

3M™
ntly
VPS





← Scientific data

Super fast setting

Material science:
Polyether setting mechanism




Setting speed



Super fast setting

Immediate setting once the working time is over is a hallmark of polyether impression material. The rapid transition from the unset to the set state is exhibited by the new 3M™ Impregum™ Super Quick Polyether Impression Materials.





← Scientific data

Super fast setting

Material science:
Polyether setting mechanism



Setting speed




Super fast setting

Material science: Polyether setting mechanism

The goal in developing the new 3M™ Impregum™ Super Quick Polyether Impression Material was to offer an advanced formula Impregum Material with all the benefits of a polyether - combined with the speed of VPS.

The goal has been reached by developing a new initiator compound. The molecule is similar to the one used in the existing Impregum™ materials: It has the same reactive group and relies on an identical reaction mechanism called cationic ring-opening polymerization. The small but crucial difference lies in the substituents of the molecule. The increased reactivity of the molecule is the basis for a faster setting reaction. Compared to Impregum Penta (Soft) Materials, the intra-oral setting time is thus reduced by approximately 50%.

New initiator - video




Setting mechanism - video





← New initiator - video

Super fast setting



Initiator

00:00:32

CLOSE



← Setting mechanism - video

Super fast setting



Polymerization

00:01:28

CLOSE



← Scientific data

Super fast setting

Material science:
Polyether setting mechanism



Setting speed



Super fast setting


Setting speed

Compared to a well known VPS brand, 3M™ Impregum™ Penta™ Super Quick Medium Body Polyether Material shows no setting and a continuous flowable viscosity during working time. A quick setting reaction starts at the end of the working time, resulting in a steep increase in viscosity to a cured stage. At this point, the material can be removed without distortions.


Method description

A controlled oscillating paddle is placed in the material sample at intra-oral temperature. As the material sets, the change in the paddle's amplitude coincides with the material's increase in viscosity during setting.

Click on image to enlarge



Setting speed



- 3M™ Impregum™ Penta™ Super Quick Medium Body Polyether Material
- Competitor VPS Impression material

Lab measurements at 32 °C, values do not necessarily represent recommended clinical working and setting times.



← Scientific Data

Precision

Clinical Evaluation Of A Novel
Super Quick Setting Polyether
Impression Material
R. Hampe et al.



Digitization of Impression Materials
R. Hampe et al.



Torque Strength of Implant Copings
in Various Impression Materials
B. Cerny et al.



Precision

3M™ Impregum™ Polyether Impression Material is known for its precision and reliability. Even in difficult clinical situations, it is able to deliver highly accurate detail reproduction, resulting in precise fitting final restorations.





← Scientific Data

Precision

Clinical Evaluation Of A Novel Super Quick Setting Polyether Impression Material
R. Hampe et al.



Digitization of Impression Materials
R. Hampe et al.



Torque Strength of Implant Copings in Various Impression Materials
B. Cerny et al.



Precision

Clinical Evaluation Of A Novel Super Quick Setting Polyether Impression Material

R. Hampe, A. Dire, B. Cerny, R. Guggenberger

Published in: J Dent Res 97 (Spec Iss B): 2522 (IADR), 2018

In an in-office evaluation, more than 6,900 impressions were taken with 3M™ Impregum™ Super Quick Polyether Impression Materials by 482 dentists. The dentists' feedback confirms the very high clinical acceptance with the tested materials. The short intra-oral setting time and high impression accuracy were specifically mentioned as very desirable.

Materials tested

- 3M™ Impregum™ Penta™ Super Quick Medium Body
- 3M™ Impregum™ Penta™ Super Quick Heavy Body
- 3M™ Impregum™ Super Quick Light Body

Method description



Customer feedback





← Method description

Precision

The research was conducted via a web-based survey which was administered after a trial period of 5 weeks. Participants were asked to use the new 3M™ Impregum™ Super Quick Polyether Impression Materials in their clinical routine and to rate different performance criteria.

Sample size: 630 dentists from Europe and the U.S. had been pre-selected to take part in this evaluation. Field dates: first testing period has been from January to March, second testing period has been from October to December 2017, 482 respondents.

Survey scope: To assess dentists' satisfaction with Impregum Super Quick Polyether Impression Materials.

CLOSE

Method description

Customer feedback




← Customer feedback

Precision

Overall, 482 participants completed the questionnaire. Quantified responses for various features and characteristics of 3M™ Impregum™ Super Quick Polyether Impression Materials, especially regarding accuracy, were analyzed using descriptive statistics.


Detail reproduction of 3M™ Impregum™ Super Quick Polyether Impression Materials

Click on image to enlarge



Satisfaction with fit of final restorations using 3M™ Impregum™ Super Quick Polyether Impression Materials


Click on image to enlarge



CLOSE X



Detail reproduction of 3M™ Impregum™ Super Quick Polyether Impression Materials



CLOSE


← Customer feedback

Precision



Satisfaction with fit of final restorations using 3M™ Impregum™ Super Quick Polyether Impression Materials

Precision



CLOSE

Customer feedback





← Scientific Data

Precision

Clinical Evaluation Of A Novel
Super Quick Setting Polyether
Impression Material
R. Hampe et al.



Digitization of Impression Materials
R. Hampe et al.



Torque Strength of Implant Copings
in Various Impression Materials
B. Cerny et al.



Precision

Digitization of Impression Materials

R. Hampe, P. Osswald, M. Krämer, A. Dire, P. Athanasiou, G.P. Morris

Published in: J Dent Res 97 (Spec Iss A): 0264 (AADR), 2018

For an easy access to the digital CAD/CAM workflow, impressions may be scanned directly. In this study the variance of the scanning accuracy was addressed. Each impression and a master model were scanned using a commercially available LED multiline scanner. The absolute mean differences between both scans for scanning accuracy was extracted.

All direct impression scans provided a high and clinically suitable accuracy with the LED multiline scanner used. If desired, impressions made with 3M™ Impregum™ Polyether Impression Material can enter the digital workflow in-office or in the lab with direct scanning of the impression with modern desktop scanners (blue light technology).

Materials tested



Method description



Detailed results





← Materials tested

Precision

Impression Material (Manufacturer)	Abbreviation
Flexitime® Fast&Scan (Kulzer)	FFS-1step
3M™ Impregum™ DuoSoft™ Quick	IPDuoQ-1step
3M™ Impregum™ Penta™ Super Quick Heavy Body / 3M™ Impregum™ Super Quick Light Body (3M)	EXPHL-1step
3M™ Impregum™ Penta™ Super Quick Medium Body / 3M™ Impregum™ Super Quick Light Body (3M)	EXPML-1step
Flexitime® Monophase Pro Scan Dynamix (Kulzer)	FM-mono
3M™ Impregum™ Penta™ Soft Quick (3M)	IPSQ-mono
3M™ Impregum™ Penta™ Super Quick Medium Body (3M)	EXP-mono

CLOSE



scanned
Each
le LED
canning


with the
gum™
the lab
light





← Method description

Precision



Accuracy of impression material scans

00:00:52

CLOSE ×



scanned
Each
le LED
canning

with the
gum™
the lab
light





As proven in a 2-sample equivalence test, all tested impressions were statistically equivalent in their accuracy. The differences to the master model were below 50 microns in all cases.

CLOSE

Materials tested

Method description

Detailes results

Precision

Clinical
Super C
Impress
R. Ham

Digitiza
R. Ham

Torque
in Variou
B. Cern

scanned
Each
le LED
canning

with the
gum™
the lab
light



← Scientific Data

Precision

Clinical Evaluation Of A Novel
Super Quick Setting Polyether
Impression Material
R. Hampe et al.



Digitization of Impression Materials
R. Hampe et al.



Torque Strength of Implant Copings
in Various Impression Materials
B. Cerny et al.



Precision

Torque Strength of Implant Copings in Various Impression Materials *B. Cerny, J. Zech, B. Kuppermann*

Published in: J Dent Res 97 (Spec Iss A): 1009 (AADR), 2018

In order to achieve accurate fabrication of an implant-retained prosthesis, the impression must record an exact registration of the position and orientation of the implant. Aim of the study was to test the suitability of different impression materials for a secure transfer of implant copings.

Even if there were significant differences in the torque strength of the VPS impression materials, all 3M polyether materials clearly showed the highest torque strength and thus the lowest chance of accidental implant displacements in comparison to all tested VPS materials.

Materials tested



Method description



Detailed results





← Materials tested

Precision

Material (Manufacturer)	Type
3M™ Impregum™ Penta™ Super Quick Medium Body (3M)	Polyether
3M™ Impregum™ Penta™ (3M)	Polyether
3M™ Impregum™ Penta™ Soft (3M)	Polyether
3M™ Impregum™ Penta™ Soft Quick (3M)	Polyether
Aquasil® Ultra+ Medium FS DECA (Dentsply Sirona)	VPS
Aquasil® Ultra+ Medium FastSet (Dentsply Sirona)	VPS
Flexitime® Monophase Cartr. (Kulzer)	VPS
Monopren® transfer (Kettenbach)	VPS
Honigum Pro-Mono (DMG)	VPS
Take 1™ Advanced™ Medium Fast Set (Kerr)	VPS
EXAFAST™ NDS Garant Mono (GC)	VPS

CLOSE



Materials

he
of the
materials


torque





← Method description

Precision



Torque test

00:00:42

CLOSE ×



Materials

the of the materials


torque





← Detailed results

Precision



Torque strength of implant copings in various impression materials

CLOSE





← Scientific Data

Improved taste

Material science:
Why the taste was improved



In-vivo study: Taste and Comfort
Evaluation of Polyether Impression
Materials
A. Syrek et al.




Improved taste

In order to further improve patient comfort, a better taste and smell was on the wish list for a new 3M™ Impregum™ Polyether Impression Material.

Improved taste and smell were achieved through the development of a new initiator compound. The compound has the same reactive group as existing Impregum materials, while the improvement comes from new substituents.


New initiator - video





← New initiator - video

Improved taste



Initiator

00:00:32

CLOSE



← Scientific Data

Improved taste

Material science:
Why the taste was improved



In-vivo study: Taste and Comfort
Evaluation of Polyether Impression
Materials
A. Syrek et al.



Improved taste

In-vivo study: Evaluation of Taste and Comfort of Polyether
Impression Materials
A. Syrek, O. Brinkmann, B. Cerny, B. Rauch

Published in: J Dent Res 97 (Spec Iss B): 3309 (IADR), 2018

Randomized, cross-over, single blinded, internal 3M study assessing the taste and comfort of two polyether impression materials. A total of 70 subjects received one impression with each material in a randomized order in dual-arch trays. After setting of the materials and removal of the impressions from the mouth, subjects were asked to give their feedback on taste and comfort of each material on a visual analog scale (VAS).

Objectives:

1. Compare the taste of the new polyether impression material with standard polyether impression material
2. Compare the comfort of both impression materials

Material tested



Results



← Scientific Data

Improved taste



← Material tested

Improved taste

3M™ Impregum™ Penta™ Super Quick Medium Body
and
3M™ Impregum™ Penta™ Polyether Impression Materials

CLOSE

Material tested

Results



The new Impregum™ Penta™ Super Quick Medium Body Polyether Impression Material revealed statistically significant higher taste and comfort ratings in this study. It had a more pleasant taste and more patient comfort compared to Impregum Penta Impression Material.

CLOSE

Material tested

Results



← Mechanical properties

Scientific Data

Mechanical properties

	3M™ Impregum™ Super Quick Light Body Polyether Impression Material*	3M™ Impregum™ Penta™ Super Quick Medium Body Polyether Impression Material**	3M™ Impregum™ Penta™ Super Quick Heavy Body Polyether Impression Material***
Consistency A+B [mm] (ISO 4823:2015)	40.5	35.0	33.0
Recovery from deformation [%] (ISO 4823:2015)	99.0	97.6	97.9
Strain in compression (ISO 4823:2015)	2.9	2.8	2.0
Detail reproduction (ISO 4823:2015)	pass	pass	pass
Linear dimensional change [%] (ISO 4823:2015)	-0.38	-0.16	-0.12
Shore hardness after 15 min/24 h (ISO 7619-1)	54/55	50/54	57/60
Tensile strength after 24 h [MPa] (3M Internal SOP)	2.33	1.85	2.16
Toughness after 24 h [J] (3M Internal SOP)	0.70	0.79	1.00
Contact angle initial, on set material [°] (3M Internal SOP)	62	66	45
Working time at 23°C [min:sec] (according to IFU)	0:45	0:45	0:45
Intra-oral setting time [min:sec] (according to IFU)	2:00	2:00	2:00

*Batch No. #VT-PENG-Abf_G-0009

**Batch No. Base #VT-PENG-B-0008, Cat #VT-PENG-K-0007

***Batch No. Base #4202097, Cat #4230412

CLOSE



Downloads

Product information:

Product brochure



Product portfolio



Competitive comparisons



Product video



Operatory guides:

Monophase impression



1-step impression




1-step dual-arch impression




Lab handling guidelines for
polyether impressions



Ordering information:



Safety Data Sheet:





Competitive comparisons

Downloads

How does your material measure up?
Discover the difference.

- | | | |
|--------------------------------------|-----------|--|
| 3M™ Impregum™ Super Quick Light Body | VS | Identium® Light Fast Impression Material |
| 3M™ Impregum™ Super Quick Light Body | VS | Panasil® Initial Contact Light VPS Impression Material |
| 3M™ Impregum™ Super Quick Light Body | VS | Flexitime® Light Flow VPS Impression Material |
| 3M™ Impregum™ Super Quick Light Body | VS | Aquasil® Ultra XLV Fast Set VPS Impression Material |
| 3M™ Impregum™ Super Quick Light Body | VS | Aquasil® Ultra LV Fast Set VPS Impression Material |

CLOSE



FAQs


General

To see the answer click on the question.

General 

Indications and compatibility with
other impression procedure products 

Storage/preparation 

Application/working and setting time/
syringing 

Material properties 

Are you planning to substitute a current material with 3M™ Impregum™ Super Quick Polyether Impression Material?

Are there specifics to be considered regarding the allergy potential of Impregum Super Quick Material?

What safety data supports Impregum Super Quick Impression Material?



FAQs

 General

FAQs

Are you planning to substitute a current material with 3M™ Impregum™ Super Quick Polyether Impression Material?

No, with the launch of Impregum Super Quick Material 3M extends the polyether impression material portfolio. All current Impregum and Impregum Quick Materials will be kept in the portfolio.

CLOSE 



FAQs

 General

FAQs

Are there specifics to be considered regarding the allergy potential of 3M™ Impregum™ Super Quick Polyether Impression Material?

Based on the biocompatibility evaluation, which included an evaluation for the sensitization potential of the product, the sensitization potential of Impregum Super Quick Material was found to be low and not different from that of other 3M impression materials.

CLOSE 



FAQs

 General

FAQs

What safety data supports 3M™ Impregum™ Super Quick Polyether Impression Material?

Impregum Super Quick Material was evaluated using international standards for biocompatibility – specifically, the evaluation was structured according to ISO 10993 (Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process) and ISO 7405 (Dentistry – Evaluation of biocompatibility of medical devices used in dentistry). These standards provide a framework of toxicological endpoints that must be included in a product evaluation based on the nature and duration of body contact.

For impression materials, the toxicological endpoints include: cytotoxicity, irritation and sensitization (ability to cause an allergic reaction). The result of the biocompatibility evaluation, which was conducted by a board-certified toxicologist, indicates that the product is safe for its intended use.

CLOSE 



General



Indications and compatibility with
other impression procedure products



Storage/preparation



Application/working and setting time/
syringing



Material properties



FAQs

Indications and compatibility with other impression procedure products

To see the answer click on the question.

Is 3M™ Impregum™ Super Quick Material indicated for implant impressions?

Can I use Impregum Super Quick Material for larger cases?

What is our recommendation for Clear Tray Aligner impressions?
Can Impregum Super Quick Material be used?

What kinds of impression trays are suitable for Impregum Super Quick Impression Material?

Does 3M™ Astringent Retraction Paste influence the setting of impression materials?

Are specific limitations to consider when using Impregum Super Quick Material (retraction, gloves)?

Is Impregum Super Quick Material compatible with other dental materials such as cements and fillings?



FAQs



Indications and compatibility with other
impression procedure products

FAQs

Is 3M™ Impregum™ Super Quick Material indicated for implant impressions?

Yes, Impregum Super Quick Material is indicated for implant impressions.
There are no constraints but it is especially suited for single-unit implant cases.

CLOSE 

General

Indications
other in

Storage

Applica
syringin

Material

on

rials?

Material

Is Impregum Super Quick Material compatible with other dental materials such as
cements and fillings?



FAQs



Indications and compatibility with other
impression procedure products

FAQs

Can I use 3M™ Impregum™ Super Quick Material for larger cases?

Impregum Super Quick Material can be used for all cases where 45 seconds working time are sufficient. There is no contraindication. 3M recommends but not limits using Impregum Super Quick Material or 1 to 2 units.

CLOSE 

General

Indications
other in

Storage

Application
syringing

Material

on


materials?

Material

Is Impregum Super Quick Material compatible with other dental materials such as cements and fillings?



FAQs

 [Indications and compatibility with other impression procedure products](#)

FAQs

What is our recommendation for Clear Tray Aligner impressions? Can 3M™ Impregum™ Super Quick Material be used?

Impregum Super Quick Material is not recommended for Clear Tray Aligner impressions, but 3M™ Imprint™ 4 VPS Impression Material is excellently suited for Clear Tray Aligner impressions. The materials should be chosen depending on the technique the dentist wants to use.

CLOSE 

General

Indications
other in

Storage

Applica
syringin

Material


on

rials?

Material

Is Impregum Super Quick Material compatible with other dental materials such as cements and fillings?



 [Indications and compatibility with other impression procedure products](#)

FAQs

What kinds of impression trays are suitable for 3M™ Impregum™ Super Quick Material?

All impression trays, full arch, quadrant and dual-arch impression trays generally used for precision impressions are suitable (metal and plastic trays, non-perforated and perforated trays, stock and custom trays). Generally, 3M recommends the use of rigid trays. Usage of a polyether tray adhesive is recommended unless using the 3M™ Impression Tray where a retention fleece makes the usage of tray adhesives unnecessary.

CLOSE 



FAQs



Indications and compatibility with other
impression procedure products

FAQs

Does 3M™ Astringent Retraction Paste influence the setting of impression materials?

As with all astringents, Astringent Retraction Paste may inhibit the setting of 3M™ Impregum™ Super Quick Material. The paste must be completely rinsed off with air-water spray and suction prior to impression taking.

CLOSE 

General

Indications
other in

Storage

Applica
syringin

Material

on

rials?

Material

is impregum super quick material compatible with other dental materials such as
cements and fillings?



FAQs

Indications and compatibility with other impression procedure products

FAQs

Are there specific limitations to consider when using 3M™ Impregum™ Super Quick Material (retraction, gloves)?

Since Impregum Super Quick Material is based on proven 3M polyether chemistry, there are no specific limitations. Retraction solutions have to be rinsed with water thoroughly before applying Impregum Super Quick Material.

CLOSE 

General

Indications
other in

Storage

Application
syringin

Material


on

rials?

Material

Is Impregum Super Quick Material compatible with other dental materials such as cements and fillings?



 **Indications and compatibility with other
impression procedure products**

FAQs

Is 3M™ Impregum™ Super Quick Material compatible with other dental materials such as cements and fillings?

There are no specific effects expected for Impregum Super Quick Material, similar to other 3M polyether or VPS impression materials. Use caution with VPS gingival mask materials, they can be used only with a separator that forms a sealed solid waxy layer (Sherasepal-U®, Shera®). Separating materials can be applied with a brush, to insulate areas that are necessary.

CLOSE 



FAQs

Storage/preparation

General



Indications and compatibility with
other impression procedure products



Storage/preparation



Application/working and setting time/
syringing



Material properties



To see the answer click on the question.

During summer, it becomes very warm in my practice, what shall I do?

What should I do in winter times when 3M™ Impregum™ Super Quick Material has been stored below 18 °C/64 °F?

Which mixing devices can be used with Impregum Penta Super Quick Impression Material?

How to use the tray adhesive?

The mixing tips cannot be attached to the Impregum Super Quick Material cartridge. I observed leakage between mixing tip and cartridge?

Can I use the standard 3M™ Impregum™ Mixing Tips?

How many applications are possible with the material in one 3M™ Intra-oral Syringe Purple?

What materials can be syringed using the Intra-oral Syringe Purple?

Is Impregum Super Quick Material performance influenced by using the Intra-oral Syringe Purple?

Paste comes out between the 50 ml cartridge and the Intra-oral Syringe Purple while loading.

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



FAQs

← Storage/preparation

FAQs

During summer, it becomes very warm in my practice, what shall I do?

In general, the working time of polyether impression materials is shortened at higher temperature. In this case the materials should be stored at a cooler place (but not below 18°C/64°F) or a regular-set version instead of a super quick setting material should be used. The setting reaction of polyether is less temperature sensitive compared to a VPS. According to the instructions for use, 3M™ Impregum™ Super Quick Material as well as all other 3M impression materials should be stored at 18°C/64°F – 25°C/77°F.

CLOSE

Paste comes out between the 50 ml cartridge and the 3M™ Intra-oral Syringe Purple while loading.

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



← Storage/preparation

FAQs

What should I do in winter times when 3M™ Impregum™ Super Quick Material has been stored below 18 °C/64 °F?

- Do not use 3M™ Penta™ Impression Materials that have been stored at temperatures below 18°C/64°F, the viscosity of the pastes will increase to such an extent that there may be mixing problems in the unit. They can be used after storing for one day at a minimum temperature of 18°C/64°F. This will re-establish the usual handling characteristics without compromising quality.
- 50 ml cartridge materials: If the material is stored at temperatures below 18°C/64°F, it needs to be discarded. Precise dosing and mixing can no longer be guaranteed even if subsequently kept at room temperature for a longer period.

CLOSE 

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



← Storage/preparation

FAQs

Which mixing devices can be used with 3M™ Impregum™ Penta™ Super Quick Impression Material?

Manufacturer	Mixing device	3M™ Pentamix™ 2 cartridge	3M™ Pentamix™ 3 cartridge
Zhermack	Modulmix	✓	✓
DMG	Mix Star emotion	✓	✓
Heraeus Kulzer	Dynamix	✓	✗
Heraeus Kulzer	Dynamix Speed	✓	✓
Kettenbach	Sympress	✓	✓
Dentsply-Sirona	Duomix	✓	✗
Dentsply-Sirona	Duomix II	✓	✓

Heraeus Kulzer, Kettenbach and Dentsply-Sirona are using the Renfert device.

- ✗ Cartridge can be inserted; however the cover cannot be closed and therefore the mixing process cannot be started.
- ✓ Cartridge is fully compatible.

CLOSE

General
Indications
other in
Storage
Applica
syringin
Material

ion Material?
purple?
nge Purple?
while loading.
only into the
into the

How should I reclose the 50 ml cartridge?



FAQs

 [Storage/preparation](#)

FAQs

How to use the tray adhesive?

For sufficient adhesion, it is important to apply a thin layer of 3M™ Polyether Adhesive, at least 1 minute before you intend to fill the tray.

CLOSE

General

Indications
other in

Storage

Application
syringin

Material

ion Material?

purple?

ange Purple?

Paste comes out between the 50 ml cartridge and the 3M™ Intra-oral Syringe Purple while loading.

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



← Storage/preparation

FAQs

The white mixing tips cannot be attached to the 3M™ Impregum™ Super Quick Material cartridge. I observed leakage between mixing tip and cartridge.

Please use the 3M™ Garant™ Mixing Tips Purple only. The white mixing tips do not fit on the Impregum Super Quick Material cartridge. Similarly, the purple tips do not fit on the 3M™ Impregum™ L DuoSoft™ cartridges.



CLOSE

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?




← Storage/preparation

FAQs

Can I use the standard 3M™ Impregum™ Mixing Tips?

Please use the 3M™ Garant™ Mixing Tips Purple only. The white mixing tips do not fit on the 3M™ Impregum™ Super Quick Material cartridge.



CLOSE

General
Indications
other in
Storage
Applica
syringin
Material

ion Material?

purple?

nge Purple?

while loading.

only into the

into the

mixing up.

How should I reclose the 50 ml cartridge?




← Storage/preparation

FAQs

How many applications are possible with the material in one 3M™ Intra-oral Syringe Purple?

According to clinical evidence and depending on the quantity syringed around each unit, the amount of material in one syringe should be sufficient for about 2 – 4 units.



CLOSE

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



← Storage/preparation

FAQs

What materials can be syringed using the 3M™ Intra-oral Syringe Purple?

The Intra-oral Syringe Purple is indicated for syringing preparations with all polyether precision impression materials from 3M in accordance with ISO 4823 Type 3 and Type 2, provided the impression materials are intended for syringing.

CLOSE 

Paste comes out between the 50 ml cartridge and the 3M™ Intra-oral Syringe Purple while loading.

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



← Storage/preparation

FAQs

Is 3M™ Impregum™ Super Quick Material performance influenced by using the 3M™ Intra-oral Syringe Purple?

All impression materials perform as usual when using the single-use Intra-oral Syringe Purple.

CLOSE

Paste comes out between the 50 ml cartridge and the 3M™ Intra-oral Syringe Purple while loading.

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?




← Storage/preparation

FAQs

Paste comes out between the 50 ml cartridge and the 3M™ Intra-oral Syringe Purple while loading.

Place the Intra-oral Syringe Purple directly onto the openings of the 50 ml cartridge and hold firmly in place while loading.



CLOSE



FAQs

 [Storage/preparation](#)

FAQs

When loading the 3M™ Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

To minimize this effect, dispense a pea-sized amount of material from the 50 ml cartridge before loading the Intra-oral Syringe Purple. Often, you cannot completely avoid uneven filling. The first small quantity of paste dispensed from the Intra-oral Syringe Purple should be discarded until the base and the catalyst pastes result in an even, homogeneous mixture.

CLOSE 

General

Indications
other in

Storage

Application
syringin

Material

ion Material?

urple?

nge Purple?

Paste comes out between the 50 ml cartridge and the 3M™ Intra-oral Syringe Purple while loading.

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



← Storage/preparation

FAQs

When loading the 3M™ Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

While loading the Intra-oral Syringe Purple, leave the mixing tip in the original angled position to prevent paste from entering the mixing tip. In an angled position overfilled paste will emerge at the side of the joint.

CLOSE 

Paste comes out between the 50 ml cartridge and the 3M™ Intra-oral Syringe Purple while loading.

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



← Storage/preparation

FAQs

How should I reclose the 50 ml cartridge?

Close the 50 ml cartridge using an original cartridge cap. Alternatively close with a mixing tip containing hardened impression material. Make sure that the catalyst paste does not come into contact with the base paste to prevent paste cross-contamination.

CLOSE 

Paste comes out between the 50 ml cartridge and the 3M™ Intra-oral Syringe Purple while loading.

When loading the Intra-oral Syringe Purple the catalyst and base pastes do not run evenly into the barrels of the syringe.

When loading the Intra-oral Syringe Purple or when placing the plungers, paste flows into the mixing tip.

How should I reclose the 50 ml cartridge?



General



Indications and compatibility with
other impression procedure products



Storage/preparation



Application/working and setting time/
syringing



Material properties



FAQs

Application/working and setting time/syringing

To see the answer click on the question.

Can 3M™ Impregum™ Super Quick Material be combined with VPS materials?

Can Impregum Super Quick Material be used with other 3M polyether materials?

Does the working time of the materials depend on the temperature?

What is the difference between the working and setting times of Impregum Super Quick, Impregum Quick and regular-setting Impregum Materials?

What is the intra-oral and what is the extra-oral working time for Impregum Super Quick Material?

Do I have to use the full working time?

The working time of Impregum Super Quick Material is too short.

When should the tray be inserted?

When should I set the timer? At start of mixing?

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE ↓



FAQs

[← Application/working and setting time/syringing](#)

FAQs

Can 3M™ Impregum™ Super Quick Material be combined with VPS materials?

No, VPS and polyether materials are chemically not compatible.

CLOSE

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE ↓

General

Indications
other in

Storage


Application
syringing

Material

ck,

ck Material?



 [Application/working and setting time/syringing](#)

FAQs


Can Impregum™ Super Quick Material be used with other 3M polyether materials?

Yes, but it is not recommended. If Impregum Super Quick Material is combined with other Impregum or 3M™ Permadyne™ Polyether Impression Materials the longer intra-oral setting time does apply.

CLOSE 

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE 



← Application/working and setting time/syringing

FAQs

Does the working time of the materials depend on the temperature?

Like all chemically setting materials the working time for 3M™ Impregum™ Super Quick Impression Material depends on the temperature. At room temperatures above 25°C/77°F, the working times for all Impregum products are shorter compared to the values given in the instructions for use. According to the IFU, Impregum Super Quick Material as well as all other 3M impression materials should be stored at 18/64°F – 25°C/77°F.

CLOSE 

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE ↓



← Application/working and setting time/syringing

FAQs

What is the difference between the working and setting times of 3M™ Impregum™ Super Quick, Impregum Quick and regular-setting Impregum Polyether Impression Materials?

Working time: The regular-setting Impregum Materials offer 1:45 min., the “Quick” products max. 1:00 min., and the new Impregum Super Quick Material offers 45 seconds maximum working time. Since 3M polyether impression materials are less affected by temperature, the extra- and intra-oral working time are the same.

Intra-oral setting time: The intra-oral setting time of Impregum Super Quick Material has been reduced significantly – to only 2 minutes. Unused working time does not need to be added to the intra-oral setting time for the Impregum Quick and Super Quick materials. In contrast to regular-setting Impregum products, 3M does not indicate a “total setting time” for Impregum Quick and Super Quick Materials.

CLOSE

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous?
Or why does the impression material not set properly?

VIEW MORE ↓

FAQs



← Application/working and setting time/syringing

FAQs

What is the intra-oral and what is the extra-oral working time for 3M™ Impregum™ Super Quick Material?

In contrast to most VPS materials, only one working time is given in the instructions for use. The setting reaction of 3M polyether materials are less temperature sensitive compared to VPS. Therefore, only one valid extra- and intra-orally working time is given.

CLOSE

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?


When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE ↓

- General
- Indications and other information
- Storage
- Application and syringing
- Material

- ck,
- ck Material?



 [Application/working and setting time/syringing](#)

FAQs

Do I have to use the full working time?

No, you can use any time up to the maximum working time of 45 seconds. The intra-oral setting time of 2 minutes remains constant, regardless of the working time used.

CLOSE 

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE 



← Application/working and setting time/syringing

FAQs

The working time of 3M™ Impregum™ Super Quick Material is too short.


Impregum Super Quick Material offers 45 seconds intra-oral syringing time and 45 seconds working time at room temperature – this time allows for syringing 1 or 2 units and inserting the tray within the 45 seconds working time. The tray should be inserted slowly – so plan in a couple of seconds for tray insertion. This is the reason why we recommend the “Super Quick” versions only for 1 or 2 unit cases. If the dentist wants to take a little more time or the 3M™ Pentamix™ Automatic Mixing Unit is not in the same operatory, 3M recommends to use a regular-setting material.

CLOSE

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous?
Or why does the impression material not set properly?

VIEW MORE ↓



 Application/working and setting time/syringing

FAQs

When should the tray be inserted?

The tray should be inserted ideally immediately after the syringing is finished, so the tray should be ready by then. The tray must be inserted within the intra-oral working time.

CLOSE 

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE ↓



FAQs

[← Application/working and setting time/syringing](#)

FAQs

When should I set the timer? At start of mixing?

Just set the timer as soon as the tray is inserted in the mouth. The tray can be removed from the mouth after the indicated intra-oral setting time which is 2:00 minutes for 3M™ Impregum™ Super Quick Material. To assure staying within the given working times you might find it helpful to use a second timer.

CLOSE

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE ↓



FAQs

[← Application/working and setting time/syringing](#)

FAQs

How does the new 3M™ Impregum™ Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

The new Impregum Super Quick Wash Material shows the same consistency according to the ISO standard measurement. Dentists might recognize a difference to the other Impregum wash materials. Impregum Super Quick Material shows a high thixotropy/structural viscosity, i.e. it flows well when syringing but stays more stable than the existing Impregum wash materials at the preparation afterwards.

CLOSE

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE ↓



FAQs

[← Application/working and setting time/syringing](#)

FAQs

When using the 3M™ Intra-oral Syringe Purple:

Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

Make sure you are using the correct syringe. Use the purple syringe for 3M™ Impregum™ Super Quick Material (2:1) and the green one for VPS (1:1) wash materials. A completely even filling is often not possible. Therefore, discard the first amount of syringed material per the instructions for use.

CLOSE

How does the new Impregum Super Quick wash material viscosity correspond to the existing Impregum wash viscosities?

When using the 3M™ Intra-oral Syringe Purple: Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

VIEW MORE ↓



FAQs

Application/working and setting time/syringing

General



Indications and compatibility with
other impression procedure products



Storage/preparation



Application/working and setting time/
syringing



Material properties



To see the answer click on the question.

When using the 3M™ Intra-oral Syringe Purple: Why is the extrusion force so high?

When using the Intra-oral Syringe Purple: Is it possible to use it with an angled mixing tip?

When using the Intra-oral Syringe Purple: During application I observed leakage at the joint.

Why is the mixing quality inhomogeneous?

Or why does the impression material not set properly?

Is it useful to air-blow the wash material onto the teeth and margins to get better detail of the impression?

I observe void formation at the preparation margin. What can I do to avoid that?

How to remove the impression from the patient's mouth?

What are the reasons for poor bonding of the wash material to the tray material?

What can be the cause of smeary or greasy impression surfaces in the preparation area?

How to disinfect 3M™ Impregum™ Super Quick Material?

Is Impregum Super Quick Material scannable?

VIEW MORE ↓

VIEW LESS ↑



FAQs

← Application/working and setting time/syringing

FAQs

When using the 3M™ Intra-oral Syringe Purple: Why is the extrusion force so high?

3M™ Impregum™ Super Quick Materials are tested for an acceptable extrusion force. Please make sure that:

- The Intra-oral Syringe Purple is completely activated and the mixing tip is in the 180° straight position.
- The material does not set in the mixing tip before or during use.

CLOSE 

[How to disinfect 3M™ Impregum™ Super Quick Material?](#)

[Is Impregum Super Quick Material scannable?](#)

VIEW MORE ↓

VIEW LESS ↑



FAQs

[← Application/working and setting time/syringing](#)

FAQs

When using the 3M™ Intra-oral Syringe Purple: Is it possible to use it with an angled mixing tip?

No, the mixing tip has to be in a 180° straight position. Otherwise, the extrusion force would be higher and a part of the paste would be emerged at the joint venting valve.

CLOSE 

[How to disinfect 3M™ Impregum™ Super Quick Material?](#)

[Is Impregum Super Quick Material scannable?](#)

VIEW MORE ↓

VIEW LESS ↑



FAQs

← Application/working and setting time/syringing

FAQs

When using the 3M™ Intra-oral Syringe Purple: During application I observed leakage at the joint.

Please be sure that the Intra-oral Syringe Purple is completely activated and the mixing tip is in the 180° straight position. Otherwise a part of the paste may emerge at the joint venting valve.

CLOSE

How to disinfect 3M™ Impregum™ Super Quick Material?

Is Impregum Super Quick Material scannable?

VIEW MORE ↓

VIEW LESS ↑



FAQs

[← Application/working and setting time/syringing](#)

FAQs

Why is the mixing quality inhomogeneous? Or why does the impression material not set properly?

Dispense a pea-sized amount of material from the 50 ml cartridge before attaching the mixing tip. The first small quantity of paste dispensed from the mixing tip should be discarded until the base and the catalyst pastes result in an even, homogeneous mixture.

CLOSE

[How to disinfect 3M™ Impregum™ Super Quick Material?](#)

[Is Impregum Super Quick Material scannable?](#)

[VIEW MORE ↓](#)

[VIEW LESS ↑](#)



FAQs

← Application/working and setting time/syringing

FAQs

Is it useful to air-blow the wash material onto the teeth and margins to get better detail of the impression?

3M recommends to carefully syringe wash materials around the prepared teeth as described above. When using the air-blow technique there is a risk of introducing bubbles into the wash material which can cause voids at the impression.

CLOSE 

How to disinfect 3M™ Impregum™ Super Quick Material?

Is Impregum Super Quick Material scannable?

VIEW MORE ↓

VIEW LESS ↑



FAQs

← Application/working and setting time/syringing

FAQs

I observe void formation at the margin of preparation. What can I do to avoid that?

Start spiral syringing around the tooth stump in the sulcus without stopping, and while ensuring that the mixing tip is not removed from the material. If the preparation has been cleaned with hydrogen peroxide, rinse carefully with water.

CLOSE 

How to disinfect 3M™ Impregum™ Super Quick Material?

Is Impregum Super Quick Material scannable?

VIEW MORE ↓

VIEW LESS ↑



FAQs

[← Application/working and setting time/syringing](#)

FAQs

How to remove the impression from the patient's mouth?

Make sure to use trays in an appropriate size (not too small). For an easier removal, block out undercuts before impression taking if information is not needed. Slightly lift the edge of the impression with one finger in order to let air creep under the impression thus overcoming the vacuum. Carefully blow air or water with the air syringe between the impression and the teeth underneath.

CLOSE

[How to disinfect 3M™ Impregum™ Super Quick Material?](#)

[Is Impregum Super Quick Material scannable?](#)

[VIEW MORE ↓](#)

[VIEW LESS ↑](#)



FAQs

← Application/working and setting time/syringing

FAQs

What are the reasons for poor bonding of the wash material to the tray material?

The main reason is that the working time has been exceeded and the wash material is already setting when the tray is seated. Be sure to use compatible materials and ideally materials from one product line with same setting speed. Do not combine VPS and polyether materials.

CLOSE 

How to disinfect 3M™ Impregum™ Super Quick Material?

Is Impregum Super Quick Material scannable?

VIEW MORE ↓

VIEW LESS ↑



← Application/working and setting time/syringing

FAQs

What can be the cause of smeary or greasy impression surfaces in the preparation area?

Such smear layers are an indication that the impression material has not set properly. One reason can be improper mixing of the material. Please be sure to discard the first amount of material. Another reason is contact with chemicals, e.g. traces of hemostatic solutions or pastes. These residues can inhibit polymerization. Also, thoroughly rinse off any retraction solution with water around the preparation. Such effects are often not seen until casting, where adhesion of plaster to the impression or the adhesion of impression material to the plaster cast can occur.

If topical anesthetics containing epinephrine are used, let the patient rinse out his mouth well before the impression is taken. If the preparation has been cleaned with hydrogen peroxide, rinse carefully with water.

CLOSE

VIEW MORE ↓

VIEW LESS ↑



← Application/working and setting time/syringing

FAQs

How to disinfect 3M™ Impregum™ Super Quick Material?

All commercial water-based disinfectants indicated for dental impression materials can be used. Glutaraldehyde solutions are recommended for disinfection. Follow the manufacturer's instructions for use, time, and level of disinfection required. After disinfection, rinse impression well with water and blow dry. The impression should be kept dry when sent to the laboratory. Temperatures exceeding 60°C/140°F will damage polyether impressions.

CLOSE 

How to disinfect 3M™ Impregum™ Super Quick Material?

Is Impregum Super Quick Material scannable?

VIEW MORE ↓

VIEW LESS ↑



FAQs

[← Application/working and setting time/syringing](#)

FAQs

Is 3M™ Impregum™ Super Quick Material scannable?

Impregum Super Quick Material can be scanned with impression scanners dedicated for impression scanning. Depending on scanner technology used applying scanning sprays or powders might be indicated.

With current LED line scanners (e.g. D2000, E3 from 3Shape), impression materials can be digitized without any surface treatment for easy access to CAD/CAM workflow.

CLOSE

[How to disinfect 3M™ Impregum™ Super Quick Material?](#)

[Is Impregum Super Quick Material scannable?](#)

VIEW MORE ↓

VIEW LESS ↑




FAQs


Application/working and setting time/syringing

To see the answer click on the question.

General 

Indications and compatibility with
other impression procedure products 

Storage/preparation 

Application/working and setting time/
syringing 

Material properties 

How to pour 3M™ Impregum™ Super Quick Material?

What needs to be considered when pouring the impression with resin based materials?

The gingival mask shows a smeary surface. What gingival mask material can be used?

Is plating of Impregum Super Quick Material possible?

The removal of the Impregum Super Quick Material impression from the model is very difficult. What shall I do?

What are the ideal storage conditions for Impregum Super Quick Material?

VIEW LESS 



FAQs

[← Application/working and setting time/syringing](#)

FAQs

How to pour 3M™ Impregum™ Super Quick Material?

The cast should be prepared from the impression with commercial gypsum stone no earlier than 30 minutes and no later than 14 days after making the impression.
Do not use surfactants (de-bubblizers).


CLOSE

General
Indications
other in
Storage
Applica
syringin
Material

l is very



FAQs

 [Application/working and setting time/syringing](#)

FAQs

What needs to be considered when pouring the impression with resin based materials?

When pouring polyether impressions with epoxy, or any urethane resin based materials, a separator must be used.

CLOSE

General

Indications
other in

Storage

Application
syringing

Material

is very



← Application/working and setting time/syringing

FAQs

The gingival mask shows a smeary surface. What gingival mask material can be used?

Choosing the right gingival mask material in combination with the right separating material is important: Polyether materials from 3M (3M™ Impregum™, 3M™ Permadyne™) can be used with common separators based on silicone, wax, or petrolatum (vaseline). C-silicones (GI-MASK®, Coltène®) can be used without a separator. Use caution with A-silicones (VPS), they can be used only with a separator that forms a sealed solid waxy layer (Sherasepal-U®, Shera®). Separating materials can be applied with a brush, to insulate areas as necessary.

CLOSE

Genera
Indicati
other in
Storage
Applica
syringin
Materia

l is very



FAQs

[← Application/working and setting time/syringing](#)

FAQs

Is plating of 3M™ Impregum™ Super Quick Material possible?

All polyether impressions can be silver-plated, but cannot be copper-plated.

CLOSE 

Genera
Indicati
other in
Storage
Applica
syringin
Materia

l is very



FAQs

← Application/working and setting time/syringing

FAQs

The removal of the 3M™ Impregum™ Super Quick Material impression from the model is very difficult. What shall I do?

Heating the impression to about 40°C/104°F after the cast has set will help in the removal of the cast. Block out, or minimize undercuts if information is not needed. Isolating or blocking out the palatal area will help in removal of the model on a maxillary impression with a deep vault. When making a custom tray, make sure that you have 2 – 3 mm material thickness per 1 mm undercut. You may section or cut the custom tray if it is “locked-on” the model.

CLOSE 



FAQs

 [Application/working and setting time/syringing](#)

FAQs

What are the ideal storage conditions for 3M™ Impregum™ Super Quick Material?

Carefully dry off impression after disinfection. Store impressions dry and below 30 °C/86 °F in the dark. Do not expose to direct sunlight and do not store when moist (e.g. storing an impression not quite dry yet in a plastic bag or together with alginate) since this may damage the impression.

CLOSE 

Genera

Indicati
other in

Storage

Applica
syringin

Material

l is very



FAQs


Material properties

To see the answer click on the question.

General 

Indications and compatibility with
other impression procedure products 

Storage/preparation 

Application/working and setting time/
syringing 

Material properties 

3M states that the taste of 3M™ Impregum™ Super Quick Material is improved. Nevertheless, some of my patients complain about unpleasant taste.

Do Impregum Polyether Impression Materials contain additional surfactants?

What is the Shore hardness of Impregum Super Quick Material?



← Application/working and setting time/syringing

FAQs

3M states that the taste of 3M™ Impregum™ Super Quick Material is improved. Nevertheless, some of my patients complain about unpleasant taste.

3M conducted a clinical study. In this study, Impregum Super Quick Material showed a more pleasant taste and increased patient comfort compared to existing Impregum Penta Impression Material. The aftertaste was also rated much better. Nevertheless, the taste sensitivity is very subjective and differs widely among individuals. Some patients will still be complaining, particularly patients who do not like the minty taste.

CLOSE

Genera

Indicati
other in


Storage

Applica
syringin

Material



FAQs

 [Application/working and setting time/syringing](#)

FAQs

Do 3M™ Impregum™ Polyether Impression Materials contain additional surfactants?

The intrinsic hydrophilicity is caused by polyether polymers. Impregum Materials do not contain typical surfactants used in VPS materials.

CLOSE

General

Indications
other in

Storage

Application
syringing

Material



FAQs

[← Application/working and setting time/syringing](#)

FAQs

What is the Shore hardness of 3M™ Impregum™ Super Quick Material?

Impregum Super Quick Material has a Shore A hardness in the range of 50 to 55.

CLOSE 

General
Indications
other in
Storage
Applica
syringin
Material

3m.com/Dental