Description

SpeedCEM® Plus is a self-adhesive, self-curing resin cement with light-curing option for the cementation of indirect restorations made of high-strength all-ceramic, metal-ceramic and metal. When SpeedCEM Plus is used, additional dentin bonding agents are not required. The automix syringe enables time-saving direct application of the luting material into the restoration.

Shades

SpeedCEM Plus is available in three shades with different degrees of translucency:

- transparent
- yellow
- white opaque

Working and setting times

The following times apply as soon as SpeedCEM Plus has been dispensed from the automix syringe (the working and setting times depend on the ambient temperature):

	At room temperature 23°C ± 1°C/ 73°F ± 2°F	Intra-orally approx. 37°C ± 1°C/ 98°F ± 2°F
Working time	approx. 2 min	approx. 1 min
Setting time (incl. working time)	approx. 6 min	approx. 3 min

Note

After SpeedCEM Plus has been dispensed from the automix syringe, the curing process can be accelerated by intensive illumination (operating light, ambient light).

Composition

The monomer matrix is composed of dimethacrylates and acidic monomers. The inorganic fillers are barium glass, ytterbium trifluoride, co-polymer and highly dispersed silicon dioxide. Additional contents are initiators, stabilizers and colour pigments (< 1 %). The primary particle size of the inorganic fillers is between 0.1 and 7 µm. The mean particle size is 5 µm. The total content of inorganic fillers is approx. 45 vol%.

Indication

- 1) Permanent cementation on natural teeth with adequate retentive features in conjunction with indirect restorations made of:
 - oxide ceramics, e.g. zirconium oxide (crowns, bridges, endodontic posts)
 - metal and metal-ceramics (inlays, onlays, crowns, bridges, endodontic posts)
 - lithium disilicate glass-ceramics (crowns, bridges)
 - fibre-reinforced composites (endodontic posts)
- 2) Permanent cementation of crowns and bridges on implant abutments that offer adequate retention and are made of the following abutment materials:
 - oxide ceramic (e.g. zirconium oxide)
 - metal (e.g. titanium)
 - lithium disilicate glass-ceramic

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Contraindication

SpeedCEM Plus is contraindicated

- in situations where the preparation/shape of the implant abutment does not provide adequate retention (e.g. veneers, short or severely tapered tooth preparations/implant abutments).
- if a patient is known to be allergic to any of the ingredients of SpeedCEM Plus.
- In general, SpeedCEM Plus should not be applied to exposed pulp or to dentin which is close to the pulp.

Side effects

Systemic side effects are not known to date. In rare cases, allergic reactions to single components have been reported.

Interactions

Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the use of products, e.g. mouth rinses and temporary cements, which contain these components should be avoided. Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process. Therefore, the preparation and syringe must not be disinfected using oxidative agents.

The syringe can be wiped e.g. with customary disinfecting wipes.

Alkaline jet media may compromise the effect of SpeedCEM Plus.

Usage

Please also refer to the relevant instructions for use of the products used in conjunction with SpeedCEM Plus for more detailed information.

1 Pre-treatment of the tooth preparation or implant abutments, try-in and isolation

1.1 Pre-treatment of the cavity or tooth preparation

Removal of the temporary restoration and thorough cleaning of the cavity

Remove possible residue of the temporary luting composite from the cavity or tooth preparation with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxyl[®] fluoride-free). Rinse with water spray. Then lightly dry with water- and oil-free air. Avoid overdrying.

Note: Cleaning with alcohol can lead to overdrying of the dentin.

Try-in of the restoration and isolation

Next, the shade, accuracy of fit and occlusion of the restoration can be checked.

Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently seated, because there is a risk of fracture. It is essential to ensure reliable isolation of the operative area with e.g. OptraGate[®], cotton rolls, saliva ejector and absorbent pads. Dental hard tissue that has become contaminated with saliva needs to be cleaned again.

1.2 Pre-treatment of the implant abutment

Removal of the temporary restoration

Remove the temporary restoration and, if present, the temporary abutment. Rinse the implant lumen and peri-implant gingiva.

Insert and check final abutment

According to the manufacturer's directions.

A retraction cord may be placed to be better able to check the accuracy of fit of the restoration and remove excess luting composite.

Try-in of the restoration and isolation

Next, the shade, accuracy of fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently seated, because there is a risk of fracture. It is essential to ensure reliable isolation of the operative area – preferably with OptraGate, cotton rolls, saliva ejector and absorbent pads – when using adhesive cementation with composites.

Pre-treatment of the abutment surface

Note: Abutments made of lithium disilicate glass-ceramic should be pre-treated extraorally using e.g. Monobond Etch & Prime[®] prior to the permanent cementation of the restoration.

- Seal the screw channel of the abutment using e.g. a temporary restorative material (e.g. Telio[®] CS Inlay). It is recommended to isolate the screw head from the temporary restorative material with cotton pellets.
- Rinse the abutment with water spray.
- Dry the abutment.

2 Cleaning and pretreatment of the restoration

Principally the manufacturer's directions must be followed.

Note: Contamination with saliva or blood must be avoided during and after the pre-treatment of the restoration (try-in). Regardless of any previous conditioning, clean contaminated restoration surfaces **using Ivoclean** as follows:

- After try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
- Cover the entire bonding surface of the restoration with a layer of Ivoclean using a microbrush or brush.
- Allow 20 seconds for Ivoclean to take effect, then thoroughly rinse with water spray and dry with oil-free air.
- Restorations made of lithium disilicate glass-ceramic (e.g. IPS e.max[®] Press/CAD) must be primed again with Monobond Etch & Prime or Monobond[®] Plus.

2.1 Restorations made of oxide ceramics - (e.g. zirconium oxide, Zenostar[®])

IMPORTANT! In order to achieve a strong bond, do not clean the surfaces with phosphoric acid.

If the restoration has been sandblasted in the laboratory:

- Apply Ivoclean, as described above, to clean the restoration.

If the restoration has not been sandblasted in the laboratory:

- Sandblast the inner surface of the restoration (use the blasting parameters stated in the manufacturer's instructions for use of the restorative material, e.g. Zenostar, max. 1 bar/14.5 psi, Al₂O₃).
- Clean the restoration e.g. in an ultrasonic unit for approx. 1 minute.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.

2.2 Metal or metal-supported restorations

IMPORTANT! In order to create a strong bond, do not clean the metal surfaces with phosphoric acid.

If the restoration has been sandblasted in the laboratory:

- Apply Ivoclean, as described above, to clean the restoration.

If the restoration has not been sandblasted in the laboratory:

- Sandblast the inner restoration surface (use the blasting parameters stated in the instructions for use of the restorative material) until an even matt surface has been achieved.
- Clean the restoration e.g. in an ultrasonic unit for approx. 1 minute.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.

If the restoration comprises precious metal, apply Monobond Plus to the pre-treated surfaces using a brush or microbrush, allow to react for 60 seconds and then disperse with a strong stream of air.

2.3 Restorations made of lithium disilicate glass-ceramic (e.g. IPS e.max® Press/CAD)

A) Procedure using Monobond Etch & Prime

Regardless of whether the restoration has been pre-treated in the laboratory:

- Following the try-in, thoroughly rinse the restoration with water spray and dry with water- and oil-free air.
- Apply Monobond Etch & Prime onto the bonding surface using a microbrush and agitate it into the surface for 20 seconds using light pressure. Allow to react for another 40 seconds.
- Then thoroughly rinse off Monobond Etch & Prime with water until the green colour is removed. If residues have remained in the micro-porosities after rinsing, the restoration may be cleaned with water in an ultrasonic cleaner for up to 5 minutes.
- Dry the restoration with a strong stream of water- and oil-free air for approx. 10 seconds.

B) Procedure using Monobond Plus

If the restoration has been pre-treated in the laboratory:

- If the restoration has already been pre-treated in the laboratory, it should be cleaned after the try-in using Ivoclean as described above.
- Subsequently, re-apply Monobond Plus onto the cleaned surface using a brush or a microbrush, allow to react for 60 seconds and then disperse with a strong stream of air.

If the restoration has not been pre-treated in the laboratory:

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic Etching Gel) for 20 seconds or according to the manufacturer's instructions for use of the restorative material.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.
- Apply Monobond Plus to the pre-treated surfaces with a brush or microbrush, leave to react for 60 seconds and then disperse with a strong stream of air.

3 Application of SpeedCEM Plus into the restoration

For each application, place a new automix tip on the double syringe. Extrude SpeedCEM Plus from the automix syringe, apply the desired quantity directly into the restoration and cover the entire bonding surface. As the luting material will cure in the used mixing tip, it may serve as a seal for the contents of the syringe until the next application.

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Manufacturer

Ivoclar Vivadent AG
FL-9494 Schaan/Liechtenstein

4 Placement of the restoration and removal of excess cement

Note: As with all composites, SpeedCEM Plus is subject to oxygen inhibition. This means that the surface layer does not polymerize during curing, as it comes in contact with atmospheric oxygen. To avoid this, cover the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after removing the excess material. After complete polymerization, the glycerine gel/air block is rinsed off with water.

4.1 solely self-curing

- Seat the restoration and retain it in place exerting uniform pressure until the excess material has been removed and the cement has completely cured.
- Remove excess material e.g. with a microbrush/brush/foam pellet/dental floss or scaler. Make sure to remove excess material within the processing time, especially in areas that are difficult to reach (proximal or gingival margins, pontics).
- To remove excess material from implant abutments, special implant scalers/curettes should be used to avoid scratching of the implant surface.

4.2 self-curing with additional light-activation to tack-cure excess material or accelerate the curing process (quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridge)

- Seat the restoration and retain it in place exerting uniform pressure until the excess material has been removed and the cement has completely cured.
- Excess cement is light-cured with a polymerization light (approx. 650 mW/cm²) from a distance of approx. 0-10 mm for 1 second per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal).
- Following this procedure, the excess material can be easily removed with a scaler.
- Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins, pontics).
- To remove excess material from implant abutments, special implant scalers/curettes should be used to avoid scratching of the implant surface.
- Subsequently, light-cure all cement joints again for 20 seconds (approx. 1,100 mW/cm²). If the light output is lower, curing takes longer.

5 Finishing the completed restoration

- Check occlusion and function and adjust as necessary.
- If necessary, rework the cement joints with finishing diamonds.
- Smooth the cement joints with finishing and polishing strips and polish with suitable polishers (e.g. OptraPol®).
- If necessary, rework the restoration margins with suitable polishers (ceramics: e.g. OptraFine®).

Special instructions for the cementation of endodontic posts:

- For the cementation of endodontic posts, carefully clean the root canal to remove any residual root canal filler material. (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite.)
- Wet the endodontic post, which has been prepared according to the instructions of the post manufacturer, with the mixed SpeedCEM Plus and apply SpeedCEM Plus directly into the root canal using endodontic tips. The material should be applied in a sufficient quantity to provide excess cement.
- Insert the endodontic post.
- Distribute excess cement on the occlusal preparation surface to achieve a thin coating. Then, light-cure SpeedCEM Plus for 20 seconds.
- Apply the core build-up material (e.g. MultiCore®) directly onto SpeedCEM Plus to serve as a bonding agent and cure according to the manufacturer's instructions (light-curing).

Special notes

SpeedCEM Plus should have room temperature when it is applied. At refrigerator temperature, the material is difficult to extrude and its processing and curing times are increased.

Do not clean the SpeedCEM Plus automix syringes with agents that have an oxidative effect.

Warning

- Avoid contact of uncured SpeedCEM Plus with the skin/mucous membrane or eyes.
- Uncured SpeedCEM Plus may cause slight irritation and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Cool storage (2–8°C/ 36–46°F) is required for SpeedCEM Plus.
- Do not use SpeedCEM Plus after the indicated expiry date.
- Do not disinfect syringes with oxidizing disinfectants.
- Store SpeedCEM Plus automix syringe with the mixing tip attached after use.
- Expiry date: see note on syringes and packages.

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

Tel. +61 3 9795 9599, Fax +61 3 9795 9645

www.ivoclarvivadent.com.au

Ivoclar Vivadent GmbH

Tech Gate Vienna, Donau-City-Strasse 1, 1220 Wien, Austria

Tel. +43 1 263 191 10, Fax: +43 1 263 191 111

www.ivoclarvivadent.at

Ivoclar Vivadent Ltda.

Alameda Caiapós, 723, Centro Empresarial Tamboré, CEP 06460-110 Barueri – SP, Brazil

Tel. +55 11 2424 7400, Fax +55 11 3466 0840

www.ivoclarvivadent.com.br

Ivoclar Vivadent Inc.

1-6600 Dixie Road, Mississauga, Ontario, L5T 2Y2, Canada

Tel. +1 905 670 8499, Fax +1 905 670 3102

www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai

Trading Co., Ltd., 2/F Building 1, 881 Wuding Road, Jing An District, 200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 3399, Fax +57 1 633 1663

www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118, F-74410 Saint-Jorioz, France

Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, D-73479 Ellwangen, Jagst, Germany

Tel. +49 7961 889 0, Fax +49 7961 6326

www.ivoclarvivadent.de

Wieland Dental + Technik GmbH & Co. KG

Lindenstrasse 2, 75175 Pforzheim, Germany

Tel. +49 7231 3705 0, Fax +49 7231 3579 59

www.wieland-dental.com

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate, Veera Desai Road, Andheri (West), Mumbai, 400 053, India

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www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale, PO Box 303011 North Harbour, Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204, P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 11 293 8345, Fax +966 11 293 8344

www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

Carretera de Fuencarral nº24, Portal 1 – Planta Baja, 28108-Alcobendas (Madrid), Spain

Telf. +34 91 375 78 20, Fax: +34 91 375 78 38

www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14, S-169 56 Solna, Sweden

Tel. +46 8 514 939 30, Fax +46 8 514 939 40

www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi, Sakayik Sokak, Nisantas' Plaza No:38/2, Kat:5 Daire:24,
34021 Sisli – Istanbul, Turkey

Tel. +90 212 343 0802, Fax +90 212 343 0842

www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building, Feldspar Close, Warrens Business Park,
Enderby, Leicester LE19 4SE, United Kingdom

Tel. +44 116 284 7880, Fax +44 116 284 7881

www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

www.ivoclarvivadent.us