FICHA TÉCNICA

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| | | |
|---|---|---|
| INSTRUMENT DISINFECTION BESTDENT | |  |
| REFERENCES | <p>60153: (1L). 60154: (5L).</p> | |
| MANUFACTURER | <p>PRISMAN GmbH Otto-Hahn-Ring 6-18 D-64653 Lorsch Phone: +49 (0)6251-866 980-0</p> | |
| CLASSIFICATION | Class IIb Medical Device, Directive 93/42/EEC. | |
| DESCRIPTION | <ul style="list-style-type: none"> - Must be diluted with water before use. - Anticorrosive. - Excellent material compatibility. - Can be used in ultrasonic baths. - Broad spectrum of activity, excellent cleaning properties. | |

TECHNICAL DATA SHEET

| | |
|-----------------------------|--|
| COMPOSITION | 15-30% nonionic surfactants, <5% cationic surfactants, <5% EDTA salts, fragrance. 100g contains 4.4g Didecyldimethylammonium chloride, 2.75g N-(3-Aminopropyl)-N-dodecyl-propan-1,3-diamine. |
| PACKAGE | 60153: Can (1L). 60154: Can (5L). |
| INTENDED USE | For cleaning and disinfection of rotary and surgical instruments of all types. |
| INSTRUCTIONS FOR USE | <p>See product instructions.</p> <p>Microbiological Effectiveness: Bactericidal, levurocidal (<i>Candida albicans</i>), tuberculocidal, “restricted viruses” (effective against membrane covered viruses, i.e. HBV, HIV, HCV).</p> <p>Effectivity: Bacteria in accordance to VAH/DGHM, EN 13727, EN 14561 (high bio load): 3% 15 min. 2% 30 min. 1% 60 min.</p> <p>Levurocidal in accordance to VAH/DGHM and EN 13624, EN 14562 (high bio load): 3% 15 min. 2% 30 min. 1% 60 min.</p> <p>TbB (Tuberculosis-Inactivation in accordance to VAH/DGHM): 3% 60 min.</p> <p>Enveloped viruses (i.e. HBV, HIV, HCV) in accordance to DVV/RKI (BVDV and Vacciniavirus EN14476 high bio load): 3% 15 min. 1% 60 min.</p> <p>Ultrasonic bath (<i>Bacteria</i> and <i>Candida albicans</i> in accordance to VAH/DGHM (high bio load): 3% 5 min.</p> |
| SAFETY WARNINGS | H314 - Causes severe skin burns and eye damage. Avoid breathing dust/fume/fume/gas/mist/mist/vapours/spray. Wear gloves and goggles. |
| INDICATIONS FOR USE | Solution is in a concentrate form which has to be diluted with water. Heavily contaminated instruments should first be cleaned before submerging into the disinfectant. After disinfection, instruments are to be thoroughly rinsed with water and dried. Do not mix with other types of disinfectant or cleansers. Depending on the level of contermination regularly exchange the solution (best every day). The instrument disinfectant solution can be used at the |

TECHNICAL DATA SHEET

| | |
|------------------------------|--|
| | concentrations indicated also in the ultrasonic bath. |
| STORAGE AND TRANSPORT | Store in a cool and dry place. Keep away from direct sunlight and heat sources. Protect from humidity. |

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