FICHA TÉCNICA

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TECHNICAL DATA SHEET

TEMPORARY CEMENT PROCLINIC EXPERT		
REFERENCES	78576	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy	
CLASSIFICATION	Medical Device, Class IIa, Directive 93/42/EEC. Notified Body: Kiwa Cermet Italia S.p.A., No. 0476	
DESCRIPTION	Temporary cement for implant retained crowns	
COMPOSITION	Multifunctional methacrylates, Urethane dimethacrylate, Catalysts	
PACKAGING	1 automix syringe 6 g 5 mixing tips	
SHELF LIFE	2 years from the date of manufacture	
INTENDED USE	Specially formulated cement for luting of provisional crowns and bridges on implant abutments.	
INSTRUCTIONS FOR USE	<ol style="list-style-type: none"> 1) Finish and polish the restoration. Check occlusion and proximal contacts and complete final adjustments. 2) Dry and isolate the area being treated. Absolute dryness of the restoration is not necessary, but contamination should be avoided. 3) Make sure that the restoration is clean and dry. It is advisable to apply a thin layer of water-based lubricant to the abutment surface to facilitate future removal of implant restorations. 4) Bleed the syringe before attaching the automix tip. Bleed the syringe again after placing Proclinic Expert Temporary Cement's automix tip to ensure correct flow of base and catalyst components. Place a thin layer of cement in the restoration. Remember that the working time for this cement is between 45-60 seconds. 5) Keep the crown in the desired position for 2-2.5 minutes by asking the patient to close on a cotton ball to ensure accurate, secure placement. Excess cement can safely be removed after 2.5 minutes. 6) Total setting time of the cement is between 4 - 5 minutes, depending on conditions. 	

SAFETY INFORMATION	Though uncommon, some individuals may be sensitive to this product. If an allergic reaction occurs, discontinue use, remove from the patient's mouth and refer patient to a doctor.
HANDLING AND STORAGE CONDITIONS	Store below 24°C (68°F). Refrigerate for maximum shelf life. It is not recommended for use beyond the expiration date marked on the packaging.