Silane

BOND ENHANCER

2 x 1,2 ml syringes
4 delivery tips

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EN SILANE ADHESION PROMOTOR

Cementation of ceramics in general, for example for:

- Preparation of porcelain crowns, veneers and inlays before definitive cementation.
- Repair of ceramic pieces or fractured metal-ceramic.

PROPERTIES

Proclinc Expert Porcelain Etch is a viscous, hydrofluoric acid for etching any type of ceramic, which creates micropores in the treated surface thereby achieving a strong bond to composite. When applied to the etched ceramic, Proclinc Expert Silane improves the restoration's bonding strength and longevity.

INSTRUCTIONS FOR USE

1. Isolate the working field intra-orally with rubber dam.
2. For fractured ceramic repairs, roughen the surface of the ceramic to be etched using a medium grain bur. Remove the shiny surface and create enough space for the composite. Rinse and dry the ceramic surface.
3. Apply Proclinc Expert Porcelain Etch to the ceramic with a brush and leave to act for 60 seconds.
4. Rinse the ceramic with copious water to remove the acid, and dry. The ceramic should have a matt, frosted appearance. If this is not the case, repeat the etching process until this appearance is achieved.
5. Apply Proclinc Expert Silane on the etched ceramic and dry well with a stream of oil-free and water-free air.

PRECAUTION!

Use with caution. Avoid contact with mucosa, eyes and skin. In case of accidental contact, wash immediately with copious amounts of water and consult a specialist. It is recommended to isolate the working field intra-orally with a rubber dam.

COMPOSITION

Proclinc Expert Silane: • Ethanol alcohol: 97% • Glycidoxy propyl methoxy silane: 3%

STORAGE Can be stored for 2 years under suitable storage conditions at room temperature of 24°C (75°F).

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