FICHA TÉCNICA

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TECHNICAL DATA SHEET

SELF-SEALING STERILIZATION POUCHES		
REFERENCES	79930: film 7x23 cm 79931: film 6x11 cm 79932: film 9x23 cm 79933: film 14x25 cm 79947: film 33x44 cm 79954: film 9x13 cm	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano, Italy	
CLASSIFICATION	Medical Device, Class I, Directive 93/42/EEC	
DESCRIPTION	Self-sealing sterilization pouches for single use	
COMPOSITION	Made of medical use laminated transparent film (PET and CPP) and medical paper.	
SPECIFICACIONS	Sterilization with steam: temperature between 121°C-134°C; the indicator colour changes from pink to brown. Sterilization with ethylene oxide: temperature between 30°C-54°C; the indicator colour changes from light blue to yellow/brown.	
PACKAGING	Box with 200 pouches	
SHELF LIFE	8 years from the date of manufacture	
INTENDED USE	Sterilization of medical devices and materials for medical use which are especially suitable for steam or ethylene oxide gas.	
HANDLING AND STORAGE	Store in a cool dry place, away from direct sunlight. Protect from moisture.	

SAFETY INFORMATION	• Before starting sterilization process, check that the pouch is filled a maximum of 3/4 of its volume, so that sealings are not obstructed at the edges. • Cover the end and the cutting surfaces of instruments with appropriate protection. • Put the instruments with the handle at the opening of the pouch. • Do not use solvent substances on the pouch. • Possible printings have to be written on a label and stuck on paper or film side; please do not write directly on the pouch, and anyway, use the maximum care to avoid contamination problem during use. • Carefully check the weight of the instruments to be sterilized in order to avoid risks of breaking.
PRESERVATION OF STERILITY	It's essential to grant the correct condition for storage and use of the material after the sterilization process. Preservation of sterility depends on type of packaging and storage conditions. The sterilized material should be stored in closed, clean and dry cupboards, protected from powder, with limited access to closed windows and doors, at a temperature between 18-22°C and humidity between 35-50%. The material should be handled as less as possible and with care. Packages broken, damaged, opened, fallen down or wet must be considered as non-sterile, and must be repacked and sterilized again. Sterility duration for steam depends on the correct and strict observance of the norms and the technical instructions for sterilization. In this context, we suggest a duration of sterility for steam in single pouch (usually used for non-invasive used on laboratory and stay in bed operation unit) of 30 days, and in double pouch (usually used with "no touch" technique for medical devices used for invasive procedures and for medical devices used in Surgery Unit) a duration of 60 days, on condition that in both cases, the conditions suitable for conservation are respected.