

POLA LIGHT NOT WORKING?

POLA
LIGHT

ADVANCED TOOTH
WHITENING SYSTEM

TRY THE FOLLOWING STEPS

STEP 1: Check that the Pola Light mouthpiece is laid flat with Pola logo facing down

STEP 2: Follow option 1 or Option 2



OPTION 1

Connect the charger horizontally - a small tick will be heard when it's connected properly and the mouthpiece will flash to indicate that it is being charged.



OPTION 2

Connect the charger roughly at a 45 angle, slightly slanted. The mouthpiece will flash to indicate that it is being charged.



YOUR SMILE. OUR VISION.



INSTAGRAM
@pola_smile



FACEBOOK
polasmiletoothwhitening

SDI Limited
3-15 Brunsdon Street,
Bayswater, Victoria 3153
Australia

www.sdi.com.au
www.polawhite.com.au