bestdent

Instructions for use

- | | |
|-------------------------------------|----------------------------------|
| de Gebrauchsanweisung | no Bruksanvisning |
| es Instrucciones de uso | cs Návod k použití |
| fr Mode d'emploi | sk Návod na použitie |
| it Istruzioni per l'uso | sl Navodila za uporabo |
| pt Instruções de uso | et Kasutusjuhised |
| nl Gebruiksaanwijzing | hu Használati utasítás |
| pl Instrukcja postępowania | hr Upute za uporabu |
| da Brugsvejledning | bg Инструкции за употреба |
| ro Instrucțiuni de utilizare | el Οδηγίες χρήσης |
| sv Bruksanvisning | fi Käyttöohjeet |

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EN - Bestdent Joya Dental

Description

It is a crystal dental jewel equipped with a base for mechanical retention, applied without any invasive procedure by means of a dental adhesive method.

Method of application

1. Isolate the field with a rubber dam, etch the enamel for 20 seconds, rinse with water and dry with air.
2. Apply a small amount of flowable composite onto the etched enamel.
3. Place the product applying slight pressure and leaving a very thin film between the tooth and the lens.
4. Light cure for 40 seconds.
5. To increase adhesion, spread a thin layer of dentin adhesive over the edge of the crystal (this does not prevent shine even if applied over the crystal).
6. If excess adhesive material occurs during application, remove it using a polishing system (for example, silicone rubbers).

Warnings

- The product must be used by professionals in the sector.
- Keep the product out of the reach of children.
- For dental use.
- Do not reuse the product.
- Non-sterile product.
- Use gloves before use.
- If the product comes into contact with oily materials, they will act as insulators; therefore the base of the Dental Jewel must be etched with orthophosphoric acid, rinsed and dried well.
- Dental Jewel can be removed at any time.
- For information on side effects and interactions of fixing products, consult the technical data sheets of the materials themselves.

Symbols

Product code

Lot number

Do not reuse

Not-sterile

IFU

Manufacturer

Distributor

Warranty

Dental World will replace any product found to be defective.

Dental World shall not be responsible for any damage or loss, direct or consequential, arising from careless use or inability to use the described product. Before using it, it is the user's responsibility to determine the suitability of the product for its intended use. The user assumes responsibility for careless use.

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