Laboratorios **Clarben S.A.**

CE 0426

OCTOPLUS®

Agujas dentales / agulhas dentais / dental needles / aiguilles dentaires / Aghi dentali

REF **01-017** OCTOPLUS agujas 30G X-CORTA 0.3 x 16 mm caja 100 u.

REF **01-020** OCTOPLUS agujas 30G CORTA 0.3 x 21 mm caja 100 u.

REF **01-016** OCTOPLUS agujas 30G CORTA 0.3 x 25 mm caja 100 u.

REF **01-018** OCTOPLUS agujas 27G LARGA 0.4 x 38 mm caja 100 u.

ESPAÑOL (ES)

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ENGLISH (EN)

Sterile single-use needle, intended for use in the dental field for anaesthesia within the oral cavity.

INDICATIONS:

Needles for dental use with threaded hub for cartridge syringes, sterile, silicone-coated, disposable, non-toxic and non-pyrogenic when in unopened and undamaged package. The welding point joining the cap to the protector ensures that the sterility of the needle is maintained.

INSTRUCTIONS FOR USE:

Positioning the needle on the syringe
1. Break the guarantee seal
2. Pull while twisting to remove the cap
3. Screw the threaded base onto the syringe, taking care not to remove the protector
4. Remove the protector with a slight twisting motion just before anaesthesia.

After use:

1. Carefully remove the needle from the syringe to avoid accidental puncture.
2. Dispose of the used needle immediately in a suitable container conform to regulations, without bending it, break it or handle it in any other way.

WARNINGS:

- Do not use if the guarantee seal is not intact.
- Do not reuse. Should the needle be accidentally reused, there is a risk of contamination through cross infection.
- Do not recap the needle by hand, as this may cause accidental puncture.
- Do not bend the needle or apply excessive pressure on it, as serious injury to you and/or your patient may occur.
- If any serious incident has occurred in relation to the device, the user should be report it to the manufacturer and the competent authority of the Member State in which the user is established.

CAUTIONS:

- Keep in original packaging until use. Store at room temperature and protect from moisture.
- Only for professional use in the dental field for the administration of anaesthetics.

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Laboratorios Clarben S.A.
Tecnofar S.p.A.