KULZER
MITSUI CHEMICALS GROUP

cara Print LEDcure

DE Betriebsanleitung
GB Operating instructions
FR Mode d'emploi
ES Instrucciones de uso
IT Istruzioni per l'uso
PT Manual de instruções
NL Gebruiksaanwijzing
SE Bruksanvisning
DK Brugsanvisning
NO Bruksanvisning
FI Käyttöohjeet
TR Kullanım kılavuzu
AE إرشادات الاستخدام لجهاز
CN 使用说明书
KR 사용 설명서

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Table of contents

1 Scope of application 2
1.1 General 2
1.2 Unit type 2
1.3 EC Declaration of Conformity 2
2 Directions for safe operation 3
2.1 Explanation of symbols 3
2.2 Transport damage 3
2.3 Operator responsibilities 3
2.4 Equipment log / operating instructions 3
2.5 Safety notes 3
3 Intended use 4
3.1 Rules for operation 5
4 Scope of delivery 5
5 Operating and indicator elements 5
6 Unit description and function 6
7 Placing the unit into operation 7
7.1 Transport 7
7.2 Setup 7
7.3 Component check before placing the unit into operation 7
7.4 Mains connection 7
7.5 Unit inspection 8
8 Operation 9
8.1 Program interface and program selection 9
8.2 Execute program 11
8.3 Coding customer-specific programs 13
8.4 Working with the unit 14
8.5 Software error messages 15
8.6 cara Cure Guide update 16
8.7 Install cara Cure Guide update 16
8.8 Errors and how to clear them 18
9 Cleaning 18
10 Maintenance / repair 19
11 Disposal 19
12 Technical details 20
12.1 Transport and storage conditions 20
12.2 Operating conditions: 20
12.3 Identification plate 20
13 Wiring diagram 21
13.1 Mains, power supply / fan 21
13.2 Controller interface 21
14 Service 22
15 Document history 22

1 Scope of application

1.1 General

cara Print LEDcure is a registered ® trademark of Kulzer GmbH.
Author f42415

These operating instructions are applicable to:

Order no.	Type and equipment	Version
66078521	cara Print LEDcure, LED light polymerization unit	2025-05

1.2 Unit type

Description	Type	Valid from serial no.
LED light polymerization unit	cara Print LEDcure	2025-05-xxxx

1.3 EC Declaration of Conformity

We, Kulzer GmbH, Leipziger Str. 2, 63450 Hanau (Germany) hereby declare, under our own responsibility, that the device described in section 1.2, based on its design and construction as well as the version we have brought to market, complies with the relevant basic health and safety requirements as defined in EC Directives and harmonized standards.

Applicable EU Directives:

Low Voltage Directive 2014/35/EU
EMC Directive 2014/30/EU

Applicable harmonized standards:

DIN EN 61010-1:2020-03
DIN EN IEC 61010-2-010:2022-12
DIN EN 61326-1:2022-11

This statement will become invalid if the device is modified without coordinating with us.

2 Directions for safe operation

2.1 Explanation of symbols

Symbols	Wordmark	Explanation
	WARNING	Observe the safety-relevant section and sub-sections of the operating instructions.
	NOTE	Note in the operating instructions concerning the optimal use of the unit.
	ELECTRIC SHOCK	Electric shock warning. Potentially life-threating consequences if instructions are not followed.
	DISCONNECT FROM MAINS	Electric shock hazard when opening unit. Unplug the power supply before opening.
	HOT SURFACE	There is a burning hazard posed by both the 3D-printed item and the inside of the chamber after the polymerization cycle has finished; wear protective gloves.
	WEAR SAFETY GOGGLES	Near-UV and blue light LED radiation hazard. Wear safety goggles during all activities while operating the unit.
	FLAMMABLE LIQUIDS	Fluid and vapor flammability hazard. Follow safety precautions and instructions.

2.2 Transport damage

After receiving the delivery, check the unit for transport damage and report it to the transportation company within 24 hours of delivery. **Do not work with a damaged unit under any circumstances.**

2.3 Operator responsibilities

In addition to complying with the manufacturer's legal provisions, the operator is required to ensure compliance and implementation of legal workplace obligations, e.g. the requirement to provide training and compliance with the Occupational Health and Safety Act as well as all other applicable laws and regulations. The provider must make written instructions available to the operators working on and with the unit in a comprehensible format and in their relevant languages that are based on the operating instructions and work to be performed.

2.4 Equipment log / operating instructions

We recommend keeping an equipment log. The equipment log is used to document inspections and all important work performed (e.g. repairs, changes). Keep these operating instructions in a safe place so that you are able to refer to safety instructions and important user information at a later date.

2.5 Safety notes

The safety instructions for protecting individuals as well as handling the unit and the material to be processed are largely dependent on the conduct of the persons operating the unit.

Before placing the unit into operation, read the operating instructions carefully and follow the information provided so as to avoid errors or potential resulting injury, particularly those related to health.

All respective national laws, provisions and directives in force must be followed in addition to the information contained in these operating instructions when setting up and operating the unit.

ELECTRIC SHOCK

Check the mains lead and plug for damage before use.
If damage is present, do not connect the unit to the mains.

WEAR SAFETY GOGGLES

The LEDs in the curing chamber emit visible light (violet and blue) at high power density. Never look into bright LEDs. Always keep the door completely closed and do not open it before the program has completed.

A glare hazard is present. If the automatic LED switch-off does not function when the door is opened, operating the unit is not permitted. Do not look into stray light from housing gaps and cooling slots for extended periods (minutes).

DISCONNECT FROM MAINS

Only **Kulzer Service** is authorised to work on the unit's electronic components when performing such work is safe (unit is de-energized).
Only approved genuine replacement parts and accessories may be used. Using other parts carries unknown risks and must be avoided in all cases.
The functionality and safety of the device are only guaranteed if the necessary inspections, maintenance, and repairs are carried out by Kulzer service partners or personnel trained by Kulzer.
Kulzer GmbH **assumes no liability** whatsoever for any damage caused to the unit resulting from improper repair, which has not been performed by **Kulzer Service Partners** or by personnel we have trained, or if genuine parts or accessories were not used when replaced.

FLAMMABLE LIQUIDS AND VAPORS

Items that contain residues or adhering material from flammable solvents or cleaning fluids or are not completely dry may not be introduced into or cured with the unit.

3 Intended use

The cara Print LEDcure light polymerization unit is used for the final polymerization of the 3D-printed objects manufactured from photosensitive fluids / resins used exclusively for dental applications.
The cara Print LEDcure light polymerization unit is manufactured in accordance with the latest state-of-the-art techniques and is safe to operate. However, this unit may present hazards if it is operated by insufficiently trained personnel or if it is not used properly or in line with its intended use.
The unit must be operated by qualified or trained members of personnel (e.g. dental laboratory technicians or dental laboratory assistants) who:

- have acquired their knowledge through specialist training or instruction,
- are able to carry out the activities he/she does commensurate with his/her qualifications,
- are able to recognise possible hazards when working with the unit.

WARNING

If certain problems occur, which are not addressed in sufficient detail in these operating instructions, please contact the supplier to ensure your own safety.

3.1 Rules for operation

WARNING

- Avoid direct visual contact with the LED light emitted. Wearing eye protection is recommended.
- Do not operate the unit with heavily soiled reflectors or without a rotating plate.
- Do not operate the unit with the housing ventilation slots closed. Ensure that the distances between the wall and the unit are observed.

WARNING

Do not use the unit for the activities described below:

- Do not use the unit to cure items that contain residues or adhering material from flammable solvents or cleaners or are not completely dry and consequently may form a flammable or explosive gas-air mixture in the polymerization chamber.
- Do not use the unit to heat food products.
- Do not use the unit on materials other than those specified for exposure in section 3.

NOTE overheating

We recommend a pause of at least 1 minute between the individual polymerization cycles to allow the unit to cool down. Continuous operation at high polymerization temperatures, high ambient temperatures or poor ventilation may result in error messages due to overheating (unit protection). Continued operation is then interrupted so the unit can cool down.


4 Scope of delivery

Check whether all components are free from defects when delivered. Please contact your supplier for any complaints.

Delivery unit:

- cara Print LEDcure
- Operating instructions on USB stick
- USB connecting cable for updating cara Cure Guide
- 2 mains cables (Europe and USA / Japan)
- Safety goggles

5 Operating and indicator elements

1) Sliding door
2) Recessed handle for opening the door
3) On/off mains switch
4) Display
5) 3 operating buttons for interactive display
6) Rotating plate
7) Mains connection w/ fuses
8) USB port for Cure Guide Update and servicing purposes

6 Unit description and function

10 LED modules in a fully mirrored chamber are used to complete the final polymerization of the material being cured. The 50 high-power LEDs emit light radiation in the visible violet range (397 nm peak) as well as in the blue range (450 nm peak), which are optimally matched to the Kulzer materials in terms of wavelength and power. Homogeneous light output is achieved at every point of the chamber through the full mirroring. The display shows all indicators for setting and controlling the unit functions. 3 interactive buttons below the display are used to operate the unit.

NO UV light (< 380 nm) with significant power is emitted. Nevertheless, safety goggles to protect the eyes against bright light radiation during operation are being provided.

The pre-set Kulzer programs are permanently stored with the name of the application and parameters (power, temperature and time). Each customer is able to create and save individual programs in parallel (power, temperature and time). The allocated customer programs are numbered consecutively - the names cannot be edited individually.
The corresponding button is pressed to start a program. Exposure starts after the rotating plate located in the curing chamber begins to turn, the heating warms the chamber and the target temperature is reached. The remaining program runtime is shown in the display.
The sliding door is not locked during the program and is monitored by a limit switch. If the door is opened during the program run, the LEDs are switched off and an error message is displayed. Operation does not resume after the door is closed. The program is cancelled and must be restarted for safety. Likewise, it is not possible to start the program if the door is open.
The display shows when the time has elapsed and is accompanied by an sound signal.
The chamber is actively warmed by a heating register and is temperature-controlled. The item temperature is also affected in parallel by the LED power setting and the absorption of the light radiation. This is very much dependent on the material, and the LED power and temperature parameters must be individually evaluated and set for customer programs.

Gloves are recommended to protect against burns caused by the potentially very warm items after polymerization. Use suitable tools (tweezers / forceps) to insert or remove the items.

Do not place or cure items in the unit that still contain residues or adhering material from flammable solvents or cleaning fluids and are not completely dry.

To clean the interior, rotating plate and the LED modules, see section 9, ‘Cleaning’.
The item(s) are placed on the rotating plate in the chamber. There are no restrictions – except for the center of the rotating plate – when selecting the position and orientation of the items. Because of the axis of rotation, there are no LEDs in the center of the rotating plate that facilitate irradiation. Therefore, items are NOT to be placed in this area.

7 Placing the unit into operation

7.1 Transport

When transporting the unit, do so carefully in a horizontal position. Avoid shaking/impact. Protect the unit from moisture, liquid spray, vapors, heat, frost and any kind of contamination (dust, soiling, etc.). For dimensions and weight, see section 12 ‘Technical details’.

7.2 Setup

The cara Print LEDcure is only to be set up indoors on a firm, non-flammable surface (laboratory table, pedestal) ensuring it is in a horizontal, safe state. The ambient temperature cannot exceed 30°C (86°F). The environment must be clean and free from excessive dust.

Setup location: Table with a load capacity of at least 15 kg.

Setup surface (W x D x H): approx. 250 x 350 x 550 mm (max. height with sliding door open)

WARNING

Following the instructions below is imperative:

- The rear side of the unit must be at least 10 cm away from the wall to avoid obstruction of exhaust air.
- Do not place paper or other objects under the unit to avoid obstruction of the fan air supply.
- There is a risk of overheating the LED modules if the air intake or output required to cool the unit is not free of obstruction. The lifecycle of the LEDs could be affected. An error message appears when the LED modules overheat.

7.3 Component check before placing the unit into operation

WARNING

The functionality of the unit is only guaranteed when the components are free from defects. Be sure to thoroughly inspect the unit components before placing the unit into operation.

Make sure that:

- all unit components are not damaged.
- there is no residue or soiling present.
- defective parts are replaced.
- there is no damage to the mains switch, the connection socket or the mains cable insulation.
- the power supply fits the wall outlet (see section 7.4 Mains connection – country-specific).
- the mains socket has a functioning ground connection.

The unit inspection must be completed at regular intervals by the operator (see section 7.6).

7.4 Mains connection

On the back of the unit there is a cold appliance plug connection socket with integrated mains fuses (2 pcs.). The device is equipped with a wide-range power supply for worldwide operation. Power outlet must have a protective earth connection. The connection is made with the supplied power cables. Any power or extension cables provided by the customer must at least correspond to the rated current of the supplied power cables and the local connection regulations.

WARNING

After setup, the power supply must be accessible to be disconnected from the mains power supply.

Mains cable included in delivery		
	Use in Europe	Voltage: 230 V
	Use in USA / Japan	Voltage USA: 115 V / voltage Japan: 100 V

ELECTRIC SHOCK

Wet conditions and moisture may lead to a short circuit.

Observe the following:

- Only operate the unit in dry areas.
- Before placing the unit into operation, check whether the rated voltage corresponds to the rated voltage of the unit specified on the identification plate.
- Only connect the unit to a mains socket with an earth connection.
- Only use the mains connection cable supplied with the unit.
- Be sure to comply with the technical specifications when the mains cable is replaced. See section 12 ‘Technical details’. If country-specific power cables are used, make sure that the cables are certified in accordance with the electrical rating details and information of the relevant standards of the respective country.

7.5 Unit inspection

The unit inspection must be completed at regular intervals by the operator (e.g. weekly). Specifically check the mains connection, mains cable, plug and fuse holder for damage.

If there are any anomalies, the unit must be taken out of operation, ensured that it is not placed back into operation and the service department notified (see section 14). For more detailed information see section 9 ‘Cleaning’ and 10 ‘Repair’.

8 Operation

Button	Explanation
SETUP	Settings
RESET	Reset
INFO	Information
TEST	Display test
user program	user program
EDIT	Edit
YES	Yes
NO	No
ESC	Cancel / back
ERR	Error
DEL	Delete
END	End

8.1 Program interface and program selection

After the power is turned on, the welcome text appears on the screen. The screen will display for approx. 3 seconds. During this time the operator can press the ‘SETUP’ (left), ‘INFO’ (middle) or ‘Test’ (right) buttons. If no button is pressed, the program selection appears on the display.

SETUP

When the ‘Setup’ button is pressed, the left ‘down’ arrow button can be used to scroll between the selectable functions in the loop. The selected function is indicated by the arrow on the left side of the display. Use the right ‘EDIT’ button (pencil) to change the setting.

NOTE

To simplify the presentation in this manual, the images of the displays are shown in ENGLISH only. On the actual unit, the display text will be shown in the language selected.

Line 1: Language selection for display text settings: DE; GB; FR; IT; ES; PT; SE; NL; JP
Line 2: Temperature unit; settings: °C / °F.
Line 3: Choice of whether each press of a button should be acknowledged with a signal tone; settings: on / off
Line 4: Abbreviated program parameter indicators above each customer program; settings: on / off
Line 5: Reset i.e. delete all customer programs. The programs are reset to zero.


Before this is executed, the following confirmation screen appears.

After confirming (YES) or cancelling (NO), you are returned to the Setup screen which you can exit by pressing 'ESC'. The settings are saved automatically.

INFO

When the 'INFO' button is pressed, the following image appears:

The screen provides information on the status of the unit (hardware) and software version (cara Cure Guide and Firmware software). As well device serial-no. and total operating hours. Press the 'ESC' button to exit the screen.


When the 'Service' button is pressed, a submenu appears in English:


This information is reserved for service personnel for service purposes (read only). Press the 'ESC' button to exit the service area.

TEST

When the right button is pressed, a display test screen appears for 10 seconds and is only relevant for service purposes, as well:


Press any button to cancel the test and return to the Main screen, program selection:

The last selected / executed program is inversely placed at the top of the list. Using the arrow keys “up” (right) and “down” (left), you can scroll through the programs. The first programs (Names and number at the time of printing this manual) are established by Kulzer and have corresponding material names. Customer-specific programs are able to be programd individually (see section 8.2) and a sequential number is applied to the program name (user program 1; user program 2, etc.).

Example program selection images:

The program options are distinguished by the icon above the middle button as follows:

Kulzer programs (non-editable): 

Coded customer programs: 

Empty customer programs: 

Kulzer programs are able to be initiated immediately. The parameters are established for the respective material. **Coded customer programs** are able to be either executed or edited. If parameters are included in the program, they are indicated in the parameter line above the program name in abbreviated form.

8.2 Execute program

The program is selected by pressing the button (Enter/Middle) quickly. The ‘EDIT’ mode is activated for customer-specific programs by pressing on the middle button for a longer (> 1.5 s) period of time (beep tone; Edit description see section 8.3). If a customer-specific program is empty, only the EDIT icon appears and the function can be activated by quickly pressing the (middle) button.

After selecting the program, the Start screen appears with the program parameters (time, temperature and LED power):

If the temperature is still too high when initiating the program, the following notification screen appears as a hint:

The actual temperature indicator with the preceding arrow signifies that the chamber temperature is > 5°C above the target program temperature. The pause symbol in the bottom line indicates that initiating the programs is not currently permitted.

NOTE

Do not place an item into the chamber when it is overheated. The item material may be damaged.
If the temperature drops to a value < 5°C above target temperature, the pause sign in the middle changes to ▶ and the program can be initiated. If the chamber temperature drops below the target temperature, the target temperature is indicated again.
If the wrong program was selected by mistake, press 'ESC' (left) to return to program selection. The selected program is initiated by pressing 'Start' (middle).

NOTE

If the program contains faulty parameters, it is not able to be initiated and only 'ERR' appears on the display. If errors are still present, the program will also not be able to be initiated. An error (ERR) is indicated in the bottom line instead of the start symbol and the symbol of the present error is displayed (e.g.: Door open). Please see section 8.7 for possible causes and workarounds.

NOTE

There are no provisions for temporary changes to program parameters; if necessary, a new program must be created.

If all parameters are correct and the program was initiated using the middle button, the flashing 'Heating' icon indicates the preheating phase. The display shows the current process temperature in the middle and the target temperature and LED power on the right.

A program run can be cancelled within the first 15 seconds by pressing the middle 'Stop' button. After 15 seconds, this option disappears and the program will be run through completely.
If the program does not have a temperature preselected (setting '- - °C'), the preheating phase is skipped.
In the next image, an animated hourglass shows the program run. If over 1 minute, the remaining runtime is indicated in minutes:seconds; if under 1 minute, it is only indicated in seconds. The actual temperature is located on the right side of the display (or alternatively '- - °C,' if no temperature is pre-set).

After the program has run successfully, the end screen appears:

Press the ‘Enter’ button (middle) to exit this screen. The End screen remains on the display until a button is touched to confirm or the door is opened. Afterwards, the program returns to the Start screen of the previously selected program.

8.3 Coding customer-specific programs

If a customer-specific program is empty, only the EDIT icon appears and the function can be activated by quickly pressing the (middle) button. If customer-specific parameters are already stored, the icon is displayed. The program is able to be selected with a quick-press or, alternatively, edited with a long-press > 1.5 seconds:

If you have inadvertently entered into the Edit function, press the ‘ESC’ button to return to program selection. If you want to delete the customer program parameters, a confirmation prompt appears after pressing the ‘DEL’ button (right):

Only after confirming by pressing the ‘YES’ button are the program parameters permanently deleted and the display returns to program selection. Press the ‘NO’ button to stay in the editor screen.

Pressing the ‘Enter’ button (middle) will toggle between the individual program parameters to modify them. Run sequence:

1. **Power:** Select 10 % – 100 % in single % increments. Set the parameter using the ‘up’ (right) and ‘down’ (left) buttons. Press the ‘Enter’ button (middle) to confirm.

2. **Temperature:** Selection from 30 – 80°C (86 – 175°F) and ‘no temperature’ (‘- -°C). Set the parameter using the ‘up’ (right) and ‘down’ (left) buttons. Press the ‘Enter’ button (middle) to confirm.

3. **Time:** Select between 1 – 30 minutes. Set the parameter using the ‘up’ (right) and ‘down’ (left) buttons. Press the ‘Enter’ button (middle) to confirm. Editing is complete, the display jumps back to the program start screen (w/ indication of the program parameters).

NOTE

- The real item temperature is very much dependent on the absorption of the LED radiation by the material. Therefore, the LED power and temperature setting must be configured individually in combination with each customer-specific material.
- The operator is solely responsible for the customer-specific parameters and processes, as well as their configuration. Always follow the specifications provided by the manufacturer of the material.
- Certain combinations of LED power and temperature are not realistic based on operational experience. The following combinations to polymerise 3D printing resins can serve as a guide:

temperature										
	80°C/176°F	√	√	√	√	√	√	√	√	√
	70°C/158°F	√	√	√	√	√	√	√	√	√
	60°C/140°F	√	√	√	√	√	√	√	√	√
	50°C/122°F	√	√	√	√	X	X	X	X	X
	40°C/104°F	√	√	√	X	X	X	X	X	X
	30°C/86°F	√	√	X	X	X	X	X	X	X
	Temp. Off	√	√	√	√	√	√	√	√	√
		20 %	30 %	40 %	50 %	60 %	70 %	80 %	90 %	100 %
		LED Power in %								

Legend: The item temperature of combinations with a red ‘X’ will probably be exceeded. When the heating is switched off (Temp. „--°C“), a temperature is only set according to the absorption properties of the material or item. **The item temperature may deviate from the air temperature indicated.**

8.4 Working with the unit

- Switch on the unit at the front.
- Open the sliding door and place the 3D-printed item on the rotating plate. Close the door.

Do not place or cure items in the unit that still contain residues or adhering material from flammable solvents or cleaning fluids and are not completely dry.

NOTE

Make sure (if possible) that the item is not placed in the unit in the center of the rotating plate, but slightly to the side. Exposure to the attached LEDs from underneath the rotating plate is supported by positioning the item outside of the center of the rotating plate.

- Edit a customer program beforehand, if necessary, then use the interactive buttons to select a polymerization program, and press the ‘Start’ button.
- The program begins with the heating phase. The actual and target temperature are indicated. The timer starts when the target temperature is achieved. If no temperature is set, the timer starts immediately after the program is initiated.

The remaining runtime and the actual temperature are indicated on the display. After the polymerization time has elapsed, a signal tone sounds and the door can be opened.
If the door is opened before the program has finished running, an error message appears. The program is cancelled, the error message must be confirmed and the program must be restarted.

NOTE

The polymerization parameters are binding as stated in the material’s instructions for use. Therefore, Kulzer programs are not able to be modified (edited) whatsoever. The provisions and processing instructions governing the manufacture of medical devices apply to some 3D-printed materials. The user is responsible for the final product when there are deviations from the processing instructions / instructions for use.

NOTE

After the program has finished running, remove the 3D-printed item(s), check the rotating plate for soiling, clean it, if necessary, and close the door.

NOTE

When working with object holders, carriers, tables or other aids, unwanted discolouration of surfaces due to LED radiation is easily possible and probable. It is the sole responsibility of the user to assess the polymerization quality and, if necessary, to carry out a post-polymerization.

HOT SURFACE

Shortly after the polymerization process finishes, the interior of the unit, especially the heating register and the set 3D item, may be very hot. This is a potential burning hazard. Observe the following:
After the polymerization process ends – especially longer processes lasting > 10 min. – for you own safety, wear protective gloves and use tools (forceps/pliers) to remove the 3D items.

8.5 Software error messages

If an error occurs during the program run, the process is terminated and an error message appears:

The icon of the triggered error is indicated on the right side of the display. Possible errors include:

LED module temperature exceeded

LED light / power supply defect

Chamber temperature sensor defect

Door open

Rotating plate blocked or motor defective

Pressing the 'Enter' button acknowledges the error message and the display jumps back to program selection. If error conditions are still present when the program is initiated, starting it is not possible → new error message. In this case, please contact the Kulzer Hotline or responsible service partner.

8.6 cara Cure Guide update

The material parameters current at the time of production are already installed on the device. The cara Cure Guide version installed on the device is displayed on the info display after the device is switched on (see Chapter 8.1). It is the operator's duty to regularly check the Kulzer website for Kulzer cara Cure Guide updates. If a newer cara Cure Guide version is available on the Kulzer website (download area), first download this file [.ZIP] to a local PC / laptop / tablet and unpack it into a common directory (no cloud).

NOTE

The 'cara Cure Guide' is an integral part of this operating manual, must always be kept up to date and observed by all operators to avoid possible risks for material processing.

8.7 Install cara Cure Guide update

The device is equipped with a USB port on the back. This is exclusively for transferring material data in the cara Cure Guide from a PC/laptop/tablet. Devices connected to the USB port must meet the requirements of IEC 62368-1.

WARNING

The USB port must not be used for any purpose other than that specified in the following description. In particular, no external chargers, power banks, other storage media or other non-intended devices may be connected via the USB port. Risk of damage! Risk of damage!

Kulzer GmbH disclaims any warranty for damage to the device or consequential damage caused by improper use of the USB port.

NOTE

Please refer to the description contained in the .ZIP file for all steps for driver installation, operation of the setup program including the update of the new cara Cure Guide! In parallel, the operating personnel must be informed immediately about the update and any parameter changes, additions, notes, etc. that may be associated with it.

NOTE

The total number of programs (Kulzer + customer programs) is limited to 165. If new Kulzer programs are added, the list of available customer programs will be reduced starting with the highest number. Therefore, always save your customer programs with an ascending numbering.

8.8 Errors and how to clear them

DISCONNECT FROM MAINS

Risk of life-threatening consequences warning
Unplug from the power supply before working on the unit.

ELECTRIC SHOCK

High voltages are life-threatening when touched. There are no components inside the unit that are serviceable or repairable by the user which require opening up the unit. Therefore, users are prohibited from opening up the unit.

Error	Indicator lights	Signal tone	Cause	How to clear
The unit does not switch on		none	Unit disconnected	Establish mains connection
		none	Fuse defect	Check fuse in the mains connection plug on the rear side of the unit and replace, if necessary. If the error reoccurs, remove the unit from operation and notify Service.
		none	Mains cable defect	Replace mains cable
		none	Incorrect mains voltage	Check permitted input voltage (see section 12 'Technical details')
Polymerization does not start or is cancelled with error message	'Door open' ERR	3 x 500 ms	Door switch triggered	Close door correctly. Re-execute polymerization, if necessary.
	'LED defect' ERR	3 x 500 ms	One or several LED module defects	Contact Service.
	'LED temperature exceeded' ERR	3 x 500 ms	LED temperatur > 65°C	Let the unit cool down; check air supply and fan function (bottom and back); check ambient conditions (if necessary, air-condition the room, ventilate)
	'Temp. sensor defect' ERR	3 x 500 ms	Chamber temperature sensor defect	Contact Service
	'Fan signal fail' ERR	3 x 500 ms	Fan defect	Contact Service
Item will be very hot	No error message	none	Air intake or output obstructed	Check the unit environment or base, keep air intake and output free of obstructions
		none	Ambient temperature too high	Air-condition the room / ventilate
		none	Check polymerization parameters for customer programs	Set of parameters according to the manufacturer specifications. Reduce LED power and temperature and extend time, if necessary.

9 Cleaning

ELECTRIC SHOCK

The LED modules operate with high voltage and current. This presents a relevant electrical hazard. If moisture enters by way of cleaning and disinfecting agents, there is a hazard of short circuit and electric shock.

Observe the following:

- Switch off the unit and unplug the power supply before cleaning.
- Disassembling the housing is not part of the operator's maintenance work.
- Do not use spray bottles – only use slightly moistened wipes.
- Only use mild cleaning and disinfecting agents composed of alcohol, allow the unit to dry completely inside and outside before placing it back into operation.
- Do not use alkaline solutions, corrosive or abrasive substances, methanol, amines, technical solvents, organic compounds, methyl methacrylate and methacrylic acid methyl ester.
- When cleaning, be sure to avoid scratching the unit – especially inside the reflector chamber.
- Do not use pointed, sharp or abrasive objects to clean the unit.
- Remove soiling completely but gently.

NOTE

Scratches or soiling on the reflecting surfaces in the chamber reduce the efficiency of the irradiation and reduce LED performance. Insufficient or inhomogeneous polymerization may result.

10 Maintenance / repair

WARNING

Unit functionality is only guaranteed with original Kulzer replacement parts. Using replacement parts of unknown origin is prohibited.

The following parts can be replaced by the end user:

- Power cable
- Fuses (see section 12 'Technical details')
- Rotating plate

WARNING

There are no maintenance or service-relevant components inside the unit for the user/customer. Therefore, opening the housing is not necessary and prohibited for the user/customer. For repairs involving opening the device, please contact the responsible country branch or the regional Kulzer service partner! See section 14 'Service'.

11 Disposal

Recyclable materials are used to manufacture the cara Print LEDcure. These do not present any particular or unusual dangers for disposal. Kulzer GmbH is naturally willing to go beyond its legal obligations and take back used devices. The device is designed for a service life of approx. 8 years.

Disposal of old equipment according to WEEE
Electrical- and Electronic Equipment Act (ElektroG 3 from 01.01.2022)

The ElektroG3 implements the European WEEE Directive 2012/19/EU into German law. It regulates the placing on the market, the take-back and the environmentally friendly disposal of electrical and electronic equipment.

For detailed information on professional disposal of disused old devices, contact your dealer or Kulzer subsidiary in your country.

CAUTION

Marked equipment must not be brought to local waste disposal centers.

12 Technical details

Value / unit	
Rated voltage	100 V – 240 V
Mains voltage tolerance	+- 10 %
Rated frequency	50 / 60 Hz
Mains power cable	EU – H05V2V2-F 3x0.75mm² 250VAC/16A US – SJT3X18AWGWB105C C13/2.5m black 125VAC/15A
Fuse	2 x T6.3H250V
Power input:	600 VA
Dimensions (W x D x H)	approx. 250 x 350 x 550 mm (max. height with sliding door open)
Weight	12.8 kg
Safety category	Safety category I
Overvoltage categorie	II
Protection category	IP 20
Degree of pollution	2
Duty cycle (DC)	100 %
Light source	10 LED modules w/ 5 LEDs with 3 watt each (=150 Watt LED power)
Wavelengths emitted (primary spectrum)	40 LEDs 397nm peak; 370 – 420 nm = 95 % spectral distribution 10 LEDs 450nm peak; 430 – 470 nm = 95 % spectral distribution
Sound pressure level	<70 dB(A)


12.1 Transport and storage conditions

Temperature	–15°C to +40°C (5°F to 104°F)
Relative humidity	Less than 100 % rel. humidity

12.2 Operating conditions:


Temperature	10°C to 30°C (50°F to 86°F)
Humidity	80 % rel. humidity up to 25°C (77°F) 50 % rel. humidity up to 40°C (104°F)
Maximum elevation above sea level	2000 m

12.3 Identification plate


13 Wiring diagram

13.1 Mains, power supply / fan

13.2 Controller interface

14 Service

www.kulzer.com

If you have any questions, comments or suggestions do not hesitate to contact us. Please use the QR code:

Contacts in other countries and further information can be found directly with the adjacent QR code on our website www.kulzer.com

Manufacturing address for technical support:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

15 Document history

2021-01	First version of cara Print LEDcure
2023-03	Update: chapt. 8.5 error messages; chapt. 12.3 identification plate; chapt. 13 circuit drawings
2025-05	Update: chapt. 1.3 EC Declaration of Conformity; chapt. 8.6, 8.7 Cara Cure guide; chapt. 11 Disposal; Editorial update of various chapters

Subject to change without notice. Dated: 2025-05

© 2025 Kulzer GmbH

All Rights reserved. No parts of this work may be reproduced in any form or by any means – graphic, electronic, or mechanical, including photocopying, recording, taping, or information storage and retrieval systems – without the written permission of the publisher.

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All rights reserved. 본 간행물의 일부 또는 전부를 발행인의 사전 서면 허가 없이