dima Print Guide & Tray

66081930/00

dima Print Guide & Tray

- DE Gebrauchsanweisung
- GB Instructions for use
- FR Mode d'emploi
- ES Instrucciones de uso
- IT Istruzioni per l'uso
- PT Instruções de uso
- NL Gebruiksaanwijzing
- KR 사용 설명서

KULZER
MITSUI CHEMICALS GROUP

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Stand: 2021-08

GB Instructions for use dima Print Guide & Tray

Medical Device – for use by dental health care professionals only.

1. Intended purpose

Photopolymer for the 3D printing of dental accessories for treatment.
This dental accessory is suitable for patients requiring dental treatment for the following indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.
Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically, or they are accessories in the rehabilitation process.

1.1 Indications
For 3D printing of surgical guides and individual impression trays

1.2 Contraindications
In case of known or suspected allergy to components of the product, its use is contraindicated. The product contains allergens: (Meth-)acrylate.

1.3 Potential side effects
Hypersensitivities to the product or its components cannot be excluded in individual cases. In case of suspicion please contact the manufacturer for ingredients information.

2. Product features

2.1 Package sizes/Delivery forms
1000 g bottle

2.2 Accessories
The material can be printed using printers with a wavelength of 405 or 385 nm (currently released devices – see table below). The cleaning takes place in the cleaning device, e.g. cara Print Clean *pro* or ultrasonic bath.
Post-curing of the material must be carried out with a suitable exposure device.

3D printer / Manufacturer		Post-curing unit / Manufacturer	
cara Print 4.0	Kulzer	• HiLite power 3D • cara Print LEDcure	Kulzer
Layer thickness:		Tray= 100 µm and 150 µm Guide= 100 µm	

Always use the latest firmware version of your 3D printer.

2.3 Composition
(Meth)acrylate-containing monomers, photo initiators

3. Product use

3.1 Processing properties
Recommended wall thickness 3 mm

Materials	Guide & Tray
Homogenization	Shake for five minutes
Cleaning solution	Isopropanol (pure)
Pre-cleaning	3 min.
Post-cleaning (fresh rinse)	2 min.
Post-curing	• HiLite power 3D 2 x 5 min (1x top, 1x bottom) • cara Print LEDcure in accordance with the exposure program

3.2 Application subitems
When filling the vat or when refilling the liquid into the bottle, it must not be mixed with other photopolymers, as otherwise negative changes in the product properties (colour, mechanical and chemical properties) may occur.

a. Shake vigorously in the original container for five minutes before use.
Inadequate shaking may result in color variations and printing errors.
b. Print the objects according to the instruction manual supplied with the 3D printer. Select the print program dima Print Guide & Tray. Processing range: 18 – 30°C (64 – 86°F), 30 – 80 % r.F.
c. Thoroughly remove remaining liquids using compressed air.
d. Pre-clean the printed objects for 3 minutes in Isopropanol using the appropriate supplies (see Item 2.2).
Note: Thoroughly rinse the inside of the structures with the spray bottle (Isopropanol).
e. Repeat step **d** using pure Isopropanol for 2 minutes.
Note:
- The objects should not remain in the Isopropanol for more than 5 minutes, otherwise they will be damaged.
- After cleaning, make sure the printed objects are dry and free from solvent residues.
For Trays: move forward with **f**
For Surgical Guides: Remove Support structures if necessary and insert drill sleeves into the drilling holes and then move forward with **f**
f. Place the printed object for final polymerization in the photo polymerization unit according to the table above.
Note: The final properties depend on the adherence to the recommended post-curing process.
For Surgical Guides: Check proper fitting of inserted sleeves after photo polymerization.
g. Carefully separate the supports with standard separating discs. The finishing is done with cross-toothed cutters and, if necessary, silicone and rubber polishers (e.g. Signum toolkit) or with the polishing motor. Polishing of the outside is mandatory for surgical guides.
h. Store the remaining printing liquid from the vat protected from light and fill it back into the bottle after 12 hours at the latest, working with a funnel and sieve to hold back any particles. Always carefully stir the liquid with a silicon spatula before resume printing.

4. General information

Please, carefully observe the following safety instructions and all processing notes in other sections. This product is intended to be used only in accordance with this instructions for use. Any use of this product inconsistent with these instructions is at the discretion and sole responsibility of the dental health care professional.

4.1 Warnings and safety information
Avoid contact with skin and eyes. May cause skin and eye irritation. Skin contact may cause sensitisation. Wear suitable protective clothing, gloves and safety goggles. In case of contact of the product with skin, immediately wash with plenty of water and soap. In case of contact with eyes, rinse with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If symptoms persist, consult an eye specialist. May cause respiratory irritation. Avoid breathing of dust/fume/gas/mist/vapours/spray. Toxic to aquatic life with long lasting effects.

4.2 Storage conditions
Do not use the material after the expiration date. Store at room temperature, not above 25°C (77°F). Keep away from sunlight. Store in a dry place.

4.3 Disposal information
Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system. European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

5. Cleaning & Disinfection (applicable for surgical guides only)

Restriction of reprocessing: The printed surgical guides are intended for single processing and single use. Please ensure that the information on cleaning & disinfection is passed on to the dental practitioner accordingly.

5.1 Point of use
No particular requirements.

5.2 Containment and transport
Transport in a closed container.

5.3 Preparation
Personal protective equipment (gloves, water repellent protective skirt, face protection mask or protective glasses and mask).

5.4 Cleaning and Disinfection: Manually
- Equipment: Non-protein-fixing VAH-listed disinfectant with proofed cleaning effect based on quaternary compounds (ID 212, Dürr Dental was validated), lint-free wipes (Braun Wipes Eco #19726 were validated), soft cleaning brush (Interlock #09098 and #09316 were validated), tank/basin for cleaner/disinfectant, tap water/flowing water (20 ± 2°C, minimum drinking water quality), medical compressed air (according to the European Pharmacopoeia).
1) The cleaning solution is prepared in accordance to the manufacturer`s instructions for the detergent (ID 212, Dürr Dental; 2% was validated).
2) Soak the medical device completely in the cleaning solution.
3) Make sure that all surfaces are completely moistened with cleaning solution.
4) Total exposure time in the cleaning solution in accordance with the manufacturer`s instructions (5 minutes was validated).
5) During the exposure time, brush the hard-to-reach areas of the immersed medical device with a soft brush for minimum 1 minute (1 minute was validated). Pay attention to the critical hard-to-reach areas and cavities where it is not possible to visually assess the cleaning effect.
6) Remove the medical device from the cleaning solution.
7) Rinse the medical device completely under tap water for minimum 1 minute (1 minute was validated) to remove the cleaning solution completely.
8) Check on cleanliness, if dirt is visible, repeat the steps above.
9) The disinfectant is prepared in accordance with the manufacturer`s instructions for the detergent (ID 212, Dürr Dental; 2% was validated).
10) Immerse the medical device completely in the disinfectant.
11) Make sure that all surfaces are completely moistened with disinfectant.
12) During the exposure time, brush the hard-to-reach areas of the immersed medical device with a soft brush for minimum 1 minute (1 minute was validated). Pay attention to the critical hard-to-reach areas and cavities where it is not possible to visually assess the cleaning effect.
13) Exposure time of the disinfectant solution in accordance with the manufacturer`s instructions (ID 212, Dürr Dental; 5 minutes was validated).
14) Remove the medical device from the disinfectant solution.
15) Rinse the medical device completely with deionized water for minimum 1 minute (1 minute was validated) to remove the disinfectant completely.
16) Dry with a lint free gauze and/or dry with medical compressed air (according to the European Pharmacopoeia).

5.5 Maintenance, Inspection and Testing
All surgical guides must be visually checked for damage. Damaged medical devices must be sorted out and not be used.

5.6 Packaging
- Equipment: Sterilization pouches (steriCLIN #3FKFS230110 and 3FKFS230112 was validated). A suitable method (sterile barrier system) must be used for packaging the products. Packaging according to DIN EN ISO 11607 (single).
- Single: Use a sterile barrier system (e.g. see-through sterilization pouches or sterilization fleece) according to DIN EN ISO 11607, which is intended by the manufacturer for steam sterilization. The products can be packed single or double. The packaging must be large enough so that the seal is not under tension (For validation it was double wrapped in see-through sterilization pouches).
Note:
- After the sealing process, the seal must be visually inspected for possible faults.
- When the package is damaged the device must be newly packaged and sealed.

5.7 Sterilization
- Equipment: Steam sterilizer according to DIN EN 285 or DIN EN 13060 with type B process. In Germany, based on the KRINKO/BfArm recommendation, a holding time of at least 5 minutes at 134°C is recommended
- Pre-vacuum process, 132°C and sterilization time at least 3 min (longer holding times are possible) (132°C, 90 seconds was validated).
1) Place the packaged medical device in the sterilizer chamber.
2) Start the sterilization program.
3) At the end of the sterilization program remove the device and let it cool down.
4) Check the packages for damages or moisture penetration. Rejected packaging must be considered non-sterile.

5.8 Storage
Storage (at least protected against dust and moisture) and storage time according to user specifications.

5.9 Additional information
The processing of medical devices should be done with validated processes.

6. Complaint notice

Please quote batch number and article number in all correspondence about the product.

According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

Additional information can be found on the download page at:
www.kulzer.com/cara-print

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