



FICHA TÉCNICA

DIQUES DE GOMA DENTALES PROCLINIC		
REFERENCIAS	0125: medios verdes 12500: medios azules 12501: finos verdes 12502: gruesos verdes 12503: medios verdes	
FABRICANTE / RESPONSABLE	Sanctuary Health Sdn Bhd No 16, Persiaran Perindustrian Kanthan 1, Kawasan Perindustrian Kanthan, 31200 Chemor, Perak, Malaysia.	
CLASIFICACIÓN	Producto Sanitario clase I, según el Reglamento (UE) 2017/745 sobre productos sanitarios.	
DESCRIPCIÓN	Diques de goma dentales.	
COMPOSICIÓN	Látex de caucho natural.	
ESPECIFICACIONES	Dimensiones: - Ref. 0125, 12500, 12501, 12502: 152 x 152 mm - Ref. 12503: 127 x 127 mm Resistencia a la tracción: 24,0 MPa (min) Alargamiento de rotura: 700% (min) Aroma: menta (verde) / normal (azul) Color: verde / azul	
ENVASE	36 unidades / caja (para 152 x 152 mm) 52 unidades / caja (para 127 x 127 mm)	
VIDA ÚTIL	5 años a partir de la fecha de fabricación	
FINALIDAD PREVISTA	Utilizar como una barrera en ortodoncia.	
INFORMACIÓN SOBRE LA SEGURIDAD	En caso de contacto con los ojos, enjuagar con abundante agua	
CONDICIONES DE TRANSPORTE Y ALMACENAMIENTO	Manipulación: Para uso dental profesional solamente. Almacenamiento: Almacenar por debajo de 30°C y un lugar fresco y seco.	

TECHNICAL DATA SHEET

PROCLINIC DENTAL DAM	
REFERENCE 0125: medium gauge, green 12500: medium gauge, blue 12501: thin gauge, green 12502: heavy gauge, green 12503: medium gauge, green	
MANUFACTURER	Sanctuary Health Sdn Bhd No 16, Persiaran Perindustrian Kanthan 1, Kawasan Perindustrian Kanthan, 31200 Chemor, Perak, Malaysia.
CLASSIFICATION	Medical Device, class I, according to Regulation (EU) 2017/745 on medical devices.
DESCRIPTION	Dental dams
COMPOSITION	Natural latex rubber
SPECIFICATIONS	Dimension: - Ref. 0125, 12500, 12501, 12502: 152 x 152 mm - Ref. 12503: 127 x 127 mm Physical strength: Tensile Strength: 24.0 MPa (min) Ultimate Elongation: 700% (min) Scent: Mint (Green) / Plain (Blue) Colour: Green / Blue
PACKAGING	36 pieces / box (for 152 x 152mm) 52 pieces / box (for 127 x 127 mm)
SHELF LIFE	5 years from the manufacturing date
INTENDED USE	Use as a barrier for orthodontic procedures
SAFETY INFORMATION	In case of eye contact, rinse with generous amount of clean water.
TRANSPORT AND STORAGE CONDITIONS	For professional dental use only. Store below 30 °C in a cool dry place.