


<b>TECHNICAL DATA SHEET</b>	
<b>Extracting Forceps "51"</b> <b>PROCLINIC</b>	
<b>REFERENCES</b>	51
<b>MANUFACTURER</b>	Carl Martin GmbH Neuenkamper Str. 80-86 42657 Solingen, Germany
<b>CLASSIFICATION</b>	Medical Devices class I, according to Directive 93/42/EWG Medical Device Regulation (EU) MDR 2017/745
<b>DESCRIPTION</b>	Extracting Forceps, upper jaw roots
<b>COMPOSITION</b>	Surgical stainless steel
<b>SPECIFICATIONS</b>	Autoclavable up to 134°C. Resistant to the chemical sterilant. Supplied in "Non-sterile" condition.
<b>PACKAGING</b>	1 unit
<b>INTENDED USE</b>	Dental forceps, root remnants and splinter forceps are used to remove (extract) teeth and roots in the upper and lower jaw.
<b>SHELF LIFE</b>	Not applicable
<b>SAFETY INFORMATION</b>	Intended for use by dental professionals only. Must be checked for fault-free operation (fractures, cracks, malfunctions etc.) and be cleaned, disinfected and sterilized before each use. New instruments are contaminated with oils and fats from the manufacturing process and must therefore be cleaned thoroughly before the first sterilisation with a suitable cleaning agent.
<b>HANDLING AND STORAGE</b>	Sterile products must be stored in a dry, clean and dust-free environment at temperatures between 5°C and 40°C.

Solingen, 03.04.2023

Verena Gerlach

  
 Carl Martin GmbH  
 Neuenkamper Str. 80-86  
 42657 Solingen  
 Tel. +49 (0)212 / 81 00 44  
 Fax +49 (0)212 / 8 73 45  
 info@carlmartin.de - www.carlmartin.de

