


TECHNICAL DATA SHEET	
Needle holder "1160/14" PROCLINIC	
REFERENCES	1160/14
MANUFACTURER	Carl Martin GmbH Neuenkamper Str. 80-86 42657 Solingen, Germany
CLASSIFICATION	Medical Devices class I, according to Directive 93/42/EWG Medical Device Regulation (EU) MDR 2017/745
DESCRIPTION	Needle holder Mathieu 14 cm
COMPOSITION	Surgical stainless steel
SPECIFICATIONS	Autoclavable up to 134°C. Resistant to the chemical sterilant. Supplied in "Non-sterile" condition.
PACKAGING	1 unit
INTENDED USE	A needle holder is used to hold (secure) needles and suture material. They are used after surgical procedures to close wounds.
SHELF LIFE	Not applicable
SAFETY INFORMATION	Intended for use by dental professionals only. Must be checked for fault-free operation (fractures, cracks, malfunctions etc.) and be cleaned, disinfected and sterilized before each use. New instruments are contaminated with oils and fats from the manufacturing process and must therefore be cleaned thoroughly before the first sterilisation with a suitable cleaning agent.
HANDLING AND STORAGE	Sterile products must be stored in a dry, clean and dust-free environment at temperatures between 5°C and 40°C.

Solingen, 08.02.2024

*O. Gellach*

Carl Martin GmbH  
Neuenkamper Str. 80-86  
D-42657 Solingen  
Tel. +49 (0)212 / 81 00 44  
Fax +49 (0)212 / 8 73 45  
info@carlmartin.de · www.carlmartin.de

