

Surgicam systems

4K **Full HD 1080P**



Features	Surgicam Lite 4K	Surgicam HD EasyRC	Surgicam 4K EasyRC	Surgicam AVC
Sensor	CMOS 1/2.3"	CMOS 1/5.8"	CMOS 1/2.3"	CMOS 1/2.3"
Video rec. resolution	3840x2160/25p	1920x1080/50i	3840x2160/25p	3840x2160/25p
Lens	10x Optical	30x Optical	10x Optical	20x Optical
Focus	Manual	Auto/Manual	Auto/Manual	Auto/Manual
Medical grade	-	-	-	Yes
Streaming	-	-	-	Yes
IP Control	-	-	-	Yes
Remote Control	Yes (LiveView)	Yes (IR)	Yes (IR)	Yes (Touchscreen)
Memory recording	Yes (microSD)	Yes (microSD)	Yes (microSD)	Yes (2 x SD)
Video Outputs	microHDMI	microHDMI	microHDMI	HDMI, 3G-SDI
Camera head mass	140 g	305 g	700 g	333 g
Cables	included	included	included	included
Colour	White	Black	Black	White
Options	Footswitch, VideoCapture HD, HVO Recorders, Arms, Monitors, Printers			



Monitors

Surgicam

Lamps & Camera Arms

Recorders

Printers

IR Remotes

Surgicam Control Unit

Twin-Arm

C-Arm

W-Arm

S-Arm




21-25

33-38

33-85

IGLU 65/*65-C/100



Lux: 50,000
Light head: 215mm

Lux: 70,000
Light head: 330mm

Lux: 85,000
Light head: 330mm

Lux: 65,000/65,000/95,000
Light head: 300mm
*Adjustable colour temp.:
3,500 - 4,000 - 4,500