


# Tooth Preparation



## MICROETCHING MADE CLEAN AND EFFICIENT

- Precise surface microetching
  - Increases surface retention
  - Insures properly prepped surface
  - Maximizes bond strength on all surfaces
- 4 hole handpiece
  - Standard airline attachment
  - Easy chairside engagement
  - Stainless steel construction
- Pre-loaded cartridge tips
  - Disposable, no-mess application
  - No aluminum oxide management
  - 50 micron aluminum oxide
  - 2 tip sizes
    - Small (6 teeth)
    - Large (18 teeth)

## CLINICAL APPLICATIONS

- Lingual bonding
- Artificial substrate bonding
- Molar prep for bonds or bands

## ETCHMASTER®

### STARTER KIT

Etchmaster Starter Kit	#EMHA
------------------------	-------

### SMALL TIPS

Regular Tips 25 Pk	#EMTRS
Lingual Tips 25 Pk	#EMTLS
Regular Tips 50 Pk	#EMTRS50
Lingual Tips 50 Pk	#EMTLS50

### LARGE TIPS

Regular Tips 25 Pk	#EMTRL
Lingual Tips 25 Pk	#EMTLL
*90 Micron 25 Pk	#EMTRL90

- \* Etchmaster Handpiece Starter kit includes
  - 4 Hole Handpiece Attachment
  - (5) Regular Small aluminum Oxide Tips
  - (5) Lingual Small aluminum Oxide Tips
- \* Recommended for removing composite from bracket bases and microetching stainless steel appliances. Not recommended for use on enamel.



SURFACE MICROETCHING



**INSTRUCTIONAL VIDEO**

Bonding to Various Artificial Substrates  
Available at RelianceU.tech