



TECHNICAL SHEET

COTTON PELLETS

I.1 Product description

Soft and consistent in shape, the highly absorbent cotton pellets are perfect for drying cavities or applying solutions.

I.2 Classification

Cotton pellets are classified as medical devices of class I according to the directive 93/42/EEC.

I.3 Product specification

Cotton pellets are provided in five different sizes:

Cotton Pellet no. 1: Ø 0.71 cm, ≈ 29 mg
Cotton Pellet no. 2: Ø 0.55 cm, ≈ 12 mg
Cotton Pellet no. 3: Ø 0.39 cm, ≈ 8 mg
Cotton Pellet no. 4: Ø 0.32 cm, ≈ 4 mg
Cotton Pellet no. 5: Ø 0.09 cm, ≈ 3 mg

Catalogue references:

N°	Refill boxes for dispenser A		Pellets in bulk for dispenser B	
	Weight	Reference	Weight	Reference
1	6 x 3 g	33251	32 g	33261
2	6 x 3 g	33252	32 g	33262
3	6 x 3 g	33253	18 g	33263
4	6 x 3 g	33254	18 g	33264
5	6 x 3 g	33255	18 g	33265

I.4 Storage conditions

Cotton pellets should be kept in a dry environment.

I.5 Shelf-life

Cotton pellets have no shelf-life.



I.6 Packaging

Cotton pellets are packaged in cardboard boxes.

I.7 Safety Data Sheet

Not applicable for this kind of product.