NO CONTACT


medel

No touch forehead thermometer



Medel **NO CONTACT** is an infrared thermometer that is able to detect the body temperature without contact. It is multifunctional and measures accurately body temperature as well as the temperature of objects and liquids.

It is fitted with a large backlight display so that the values can be easily read during day and night. It detects the temperature in one second with a sound indicator at the end of the reading, which can be disabled in night mode. Class IIa Medical Device.



Temperature values

Red LED light for fever alarm



TECHNICAL SPECIFICATIONS	NO CONTACT body	NO CONTACT objects/liquids
Measuring range	34°C \sim 43°C (93,2°F \sim 109,4°F)	0°C \sim 100°C (32,0°F \sim 212,0°F)
Calibration accuracy	36°C \sim 39°C : +0,2°C else: +0,3°C	$<20^{\circ}C:+1^{\circ}C>20^{\circ}C:+5\%$
Weight	80g (with batteries)	
Size	128,5mm X 48,83mm X 38,85mm	
Working distance	0,5 ~ 3 cm	
Low battery indicator	\checkmark	
Body and objects temperature	\checkmark	
Date and time		
Fever indicator	\checkmark	
Memory	30 memories	



Also measures objects/liquids temperature



