

	TECHNICAL DATA SHEET	
		DATE: 01/06/18
	Proclinic	20170509

<b>PRODUCT</b>	DISPOSABLE BRUSHES FOR ADHESIVE-100u	
<b>REFERENCE</b>	2145	
<b>COMPANY/Manufacturer</b>	HUANGHUA PROMISEE DENTAL CO., LTD MOULD INDUSTRY GARDEN EMPOLDER SECTION HUANGHUA CITY HEBEI PROVINCE CHINA	
<b>CLASSIFICATION</b>	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
<b>PRODUCT DESCRIPTION</b>	DISPOSABLE PENCILS FOR ADHESIVE-100u Replaceable brush tips, adopt imported Dupont hair, large size and small size for your choice according to your requirement.	
<b>SPECIFICATIONS</b>	(stick) 3.2*14mm Hair:(Large)8-9mm, Tip color: Black	
<b>PACKAGING</b>	100pcs per bag	
<b>SHELF LIFE</b>	3 YEARS	
<b>INTENDED OF USE</b>	brushes tips is used to adhesive. It can bend in all angles, the dose can be controlled by it and it can reach at any narrow place and which it can be used for etch gel.	
<b>SECURITY</b>	This brushes tips needn't sterilization before use. Gingival brushes is a single use product. Pls do not re-use it. When disposal, please follow the related regulation in your region	
<b>TRANSPORT AND STORAGE CONDITIONS</b>	The brush holder r shall be storage in a temperature range: -20C°---70C°. Under the specified storage condition, the expiration is 3 years.	



和清逸  
2017.05.09