


REFERENCIA	16113500	CÓDIGO EAN	842965900608
PRODUCTO	Equipo dentista (implantación)		

## COMPONENTES

<p>1 Paño quirúrgico (SMS) 60x60cm</p> <p>3 Bata Quirúrgica Protectora p.p. 28/30 grs.con puños. No estéril.</p> <p>1 Guantes Látex c/carpetilla talla "G" (par)</p> <p>2 Guantes Látex c/carpetilla talla "M" (par)</p> <p>3 Mascarilla quirúrgica con cintas</p> <p>2 Gorro enfermera elástico posterior</p> <p>1 Gorro cirujano pp.cintas</p> <p>2 Paño quirúrgico (L) 50x60cm</p> <p>1 Paño quirúrgico (L) 75x100cm</p> <p>1 Paño quirúrgico fenestrado sin adhesivo (IS) 50x60cm</p>	
---	--

## PRESENTACIÓN

**Unidades embalaje 12**

**Presentación Estéril por óxido de Etileno**

**Alta 1/01/2000**

**Revisión 8/07/2016**

## NORMATIVA

Todos los productos de Cardiva Soluciones Integrales disponen de marcado CE, cumpliendo con la Directiva Europea 93/42/CE relativa a productos sanitarios.


Cardiva Soluciones Integrales tiene implantado y certificado un sistema de gestión de calidad y medio ambiente acorde a las normas ISO 9001; ISO 13485 e ISO 14001.

## COMPOSICIÓN

### 1 Paño quirúrgico (SMS) 60x60cm

Ancho 60cm; Largo 60cm.

- Tejido SMS (45 g/m<sup>2</sup>): Polipropileno formado por tres capas (Spunbonded-Metblown-Spundonded).



### 3 Bata Quirúrgica Protectora p.p. 28/30 grs.con puños. No estéril.

Largo 116cm; Ancho 142cm; Mangas 56cm; Puños 6cm.

- Tela sin tejer polipropileno 100% de 28/30 g/m<sup>2</sup>.
- Puños: poliéster.




### 1 Guantes Látex c/carpetilla talla "G" (par)

### 2 Guantes Látex c/carpetilla talla "M" (par)

### 3 Mascarilla quirúrgica con cintas


Medida de las cintas 42cm; Largo mascarilla 17,5cm; Ancho mascarilla 9,5cm.



### 2 Gorro enfermera elástico posterior

Largo de banda 65cm; Ancho de banda 15cm; Diámetro corona 20,5cm; Elástico relajado 10cm.


- Polipropileno 100%



## 1 Gorro cirujano pp.cintas

Largo banda 61cm; Ancho banda 11cm; Largo de tapa 23,5cm; Ancho de tapa 16cm


- 100% polipropileno tipo spunbonded 20 g.



## 2 Paño quirúrgico (L) 50x60cm

Ancho 50cm; Largo 60cm.


- Tejido Rayón de 24 grs/m<sup>2</sup>. (100% Viscosa).



## 1 Paño quirúrgico (L) 75x100cm

Ancho 100cm; Largo 75cm.


- Tejido Rayón de 24 grs/m<sup>2</sup>. (100% Viscosa).



## 1 Paño quirúrgico fenestrado sin adhesivo (IS) 50x60cm


Ancho 60cm; Largo 50cm; Diámetro fenestra 8cm.

- Tejido IS (56 g/m<sup>2</sup>): PP 30g laminado PE 25μ



REFERENCE	16113500	BAR CODE	842965900608
PRODUCT	Dental Implant Set		

## COMPONENTS

<ul style="list-style-type: none"> <li>1 Surgical drape (SMS) 60x60cm</li> <li>3 Protector Gown p.p. 28/30 gr.w/cuffs Not sterile.</li> <li>1 Guantes Látex c/carpetilla talla "G" (par)</li> <li>2 Guantes Látex c/carpetilla talla "M" (par)</li> <li>3 Surgical mask with tapes</li> <li>2 Nurse Cap Elastic back</li> <li>1 Surgeon Cap pp. with laces</li> <li>2 Surgical Drape (L) 50x60cm</li> <li>1 Surgical Drape (L) 75x100cm</li> <li>1 Fenestrated Surgical Drape without adhesive (IS) 50x60cm</li> </ul>	
--	---

## PRESENTATION

**Packing Units** 12

**Presentation** Oxide Etylene

**Entry date** 1/01/2000

**Review** 8/07/2016

## APPLICABLE STANDARD

All Cardiva Integral Solutions' products have CE mark according to European Directive 93/42/CE concerning Medical Devices.


Cardiva Integral Solutions has implemented and certified a quality and environment management system according to the standards ISO 9001; ISO 13485 and ISO 14001.

## COMPOSITION

1 Surgical drape (SMS) 60x60cm

Width 60cm; Length 60cm.

- SMS Fabric (45 g/m<sup>2</sup>): Polypropylene made up of three layers (Spunbonded-Metblown-Spundonded).



3 Protector Gown p.p. 28/30 gr.w/cuffs Not sterile.

Length 116 cm; Width 142 cm; Sleeve 56 cm; Cuff 6cm.

- Polypropylene 100% non-woven (28/30g/m<sup>2</sup>)
- Cuffs: Polyester.




1 Guantes Látex c/carpetilla talla "G" (par)

2 Guantes Látex c/carpetilla talla "M" (par)

3 Surgical mask with tapes


Tapes 18.5cm; Length 17.5cm; Width 9.5cm.



2 Nurse Cap Elastic back

Band long 65cm; Band width 5cm; Top diameter 20,5cm; Relaxed elastic 10cm.


- Polypropylene 100%



## 1 Surgeon Cap pp. with laces

Long of the band 61cm; Width of the band 11cm; Top 23,5x16cm.


- 100% polypropylene



## 2 Surgical Drape (L) 50x60cm

Width 60cm; Length 60cm.


- Rayon Fabric 24 grs/m<sup>2</sup>. (100% Viscose).



## 1 Surgical Drape (L) 75x100cm

Width 100cm; Length 75cm.

- Rayon Fabric 24 grs/m<sup>2</sup>. (100% Viscose).



## 1 Fenestrated Surgical Drape without adhesive (IS) 50x60cm

Width 60cm; Length 50cm; Aperture diameter: 8cm

- IS Fabric (56g/m<sup>2</sup>): PP 30g/m<sup>2</sup> laminated PE 25μ.

