




	TECHNICAL SHEET	DATE
		14-05-18

DISPOSABLE DENTAL TRAY 28x18 COLORED - ROUGH BOTTOM		IMAGE :	
REFERENCE	BD4182.XX	CLIENTE:	
PACKAGING:	CONVENCIONAL	Customer ref:	
DRAWING:	☒ ADJUNTO		
			
MATERIAL :			
P.E.T.			
THICKNESS :			180 +/- 5% MICRONS
COLOR :			
			DARK BLUE, CYAN, ORANGE, YELLOW
AVERAGE UNIT WEIGH:		14,50	grs.
PACKAGING:		TECHNIQUE: :	
		THERMOFORMED	
EMBALAJE :	LOGO :	PROPIO CLIENTE	DISTRIBUTION :
	BOX MODEL:	14 click coin blanca	
	SIZE:	400X300X110 mm	
QUANTITY x BOX :	400 UNITS		2 piles of 200 units each, in a unique polyethylene bag
PALETIZATION:			
PALLET:	EUROPALET 800 X 1200		
MAXIMUM HEIGH:	1.780,00 mm		
DISTRIBUTION:	18 boxes x 5 stages + 6 boxes		
QUANTITY x PALLET:	38.400 UNITS		
FINAL WRAPPING:	Shrink-Wrapping		
Packaging observations:			
REGULATION:			
Medical device class 1, European directive 93/42 CEE			
QUALITY DEPT.	CLICK COIN, S.L.	CUSTOMER COMPLIANCE	
QUALITY MANAGER	Gemma Bochaca Royo	RESPONSIBLE:	
DATE:	14-05-18	DATE:	
SIGNATURE:		SIGNATURE:	(If needed)