



DISPOSABLE DENTAL TRAY 18x14 COLORED - ROUGH BOTTOM

REFERENCE: VBD4278.BL CLIENT:

PACKAGING: CONVENCIONAL Customer ref:

DRAWING: ADJUNTO



MATERIAL : P.E.T.

THICKNESS : 180 +/- 5% MICRONS

COLOR :

DARK BLUE, CYAN, ORANGE, YELLOW

AVERAGE UNIT WEIGH: 7,25 grs. TECHNIQUE: : THERMOFORMED

PACKAGING:

EMBALAJE :	LOGO :	PROPIO CLIENTE <input type="text"/>	DISTRIBUTION :
	BOX MODEL:	17 click coin blanca <input type="text"/>	
	SIZE:	305x200x110 mm	

QUANTITY x BOX : 400 UNITS

2 piles of 200 units each, in a unique polyethylene bag

PALETIZATION:

PALLET:	EUROPALET 800 X 1200 <input type="text"/>
MAXIMUM HEIGH:	1.550,00 mm
DISTRIBUTION:	20 boxes x 7 stages
QUANTITY x PALLET:	56.000 UNITS
FINAL WRAPPING:	Shrink-Wrapping
Packaging observations:	<input type="text"/>

REGULATION: Medical device class 1, European directive 93/42 CEE

QUALITY DEPT.	CLICK COIN, S.L.	CUSTOMER COMPLIANCE	
QUALITY MANAGER	Gemma Bochaca Royo	RESPONSIBLE:	<input type="text"/>
DATE:	14-05-18	DATE:	<input type="text"/>
SIGNATURE:		SIGNATURE:	<input type="text"/>

(If needed)