



AKZENTA ONE SOFT COVER TISSUE-PE

DISPOSABLE HEADREST COVER TO PROTECT THE PATIENT



CONFIGURATION OVERVIEW

Type: Non-Sterile
Latex free

Features: The special features of the materials offer high absorbency and waterproofing, so to assure a high comfort to the patient.

TECHNICAL SPECIFICATIONS

1 ply of pure wadding of long fiber cellulose (TISSUE) partially laminated with a second ply of high density polyethylene (21gsm tissue with 17 gsm PE laminated)

Dimensions: 33x34 cm (tolerance \pm 1%)

PRODUCT VARIATIONS

Art-No.	HRCBH26-AKZ-LU-C86	HRCBH26-AKZ-MC-C86	HRCBH26-AKZ-WH-C86	HRCBH26-AKZ-LI-C86
Color	Laguna Blue	Magic Blue	White	Lilac
	●	●	○	●

PACKAGING

Size: 33x34 cm
Pieces/carton: 4x75 pcs
Carton/size: 37,5x19x22 cm
Carton weight: 2.50 kg
Carton/pallet: 108
Carton/layer: 12
Layers/pallet: 9
Pallet size: 80x120x169 cm
Pallet weight: 231kg

STORAGE

In a dry place at ambient temperature, away from heat sources, flames or sparks