

**LARIDENT SRL**  
VIA LAMANEIGRA, 12R  
16030 TRIBOGNA (GE) ITALY  
TEL.00390185938068 FAX. 00390185938595 info@larident.it

# TECHNICAL SHEET

OF PRODUCT

<b>PRODUCT NAME</b>	ENDO ORGANIZER SET
<b>CODE</b>	Q 65
<b>MATERIAL</b>	POLYAMIDE 6 WITH GLASS FIBRE
<b>THERMAL PROPERTIES</b>	MATERIAL AUTOCLAVABLE UP TO 134°C.
<b>MECHANICAL, PHYSICAL PROPERTIES OR CHARACTERISTICS</b>	STERY HOLDER FOR 16 ENDODONTIC INSTRUMENTS WITH ENDOMETER MEASURER AND SPONGES FOR CLEANING AVAILABLE IN : BLUE YELLOW AND GREEN  Q65 WITH 25 SPONGES : CODE Q67
<b>CHEMICAL PROPERTIES</b>	YOU CAN CLEAN IT WITH ALL THE DISINFECTANT PRODUCTS ON THE MARKET
<b>INSTRUCTIONS FOR USE</b>	IT IS USED BY HIGHLY TRAINED OR MEDICAL SPECIALIST IT IS SUPPLIED IN NON STERIL CONDITION , SO SHOULD BE METICULOSLY CLEANED AND STERILIZED BEFORE USING
<b>CLASSIFICATION</b> 93/42 EEC	NO MEDICAL DEVICE

TRIBOGNA , 20/04/2016