



Created by: **jth** 09.09.2015
Checked by: **wah** 10.09.2015

Type description:

S-11

REF

30056

Category:

Hand- und Winkelstücke

Manufacturer:

W&H

		Instructions for use	Leaflet
Valid from SN:	01001		
Intended use:	D001465	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Transmission ratio:	1:1	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Colour coding:	blue	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Maximum motor speed:	50.000 rpm	<input type="checkbox"/>	<input type="checkbox"/>
Maximum bur speed:	50.000 rpm	<input type="checkbox"/>	<input type="checkbox"/>
Dimension:			
• Head diameter:	10,6 mm	<input type="checkbox"/>	<input type="checkbox"/>
• Addendum without bur (19mm):	mm	<input type="checkbox"/>	<input type="checkbox"/>
• Addendum with bur (19mm):	mm	<input type="checkbox"/>	<input type="checkbox"/>
• Length	110,9 mm	<input type="checkbox"/>	<input type="checkbox"/>
Weight:	95,6 g	<input type="checkbox"/>	<input type="checkbox"/>
Noise level:	< 70 dBA	<input type="checkbox"/>	<input type="checkbox"/>
Connection motor - instrument according to ISO 3964:		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Chucking system:			
• Type _____	lever chuck	<input checked="" type="checkbox"/>	<input type="checkbox"/>
• Bur diameter accord. to ISO 1797-1/2:	2,35 mm	<input checked="" type="checkbox"/>	<input type="checkbox"/>
• Maximum permitted bur length:	45 mm	<input checked="" type="checkbox"/>	<input type="checkbox"/>
• Maximum permitted bur diameter (drill head):	mm	<input type="checkbox"/>	<input type="checkbox"/>
• operating- or insert force at chuck system:	N	<input type="checkbox"/>	<input type="checkbox"/>
• Maximum permitted torque (bur shaft, soft)	Ncm	<input type="checkbox"/>	<input type="checkbox"/>
• Maximum permitted torque (bur shaft, hardness > 50 HRC, >520 HV)	Ncm	<input type="checkbox"/>	<input type="checkbox"/>
• Retention force	N	<input type="checkbox"/>	<input type="checkbox"/>
Light system:	no light	<input type="checkbox"/>	<input type="checkbox"/>
Maximum light intensity:	lx	<input type="checkbox"/>	<input type="checkbox"/>
Cooling system:			
• Type _____	external	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
• Minimum amount of spray water:	50 ml/min	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Document type

Product Data Sheet

Number

00188-AEN

Index

00

PEOPLE HAVE PRIORITY



Created by:

jth 09.09.2015

Checked by:

wah 10.09.2015

Type description:

S-11

REF

30056

Category:

Hand- und Winkelstücke

Manufacturer:

W&H

		Instructions for use	Leaflet
• Type of ball bearings (steel/ceramic):	steel	<input type="checkbox"/>	<input type="checkbox"/>
Permitted hygiene methods:		<input checked="" type="checkbox"/>	<input type="checkbox"/>
• Sterilizer:	<input checked="" type="checkbox"/>		
• Thermo washer disinfectant:	<input checked="" type="checkbox"/>		
• Chemical autoclave:	<input type="checkbox"/>		
• Ultra sonic bath:	<input type="checkbox"/>		
• Disinfection solutions which are released for the product:	<input checked="" type="checkbox"/>		
Permitted lubrication methods:		<input type="checkbox"/>	<input type="checkbox"/>
• W&H serviceoil	<input checked="" type="checkbox"/>		
• None	<input type="checkbox"/>		
Warranty:	12 Monate	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Operation info:		<input type="checkbox"/>	<input type="checkbox"/>