

Meisinger Qualitäts- Management	Certificate for application and function	
Hersteller Manufacturer	HAGER & MEISINGER GMBH Hansemannstraße 10 41468 Neuss	

Fabrikat Article	Diamond burs (article no. 5100801314014, 5100801314018)
---------------------	--

TECHNISCHE DATEN / TECHNICAL DATA

Form Form	acc. to ISO 6360-2
Kopf Ø Head Ø	round
Schaft Ø Shank Ø	1.60 mm (FG)
Schaftausführung Design of shank	acc. to ISO 1797-1
Toleranzen Tolerance	acc. to ISO 7711-1

Material Material	Stainless Steel 1.4034, Natural Diamond
Härte Hardness	10.060 HV2
Anwendung Application	Dental and surgical application

Durchgeführte Prüfungen Performed tests	Checks regarding dimensional accuracy, visual and functional parameters at a) first article inspection for every new lot in serial production b) in-process control c) final inspection
--	--

Weitere angewandte Normen Other applied standards	EN 1639, EN 1641, ISO 14971
--	-----------------------------

Neuss, 13.04.2017
Ort, Datum / Place, Date


i.A. Dr. Melanie May

HAGER & MEISINGER GMBH
HANSEMANNSTR. 10
41468 NEUSS
TELEFON (02131) 2 01 20 FAX 2 01 22 22