



SOFT COVER TISSUE-PE

DISPOSABLE HEADREST COVER TO PROTECT THE PATIENT



CONFIGURATION OVERVIEW

Type: Non-Sterile
Latex free

Features: The special features of the materials offer high absorbency and waterproofing, so to assure a high comfort to the patient.

TECHNICAL SPECIFICATIONS

1 ply of pure wadding of long fiber cellulose (TISSUE) partially laminated with a second ply of high density polyethylene. (21 gsm Tissue with 17 gsm PE laminated)

Dimensions: 25x33 cm (tolerance ± 1%) / 25x25 (tolerance ± 1%)

Article reference number: - HRCAH26-AKZ-...-C28 (25x33 cm)
-HRCAH26-AKZ-...-C87 (25X25 cm)

PRODUCT VARIATIONS

Art-No.	HRCAH26-AKZ-LU	HRCAH26-AKZ-MC	HRCAH26-AKZ-WH	HRCAH26-AKZ-LI	HRCAH26-AKZ-PI	HRCAH26-AKZ-FU	HRCAH26-AKZ-FG	HRCAH26-AKZ-EG	HRCAH26-AKZ-IB
Color	Laguna Blue	Magic Blue	White	Lilac	Pink panther	Flamingo fuchsia	Fresh green	Emerald green	Ice Blue

PACKAGING

Size: 25x33 cm / 25x25 cm
Pieces/carton: 500
Carton/size: 38x29x16 cm
Carton weight: 3.90 kg
Carton/pallet: 88
Carton/layer: 8
Layers/pallet: 11
Pallet size: 80x120x191 cm
Pallet weight: 363 kg

STORAGE

In a dry place at ambient temperature, away from heat sources, flames, or sparks.