

## Space requirements




Fig. 1: View from the front, the right and above

Dimensions		Vacuklav 23 B+	Vacuklav 31 B+
Width	A	42.5 cm	
Height	B	49 cm	
Depth, total	C	74 cm	62 cm
Clearance between the device feet	C <sub>1</sub>	45 cm	37 cm
Clearance from rear device foot up to the rear panel	C <sub>2</sub>	7 cm	8 cm
Min. clearance to the side	D <sub>1</sub>	5 cm	
Min. clearance to the side of the door hinge	D <sub>2</sub>	10 cm	
Min. clearance to the rear	E	10 cm	
Tank lid (width)	F	23 cm	
Tank lid (depth)	F <sub>1</sub>	19.5 cm	
Tank lid (clearance to the rear)	F <sub>2</sub>	23 cm	
Free area with a fully-opened door	G	66 cm	54 cm
Min. clearance to the top (can be pulled out / with exhaust shaft)	--	5/20 cm	

The area above the steam sterilizer should be freely accessible in order to enable easy filling of the storage tank and to ensure good ventilation.