


Space requirements



Dimensions		Vacuklav 40 B+	Vacuklav 44 B+
Width	A	46 cm	46 cm
Height	B	56.5 cm	56.5 cm
Depth, total ^{*)}	C	58 cm	69 cm
Clearance between the device feet	C1	24 cm	35 cm
Clearance from rear device foot up to the rear panel	C2	17.6 cm	17.6 cm
Min. clearance to the side	D1	5 cm	5 cm
Min. clearance to the side of the door hinge	D2	14 cm	14 cm
Min. clearance to the rear	E	5 cm	5 cm
Free area with a fully-opened door	F	52 cm	63 cm
Clearance above	G	5 cm	5 cm

^{*)} Includes a safety combination in accordance with EN 1717

Maintain the side clearances when installing the steam sterilizer.