

Hygiene sleeves

Hygiene sleeves are a cover used in dental industry to avoid cross contamination that are soft, comfortable and tasteless.

Ref	Size	Measure	weight/pcs.	Box of pcs.	Box measure
70302	51x21 mm	51x21x0,065 mm	0,3g	500	72x61x60 mm
70303	80x40 mm	80x40x0,065 mm	0,3g	500	98x66x65 mm
70305	80x26 mm	80x26x0,09 mm	0,2g	500	96x61x59 mm
70306	110x55 mm	110x55x0,065 mm	0,7g	500	125x75x75 mm
70307	210(140)x55 mm	210(140)x55x0,13 mm	1,1g	500	240x75x75 mm
70308	40x21 mm	40x21x0,065 mm	0,3g	200	65x40x40 mm
70325	56x243 mm	56x243x0,05 mm	0,5g	200	255x65x25 mm

Contents

Hygienic sleeves
200 or 500 pieces
Packed in roll or individually
1 sticker

Features

Hygiene Sleeves are made out of latex-free, odorless plastics. The material has been tested not to provoke any skin irritation. The sleeves serve as a support for dental hygiene measures but do not replace the manufacturer's instructions for the disinfection or sterilization. The Hygiene Sleeves are available in different sizes and as such compatible with a variety of bite blocks, lights and other tools in the market.

Material Specification

Product: LDPE-foil / different thicknesses
transparent
Package: carton box

Information

Hygiene Sleeves are single-use, transparent plastic protection sleeves to prevent cross contamination between patients and to protect the tool from dust, dirt or physical damage during digital intraoral dental x ray procedures. Protective Hygienic sleeves are single use only and should never be re-used. Dispose protective hygienic sleeves after each patient use. Never use a damaged protective hygienic sleeve.

Medical device

Class 1 medical device. Classification rule 5. The product meets the applicable provisions of the DIRECTIVE 93/42/EEC amended by the DIRECTIVE 07/47/EEC concerning medical devices.

For more information: see operation instructions
phone 0049-415181011